

Ab✓Chek[®] NuLife



HIGH-QUALITY ELISA KITS

Discover Precision in
Immunoassay



Advanced **ELISA** solutions
engineered for accuracy, reliability
and innovation.

EXPLORE THE
ENTIRE
CATALOGUE

AbChek offers cutting-edge technologies designed to support laboratories, researchers and healthcare professionals with confidence reliability and consistency.



Complete
Species Coverage



Affordable



Stringent Quality
Inspection



Expert Technical
Assistance

CONTENTS

01-22 ABCHEK HUMAN ELISA KIT



23-34 ABCHEK MOUSE ELISA KIT



35-45 ABCHEK RAT ELISA KIT



46 ABCHEK UNIVERSAL ELISA KIT



47-48 ABCHEK MF SERIES ELISA KIT



SR. NO.	AbChek Cat No.	Product Name
1	abk1001H	Human ANG2(Angiopoietin 2) ELISA Kit
2	abk2001H	Human BAFF/CD257(B-Cell Activating Factor) ELISA Kit
3	abk3001H	Human BMP-2(Bone Morphogenetic Protein 2) ELISA Kit
4	abk4001H	Human BMP-4(Bone Morphogenetic Protein 4) ELISA Kit
5	abk5001H	Human BMP-7(Bone Morphogenetic Protein 7) ELISA Kit
6	abk6001H	Human E-Cad(E-Cadherin) ELISA Kit
7	abk7001H	Human CA9(Carbonic Anhydrase IX) ELISA Kit
8	abk8001H	Human CTSV(Cathepsin V) ELISA Kit
9	abk9001H	Human MIP-1 α (Macrophage Inflammatory Protein 1 Alpha) ELISA Kit
10	abk0101H	Human MIP-1 β (Macrophage Inflammatory Protein 1 Beta) ELISA Kit
11	abk1101H	Human MCP-3(Monocyte Chemotactic Protein 3) ELISA Kit
12	abk2101H	Human ECF/CCL11(Eosinophil Chemotactic Factor) ELISA Kit
13	abk3101H	Human TARC(Thymus Activation Regulated Chemokine) ELISA Kit
14	abk4101H	Human SLC(Secondary Lymphoid Tissue Chemokine) ELISA Kit
15	abk5101H	Human MPIF2(Myeloid Progenitor Inhibitory Factor 2) ELISA Kit
16	abk6101H	Human MIP4 α (Macrophage Inflammatory Protein 4 Alpha) ELISA Kit
17	abk7101H	Human CTACK(Cutaneous T-cell Attracting Chemokine) ELISA Kit
18	abk8101H	Human Fc ϵ R II /CD23(Receptor II for the Fc Region of Immunoglobulin E) ELISA Kit
19	abk9101H	Human sCD40L(Soluble Cluster of Differentiation 40 Ligand) ELISA Kit
20	abk0201H	Human sCD163(Soluble Cluster of Differentiation 163) ELISA Kit
21	abk1201H	Human GP39(Glycoprotein 39, Cartilage) ELISA Kit
22	abk2201H	Human CNTF(Ciliary Neurotrophic Factor) ELISA Kit
23	abk3201H	Human F3/TF(Tissue Factor) ELISA Kit
24	abk4201H	Human CRN(Corin) ELISA Kit
25	abk5201H	Human CX3CL1(Chemokine C-X3-C-Motif Ligand 1) ELISA Kit
26	abk7201H	Human GCP-2(Granulocyte Chemotactic Protein 2) ELISA Kit
27	abk0301H	Human BLC(B-Lymphocyte Chemoattractant) ELISA Kit
28	abk1301H	Human Cyt-C(Cytochrome C) ELISA Kit
29	abk2301H	Human DKK1(Dickkopf Relat Protein 1) ELISA Kit
30	abk3301H	Human EGF(Epidermal Growth Factor) ELISA Kit
31	abk4301H	Human EGFR(Epidermal Growth Factor Receptor) ELISA Kit
32	abk5301H	Human EG-VEGF(Endocrine Gland Derived Vascular Endothelial Growth Factor) ELISA Kit
33	abk6301H	Human ET-1(Endothelin 1) ELISA Kit
34	abk7301H	Human EPCR(Endothelial Protein C receptor) ELISA Kit
35	abk8301H	Human FAS/CD95(Factor Related Apoptosis) ELISA Kit
36	abk9301H	Human FASL/TNFSF6(Factor Related Apoptosis Ligand) ELISA Kit
37	abk0401H	Human AFGF/FGF1(Acidic Fibroblast Growth Factor 1) ELISA Kit
38	abk1401H	Human FGF19(Fibroblast Growth Factor 19) ELISA Kit
39	abk2401H	Human FGF21(Fibroblast Growth Factor 21) ELISA Kit
40	abk3401H	Human FST(Follistatin) ELISA Kit
41	abk4401H	Human GDF15(Growth Differentiation Factor 15) ELISA Kit
42	abk5401H	Human HGF(Hepatocyte Growth Factor) ELISA Kit
43	abk6401H	Human IFN- β (Interferon Beta) ELISA Kit
44	abk7401H	Human IL-1 α (Interleukin 1 Alpha) ELISA Kit
45	abk8401H	Human IL-1ra/IL-1F3(Interleukin 1 Receptor Antagonist) ELISA Kit
46	abk9401H	Human IL-1R2(Interleukin 1 Receptor Type II) ELISA Kit
47	abk0501H	Human LIF(Leukemia Inhibitory Factor) ELISA Kit
48	abk1501H	Human TNFSF14(Tumor Necrosis Factor Ligand Superfamily, Member 14) ELISA Kit
49	abk2501H	Human MCSF(Macrophage Colony Stimulating Factor 1) ELISA Kit
50	abk3501H	Human Pro-MMP-1(pro-Matrix Metalloproteinase 1) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
51	abk4501H	Human IL-3(Interleukin 3) ELISA Kit
52	abk5501H	Human IL-4(Interleukin 4) ELISA Kit
53	abk6501H	Human IL-13(Interleukin 13) ELISA Kit
54	abk7501H	Human IL-22(Interleukin 22) ELISA Kit
55	abk8501H	Human IL-23(Interleukin 23) ELISA Kit
56	abk1601H	Human VEGF-A(Vascular Endothelial Cell Growth Factor A) ELISA Kit
57	abk2601H	Human EREG(Epiregulin) ELISA Kit
58	abk3601H	Human Bcl-2(B-cell Leukemia/Lymphoma 2) ELISA Kit
59	abk4601H	Human REN(Renin) ELISA Kit
60	abk5601H	Human CTSC(Cathepsin C) ELISA Kit
61	abk7601H	Human MMP-25(Matrix Metalloproteinase 25) ELISA Kit
62	abk8601H	Human MMP-12(Matrix Metalloproteinase 12) ELISA Kit
63	abk9601H	Human CT(Calcitonin) ELISA Kit
64	abk2701H	Human PAPPA(Pregnancy Associated Plasma Protein A) ELISA Kit
65	abk3701H	Human CSTA(Cystatin A) ELISA Kit
66	abk5701H	Human IL-12(Interleukin 12) ELISA Kit
67	abk6701H	Human IL-12 p40(Interleukin 12 p40) ELISA Kit
68	abk8701H	Human CD42(Cluster of Differentiation 42) ELISA Kit
69	abk9701H	Human HLA-B27(Human Leukocyte Antigen B27) ELISA Kit
70	abk0801H	Human FGFR2(Fibroblast Growth Factor Receptor 2) ELISA Kit
71	abk1801H	Human MBP(Myelin Basic Protein) ELISA Kit
72	abk2801H	Human TE(Telomerase) ELISA Kit
73	abk3801H	Human TFPI(Tissue Factor Pathway Inhibitor) ELISA Kit
74	abk6353H	Human TM(Thrombomodulin) ELISA Kit
75	abk5801H	Human FSTL1(Follistatin Like Protein 1) ELISA Kit
76	abk6801H	Human LOX(Lysyl Oxidase) ELISA Kit
77	abk7801H	Human GH(Growth Hormone) ELISA Kit
78	abk8801H	Human COL4(Collagen Type IV) ELISA Kit
79	abk9801H	Human FN(Fibronectin) ELISA Kit
80	abk0901H	Human IL-9(Interleukin 9) ELISA Kit
81	abk1901H	Human PC I (Procollagen I) ELISA Kit
82	abk2901H	Human PCIII(Procollagen III) ELISA Kit
83	abk3901H	Human PIINP(N-Terminal Procollagen III Propeptide) ELISA Kit
84	abk5901H	Human ADCY2(Adenylate Cyclase 2) ELISA Kit
85	abk6901H	Human AGGF1(Angiogenic Factor with G Patch and FHA Domains 1) ELISA Kit
86	abk7901H	Human C4a(Complement Component 4a) ELISA Kit
87	abk8901H	Human C5a(Complement Component 5a) ELISA Kit
88	abk9901H	Human IL-5(Interleukin 5) ELISA Kit
89	abk0011H	Human IL-6R(Interleukin 6 Receptor) ELISA Kit
90	abk1011H	Human Fc γ R I (Receptor I for the Fc Region of Immunoglobulin G) ELISA Kit
91	abk2011H	Human NCAD(Neural Cadherin) ELISA Kit
92	abk3011H	Human PDGFRL(Platelet Derived Growth Factor Receptor Like Protein) ELISA Kit
93	abk4011H	Human LUM(Lumican) ELISA Kit
94	abk5011H	Human PRR4(Proline Rich Protein 4, Lacrimal) ELISA Kit
95	abk6011H	Human PSAP(Prosaposin) ELISA Kit
96	abk7011H	Human PTGES(Prostaglandin E Synthase, Microsomal) ELISA Kit
97	abk8011H	Human ISG15(Interferon Stimulated Gene 15) ELISA Kit
98	abk9011H	Human AIMP1(Aminoacyl tRNA Synthetase Complex Interacting Multifunctional Protein 1) ELISA Kit
99	abk1111H	Human 15-LO(Arachidonate 15-Lipoxygenase) ELISA Kit
100	abk2111H	Human BACE1(Beta-Site APP Cleaving Enzyme 1) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
101	abk3111H	Human IL-32(Interleukin 32) ELISA Kit
102	abk4111H	Human TNFRSF1A(Tumor Necrosis Factor Receptor Superfamily, Member 1A) ELISA Kit
103	abk5111H	Human IL-15(Interleukin 15) ELISA Kit
104	abk6111H	Human 5-LO(Arachidonate 5-Lipoxygenase) ELISA Kit
105	abk7111H	Human 5-NT(5-Nucleotidase) ELISA Kit
106	abk8111H	Human α 1ACT(Alpha 1 Antichymotrypsin) ELISA Kit
107	abk9111H	Human IL-16(Interleukin 16) ELISA Kit
108	abk6721H	Human AR(Amphiregulin) ELISA Kit
109	abk1211H	Human ARO(Aromatase) ELISA Kit
110	abk2211H	Human BLVRB(Biliverdin Reductase B) ELISA Kit
111	abk3211H	Human IL-18(Interleukin 18) ELISA Kit
112	abk4211H	Human IL-19(Interleukin 19) ELISA Kit
113	abk5211H	Human ACTG1(Actin Gamma 1) ELISA Kit
114	abk6211H	Human ACTN2(Actinin Alpha 2) ELISA Kit
115	abk7211H	Human ACTN3(Actinin Alpha 3) ELISA Kit
116	abk8211H	Human ADA(Adenosine Deaminase) ELISA Kit
117	abk9211H	Human ADAM10(A Disintegrin And Metalloprotease 10) ELISA Kit
118	abk0311H	Human ADAM8(A Disintegrin And Metalloprotease 8) ELISA Kit
119	abk1311H	Human ADAM9(A Disintegrin And Metalloprotease 9) ELISA Kit
120	abk2311H	Human ADAMTS4(ADAM with Thrombospondin Type 1 Motif 4) ELISA Kit
121	abk3311H	Human ADCY1(Adenylate Cyclase 1, Brain) ELISA Kit
122	abk4311H	Human ADH1(Alcohol Dehydrogenase 1) ELISA Kit
123	abk6311H	Human ADIPOR1(Adiponectin Receptor 1) ELISA Kit
124	abk7311H	Human ADIPOR2(Adiponectin Receptor 2) ELISA Kit
125	abk8311H	Human ADM(Adrenomedullin) ELISA Kit
126	abk9311H	Human ADRA1A(Adrenergic Receptor alpha-1A) ELISA Kit
127	abk0411H	Human ADRP(Adipose Differentiation Related Protein) ELISA Kit
128	abk1411H	Human IL-20(Interleukin 20) ELISA Kit
129	abk2411H	Human DLK-1(Delta Like 1 Homolog) ELISA Kit
130	abk3411H	Human ACE2(Angiotensin I Converting Enzyme 2) ELISA Kit
131	abk4411H	Human BMP-3(Bone Morphogenetic Protein 3) ELISA Kit
132	abk6411H	Human CAPS(Calcyphosine) ELISA Kit
133	abk7411H	Human FABP4(Fatty Acid Binding Protein 4, Adipocyte) ELISA Kit
134	abk8411H	Human AFP-L3(Alpha-Fetoprotein Lens Culinaris Agglutinin 3) ELISA Kit
135	abk9411H	Human FUCA(Alpha-L-Fucosidase, Tissue) ELISA Kit
136	abk0511H	Human CCNB1(Cyclin-B1) ELISA Kit
137	abk1511H	Human AGC(Aggregan) ELISA Kit
138	abk2511H	Human AGER(Total Advanced Glycosylation End Product Specific Receptor) ELISA Kit
139	abk3511H	Human CD3d(Cluster of Differentiation 3d) ELISA Kit
140	abk4511H	Human AGR2(Anterior Gradient Protein 2) ELISA Kit
141	abk5511H	Human AGRN(Agrin) ELISA Kit
142	abk6511H	Human AGT(Angiotensinogen) ELISA Kit
143	abk7511H	Human AIF1(Allograft Inflammatory Factor 1) ELISA Kit
144	abk8511H	Human ALDM(Aldehyde Dehydrogenase, Mitochondrial) ELISA Kit
145	abk9511H	Human Alk-Smase(Alkaline Sphingomyelinase) ELISA Kit
146	abk0611H	Human INHB(Inhibin B) ELISA Kit
147	abk1611H	Human CDK6(Cyclin Dependent Kinase 6) ELISA Kit
148	abk2611H	Human AMBP(Alpha-1-Microglobulin/Bikunin Precursor) ELISA Kit
149	abk3611H	Human CDK8(Cyclin Dependent Kinase 8) ELISA Kit
150	abk4611H	Human AMY1(Amylase Alpha 1, Salivary) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
151	abk6611H	Human CDKN3(Cyclin Dependent Kinase Inhibitor 3) ELISA Kit
152	abk7353H	Human ANG(Angiostatin) ELISA Kit
153	abk8611H	Human ANG II R1(Angiotensin II Receptor 1) ELISA Kit
154	abk0711H	Human ANG II R2(Angiotensin II Receptor 2) ELISA Kit
155	abk1711H	Human ANGPT4(Angiopietin 4) ELISA Kit
156	abk2711H	Human ANGPTL1(Angiopietin Like Protein 1) ELISA Kit
157	abk3711H	Human ANGPTL4(Angiopietin Like Protein 4) ELISA Kit
158	abk4711H	Human ANGPTL6(Angiopietin Like Protein 6) ELISA Kit
159	abk5711H	Human ANG-R-Tie1(Angiopietin Receptor Tie1) ELISA Kit
160	abk6711H	Human ANG-R-Tie2(Angiopietin Receptor Tie2) ELISA Kit
161	abk7711H	Human Ank-B(Ankyrin B) ELISA Kit
162	abk8711H	Human MAG(Myelin Associated Glycoprotein) ELISA Kit
163	abk9711H	Human CDNF(Cerebral Dopamine Neurotrophic Factor) ELISA Kit
164	abk0811H	Human FGFR1(Fibroblast Growth Factor Receptor 1) ELISA Kit
165	abk1811H	Human CNN1(Calponin 1, Basic, Smooth Muscle) ELISA Kit
166	abk2811H	Human CNTNAP2(Contactin Associated Protein Like Protein 2) ELISA Kit
167	abk3811H	Human ACHAP(Acetylcholinesterase Associated Protein) ELISA Kit
168	abk4811H	Human ACO1(Aconitase 1) ELISA Kit
169	abk5811H	Human CPE(Carboxypeptidase E) ELISA Kit
170	abk7811H	Human CG α (Chorionic Gonadotropin Alpha Polypeptide) ELISA Kit
171	abk8811H	Human FGFR3(Fibroblast Growth Factor Receptor 3) ELISA Kit
172	abk9811H	Human DPP10(Dipeptidyl Peptidase X) ELISA Kit
173	abk0911H	Human CPR(Cytochrome P450 Reductase) ELISA Kit
174	abk1911H	Human CRLF1(Cytokine Receptor Like Factor 1) ELISA Kit
175	abk2911H	Human CSF2R β (Colony Stimulating Factor 2 Receptor Beta) ELISA Kit
176	abk3911H	Human CA3(Carbonic Anhydrase III, Muscle Specific) ELISA Kit
177	abk5911H	Human CYP1A2(Cytochrome P450, family 1, subfamily A, polypeptide 2) ELISA Kit
178	abk7911H	Human DCLK1(Doublecortin Like Kinase 1) ELISA Kit
179	abk8911H	Human DPYD(Dihydropyrimidine Dehydrogenase) ELISA Kit
180	abk9911H	Human DSE(Dermatan Sulfate Epimerase) ELISA Kit
181	abk0021H	Human EPYC(Epiphycan) ELISA Kit
182	abk1021H	Human FADS2(Fatty Acid Desaturase 2) ELISA Kit
183	abk2021H	Human FAIM3(Fas Apoptotic Inhibitory Molecule 3) ELISA Kit
184	abk3021H	Human FASN(Fatty Acid Synthase) ELISA Kit
185	abk4021H	Human FBLN4(Fibulin 4) ELISA Kit
186	abk5021H	Human CYP24A1(Cytochrome P450, family 24, subfamily A, polypeptide 1) ELISA Kit
187	abk6021H	Human FCN1(Ficolin 1) ELISA Kit
188	abk8021H	Human Fc γ BP(Fc Fragment of IgG Binding Protein) ELISA Kit
189	abk9021H	Human Fc γ R3A(Fc Fragment of IgG Low Affinity IIIa Receptor) ELISA Kit
190	abk0121H	Human CYP1B1(Cytochrome P450, family 1, subfamily B, polypeptide 1) ELISA Kit
191	abk1121H	Human FETUA(Fetuin A) ELISA Kit
192	abk2121H	Human FG γ (Fibrinogen Gamma chain) ELISA Kit
193	abk3121H	Human FLOT2(Flotillin 2) ELISA Kit
194	abk9032H	Human GC(Glycocalicin) ELISA Kit
195	abk5121H	Human GLRX3(Glutaredoxin 3) ELISA Kit
196	abk6121H	Human GIF(Gastric Intrinsic Factor) ELISA Kit
197	abk7121H	Human C4BP β (C4 Binding Protein Beta) ELISA Kit
198	abk8121H	Human CA4(Carbonic Anhydrase IV) ELISA Kit
199	abk9121H	Human C4BP α (C4 Binding Protein Alpha) ELISA Kit
200	abk0221H	Human GPLD1(Glycosylphosphatidylinositol Specific Phospholipase D1) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
201	abk1221H	Human PRELP(Proline Arginine Rich End Leucine Rich Repeat Protein) ELISA Kit
202	abk3221H	Human TARS(Threonyl tRNA Synthetase, cytoplasmic) ELISA Kit
203	abk4221H	Human GTF2A1(General Transcription Factor II A, Polypeptide 1) ELISA Kit
204	abk5221H	Human GUCA2A(Guanylate Cyclase Activator 2A) ELISA Kit
205	abk6221H	Human INSL3(Insulin Like Protein 3) ELISA Kit
206	abk7221H	Human INSL5(Insulin Like Protein 5) ELISA Kit
207	abk8221H	Human GUCA2B(Guanylate Cyclase Activator 2B) ELISA Kit
208	abk9221H	Human GZMK(Granzyme K) ELISA Kit
209	abk0321H	Human ALOXE3(Arachidonate Lipoxygenase 3) ELISA Kit
210	abk1321H	Human CYP26A1(Cytochrome P450, family 26, subfamily A, polypeptide 1) ELISA Kit
211	abk2321H	Human HCCS(Holocytochrome C Synthase) ELISA Kit
212	abk3321H	Human CAMK2α(Calcium/Calmodulin-Dependent Protein Kinase II Alpha) ELISA Kit
213	abk4321H	Human BAG3(Bcl2 Associated Athanogene 3) ELISA Kit
214	abk5321H	Human ARSB(Arylsulfatase B) ELISA Kit
215	abk6321H	Human ANXA5(Annexin A5) ELISA Kit
216	abk7321H	Human HTRA1(HtrA Serine Peptidase 1) ELISA Kit
217	abk8321H	Human IGF1R(Insulin Like Growth Factor 1 Receptor) ELISA Kit
218	abk9321H	Human HYOU1(Hypoxia Up Regulated 1) ELISA Kit
219	abk0421H	Human ICAM4(Intercellular Adhesion Molecule 4) ELISA Kit
220	abk1421H	Human IDE(Insulin Degrading Enzyme) ELISA Kit
221	abk3421H	Human LIPH(Lipase H) ELISA Kit
222	abk4421H	Human LTβR(Lymphotoxin Beta Receptor) ELISA Kit
223	abk6421H	Human IGF2BP3(Insulin Like Growth Factor 2 mRNA Binding Protein 3) ELISA Kit
224	abk7421H	Human TTN(Titin) ELISA Kit
225	abk8421H	Human IGF2R(Insulin Like Growth Factor 2 Receptor) ELISA Kit
226	abk9421H	Human IGFALS(Insulin Like Growth Factor Binding Protein, Acid Labile Subunit) ELISA Kit
227	abk0521H	Human IGFBP-7(Insulin Like Growth Factor Binding Protein 7) ELISA Kit
228	abk1521H	Human ANXA2(Annexin A2) ELISA Kit
229	abk2521H	Human INSR(Insulin Receptor) ELISA Kit
230	abk5521H	Human LCN1(Lipocalin 1) ELISA Kit
231	abk6521H	Human ApoA4(Apolipoprotein A4) ELISA Kit
232	abk7521H	Human ApoB(Apolipoprotein B) ELISA Kit
233	abk8521H	Human ApoC1(Apolipoprotein C1) ELISA Kit
234	abk9521H	Human ApoC2(Apolipoprotein C2) ELISA Kit
235	abk0621H	Human ApoC4(Apolipoprotein C4) ELISA Kit
236	abk1621H	Human ApoD(Apolipoprotein D) ELISA Kit
237	abk2621H	Human ApoE(Apolipoprotein E) ELISA Kit
238	abk3621H	Human ApoH(Apolipoprotein H) ELISA Kit
239	abk4621H	Human ApoM(Apolipoprotein M) ELISA Kit
240	abk5621H	Human CA6(Carbonic Anhydrase VI) ELISA Kit
241	abk6621H	Human AATF(Apoptosis Antagonizing Transcription Factor) ELISA Kit
242	abk7621H	Human LIPC(Lipase, Hepatic) ELISA Kit
243	abk8621H	Human AIF(Apoptosis Inducing Factor) ELISA Kit
244	abk9621H	Human LIPF(Lipase, Gastric) ELISA Kit
245	abk0721H	Human APAF1(Apoptosis Protease Activating Factor 1) ELISA Kit
246	abk1721H	Human ASK1(Apoptosis Signal Regulating Kinase 1) ELISA Kit
247	abk2721H	Human AQP-1(Aquaporin 1) ELISA Kit
248	abk3721H	Human AQP-2(Aquaporin 2) ELISA Kit
249	abk4721H	Human AQP-3(Aquaporin 3) ELISA Kit
250	abk5721H	Human AQP-4(Aquaporin 4) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
251	abk0211H	Human AR(Androgen Receptor) ELISA Kit
252	abk7721H	Human AVPR1A(Arginine Vasopressin Receptor 1A) ELISA Kit
253	abk8721H	Human KLK6(Kallikrein 6) ELISA Kit
254	abk9721H	Human KLK7(Kallikrein 7) ELISA Kit
255	abk0821H	Human ARHGDlα(Rho GDP Dissociation Inhibitor Alpha) ELISA Kit
256	abk1821H	Human ARRβ1(Arrestin Beta 1) ELISA Kit
257	abk2821H	Human ARRβ2(Arrestin Beta 2) ELISA Kit
258	abk3821H	Human MANF(Mesencephalic Astrocyte Derived Neurotrophic Factor) ELISA Kit
259	abk4821H	Human MAP2K1(Mitogen Activated Protein Kinase Kinase 1) ELISA Kit
260	abk5821H	Human AhR(Aryl Hydrocarbon Receptor) ELISA Kit
261	abk6821H	Human ARSA(Arylsulfatase A) ELISA Kit
262	abk7821H	Human ASGR2(Asialoglycoprotein Receptor 2) ELISA Kit
263	abk8821H	Human ASGR1(Asialoglycoprotein Receptor 1) ELISA Kit
264	abk9821H	Human ASM(Acid Sphingomyelinase) ELISA Kit
265	abk0921H	Human MAPK14(Mitogen Activated Protein Kinase 14) ELISA Kit
266	abk1921H	Human AST(Aspartate Aminotransferase) ELISA Kit
267	abk2921H	Human ASPN(Asporin) ELISA Kit
268	abk3921H	Human ATA(Angiotensinase A) ELISA Kit
269	abk4921H	Human ATM(Ataxia Telangiectasia Mutated) ELISA Kit
270	abk5921H	Human ATXN1(Ataxin 1) ELISA Kit
271	abk6921H	Human ATF1(Activating Transcription Factor 1) ELISA Kit
272	abk7921H	Human ATF2(Activating Transcription Factor 2) ELISA Kit
273	abk8921H	Human ATF3(Activating Transcription Factor 3) ELISA Kit
274	abk9921H	Human ATF4(Activating Transcription Factor 4) ELISA Kit
275	abk0031H	Human ATF6(Activating Transcription Factor 6) ELISA Kit
276	abk2031H	Human MGL(Monoacylglycerol Lipase) ELISA Kit
277	abk3031H	Human AIRE(Autoimmune Regulator) ELISA Kit
278	abk4031H	Human AXIN2(Axis Inhibition Protein 2) ELISA Kit
279	abk6031H	Human MCRS1(Microspherule Protein 1) ELISA Kit
280	abk7031H	Human BAFFR(B-Cell Activation Factor Receptor) ELISA Kit
281	abk8031H	Human BCL2L1/BCL-X(Bcl-2 Like Protein 1) ELISA Kit
282	abk9031H	Human BPI(Bactericidal/Permeability Increasing Protein) ELISA Kit
283	abk0131H	Human FGF6(Fibroblast Growth Factor 6) ELISA Kit
284	abk1131H	Human LDHA(Lactate Dehydrogenase A) ELISA Kit
285	abk2131H	Human Bcl-10(B-cell Lymphoma/Leukemia 10) ELIS Kit ELISA Kit
286	abk3131H	Human BAK1(Bcl2 Antagonist/Killer 1) ELISA Kit
287	abk4131H	Human BAD(Bcl2 Associated Death Promoter) ELISA Kit
288	abk5131H	Human BAX(Bcl-2 Associated X Protein) ELISA Kit
289	abk6131H	Human BMF(Bcl-2 Modifying Factor) ELISA Kit
290	abk7131H	Human BECN1(Beclin 1) ELISA Kit
291	abk8131H	Human SEMA3B(Semaphorin 3B) ELISA Kit
292	abk0231H	Human BACE2(Beta-Site APP Cleaving Enzyme 2) ELISA Kit
293	abk1231H	Human BID(BH3 Interacting Domain Death Agonist) ELISA Kit
294	abk2231H	Human MMRN1(Multimerin 1) ELISA Kit
295	abk3231H	Human Bcl2L11(Bcl-2 Like Protein 11) ELISA Kit
296	abk4231H	Human BMP-10(Bone Morphogenetic Protein 10) ELISA Kit
297	abk5231H	Human BMPR2(Bone Morphogenetic Protein Receptor II) ELISA Kit
298	abk6231H	Human MOG(Myelin Oligodendrocyte Glycoprotein) ELISA Kit
299	abk7231H	Human MT1E(Metallothionein 1E) ELISA Kit
300	abk9231H	Human MSMβ(Microseminoprotein Beta) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
301	abk0331H	Human BRAK/CXCL14(Breast and Kidney Expressed Chemokine) ELISA Kit
302	abk1331H	Human BCAN(Brevican) ELISA Kit
303	abk2331H	Human MT1M(Metallothionein 1M) ELISA Kit
304	abk3331H	Human MT3(Metallothionein 3) ELISA Kit
305	abk6331H	Human MUC13(Mucin 13) ELISA Kit
306	abk7331H	Human CA724(Tumor Marker CA724) ELISA Kit
307	abk8331H	Human CNR-1(Cadherin Related Neuronal Receptor 1) ELISA Kit
308	abk9331H	Human MUC1(Mucin 1) ELISA Kit
309	abk0431H	Human CABIN1(Calcineurin Binding Protein 1) ELISA Kit
310	abk1431H	Human CTR(Calcitonin Receptor) ELISA Kit
311	abk2431H	Human CACYBP(Calcyclin Binding Protein) ELISA Kit
312	abk3431H	Human CALD(Caldesmon) ELISA Kit
313	abk4431H	Human MSR1(Macrophage Scavenger Receptor 1) ELISA Kit
314	abk5431H	Human CAM(Calmodulin) ELISA Kit
315	abk6431H	Human CNX(Calnexin) ELISA Kit
316	abk7431H	Human CRT(Calreticulin) ELISA Kit
317	abk8431H	Human CR(Calretinin) ELISA Kit
318	abk9431H	Human CALU(Calumenin) ELISA Kit
319	abk1531H	Human MUC17(Mucin 17) ELISA Kit
320	abk2531H	Human MUC2(Mucin 2) ELISA Kit
321	abk3531H	Human MUC20(Mucin 20) ELISA Kit
322	abk4531H	Human CNR1(Cannabinoid Receptor 1, Brain) ELISA Kit
323	abk5531H	Human CA1(Carbonic Anhydrase I) ELISA Kit
324	abk6531H	Human CA2(Carbonic Anhydrase II) ELISA Kit
325	abk7531H	Human ANKRD1(Ankyrin Repeat Domain 1) ELISA Kit
326	abk0631H	Human CLS(Cardiopilin Synthase) ELISA Kit
327	abk1631H	Human CLCF1(Cardiotrophin Like Cytokine Factor 1) ELISA Kit
328	abk2631H	Human PAR1(Protease Activated Receptor 1) ELISA Kit
329	abk3631H	Human COMP(Cartilage Oligomeric Matrix Protein) ELISA Kit
330	abk4631H	Human CASP14(Caspase 14) ELISA Kit
331	abk5631H	Human CASP8(Caspase 8) ELISA Kit
332	abk6631H	Human CASP4(Caspase 4) ELISA Kit
333	abk7631H	Human CASP7(Caspase 7) ELISA Kit
334	abk8631H	Human CAD(Caspase Activated Deoxyribonuclease) ELISA Kit
335	abk9631H	Human CASP 9(Caspase 9) ELISA Kit
336	abk0731H	Human CTNNA1(Catenin, Alpha 1) ELISA Kit
337	abk1731H	Human CTNNB1(Catenin, Beta 1) ELISA Kit
338	abk2731H	Human CTSD(Cathepsin D) ELISA Kit
339	abk3731H	Human CTSK(Cathepsin K) ELISA Kit
340	abk4731H	Human CTSL(Cathepsin L) ELISA Kit
341	abk5731H	Human CDX2(Caudal Type Homeobox Transcription Factor 2) ELISA Kit
342	abk6731H	Human C/EBPα(CCAAT/Enhancer Binding Protein Alpha) ELISA Kit
343	abk7731H	Human C/EBPβ(CCAAT/Enhancer Binding Protein Beta) ELISA Kit
344	abk8731H	Human C/EBPδ(CCAAT/Enhancer Binding Protein, Delta) ELISA Kit
345	abk9731H	Human C/EBPε(CCAAT/Enhancer Binding Protein, Epsilon) ELISA Kit
346	abk0831H	Human C/EBPγ(CCAAT/Enhancer Binding Protein Gamma) ELISA Kit
347	abk1831H	Human CD59(Protectin) ELISA Kit
348	abk2831H	Human CENPA(Centromere Protein A) ELISA Kit
349	abk3831H	Human CA5B(Carbonic Anhydrase VB, Mitochondrial) ELISA Kit
350	abk4831H	Human CFL1(Cofilin 1, Non-Muscle) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
351	abk5831H	Human FOS(V-Fos FBJ Murine Osteosarcoma Viral Oncogene Homolog) ELISA Kit
352	abk6831H	Human MUC21(Mucin 21) ELISA Kit
353	abk7831H	Human CHEM(Chemerin) ELISA Kit
354	abk8831H	Human CCL14(Cheokine C-C-Motif Ligand 14) ELISA Kit
355	abk9831H	Human CCL16(Cheokine C-C-Motif Ligand 16) ELISA Kit
356	abk0931H	Human CCL3L1(Cheokine C-C-Motif Ligand 3 Like Protein 1) ELISA Kit
357	abk1931H	Human CCR8(Cheokine C-C-Motif Receptor 8) ELISA Kit
358	abk2931H	Human TERT(Telomerase Reverse Transcriptase) ELISA Kit
359	abk3931H	Human CKLFSF4(Cheokine Like Factor Superfamily 4) ELISA Kit
360	abk4931H	Human CKLFSF5(Cheokine Like Factor Superfamily 5) ELISA Kit
361	abk5931H	Human CKLFSF6(Cheokine Like Factor Superfamily 6) ELISA Kit
362	abk6931H	Human MUC4(Mucin 4) ELISA Kit
363	abk7931H	Human MUC6(Mucin 6) ELISA Kit
364	abk8931H	Human MYBPC3(Myosin Binding Protein C, Cardiac) ELISA Kit
365	abk9931H	Human MYH10(Myosin Heavy Chain 10, Non Muscle) ELISA Kit
366	abk0041H	Human MYH9(Myosin Heavy Chain 9, Non Muscle) ELISA Kit
367	abk1041H	Human MYL9(Myosin Light Chain 9, Regulatory) ELISA Kit
368	abk2041H	Human CCKAR(Cholecystokinin A Receptor) ELISA Kit
369	abk3041H	Human CETP(Cholesteryl Ester Transfer Protein) ELISA Kit
370	abk4041H	Human CENPF(Centromere Protein F) ELISA Kit
371	abk5041H	Human CHRM2(Cholinergic Receptor, Muscarinic 2) ELISA Kit
372	abk6041H	Human CA12(Carbonic Anhydrase XII) ELISA Kit
373	abk7041H	Human CBX3(Chromobox Homolog 3) ELISA Kit
374	abk8041H	Human CHD5(Chromodomain Helicase DNA Binding Protein 5) ELISA Kit
375	abk9041H	Human NEF3(Neurofilament 3) ELISA Kit
376	abk0141H	Human NEFL(Neurofilament, Light Polypeptide) ELISA Kit
377	abk1141H	Human NOS1/nNOS(Nitric Oxide Synthase 1, Neuronal) ELISA Kit
378	abk2141H	Human C-jun(Proto-oncogene c-Jun) ELISA Kit
379	abk3141H	Human CC16(Clara Cell 16kD protein) ELISA Kit
380	abk4141H	Human CLDN1(Claudin 1) ELISA Kit
381	abk5141H	Human CD30(Cluster of Differentiation 30) ELISA Kit
382	abk6141H	Human CD97(Cluster of Differentiation 97) ELISA Kit
383	abk7141H	Human CD99(Cluster of Differentiation 99) ELISA Kit
384	abk8141H	Human CD4(Cluster of Differentiation 4) ELISA Kit
385	abk9141H	Human NOS2/iNOS(Nitric Oxide Synthase 2, Inducible) ELISA Kit
386	abk0241H	Human CLDN3(Claudin 3) ELISA Kit
387	abk2241H	Human F7(Coagulation Factor VII) ELISA Kit
388	abk3241H	Human F9(Coagulation Factor IX) ELISA Kit
389	abk4241H	Human F10(Coagulation Factor X) ELISA Kit
390	abk5241H	Human F13(Coagulation Factor XIII) ELISA Kit
391	abk6241H	Human F12(Coagulation Factor XII) ELISA Kit
392	abk7241H	Human F5(Coagulation Factor V) ELISA Kit
393	abk8241H	Human F11(Coagulation Factor XI) ELISA Kit
394	abk9241H	Human F13A1(Coagulation Factor XIII A1 Polypeptide) ELISA Kit
395	abk0341H	Human F13B(Coagulation Factor XIII B Polypeptide) ELISA Kit
396	abk1341H	Human CFL2(Cofilin 2, Muscle) ELISA Kit
397	abk2341H	Human COL15(Collagen Type XV) ELISA Kit
398	abk3341H	Human COL6(Collagen Type VI) ELISA Kit
399	abk4341H	Human COL1α2(Collagen Type I Alpha 2) ELISA Kit
400	abk5341H	Human COL2α1(Collagen Type II Alpha 1) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
401	abk6341H	Human COL3α1(Collagen Type III Alpha 1) ELISA Kit
402	abk7341H	Human COL4α1(Collagen Type IV Alpha 1) ELISA Kit
403	abk8341H	Human COL6α1(Collagen Type VI Alpha 1) ELISA Kit
404	abk9341H	Human COL7(Collagen Type VII) ELISA Kit
405	abk0441H	Human COL8α1(Collagen Type VIII Alpha 1) ELISA Kit
406	abk1441H	Human COL11α1(Collagen Type XI Alpha 1) ELISA Kit
407	abk2441H	Human COL17(Collagen Type X VII) ELISA Kit
408	abk3441H	Human Collagenase I ELISA Kit
409	abk4441H	Human NPYR1(Neuropeptide Y Receptor Y1) ELISA Kit
410	abk5441H	Human CLL1(Collectin Liver 1) ELISA Kit
411	abk6441H	Human CCSA-4(Colon Cancer Specific Antigen-4) ELISA Kit
412	abk7441H	Human CLDN9(Claudin 9) ELISA Kit
413	abk8441H	Human NEP(Neprilysin) ELISA Kit
414	abk9441H	Human CLDN4(Claudin 4) ELISA Kit
415	abk0541H	Human C7(Complement Component 7) ELISA Kit
416	abk1541H	Human C9(Complement Component 9) ELISA Kit
417	abk2541H	Human CFP(Complement Factor P) ELISA Kit
418	abk3541H	Human C3a(Complement Component 3a) ELISA Kit
419	abk4541H	Human C3b(Complement Fragment 3b) ELISA Kit
420	abk5541H	Human C3c(Complement Fragment 3c) ELISA Kit
421	abk6541H	Human NPYR2(Neuropeptide Y Receptor Y2) ELISA Kit
422	abk7541H	Human C5b(Complement Fragment 5b) ELISA Kit
423	abk8541H	Human NPYR5(Neuropeptide Y Receptor Y5) ELISA Kit
424	abk9541H	Human CTGF(Connective Tissue Growth Factor) ELISA Kit
425	abk0641H	Human CX26(Connexin 26) ELISA Kit
426	abk1641H	Human CHGB(Chromogranin B) ELISA Kit
427	abk2641H	Human DEFα1(Defensin Alpha 1, Neutrophil) ELISA Kit
428	abk3641H	Human CTX I (Cross Linked C-telopeptide of Type I Collagen) ELISA Kit
429	abk4641H	Human NTX I (Cross Linked N-telopeptide of Type I Collagen) ELISA Kit
430	abk5641H	Human CTX II (Cross Linked C-telopeptide of Type II Collagen) ELISA Kit
431	abk6641H	Human LDHB(Lactate Dehydrogenase B) ELISA Kit
432	abk7641H	Human HMGCS1(3-Hydroxy-3-Methylglutaryl-Coenzyme A Synthase 1) ELISA Kit
433	abk8641H	Human PIIICP(Procollagen III C-Terminal ProPeptide) ELISA Kit
434	abk9641H	Human DKK4(Dickkopf Related Protein 4) ELISA Kit
435	abk1741H	Human CLDN11(Claudin 11) ELISA Kit
436	abk2741H	Human CXCR1(CXC-Chemokine Receptor 1) ELISA Kit
437	abk3741H	Human SCG2(Secretogranin II) ELISA Kit
438	abk4741H	Human MYCBP(C-Myc Binding Protein) ELISA Kit
439	abk8353H	Human MYC(V-Myc Myelocytomatosis Viral Oncogene Homolog) ELISA Kit
440	abk6741H	Human CLEC4E(C-Type Lectin Domain Family 4 Member E) ELISA Kit
441	abk7741H	Human LDHC(Lactate Dehydrogenase C) ELISA Kit
442	abk8741H	Human EphA2(Tyrosine protein kinase receptor A2) ELISA Kit
443	abk9741H	Human DKK2(Dickkopf Related Protein 2) ELISA Kit
444	abk0841H	Human D-LDH(D-Lactate Dehydrogenase) ELISA Kit
445	abk1841H	Human PA2G4(Proliferation Associated Protein 2G4) ELISA Kit
446	abk2841H	Human COL1α1(Collagen Type I Alpha 1) ELISA Kit
447	abk3841H	Human COL6α3(Collagen Type VI Alpha 3) ELISA Kit
448	abk4841H	Human COL13(Collagen Type X III) ELISA Kit
449	abk5841H	Human COL14(Collagen Type X IV) ELISA Kit
450	abk8841H	Human EFNA1(Ephrin A1) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
451	abk9841H	Human GATA3(GATA Binding Protein 3) ELISA Kit
452	abk0941H	Human GATA4(GATA Binding Protein 4) ELISA Kit
453	abk1941H	Human GNβ1(G Protein Beta 1) ELISA Kit
454	abk2941H	Human HHIP(Hedgehog Interacting Protein) ELISA Kit
455	abk3941H	Human H-Cadherin(Cadherin, Heart) ELISA Kit
456	abk4941H	Human NRG-4(Neuregulin 4) ELISA Kit
457	abk5941H	Human PTPRK(Protein Tyrosine Phosphatase Receptor Type K) ELISA Kit
458	abk6941H	Human SELL(L-Selectin) ELISA Kit
459	abk7941H	Human MXD1(MAX Dimerization Protein 1) ELISA Kit
460	abk8941H	Human HBμ(Hemoglobin, Mu) ELISA Kit
461	abk9941H	Human NAGase(N-Acetyl Beta-D-Glucosaminidase) ELISA Kit
462	abk0051H	Human NRG-2(Neuregulin 2) ELISA Kit
463	abk1051H	Human NRG-3(Neuregulin 3) ELISA Kit
464	abk2051H	Human N-MID-OT(N-MID Osteocalcin) ELISA Kit
465	abk3051H	Human PTPRN(Protein Tyrosine Phosphatase Receptor Type N) ELISA Kit
466	abk4051H	Human GAS(N-Acetylgalactosamine 6-Sulfatase) ELISA Kit
467	abk5051H	Human PTPRO(Protein Tyrosine Phosphatase Receptor Type O) ELISA Kit
468	abk6051H	Human P27(P27 Protein) ELISA Kit
469	abk7051H	Human TP53(Tumor Protein p53) ELISA Kit
470	abk8051H	Human PDAP1(Pdgfa Associated Protein 1) ELISA Kit
471	abk9051H	Human PDCN(Podocin) ELISA Kit
472	abk0151H	Human Preptin ELISA Kit
473	abk1151H	Human P-Cadherin(Cadherin, Placental) ELISA Kit
474	abk2151H	Human P-GP(Permeability Glycoprotein) ELISA Kit
475	abk3151H	Human RAPGEF1(Rap Guanine Nucleotide Exchange Factor 1) ELISA Kit
476	abk4151H	Human Rac-GAP1(Rac-GTPase Activating Protein 1) ELISA Kit
477	abk5151H	Human RHOB(Ras Homolog Gene Family, Member B) ELISA Kit
478	abk6151H	Human RCOR3(Rest Corepressor 3) ELISA Kit
479	abk7151H	Human RHOA(Ras Homolog Gene Family, Member A) ELISA Kit
480	abk8151H	Human Rock-2(Rho Associated Coiled Coil Containing Protein Kinase 2) ELISA Kit
481	abk9151H	Human RND1(Rho Family GTPase 1) ELISA Kit
482	abk0251H	Human RUNX2(Runt-Related Transcription Factor 2) ELISA Kit
483	abk1251H	Human COL9α3(Collagen Type IX Alpha 3) ELISA Kit
484	abk3251H	Human SHC/SLP-76(SH2-Containing Protein) ELISA Kit
485	abk4251H	Human Slit2(Slit Homolog 2) ELISA Kit
486	abk5251H	Human Smad1/MADH1(Mothers Against Decapentaplegic Homolog 1) ELISA Kit
487	abk6251H	Human Smad4/MADH4(Mothers Against Decapentaplegic Homolog 4) ELISA Kit
488	abk7251H	Human Smad7/MADH7(Mothers Against Decapentaplegic Homolog 7) ELISA Kit
489	abk8251H	Human Sp100(Sp100 Nuclear Antigen) ELISA Kit
490	abk9251H	Human PIAS1(Protein Inhibitor of Activated STAT 1) ELISA Kit
491	abk0351H	Human PIAS3(Protein Inhibitor of Activated STAT3) ELISA Kit
492	abk1351H	Human PCDHβ16(Protocadherin Beta 16) ELISA Kit
493	abk2351H	Human COL5α2(Collagen Type V Alpha 2) ELISA Kit
494	abk3351H	Human SAG(S Antigen) ELISA Kit
495	abk4351H	Human THP(Tamm-Horsfall Glycoprotein) ELISA Kit
496	abk5351H	Human TBP(TATA Binding Protein) ELISA Kit
497	abk6351H	Human MAPτ(Microtubule Associated Protein Tau/Tau Protein) ELISA Kit
498	abk7351H	Human TIEG1(TGF Beta Inducible Early Response Gene 1) ELISA Kit
499	abk8351H	Human TLR9(Toll-Like Receptor 9) ELISA Kit
500	abk9351H	Human PAR3(Protease Activated Receptor 3) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
501	abk1451H	Human XCR1(XC-Chemokine Receptor 1) ELISA Kit
502	abk2451H	Human PCDH15(Protocadherin 15) ELISA Kit
503	abk3451H	Human α 1BG(Alpha-1-B-Glycoprotein) ELISA Kit
504	abk4451H	Human α 1-AT(Alpha 1-Antitrypsin) ELISA Kit
505	abk5451H	Human PLC γ 2(Phospholipase C Gamma 2, Phosphatidylinositol Specific) ELISA Kit
506	abk6451H	Human α 2-AP(α 2-Antiplasmin) ELISA Kit
507	abk7451H	Human IDU α (Alpha-L-Iduronidase) ELISA Kit
508	abk8451H	Human NAGLU(N-Acetyl Alpha-D-Glucosaminidase) ELISA Kit
509	abk9451H	Human GL α (Galactosidase Alpha) ELISA Kit
510	abk0551H	Human α Actin(Alpha-Actin) ELISA Kit
511	abk1551H	Human α -La(Alpha-Lactalbumin) ELISA Kit
512	abk2551H	Human SNC α (Synuclein Alpha) ELISA Kit
513	abk3551H	Human PDGFR β (Platelet Derived Growth Factor Receptor Beta) ELISA Kit
514	abk4551H	Human P II CP(Procollagen II C-Terminal ProPeptide) ELISA Kit
515	abk5551H	Human PCDH γ A2(Protocadherin Gamma A2) ELISA Kit
516	abk6551H	Human MAN1A1(Mannosidase Alpha Class 1A Member 1) ELISA Kit
517	abk7551H	Human β -LPH(β -Lipotropic Hormone) ELISA Kit
518	abk8551H	Human DEFB1(Defensin Beta 1) ELISA Kit
519	abk9551H	Human DEF β 2/DEFB2(Defensin Beta 2) ELISA Kit
520	abk0651H	Human β Manase(β Mannosidase) ELISA Kit
521	abk1651H	Human CYSLTR1(Cysteiny Leukotriene Receptor 1) ELISA Kit
522	abk2651H	Human IRAK3(Interleukin 1 Receptor Associated Kinase 3) ELISA Kit
523	abk3651H	Human GUS β (β -glucuronidase) ELISA Kit
524	abk4651H	Human SYCN(Syncollin) ELISA Kit
525	abk5651H	Human TUB β (β -Tubulin) ELISA Kit
526	abk6651H	Human γ -ECS(γ -Glutamyl Systeine Synthetase) ELISA Kit
527	abk7651H	Human PLA2G10(Phospholipase A2, Group X) ELISA Kit
528	abk8651H	Human γ GT1(Gamma Glutamyltransferase 1) ELISA Kit
529	abk9651H	Human POMC(Pro-Opiomelanocortin) ELISA Kit
530	abk1751H	Human PIK3C β (Phosphoinositide-3-Kinase Catalytic Beta Polypeptide) ELISA Kit
531	abk2751H	Human OIT3(Oncoprotein Induced Transcript 3) ELISA Kit
532	abk3751H	Human PIK3C α (Phosphoinositide-3-Kinase Catalytic Delta Polypeptide) ELISA Kit
533	abk4751H	Human OCT2(Octamer-binding protein 2) ELISA Kit
534	abk5751H	Human OCT4(Octamer-binding protein 4) ELISA Kit
535	abk6751H	Human PLA2G12B(Phospholipase A2, Group XII B) ELISA Kit
536	abk7751H	Human IL-2R α (Interleukin 2 Receptor Alpha) ELISA Kit
537	abk8751H	Human IL-2R β (Interleukin 2 Receptor Beta) ELISA Kit
538	abk9751H	Human NFIL3(Nuclear Factor, Interleukin 3 Regulated) ELISA Kit
539	abk0851H	Human LILRB4(Leukocyte Immunoglobulin Like Receptor Subfamily B Member 4) ELISA Kit
540	abk1851H	Human LIFR(Leukemia Inhibitory Factor Receptor) ELISA Kit
541	abk2851H	Human GARS(Glycyl tRNA Synthetase) ELISA Kit
542	abk3851H	Human EPHA1(Ephrin Type A Receptor 1) ELISA Kit
543	abk4851H	Human SPOP(Speckle Type POZ Protein) ELISA Kit
544	abk5851H	Human GP4/CD36(Platelet Membrane Glycoprotein IV) ELISA Kit
545	abk6851H	Human PLA2G4D(Phospholipase A2, Group IV D) ELISA Kit
546	abk7851H	Human CYSLTR2(Cysteiny Leukotriene Receptor 2) ELISA Kit
547	abk8851H	Human GAL1(Galectin 1) ELISA Kit
548	abk9851H	Human GAL12(Galectin 12) ELISA Kit
549	abk0951H	Human GAL14(Galectin 14) ELISA Kit
550	abk1951H	Human GAL2(Galectin 2) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
551	abk2951H	Human EPHA3(Ephrin Type A Receptor 3) ELISA Kit
552	abk3951H	Human GAL4(Galectin 4) ELISA Kit
553	abk4951H	Human GAL8(Galectin 8) ELISA Kit
554	abk5951H	Human GAL7(Galectin 7) ELISA Kit
555	abk6951H	Human GAL9(Galectin 9) ELISA Kit
556	abk7951H	Human UBC(Urinary Bladder Cancer Antigen) ELISA Kit
557	abk8951H	Human PLA2R1(Phospholipase A2 Receptor 1) ELISA Kit
558	abk9951H	Human PLB(Phospholipase B) ELISA Kit
559	abk1061H	Human EFNA3(Ephrin A3) ELISA Kit
560	abk2061H	Human PLD2(Phospholipase D2) ELISA Kit
561	abk3061H	Human PDXK(Pyridoxal Kinase) ELISA Kit
562	abk4061H	Human OCLN(Occludin) ELISA Kit
563	abk5061H	Human SEMA7A(Semaphorin 7A) ELISA Kit
564	abk6061H	Human PLVAP(Plasmalemma Vesicle Associated Protein) ELISA Kit
565	abk7061H	Human HAGH(Hydroxyacylglutathione Hydrolase) ELISA Kit
566	abk8061H	Human SPJ(Surfactant Associated Protein J) ELISA Kit
567	abk9061H	Human GPRC5A(G Protein Coupled Receptor, Family C, Group 5, Member A) ELISA Kit
568	abk0161H	Human FABP5(Fatty Acid Binding Protein 5, Epidermal) ELISA Kit
569	abk1161H	Human PDK4(Pyruvate Dehydrogenase Kinase Isozyme 4) ELISA Kit
570	abk2161H	Human PDP(Pyruvate Dehydrogenase Phosphatase) ELISA Kit
571	abk3161H	Human PPOX(Protoporphyrinogen Oxidase) ELISA Kit
572	abk4161H	Human PRB1(Basic Salivary Proline Rich Protein 1) ELISA Kit
573	abk5161H	Human FOXC1(Forkhead Box Protein C1) ELISA Kit
574	abk6161H	Human FOXO1(Forkhead Box Protein O1) ELISA Kit
575	abk7161H	Human FOXO3(Forkhead Box Protein O3) ELISA Kit
576	abk8161H	Human FOXO4(Forkhead Box Protein O4) ELISA Kit
577	abk9161H	Human FOXP3(Forkhead Box Protein P3) ELISA Kit
578	abk0261H	Human LATS(Long Acting Thyroid Stimulator) ELISA Kit
579	abk1261H	Human EFNA4(Ephrin A4) ELISA Kit
580	abk2261H	Human TFF3(Trefoil Factor 3, Intestinal) ELISA Kit
581	abk3261H	Human TAGLN(Transgelin) ELISA Kit
582	abk4261H	Human TEP1(Telomerase Associated Protein 1) ELISA Kit
583	abk5261H	Human PRB2(Basic Salivary Proline Rich Protein 2) ELISA Kit
584	abk8261H	Human FGF23(Fibroblast Growth Factor 23) ELISA Kit
585	abk9261H	Human FGF9(Fibroblast Growth Factor 9) ELISA Kit
586	abk0361H	Human FGFBP1(Fibroblast Growth Factor Binding Protein 1) ELISA Kit
587	abk1361H	Human FSHR(Follicle Stimulating Hormone Receptor) ELISA Kit
588	abk2361H	Human TGM1(Transglutaminase 1, Keratinocyte) ELISA Kit
589	abk3361H	Human PACAPR(Pituitary Adenylate Cyclase Activating Polypeptide Receptor 1) ELISA Kit
590	abk4361H	Human EPHA10(Ephrin Type A Receptor 10) ELISA Kit
591	abk5361H	Human ER β (Estrogen Receptor Beta) ELISA Kit
592	abk6361H	Human ER α (Estrogen Receptor Alpha) ELISA Kit
593	abk7361H	Human EFNA5(Ephrin A5) ELISA Kit
594	abk8361H	Human MAPKAPK3(Mitogen Activated Protein Kinase Activated Protein Kinase 3) ELISA Kit
595	abk0461H	Human TEF(Thyrotroph Embryonic Factor) ELISA Kit
596	abk2461H	Human GHRH(Growth Hormone Releasing Hormone) ELISA Kit
597	abk4461H	Human SCTR(Secretin Receptor) ELISA Kit
598	abk6461H	Human FABP7(Fatty Acid Binding Protein 7, Brain) ELISA Kit
599	abk7461H	Human MCP-4/CCL13(Monocyte Chemotactic Protein 4) ELISA Kit
600	abk8461H	Human GPER1(G Protein Coupled Estrogen Receptor 1) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
601	abk9461H	Human ELN(Elastin) ELISA Kit
602	abk0561H	Human ELA4(Elastase 4) ELISA Kit
603	abk2561H	Human PDIA2(Protein Disulfide Isomerase A2) ELISA Kit
604	abk3561H	Human PDIA3(Protein Disulfide Isomerase A3) ELISA Kit
605	abk4561H	Human PRSS2(Protease, Serine, 2) ELISA Kit
606	abk5561H	Human PKBG(Protein Kinase B Gamma) ELISA Kit
607	abk6561H	Human PICK1(Protein Interacting With Protein Kinase C, Alpha 1) ELISA Kit
608	abk7561H	Human PKCβ1(Protein Kinase C Beta 1) ELISA Kit
609	abk8561H	Human PKCγ(Protein Kinase C Gamma) ELISA Kit
610	abk9561H	Human PKCδ1(Protein Kinase C Delta) ELISA Kit
611	abk0661H	Human PKCε(Protein Kinase C Epsilon) ELISA Kit
612	abk1661H	Human PKCζ(Protein Kinase C Zeta) ELISA Kit
613	abk2661H	Human PKCθ(Protein Kinase C Theta) ELISA Kit
614	abk3661H	Human PRSS23(Protease, Serine 23) ELISA Kit
615	abk4661H	Human PKD1(Protein Kinase D1) ELISA Kit
616	abk5661H	Human PKN1(Protein Kinase N1) ELISA Kit
617	abk6661H	Human PKIα(Protein Kinase Inhibitor Alpha) ELISA Kit
618	abk7661H	Human PKIβ(Protein Kinase Inhibitor Beta) ELISA Kit
619	abk8661H	Human PKIγ(Protein Kinase Inhibitor Gamma) ELISA Kit
620	abk9661H	Human PTPRJ(Protein Tyrosine Phosphatase Receptor Type J) ELISA Kit
621	abk0761H	Human PPP1R1A(Protein Phosphatase 1 Regulatory Subunit A) ELISA Kit
622	abk1761H	Human PRSS1(Protease, Serine, 1) ELISA Kit
623	abk2761H	Human PRSS12(Protease, Serine 12) ELISA Kit
624	abk4761H	Human PAR2(Protease Activated Receptor 2) ELISA Kit
625	abk5761H	Human PCI(Protein C Inhibitor) ELISA Kit
626	abk6761H	Human PDI(Protein Disulfide Isomerase) ELISA Kit
627	abk7761H	Human PDIA5(Protein Disulfide Isomerase A5) ELISA Kit
628	abk8761H	Human NGF(Nerve Growth Factor) ELISA Kit
629	abk9761H	Human PRSS33(Protease, Serine 33) ELISA Kit
630	abk0861H	Human LDL(Low Density Lipoprotein) ELISA Kit
631	abk1861H	Human LDLR(Low Density Lipoprotein Receptor) ELISA Kit
632	abk2861H	Human LRP6(Low Density Lipoprotein Receptor Related Protein) ELISA Kit
633	abk3861H	Human RETN(Resistin) ELISA Kit
634	abk4861H	Human RPS6Kβ1(Ribosomal Protein S6 Kinase Beta 1) ELISA Kit
635	abk5861H	Human ETFβ(Electron Transfer Flavoprotein Beta Polypeptide) ELISA Kit
636	abk6861H	Human APP(Amyloid Precursor Protein) ELISA Kit
637	abk7861H	Human M-AChR M1(Muscarinic Acetylcholine Receptor M1) ELISA Kit
638	abk8861H	Human M-AChR M3(Muscarinic Acetylcholine Receptor M3) ELISA Kit
639	abk9861H	Human D1R/DRD1(Dopamine Receptor D1) ELISA Kit
640	abk0961H	Human D2R/DRD2(Dopamine Receptor D2) ELISA Kit
641	abk1961H	Human D3R/DRD3(Dopamine Receptor D3) ELISA Kit
642	abk2961H	Human D4R/DRD4(Dopamine Receptor D4) ELISA Kit
643	abk3961H	Human D5R/DRD5(Dopamine Receptor D5) ELISA Kit
644	abk4961H	Human DDC(Dopamine Decarboxylase) ELISA Kit
645	abk5961H	Human DAT(Dopamine Transporter) ELISA Kit
646	abk6961H	Human DT(Dopachrome Tautomerase) ELISA Kit
647	abk7961H	Human PIGR/SC(Polymeric Immunoglobulin Receptor/Membrane Secretory Component) ELISA Kit
648	abk8961H	Human PTN(Pleiotrophin) ELISA Kit
649	abk9961H	Human DAO(Diamine Oxidase) ELISA Kit
650	abk0071H	Human UGT2B4(Uridine Diphosphate Glucuronosyltransferase 2 Family Polypeptide B4) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
651	abk1071H	Human DLD(Dihydrolipoyl Dehydrogenase) ELISA Kit
652	abk2071H	Human DLAT(Dihydrolipoyl Transacetylase) ELISA Kit
653	abk3071H	Human DPYSL2(Dihydropyrimidinase Like Protein 2) ELISA Kit
654	abk4071H	Human DNM1(Dynamin 1) ELISA Kit
655	abk5071H	Human DNM2(Dynamin 2) ELISA Kit
656	abk6071H	Human Ub(Ubiquitin) ELISA Kit
657	abk7071H	Human UBD(Ubiquitin D) ELISA Kit
658	abk8071H	Human UBE1L(Ubiquitin Activating Enzyme E1 Like Protein) ELISA Kit
659	abk9071H	Human UBE2C(Ubiquitin Conjugating Enzyme E2C) ELISA Kit
660	abk0171H	Human UBE3A(Ubiquitin Protein Ligase E3A) ELISA Kit
661	abk1171H	Human NNE(Non-Neuronal Enolase) ELISA Kit
662	abk2171H	Human SDC3(Syndecan 3) ELISA Kit
663	abk3171H	Human FBLN5(Fibulin 5) ELISA Kit
664	abk4171H	Human REG1β(Regenerating Islet Derived Protein 1 Beta) ELISA Kit
665	abk5171H	Human SPC(Pulmonary Surfactant Associated Protein C) ELISA Kit
666	abk6171H	Human PARC/CCL18(Pulmonary Activation Regulated Chemokine) ELISA Kit
667	abk7171H	Human LRP(Lung Resistance-Related Protein) ELISA Kit
668	abk8171H	Human PRSS8(Protease, Serine 8) ELISA Kit
669	abk9171H	Human SPARCL1(SPARC Like Protein 1) ELISA Kit
670	abk1271H	Human sIgA(Secretory Immunoglobulin A) ELISA Kit
671	abk2271H	Human SLPI(Secretory Leukocyte Protease Inhibitor) ELISA Kit
672	abk3271H	Human SAP(Serum Amyloid P Component) ELISA Kit
673	abk4271H	Human PAR4(Protease Activated Receptor 4) ELISA Kit
674	abk5271H	Human GPR120(G Protein Coupled Receptor 120) ELISA Kit
675	abk6271H	Human REG3α(Regenerating Islet Derived Protein 3 Alpha) ELISA Kit
676	abk7271H	Human Cyr61(Cysteine Rich Protein, Angiogenic Inducer 61) ELISA Kit
677	abk8271H	Human REG3γ(Regenerating Islet Derived Protein 3 Gamma) ELISA Kit
678	abk9271H	Human S100A8(S100 Calcium Binding Protein A8) ELISA Kit
679	abk0371H	Human S100A9(S100 Calcium Binding Protein A9) ELISA Kit
680	abk1371H	Human GPR37(G Protein Coupled Receptor 37) ELISA Kit
681	abk2371H	Human S100A11(S100 Calcium Binding Protein A11) ELISA Kit
682	abk3371H	Human S100A12(S100 Calcium Binding Protein A12) ELISA Kit
683	abk4371H	Human S100A13(S100 Calcium Binding Protein A13) ELISA Kit
684	abk5371H	Human S100A6(S100 Calcium Binding Protein A6) ELISA Kit
685	abk6371H	Human S100A7(S100 Calcium Binding Protein A7) ELISA Kit
686	abk8371H	Human SDC1(Syndecan 1) ELISA Kit
687	abk9371H	Human SCF(Stem Cell Factor) ELISA Kit
688	abk0471H	Human SCFR(Stem Cell Growth Factor Receptor) ELISA Kit
689	abk1471H	Human GAL(Galanin) ELISA Kit
690	abk2471H	Human GALR1(Galanin Receptor 1) ELISA Kit
691	abk3471H	Human PGAM1(Phosphoglycerate Mutase 1) ELISA Kit
692	abk4471H	Human PK(Pyruvate Kinase, Liver And RBC) ELISA Kit
693	abk5471H	Human PYGL(Glycogen Phosphorylase, Liver) ELISA Kit
694	abk6471H	Human SFRP2(Secreted Frizzled Related Protein 2) ELISA Kit
695	abk7471H	Human GP-73(Golgi Protein 73) ELISA Kit
696	abk8471H	Human SIAE(Sialic Acid Acetylesterase) ELISA Kit
697	abk9471H	Human ST2(Syntenin 2) ELISA Kit
698	abk0571H	Human SIGLEC1(Sialic Acid Binding Immunoglobulin Like Lectin 1) ELISA Kit
699	abk1571H	Human SIGLEC10(Sialic Acid Binding Immunoglobulin Like Lectin 10) ELISA Kit
700	abk2571H	Human SIGLEC2(Sialic Acid Binding Immunoglobulin Like Lectin 2) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
701	abk3571H	Human PADI4(Peptidyl Arginine Deiminase Type IV) ELISA Kit
702	abk4571H	Human SIGLEC8(Sialic Acid Binding Immunoglobulin Like Lectin 8) ELISA Kit
703	abk5571H	Human Smad3/MADH3(Mothers Against Decapentaplegic Homolog 3) ELISA Kit
704	abk6571H	Human OFQ/N(Orphanin FQ/Nociceptin) ELISA Kit
705	abk7571H	Human GCLM(Glutamate Cysteine Ligase, Modifier Subunit) ELISA Kit
706	abk8571H	Human GSTθ1/GSTt1(Glutathione S Transferase Theta 1) ELISA Kit
707	abk9571H	Human GSTθ2/GSTt2(Glutathione S Transferase Theta 2) ELISA Kit
708	abk0671H	Human GSTω1(Glutathione S Transferase Omega 1) ELISA Kit
709	abk1671H	Human HDAC2(Histone Deacetylase 2) ELISA Kit
710	abk2671H	Human C1QTNF1(C1q and Tumor Necrosis Factor Related Protein 1) ELISA Kit
711	abk4453H	Human GR(Glutathione Reductase) ELISA Kit
712	abk4671H	Human GSTpi(Glutathione S Transferases Pi) ELISA Kit
713	abk5671H	Human OPG(Osteoprotegerin) ELISA Kit
714	abk6671H	Human OC/BGP(Osteocalcin) ELISA Kit
715	abk7671H	Human ELA3B(Elastase 3B) ELISA Kit
716	abk8671H	Human MUSK(Muscle Skeletal Receptor Tyrosine Kinase) ELISA Kit
717	abk9671H	Human GSTα5(Glutathione S Transferase Alpha 5) ELISA Kit
718	abk0771H	Human OPN(Osteopontin) ELISA Kit
719	abk1771H	Human MPIF-1/CCL23(Myeloid Progenitor Inhibitory Factor 1) ELISA Kit
720	abk2771H	Human SOX9(Sex Determining Region Y Box Protein 9) ELISA Kit
721	abk3771H	Human GREM1(Gremlin 1) ELISA Kit
722	abk4771H	Human ON(Osteonectin) ELISA Kit
723	abk5771H	Human SF1(Steroidogenic Factor 1) ELISA Kit
724	abk6771H	Human PAPPA2(Pappalysin 2) ELISA Kit
725	abk7771H	Human P-CK(Pan-Cytokeratin) ELISA Kit
726	abk8771H	Human ALDOC(Aldolase C, Fructose Bisphosphate) ELISA Kit
727	abk9771H	Human PMP3(Peroxisomal Membrane Protein 3) ELISA Kit
728	abk0871H	Human PPARγC1α(Peroxisome Proliferator Activated Receptor Gamma Coactivator 1 Alpha) ELISA Kit
729	abk1871H	Human PPAR-α(Peroxisome Proliferator Activated Receptor Alpha) ELISA Kit
730	abk2871H	Human PPAR-γ(Peroxisome Proliferator Activated Receptor Gamma) ELISA Kit
731	abk3871H	Human NXF1(Nuclear RNA Export Factor 1) ELISA Kit
732	abk4871H	Human NMP-22(Nuclear Matrix Protein 22) ELISA Kit
733	abk5871H	Human NUP107(Nucleoporin 107kDa) ELISA Kit
734	abk6871H	Human NUP133(Nucleoporin 133kDa) ELISA Kit
735	abk7871H	Human NUP155(Nucleoporin 155kDa) ELISA Kit
736	abk8871H	Human NUP188(Nucleoporin 188kDa) ELISA Kit
737	abk9871H	Human NUP205(Nucleoporin 205kDa) ELISA Kit
738	abk0971H	Human NUP214(Nucleoporin 214kDa) ELISA Kit
739	abk1971H	Human NUP35(Nucleoporin 35kDa) ELISA Kit
740	abk2971H	Human NUP37(Nucleoporin 37kDa) ELISA Kit
741	abk3971H	Human NUP50(Nucleoporin 50kDa) ELISA Kit
742	abk4971H	Human NUP85(Nucleoporin 85kDa) ELISA Kit
743	abk5971H	Human NUP98(Nucleoporin 98kDa) ELISA Kit
744	abk6971H	Human NRIP1(Nuclear Receptor Interacting Protein 1) ELISA Kit
745	abk7971H	Human RNASE3/ECP(Ribonuclease A3/Eosinophil Cationic Protein) ELISA Kit
746	abk8971H	Human RI(Ribonuclease Inhibitor) ELISA Kit
747	abk9971H	Human RNASE(Ribonuclease) ELISA Kit
748	abk0081H	Human RPL10(Ribosomal Protein L10) ELISA Kit
749	abk1081H	Human RPL26(Ribosomal Protein L26) ELISA Kit
750	abk2081H	Human MART/Melan-A(Melanoma Marker A) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
751	abk3081H	Human MAGE(Melanoma-associated Antigen) ELISA Kit
752	abk4081H	Human MLZE(Melanoma Derived Leucine Zipper Extra Nuclear Factor) ELISA Kit
753	abk5081H	Human THRSP(Thyroid Hormone Responsive) ELISA Kit
754	abk6081H	Human EPOR(Erythropoietin Receptor) ELISA Kit
755	abk7081H	Human EXT1(Exostoses 1) ELISA Kit
756	abk8081H	Human GSTκ1(Glutathione S Transferase Kappa 1) ELISA Kit
757	abk9081H	Human LTA4H(Leukotriene A4 Hydrolase) ELISA Kit
758	abk0181H	Human GSTα4(Glutathione S Transferase Alpha 4) ELISA Kit
759	abk1181H	Human SNCγ(Synuclein Gamma) ELISA Kit
760	abk2181H	Human TXLNa(Taxilin Alpha) ELISA Kit
761	abk3181H	Human TXLNb(Taxilin Beta) ELISA Kit
762	abk4181H	Human Sec23b(Protein transport protein Sec23B) ELISA Kit
763	abk5181H	Human RNF4(Ring Finger Protein 4) ELISA Kit
764	abk6181H	Human FABP6(Fatty Acid Binding Protein 6, Ileal) ELISA Kit
765	abk7181H	Human TMPPRS2(Transmembrane Protease, Serine 2) ELISA Kit
766	abk8181H	Human TNNI1(Troponin I Type 1, Slow Skeletal) ELISA Kit
767	abk9181H	Human FSCN(Fascin) ELISA Kit
768	abk0281H	Human TNNC2(Troponin C Type 2) ELISA Kit
769	abk1281H	Human TNC(Tenascin C) ELISA Kit
770	abk2281H	Human TNR(Tenascin R) ELISA Kit
771	abk3281H	Human TNX(Tenascin X) ELISA Kit
772	abk4281H	Human IgA1(Immunoglobulin A1) ELISA Kit
773	abk5281H	Human MYL1(Myosin Light Chain 1, Alkali, Fast Skeletal) ELISA Kit
774	abk6281H	Human MLCK/MYLK(Myosin Light Chain Kinase) ELISA Kit
775	abk7281H	Human MYH1(Myosin Heavy Chain 1, Skeletal Muscle, Adult) ELISA Kit
776	abk8281H	Human FABP3(Fatty Acid Binding Protein 3, Muscle and Heart) ELISA Kit
777	abk9281H	Human PYGM(Glycogen Phosphorylase, Muscle) ELISA Kit
778	abk0381H	Human CKM(Creatine Kinase, Muscle) ELISA Kit
779	abk1381H	Human MYOG(Myogenin) ELISA Kit
780	abk2381H	Human DMD(Dystrophin) ELISA Kit
781	abk3381H	Human MSTN(Myostatin) ELISA Kit
782	abk4381H	Human PUMA(P53 Upregulated Modulator of Apoptosis) ELISA Kit
783	abk5381H	Human MGP(Matrix Gla Protein) ELISA Kit
784	abk6381H	Human MMP-11(Matrix Metalloproteinase 11) ELISA Kit
785	abk7381H	Human MGST1(Microsomal Glutathione S Transferase 1) ELISA Kit
786	abk8381H	Human MMP-2(Matrix Metalloproteinase 2) ELISA Kit
787	abk9381H	Human MMP-3(Matrix Metalloproteinase 3) ELISA Kit
788	abk0481H	Human MMP-17(Matrix Metalloproteinase 17) ELISA Kit
789	abk1481H	Human MMP-24(Matrix Metalloproteinase 24) ELISA Kit
790	abk2481H	Human MMP-7(Matrix Metalloproteinase 7) ELISA Kit
791	abk4481H	Human ELA2B(Elastase 2B) ELISA Kit
792	abk5481H	Human TIMP-2(Tissue Inhibitors of Metalloproteinase 2) ELISA Kit
793	abk6481H	Human TIMP-3(Tissue Inhibitors of Metalloproteinase 3) ELISA Kit
794	abk7481H	Human TIMP-4(Tissue Inhibitors of Metalloproteinase 4) ELISA Kit
795	abk8481H	Human LGALS3BP(Lectin Galactoside Binding, Soluble 3 Binding Protein) ELISA Kit
796	abk9481H	Human TNNI2(Troponin I Type 2, Fast Skeletal) ELISA Kit
797	abk0581H	Human TP/TXA2R(Thromboxane Receptor) ELISA Kit
798	abk1581H	Human SDF-2(Stromal Cell Derived Factor 2) ELISA Kit
799	abk2581H	Human SDF2L1(Stromal Cell Derived Factor 2 Like Protein 1) ELISA Kit
800	abk3581H	Human SDF-4(Stromal Cell Derived Factor 4) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
801	abk4581H	Human TDGF1(Teratocarcinoma Derived Growth Factor 1) ELISA Kit
802	abk5581H	Human KLK8(Kallikrein 8) ELISA Kit
803	abk6581H	Human VLDLR(Very Low Density Lipoprotein Receptor) ELISA Kit
804	abk7581H	Human GAL3(Galectin 3) ELISA Kit
805	abk8581H	Human FPR2(Formyl Peptide Receptor 2) ELISA Kit
806	abk9581H	Human PTHR2(Parathyroid Hormone Receptor 2) ELISA Kit
807	abk0681H	Human PTHrP(Parathyroid Hormone Related Protein) ELISA Kit
808	abk1681H	Human TPP1(Tripeptidyl Peptidase 1) ELISA Kit
809	abk2681H	Human CD276(Cluster of Differentiation 276) ELISA Kit
810	abk3681H	Human TBG(Thyroxine Binding Globulin) ELISA Kit
811	abk4681H	Human PAH(Phenylalanine Hydroxylase) ELISA Kit
812	abk5681H	Human GJβ1(Gap Junction Protein Beta 1) ELISA Kit
813	abk7681H	Human PCPE1(Procollagen C Proteinase Enhancer 1) ELISA Kit
814	abk8681H	Human GDNF(Glial Cell Line Derived Neurotrophic Factor) ELISA Kit
815	abk9681H	Human CDSN(Corneodesmosin) ELISA Kit
816	abk0781H	Human CNTN1(Contactin 1) ELISA Kit
817	abk1781H	Human CNTN2(Contactin 2) ELISA Kit
818	abk2781H	Human CNTN3(Contactin 3) ELISA Kit
819	abk3781H	Human CNTNAP1(Contactin Associated Protein 1) ELISA Kit
820	abk4781H	Human CNTNAP4(Contactin Associated Protein 4) ELISA Kit
821	abk5781H	Human ELMO2(Engulfment And Cell Motility 2) ELISA Kit
822	abk6781H	Human Des(Desmin) ELISA Kit
823	abk7781H	Human SSRP1(Structure Specific Recognition Protein 1) ELISA Kit
824	abk8781H	Human HPR(Haptoglobin Related Protein) ELISA Kit
825	abk9781H	Human TNNT1(Troponin T Type 1, Slow Skeletal) ELISA Kit
826	abk0881H	Human TNPO2(Transportin 2) ELISA Kit
827	abk1881H	Human UCP-2(Uncoupling Protein 2, Mitochondrial) ELISA Kit
828	abk2881H	Human HLTF(Helicase Like Transcription Factor) ELISA Kit
829	abk3881H	Human MT2(Metalllothionein 2) ELISA Kit
830	abk4881H	Human TJP1(Tight junction protein ZO-1) ELISA Kit
831	abk5881H	Human CLC(Charcot Leyden Crystal Protein) ELISA Kit
832	abk6881H	Human ACVB(Activin B) ELISA Kit
833	abk7881H	Human CRYαB(Crystallin Alpha B) ELISA Kit
834	abk8881H	Human Sp17(Sperm Protein 17) ELISA Kit
835	abk9881H	Human MSP(Macrophage Stimulating Protein) ELISA Kit
836	abk0981H	Human MCSFR(Macrophage Colony Stimulating Factor Receptor) ELISA Kit
837	abk1981H	Human MIP-5(Macrophage Inflammatory Protein 5) ELISA Kit
838	abk2981H	Human FZD8(Frizzled Homolog 8) ELISA Kit
839	abk3981H	Human HAAO(3-Hydroxyanthranilate 3,4-Dioxygenase) ELISA Kit
840	abk4981H	Human CRHBP(Corticotropin Releasing Hormone Binding Protein) ELISA Kit
841	abk5981H	Human αHSP(Alpha-Hemoglobin Stabilizing Protein) ELISA Kit
842	abk6981H	Human MMP-28(Matrix Metalloproteinase 28) ELISA Kit
843	abk7981H	Human SPTAN1(Alpha-Fodrin) ELISA Kit
844	abk8981H	Human AZGP1(Zinc-alpha-2-glycoprotein) ELISA Kit
845	abk9981H	Human MSLN(Mesothelin) ELISA Kit
846	abk0091H	Human SOST(Sclerostin) ELISA Kit
847	abk1091H	Human IRF5(Interferon Regulatory Factor 5) ELISA Kit
848	abk2091H	Human SIRT1(Sirtuin 1) ELISA Kit
849	abk3091H	Human VWFCP/ADAMTS13(Von Willebrand Factor Cleaving Protease) ELISA Kit
850	abk4091H	Human MTR(5-Methyltetrahydrofolate Homocysteine Methyltransferase) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
851	abk5091H	Human PGF(Placental Growth Factor) ELISA Kit
852	abk6091H	Human FOXC2(Forkhead Box Protein C2) ELISA Kit
853	abk7091H	Human ESM1(Endothelial Cell Specific Molecule 1) ELISA Kit
854	abk8091H	Human Rac1(Ras Related C3 Botulinum Toxin Substrate 1) ELISA Kit
855	abk9091H	Human LILRB2(Leukocyte Immunoglobulin Like Receptor Subfamily B Member 2) ELISA Kit
856	abk0191H	Human OTR(Oxytocin Receptor) ELISA Kit
857	abk1191H	Human PDIA6(Protein Disulfide Isomerase A6) ELISA Kit
858	abk2191H	Human HDL(High Density Lipoprotein) ELISA Kit
859	abk3191H	Human Pro-MMP-10(Pro-Matrix Metalloproteinase 10) ELISA Kit
860	abk4191H	Human Pro-MMP-13(Pro-Matrix Metalloproteinase 13) ELISA Kit
861	abk5191H	Human SIRPA(Signal Regulatory Protein Alpha) ELISA Kit
862	abk6191H	Human PDGF-AB(Platelet Derived Growth Factor AB) ELISA Kit
863	abk7191H	Human PGRN(Progranulin) ELISA Kit
864	abk8191H	Human RLN2(Relaxin-2) ELISA Kit
865	abk9191H	Human RBP1(Retinol Binding Protein 1, Cellular) ELISA Kit
866	abk0291H	Human Surv(Survivin) ELISA Kit
867	abk1291H	Human RBP3(Retinol Binding Protein 3, Interstitial) ELISA Kit
868	abk6453H	Human TPO(Thrombopoietin) ELISA Kit
869	abk3291H	Human TSP-1(Thrombospondin-1) ELISA Kit
870	abk4291H	Human TSP-2(Thrombospondin-2) ELISA Kit
871	abk5291H	Human IL7R(Interleukin 7 Receptor) ELISA Kit
872	abk6291H	Human TRAIL/TNFSF10(Tumor Necrosis Factor Related Apoptosis Inducing Ligand) ELISA Kit
873	abk7291H	Human PLD1(Phospholipase D1) ELISA Kit
874	abk8291H	Human TREM-1(Triggering Receptor Expressed on Myeloid Cells-1) ELISA Kit
875	abk9291H	Human TSLP(Thymic Stromal Lymphopoietin) ELISA Kit
876	abk0391H	Human PLAUR/uPAR(Plasminogen Activator, Urokinase Receptor) ELISA Kit
877	abk1391H	Human VEGF-C(Vascular Endothelial Growth Factor C) ELISA Kit
878	abk2391H	Human VEGF-D(Vascular Endothelial Growth Factor D) ELISA Kit
879	abk3391H	Human VEGFR-2/KDR(Vascular Endothelial Growth Factor Receptor 2) ELISA Kit
880	abk4391H	Human CD34(Cluster Of Differentiation 34) ELISA Kit
881	abk5391H	Human CYB5A(Cytochrome b5, type A) ELISA Kit
882	abk6391H	Human sOX40L(soluble OX40L) ELISA Kit
883	abk7391H	Human BRAF(B-Raf Proto Oncogene Serine/Threonine Protein Kinase) ELISA Kit
884	abk8391H	Human GzmA(Granzyme A) ELISA Kit
885	abk9391H	Human GzmB(Granzyme B) ELISA Kit
886	abk0491H	Human GNLY(Granulysin) ELISA Kit
887	abk1491H	Human S100P(S100 Calcium Binding Protein P) ELISA Kit
888	abk2491H	Human CR2/CD21(Complement Receptor 2) ELISA Kit
889	abk3491H	Human cADPRH/CD38(Cyclic ADP Ribose Hydrolase) ELISA Kit
890	abk4491H	Human CLDN5(Claudin 5) ELISA Kit
891	abk5491H	Human GRIN2B(Glutamate Receptor, Ionotropic, N-Methyl-D-Aspartate 2B) ELISA Kit
892	abk6491H	Human GRIN1(Glutamate Receptor, Ionotropic, N-Methyl-D-Aspartate 1) ELISA Kit
893	abk7491H	Human GRIN2A(Glutamate Receptor, Ionotropic, N-Methyl-D-Aspartate 2A) ELISA Kit
894	abk8491H	Human FXN(Frataxin) ELISA Kit
895	abk9491H	Human LOX-1(Lectin Like Oxidized Low Density Lipoprotein Receptor 1) ELISA Kit
896	abk0591H	Human ACTN4(Actinin Alpha 4) ELISA Kit
897	abk1591H	Human sAC(Soluble Adenylate Cyclase) ELISA Kit
898	abk2591H	Human sFMC(Soluble Fibrin Monomer Complex) ELISA Kit
899	abk3591H	Human PECAM1/CD31(Platelet/Endothelial Cell Adhesion Molecule 1) ELISA Kit
900	abk4591H	Human COL18α1(Collagen Type XVIII Alpha 1) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
901	abk5591H	Human MITF(Microphthalmia Associated Transcription Factor) ELISA Kit
902	abk6591H	Human BMP-5(Bone Morphogenetic Protein 5) ELISA Kit
903	abk7591H	Human IL-12 p35(Interleukin 12 p35) ELISA Kit
904	abk8591H	Human IL-25(Interleukin 25) ELISA Kit
905	abk9591H	Human IL-34(Interleukin 34) ELISA Kit
906	abk0691H	Human RLN1(Relaxin-1) ELISA Kit
907	abk1691H	Human RLN3(Relaxin-3) ELISA Kit
908	abk2691H	Human Tyk-2(Tyrosine Kinase 2) ELISA Kit
909	abk3691H	Human TH(Tyrosine Hydroxylase) ELISA Kit
910	abk4691H	Human mTOR(Mammalian Target of Rapamycin) ELISA Kit
911	abk5691H	Human PSPN(Persephin) ELISA Kit
912	abk6691H	Human TGFβ1(Transforming Growth Factor Beta Induced Protein) ELISA Kit
913	abk7691H	Human TSC22(Transforming Growth Factor Beta Stimulated Protein Clone 22) ELISA Kit
914	abk8691H	Human HLA-A(Leukocyte Antigen A) ELISA Kit
915	abk9691H	Human PKCη(Protein Kinase C Eta) ELISA Kit
916	abk0791H	Human UCP-1(Uncoupling Protein 1, Mitochondrial) ELISA Kit
917	abk1791H	Human IL18BP(Interleukin 18 Binding Protein) ELISA Kit
918	abk2791H	Human TYR(Tyrosinase) ELISA Kit
919	abk3791H	Human NFκB2(Nuclear Factor Kappa B2) ELISA Kit
920	abk4791H	Human FGL1(Fibrinogen Like Protein 1) ELISA Kit
921	abk5791H	Human FBLN2(Fibulin 2) ELISA Kit
922	abk6791H	Human FBLN3(Fibulin 3) ELISA Kit
923	abk8791H	Human PRG4(Proteoglycan 4) ELISA Kit
924	abk9791H	Human FBLN1(Fibulin 1) ELISA Kit
925	abk0891H	Human CX40(Connexin 40) ELISA Kit
926	abk1891H	Human CX43(Connexin 43) ELISA Kit
927	abk2891H	Human SDC4(Syndecan 4) ELISA Kit
928	abk3891H	Human LAP(Leucine Aminopeptidase) ELISA Kit
929	abk4891H	Human LTβ(Lymphotoxin Beta) ELISA Kit
930	abk5891H	Human LFA-1(Lymphocyte Function Associated Antigen 1) ELISA Kit
931	abk6891H	Human LFA-2/CD2(Lymphocyte Function Associated Antigen 2) ELISA Kit
932	abk7891H	Human Lptn/XCL1(Lymphotoactin) ELISA Kit
933	abk8891H	Human PGAM2(Phosphoglycerate Mutase 2) ELISA Kit
934	abk9891H	Human PFKM(Phosphofructokinase, Muscle) ELISA Kit
935	abk0991H	Human plkBα(Phospho-Inhibitory Subunit of NF Kappa B Alpha) ELISA Kit
936	abk1991H	Human PGM(Phosphoglucomutase) ELISA Kit
937	abk2991H	Human PCK1(Phosphoenolpyruvate Carboxykinase 1) ELISA Kit
938	abk9353H	Human GLO1(Glyoxalase I) ELISA Kit
939	abk4991H	Human PLSCR1(Phospholipid Scramblase 1) ELISA Kit
940	abk5991H	Human PLSCR2(Phospholipid Scramblase 2) ELISA Kit
941	abk6991H	Human PI3K2α(Phosphoinositide-3-Kinase Class-2-Alpha Polypeptide) ELISA Kit
942	abk7991H	Human GPC1(Glypican 1) ELISA Kit
943	abk8991H	Human GPC2(Glypican 2) ELISA Kit
944	abk9991H	Human GPC3(Glypican 3) ELISA Kit
945	abk2002H	Human PLTP(Phospholipid Transfer Protein) ELISA Kit
946	abk3002H	Human PLCγ1(Phospholipase C Gamma 1) ELISA Kit
947	abk4002H	Human FUR(Furin) ELISA Kit
948	abk6002H	Human SCCA2(Squamous Cell Carcinoma Antigen 2) ELISA Kit
949	abk7002H	Human pPLA2(Phospholipase A2, Pancreas) ELISA Kit
950	abk8002H	Human DS(Dermatan Sulfate) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
951	abk9002H	Human Trx(Thioredoxin) ELISA Kit
952	abk0102H	Human TrxR1(Thioredoxin Reductase 1) ELISA Kit
953	abk1102H	Human TBP2(Thioredoxin Binding Protein 2) ELISA Kit
954	abk0453H	Human IAP(immunosuppressive acidic protein) ELISA Kit
955	abk3102H	Human MCL1(Myeloid Cell Leukemia Sequence 1) ELISA Kit
956	abk5102H	Human FcαR/CD89(Receptor I for the Fc Fragment of IgA) ELISA Kit
957	abk6102H	Human FcεR I (Receptor I for the Fc Fragment of IgE) ELISA Kit
958	abk7102H	Human FcγR II/CD32(Receptor II for the Fc Fragment of IgG) ELISA Kit
959	abk8102H	Human FcγR3B(Fc Fragment of IgG Low Affinity IIIb Receptor) ELISA Kit
960	abk9102H	Human Fcγ(Fc Fragment of IgG) ELISA Kit
961	abk0202H	Human MCP/CD46(Membrane Cofactor Protein) ELISA Kit
962	abk1202H	Human TSP-4(Thrombospondin-4) ELISA Kit
963	abk2202H	Human TDT(Terminal Deoxynucleotidyl Transferase) ELISA Kit
964	abk3202H	Human SMN1(Survival of Motor Neuron 1, Telomeric) ELISA Kit
965	abk4202H	Human ETRA(Endothelin Receptor A) ELISA Kit
966	abk5202H	Human ECE1(Endothelin Converting Enzyme 1) ELISA Kit
967	abk6202H	Human REPIN1(Replication Initiator 1) ELISA Kit
968	abk7202H	Human EL/LIPG(Endothelial Lipase) ELISA Kit
969	abk8202H	Human VASPIN(Visceral Adipose Specific Serine Protease Inhibitor) ELISA Kit
970	abk9202H	Human VF(Visfatin) ELISA Kit
971	abk0302H	Human ERN1(Endoplasmic Reticulum to Nucleus Signalling 1) ELISA Kit
972	abk1302H	Human CFTR(Cystic Fibrosis Transmembrane Conductance Regulator) ELISA Kit
973	abk2302H	Human NGB(Neuroglobin) ELISA Kit
974	abk3302H	Human MAdCAM1(Mucosal Addressin Cell Adhesion Molecule 1) ELISA Kit
975	abk4302H	Human FAK(Focal Adhesion Kinase) ELISA Kit
976	abk5302H	Human OTC(Ornithine Transcarbamylase) ELISA Kit
977	abk6302H	Human ODC(Ornithine Decarboxylase) ELISA Kit
978	abk7302H	Human GBP4(Guanylate Binding Protein 4) ELISA Kit
979	abk8302H	Human UGCG(UDP-glucose ceramide glucosyltransferase) ELISA Kit
980	abk9302H	Human PLAU/uPA(Urokinase-Type Plasminogen Activator) ELISA Kit
981	abk0402H	Human UPP1(Uridine Phosphorylase 1) ELISA Kit
982	abk1402H	Human ETRB(Endothelin Receptor B) ELISA Kit
983	abk3453H	Human GS(Gelsolin) ELISA Kit
984	abk4801H	Human TM(Thrombin) ELISA Kit
985	abk4402H	Human TAT(Thrombin/Antithrombin Complex) ELISA Kit
986	abk5402H	Human PT(Prothrombin) ELISA Kit
987	abk6402H	Human F1+2(Prothrombin Fragment 1+2) ELISA Kit
988	abk7402H	Human KERA(Keratocan) ELISA Kit
989	abk8402H	Human VCL(Vinculin) ELISA Kit
990	abk9402H	Human EED(Embryonic Ectoderm Development) ELISA Kit
991	abk0502H	Human DPT(Dermatopontin) ELISA Kit
992	abk1502H	Human SYK(Spleen Tyrosine Kinase) ELISA Kit
993	abk2502H	Human PDPN(Podoplanin) ELISA Kit
994	abk3502H	Human OSF(Osteoclast Stimulating Factor 1) ELISA Kit
995	abk4502H	Human OSCAR(Osteoclast Associated Receptor) ELISA Kit
996	abk5502H	Human G6PD(Glucose 6 Phosphate Dehydrogenase) ELISA Kit
997	abk6502H	Human GPI(Glucose 6 Phosphate Isomerase) ELISA Kit
998	abk7502H	Human GCK(Glucokinase) ELISA Kit
999	abk8502H	Human GKRP(Glucokinase Regulatory Protein) ELISA Kit
1000	abk9502H	Human GLUT1(Glucose Transporter 1) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
1001	abk0602H	Human GLUT2(Glucose Transporter 2) ELISA Kit
1002	abk1602H	Human GLUT3(Glucose Transporter 3) ELISA Kit
1003	abk2602H	Human GLUT4(Glucose Transporter 4) ELISA Kit
1004	abk3602H	Human GTF2B(General Transcription Factor II B) ELISA Kit
1005	abk4602H	Human GTF2F1(General Transcription Factor II F, Polypeptide 1) ELISA Kit
1006	abk5602H	Human PK2(Prokineticin 2) ELISA Kit
1007	abk6602H	Human PKR1(Prokineticin Receptor 1) ELISA Kit
1008	abk7602H	Human PKR2(Prokineticin Receptor 2) ELISA Kit
1009	abk8602H	Human EP3(Prostaglandin E Receptor 3) ELISA Kit
1010	abk0702H	Human FP(Prostaglandin F Receptor) ELISA Kit
1011	abk1702H	Human PTGIS(Prostaglandin I Synthase) ELISA Kit
1012	abk2702H	Human PTGES2(Prostaglandin E Synthase 2) ELISA Kit
1013	abk3702H	Human PTGS1(Prostaglandin Endoperoxide Synthase 1) ELISA Kit
1014	abk4702H	Human PTGS2/COX-2(Prostaglandin Endoperoxide Synthase 2) ELISA Kit
1015	abk6702H	Human DSP(Desmoplakin) ELISA Kit
1016	abk7702H	Human ERCC1(Excision Repair Cross Complementing Rodent Repair Deficiency Complementation 1) ELISA Kit
1017	abk8702H	Human KLK1(Kallikrein 1) ELISA Kit
1018	abk9702H	Human LXRA(Liver X Receptor Alpha) ELISA Kit
1019	abk0802H	Human HSPβ6/HSP-20(Heat Shock Protein Beta 6) ELISA Kit
1020	abk1802H	Human HSP-27/HSPB1(Heat Shock Protein 27) ELISA Kit
1021	abk2802H	Human HSP-40(Heat Shock Protein 40) ELISA Kit
1022	abk3802H	Human HSP-60/HSPD1(Heat Shock Protein 60) ELISA Kit
1023	abk4802H	Human HSP-90(Heat Shock Protein 90) ELISA Kit
1024	abk5802H	Human HSP gp96(Heat Shock Protein glycoprotein 96) ELISA Kit
1025	abk6802H	Human HSF2(Heat Shock Transcription Factor 2) ELISA Kit
1026	abk7802H	Human PSG1/SP1(Pregnancy Specific Beta-1-Glycoprotein 1) ELISA Kit
1027	abk9802H	Human TERF2(Telomeric Repeat Binding Factor 2) ELISA Kit
1028	abk0902H	Human TFF1(Trefoil Factor 1) ELISA Kit
1029	abk1902H	Human TFF2(Trefoil Factor 2) ELISA Kit
1030	abk2902H	Human PEDF(Pigment Epithelium Derived Factor) ELISA Kit
1031	abk3902H	Human SDH(Sorbitol Dehydrogenase) ELISA Kit
1032	abk4902H	Human CCR2(Chemokine C-C-Motif Receptor 2) ELISA Kit
1033	abk5902H	Human Ep-CAM/CD326(Epithelial Cell Adhesion Molecule) ELISA Kit
1034	abk8902H	Human GMFB(Glia Maturation Factor, Beta) ELISA Kit
1035	abk0012H	Human HK1(Hexokinase 1) ELISA Kit
1036	abk1012H	Human P2(Myelin Protein 2) ELISA Kit
1037	abk3012H	Human NCAM/CD56(Neural Cell Adhesion Molecule) ELISA Kit
1038	abk4012H	Human L1CAM/CD171(L1-Cell Adhesion Molecule) ELISA Kit
1039	abk5012H	Human NTRK2(Neurotrophic Tyrosine Kinase Receptor Type 2) ELISA Kit
1040	abk6012H	Human NTRK3(Neurotrophic Tyrosine Kinase Receptor Type 3) ELISA Kit
1041	abk7012H	Human NAIP(Neural Apoptosis Inhibitory Protein) ELISA Kit
1042	abk8012H	Human NPHN(Nephrin) ELISA Kit
1043	abk9012H	Human GTEAGE(Glomerular Tissue Extract Advanced Glycosylation End Products) ELISA Kit
1044	abk0112H	Human GROγ/CXCL3(Growth Regulated Oncogene Gamma) ELISA Kit
1045	abk1112H	Human GDF1(Growth Differentiation Factor 1) ELISA Kit
1046	abk2112H	Human GDF10(Growth Differentiation Factor 10) ELISA Kit
1047	abk3112H	Human GDF11(Growth Differentiation Factor 11) ELISA Kit
1048	abk4112H	Human GDF3(Growth Differentiation Factor 3) ELISA Kit
1049	abk5112H	Human GDF6(Growth Differentiation Factor 6) ELISA Kit
1050	abk6112H	Human GDF7(Growth Differentiation Factor 7) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
1051	abk7112H	Human GDF9(Growth Differentiation Factor 9) ELISA Kit
1052	abk8112H	Human HK2(Hexokinase 2) ELISA Kit
1053	abk9112H	Human GH2(Growth Hormone 2) ELISA Kit
1054	abk0212H	Human OSBPL1A(Oxysterol Binding Protein Like Protein 1A) ELISA Kit
1055	abk2212H	Human GHITM(Growth Hormone Inducible Transmembrane Protein) ELISA Kit
1056	abk3212H	Human GAP43(Growth Associated Protein 43) ELISA Kit
1057	abk5212H	Human SSTR1(Somatostatin Receptor 1) ELISA Kit
1058	abk6212H	Human Grb10(Growth Factor Receptor Bound Protein 10) ELISA Kit
1059	abk7212H	Human Grb14(Growth Factor Receptor Bound Protein 14) ELISA Kit
1060	abk8212H	Human Grb2(Growth Factor Receptor Bound Protein 2) ELISA Kit
1061	abk9212H	Human Grb7(Growth Factor Receptor Bound Protein 7) ELISA Kit
1062	abk1312H	Human RB1(Retinoblastoma Protein 1) ELISA Kit
1063	abk2312H	Human CYPB(Cyclophilin B) ELISA Kit
1064	abk3312H	Human CYPD(Cyclophilin D) ELISA Kit
1065	abk4312H	Human NRP2(Neuropilin 2) ELISA Kit
1066	abk5312H	Human CS(Citrate Synthase) ELISA Kit
1067	abk6312H	Human OGG1(8-Oxoguanine Glycosylase 1) ELISA Kit
1068	abk7312H	Human ECFR(Eosinophil Chemotactic Factor Receptor) ELISA Kit
1069	abk8312H	Human EAR1(Eosinophil Associated Ribonuclease A family Member 1) ELISA Kit
1070	abk9312H	Human EAR2(Eosinophil Associated Ribonuclease A family Member 2) ELISA Kit
1071	abk0412H	Human EAR3(Eosinophil Associated Ribonuclease A family Member 3) ELISA Kit
1072	abk1412H	Human NOG(Noggin) ELISA Kit
1073	abk2412H	Human PLN(Phospholamban) ELISA Kit
1074	abk3412H	Human RIPK2(Receptor TNFRSF Interacting Serine Threonine Kinase 2) ELISA Kit
1075	abk4412H	Human TNK2(Tyrosine Kinase, Non Receptor 2) ELISA Kit
1076	abk5412H	Human DAF(Decay Accelerating Factor) ELISA Kit
1077	abk6412H	Human TNNI3K(Serine/threonine-protein kinase TNNI3K) ELISA Kit
1078	abk7412H	Human MAP2K6/MKK6(Mitogen Activated Protein Kinase Kinase 6) ELISA Kit
1079	abk8412H	Human RXFP1(Relaxin/Insulin Like Family Peptide Receptor 1) ELISA Kit
1080	abk9412H	Human TACR1(Tachykinin Receptor 1) ELISA Kit
1081	abk0512H	Human TACR2(Tachykinin Receptor 2) ELISA Kit
1082	abk1512H	Human TACR3(Tachykinin Receptor 3) ELISA Kit
1083	abk2512H	Human PLP1(Proteolipid Protein 1, Myelin) ELISA Kit
1084	abk3512H	Human MyD88(Myeloid Differentiation Factor 88) ELISA Kit
1085	abk4512H	Human POT1(Protection of Telomeres Protein 1) ELISA Kit
1086	abk5512H	Human HBF(Fetal Hemoglobin) ELISA Kit
1087	abk6512H	Human PR3(Proteinase 3) ELISA Kit
1088	abk7512H	Human PL(Placental Lactogen) ELISA Kit
1089	abk8512H	Human PP13(Placental Protein 13) ELISA Kit
1090	abk9512H	Human PLAP/ALPP(Placental Alkaline Phosphatase) ELISA Kit
1091	abk0612H	Human CPT1A(Carnitine Palmitoyltransferase 1A, Liver) ELISA Kit
1092	abk1612H	Human CPT1B(Carnitine Palmitoyltransferase 1B, Muscle) ELISA Kit
1093	abk2612H	Human OSBPL8(Oxysterol Binding Protein Like Protein 8) ELISA Kit
1094	abk3612H	Human CFHR3(Complement Factor H Related Protein 3) ELISA Kit
1095	abk4612H	Human PGRP-S(Peptidoglycan Recognition Protein-S) ELISA Kit
1096	abk5612H	Human PPIL1(Peptidylprolyl Isomerase/Cyclophilin Like 1) ELISA Kit
1097	abk6612H	Human RAP1A(Ras-related protein Rap-1A) ELISA Kit
1098	abk7612H	Human GP5(Glycoprotein V, Platelet) ELISA Kit
1099	abk8612H	Human GP6(Glycoprotein VI, Platelet) ELISA Kit
1100	abk9612H	Human GlyCAM1(Glycosylation Dependent Cell Adhesion Molecule 1) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
1101	abk0712H	Human GR α (Glucocorticoid Receptor Alpha) ELISA Kit
1102	abk1712H	Human GR β (Glucocorticoid Receptor Beta) ELISA Kit
1103	abk2712H	Human GSK3 β (Glycogen Synthase Kinase 3 Beta) ELISA Kit
1104	abk3712H	Human GSK3 α (Glycogen Synthase Kinase 3 Alpha) ELISA Kit
1105	abk4712H	Human A-GHRL(Acylated Ghrelin) ELISA Kit
1106	abk5712H	Human GLTP(Glycolipid Transfer Protein) ELISA Kit
1107	abk6712H	Human CFHR5(Complement Factor H Related Protein 5) ELISA Kit
1108	abk7712H	Human FTL(Ferritin, Light Polypeptide) ELISA Kit
1109	abk8712H	Human FTH(Ferritin, Heavy Polypeptide) ELISA Kit
1110	abk9712H	Human PIR(Pirin) ELISA Kit
1111	abk0812H	Human MRP1(Multidrug Resistance Associated Protein 1) ELISA Kit
1112	abk1812H	Human SYNPO(Synaptopodin) ELISA Kit
1113	abk2812H	Human SYP(Synaptophysin) ELISA Kit
1114	abk3812H	Human SAP25(Synaptosome Associated Protein 25) ELISA Kit
1115	abk7812H	Human DNase- I (Deoxyribonuclease I) ELISA Kit
1116	abk8812H	Human TERF1(Telomeric Repeat Binding Factor 1) ELISA Kit
1117	abk0912H	Human SIGLEC12(Sialic Acid Binding Ig Like Lectin 12) ELISA Kit
1118	abk1912H	Human SIGLEC3(Sialic Acid Binding Ig Like Lectin 3) ELISA Kit
1119	abk2912H	Human ITLN1(Intelectin 1/Omentin) ELISA Kit
1120	abk3912H	Human TUB α 3C(Tubulin Alpha 3c) ELISA Kit
1121	abk4912H	Human CFHR4(Complement Factor H Related Protein 4) ELISA Kit
1122	abk5912H	Human STMN1(Stathmin 1) ELISA Kit
1123	abk6912H	Human DCX(Doublecortin) ELISA Kit
1124	abk7912H	Human VDR(Vitamin D Receptor) ELISA Kit
1125	abk8912H	Human CFHR2(Complement Factor H Related Protein 2) ELISA Kit
1126	abk0022H	Human CFHR1(Complement Factor H Related Protein 1) ELISA Kit
1127	abk1022H	Human PP(Pepsin) ELISA Kit
1128	abk4022H	Human GKN1(Gastrokine 1) ELISA Kit
1129	abk5022H	Human GKN2(Gastrokine 2) ELISA Kit
1130	abk7022H	Human GT-34(Gastrin 34) ELISA Kit
1131	abk8022H	Human ProGRP(Pro-Gastrin Releasing Peptide) ELISA Kit
1132	abk0122H	Human GIP(Gastric Inhibitory Polypeptide) ELISA Kit
1133	abk1122H	Human MFGE8(Milk Fat Globule EGF Factor 8) ELISA Kit
1134	abk2122H	Human MMP23B(Matrix Metalloproteinase 23B) ELISA Kit
1135	abk3122H	Human CK-13/KRT13(Cytokeratin 13) ELISA Kit
1136	abk4122H	Human CK-19/KRT19(Cytokeratin 19) ELISA Kit
1137	abk5122H	Human CK-20/KRT20(Cytokeratin 20) ELISA Kit
1138	abk6122H	Human CYFRA21-1(Cytokeratin Fragment Antigen 21-1) ELISA Kit
1139	abk7122H	Human Bcl11A(B-Cell Lymphoma 11A) ELISA Kit
1140	abk8122H	Human CVL/EZR(Cytovillin/Ezrin) ELISA Kit
1141	abk9122H	Human CYB561(Cytochrome b-561) ELISA Kit
1142	abk0222H	Human CYP21B(Cytochrome P450 21B) ELISA Kit
1143	abk1222H	Human ECM1(Extracellular Matrix Protein 1) ELISA Kit
1144	abk2222H	Human MEPE(Matrix Extracellular Phosphoglycoprotein) ELISA Kit
1145	abk3222H	Human ERK1(Extracellular Signal Regulated Kinase 1) ELISA Kit
1146	abk4222H	Human ERK2(Extracellular Signal Regulated Kinase 2) ELISA Kit
1147	abk5222H	Human ERK3(Extracellular Signal Regulated Kinase 3) ELISA Kit
1148	abk6222H	Human SOCS3(Suppressors of Cytokine Signaling 3) ELISA Kit
1149	abk7222H	Human CYTH1(Cytohesin 1) ELISA Kit
1150	abk8222H	Human CYTIP(Cytohesin Interacting Protein) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
1151	abk9222H	Human FLNB(Filamin B, Beta) ELISA Kit
1152	abk0322H	Human FLNC(Filamin C, Gamma) ELISA Kit
1153	abk1322H	Human FMOD(Fibromodulin) ELISA Kit
1154	abk2322H	Human FLRT2(Fibronectin Leucine Rich Transmembrane Protein 2) ELISA Kit
1155	abk3322H	Human PAP(Plasmin-Antiplasmin Complex) ELISA Kit
1156	abk4322H	Human Plg(Plasminogen) ELISA Kit
1157	abk5322H	Human PAI1(Plasminogen Activator Inhibitor 1) ELISA Kit
1158	abk6322H	Human PAI2(Plasminogen Activator Inhibitor 2) ELISA Kit
1159	abk7322H	Human tPA(Tissue-type Plasminogen Activator) ELISA Kit
1160	abk8322H	Human TAFI(Thrombin Activatable Fibrinolysis Inhibitor) ELISA Kit
1161	abk9322H	Human FBL(Fibrillarin) ELISA Kit
1162	abk0422H	Human GHR(Growth Hormone Receptor) ELISA Kit
1163	abk1422H	Human FG α (Fibrinogen Alpha) ELISA Kit
1164	abk2422H	Human FG β (Fibrinogen Beta) ELISA Kit
1165	abk3422H	Human FGL2(Fibrinogen Like Protein 2) ELISA Kit
1166	abk5422H	Human SNRPC(Small Nuclear Ribonucleoprotein Polypeptide C) ELISA Kit
1167	abk6422H	Human SBEM(Small Breast Epithelial Mucin) ELISA Kit
1168	abk7422H	Human TINF2(TERF1 Interacting Nuclear Factor 2) ELISA Kit
1169	abk9422H	Human SRP9(Signal Recognition Particle 9kDa) ELISA Kit
1170	abk0522H	Human NPTX1(Neuronal Pentraxin 1) ELISA Kit
1171	abk1522H	Human STAT1(Signal Transducer and Activator of Transcription 1) ELISA Kit
1172	abk2522H	Human STAT2(Signal Transducer and Activator of Transcription 2) ELISA Kit
1173	abk3522H	Human STAT4(Signal Transducer and Activator of Transcription 4) ELISA Kit
1174	abk4522H	Human STAT6(Signal Transducer and Activator of Transcription 6) ELISA Kit
1175	abk5522H	Human SOX1(Sex Determining Region Y Box Protein 1) ELISA Kit
1176	abk6522H	Human SOX18(Sex Determining Region Y Box Protein 18) ELISA Kit
1177	abk7522H	Human SOX3(Sex Determining Region Y Box Protein 3) ELISA Kit
1178	abk8522H	Human TS/TYMS(Thymidylate Synthetase) ELISA Kit
1179	abk9522H	Human TLA(Thymus Leukemia Antigen) ELISA Kit
1180	abk0622H	Human TP(Thymidine Phosphorylase) ELISA Kit
1181	abk3622H	Human NPTX2(Neuronal Pentraxin 2) ELISA Kit
1182	abk5622H	Human VEGF165(Vascular Endothelial Growth Factor 165) ELISA Kit
1183	abk6622H	Human GST μ 2/GSTm2(Glutathione S Transferase Mu 2, Muscle) ELISA Kit
1184	abk7622H	Human SUMO2(Small Ubiquitin Related Modifier 2) ELISA Kit
1185	abk8622H	Human IDO(Indoleamine-2,3-Dioxygenase) ELISA Kit
1186	abk9622H	Human OGN(Osteoglycin) ELISA Kit
1187	abk0722H	Human VEGF-B(Vascular Endothelial Cell Growth Factor B) ELISA Kit
1188	abk1722H	Human VEGFR3/Flt4(Vascular Endothelial Cell Growth Factor Receptor 3) ELISA Kit
1189	abk3722H	Human VWF(Von Willebrand Factor) ELISA Kit
1190	abk4722H	Human HbC(Hemoglobin C) ELISA Kit
1191	abk5722H	Human HB β (Hemoglobin Beta) ELISA Kit
1192	abk6722H	Human HO1(Heme Oxygenase 1) ELISA Kit
1193	abk7722H	Human HO2(Heme Oxygenase 2) ELISA Kit
1194	abk8722H	Human SGK2(Serum/Glucocorticoid Regulated Kinase 2) ELISA Kit
1195	abk9722H	Human KRAS(V-Ki-Ras2 Kirsten Rat Sarcoma Viral Oncogene Homolog) ELISA Kit
1196	abk0822H	Human SEMA4C(Semaphorin 4C) ELISA Kit
1197	abk1822H	Human SERT(Serotonin Transporter) ELISA Kit
1198	abk2822H	Human TPP(Thrombus Precursor Protein) ELISA Kit
1199	abk4822H	Human STC2(Stanniocalcin 2) ELISA Kit
1200	abk5822H	Human FDP(Fibrinogen Degradation Product) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
1201	abk6822H	Human FPA(Fibrinopeptide A) ELISA Kit
1202	abk7822H	Human FPB(Fibrinopeptide B) ELISA Kit
1203	abk8822H	Human SSTR2(Somatostatin Receptor 2) ELISA Kit
1204	abk9822H	Human PHLDA2(Pleckstrin Homology Like Domain Family A, Member 2) ELISA Kit
1205	abk0922H	Human PAFAHIB3(Platelet Activating Factor Acetylhydrolase Ib3) ELISA Kit
1206	abk1922H	Human HAase(Hyaluronidase) ELISA Kit
1207	abk2922H	Human GP II b/IIIa(Platelet Membrane Glycoprotein II b/IIIa complex) ELISA Kit
1208	abk3922H	Human ITGB3(Integrin Beta 3) ELISA Kit
1209	abk4922H	Human PKP1(Plakophilin 1) ELISA Kit
1210	abk5922H	Human PKP2(Plakophilin 2) ELISA Kit
1211	abk6922H	Human ANGPTL8(Angiopoietin Like Protein 8) ELISA Kit
1212	abk7922H	Human PFKP(Phosphofructokinase, Platelet) ELISA Kit
1213	abk8922H	Human PDGFRα(Platelet Derived Growth Factor Receptor Alpha) ELISA Kit
1214	abk9922H	Human PF4V1(Platelet Factor 4 Variant 1) ELISA Kit
1215	abk0032H	Human GYPE(Glycophorin E) ELISA Kit
1216	abk1032H	Human SSTR4(Somatostatin Receptor 4) ELISA Kit
1217	abk2032H	Human NOX1(Nicotinamide Adenine Dinucleotide Phosphate Oxidase 1) ELISA Kit
1218	abk3032H	Human EloA(Elongin A) ELISA Kit
1219	abk5032H	Human SSTR3(Somatostatin Receptor 3) ELISA Kit
1220	abk6032H	Human PI(Proinsulin) ELISA Kit
1221	abk7032H	Human AMY2(Amylase Alpha 2, Pancreatic) ELISA Kit
1222	abk2453H	Human PP(Pancreatic Polypeptide) ELISA Kit
1223	abk1132H	Human JAK2(Janus Kinase 2) ELISA Kit
1224	abk2132H	Human PTF1α(Pancreas Specific Transcription Factor 1α) ELISA Kit
1225	abk3132H	Human REG1α(Regenerating Islet Derived Protein 1 Alpha) ELISA Kit
1226	abk4132H	Human REG4(Regenerating Islet Derived Protein 4) ELISA Kit
1227	abk5132H	Human SULT1A1(Sulfotransferase 1A1) ELISA Kit
1228	abk1453H	Human PL(Pancreatic Lipase) ELISA Kit
1229	abk7132H	Human PTPN11(Protein Tyrosine Phosphatase, Non Receptor Type11) ELISA Kit
1230	abk8132H	Human OSM(Oncostatin M) ELISA Kit
1231	abk9132H	Human DCN(Decorin) ELISA Kit
1232	abk0232H	Human PHB(Prohibitin) ELISA Kit
1233	abk1232H	Human CRY1(Cryptochrome 1) ELISA Kit
1234	abk2232H	Human SEMA4D(Semaphorin 4D) ELISA Kit
1235	abk3232H	Human SSTR5(Somatostatin Receptor 5) ELISA Kit
1236	abk4232H	Human KNG1(Kininogen 1) ELISA Kit
1237	abk5232H	Human S100A2(S100 Calcium Binding Protein A2) ELISA Kit
1238	abk6232H	Human VAP-1(Vascular Adhesion Protein 1) ELISA Kit
1239	abk7232H	Human PTPN6(Protein Tyrosine Phosphatase, Non Receptor Type 6) ELISA Kit
1240	abk8232H	Human MGMT(O-6-Methylguanine DNA Methyltransferase) ELISA Kit
1241	abk9232H	Human PCDH1(Protocadherin 1) ELISA Kit
1242	abk0332H	Human PCDHβ2(Protocadherin Beta 2) ELISA Kit
1243	abk1332H	Human FBN1(Fibrillin 1) ELISA Kit
1244	abk2332H	Human FBN2(Fibrillin 2) ELISA Kit
1245	abk3332H	Human FBN3(Fibrillin 3) ELISA Kit
1246	abk4332H	Human PAQR7(Progestin and AdipoQ Receptor Family Member VII) ELISA Kit
1247	abk5332H	Human PGR(Progesterone Receptor) ELISA Kit
1248	abk6332H	Human sCD44(Soluble Cluster of Differentiation 44) ELISA Kit
1249	abk7332H	Human PS2(Presenilin 2) ELISA Kit
1250	abk8332H	Human EGR1(Early Growth Response Protein 1) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
1251	abk9332H	Human EGR2(Early Growth Response Protein 2) ELISA Kit
1252	abk0432H	Human EGR3(Early Growth Response Protein 3) ELISA Kit
1253	abk1432H	Human EGR4(Early Growth Response Protein 4) ELISA Kit
1254	abk2432H	Human NFYC(Nuclear Transcription Factor Y Gamma) ELISA Kit
1255	abk3432H	Human MUC5AC(Mucin 5 Subtype AC) ELISA Kit
1256	abk4432H	Human MUC5B(Mucin 5 Subtype B) ELISA Kit
1257	abk5432H	Human MUC7(Mucin 7) ELISA Kit
1258	abk6432H	Human eIF3a(Eukaryotic Translation Initiation Factor 3a) ELISA Kit
1259	abk7432H	Human PC4(Positive Cofactor 4) ELISA Kit
1260	abk8432H	Human IL4I1(Interleukin 4 Induced Protein 1) ELISA Kit
1261	abk9432H	Human GYS1(Glycogen Synthase 1, Muscle) ELISA Kit
1262	abk0532H	Human GYS2(Glycogen Synthase 2, Liver) ELISA Kit
1263	abk1532H	Human OSMR(Oncostatin M Receptor) ELISA Kit
1264	abk2532H	Human MK(Midkine) ELISA Kit
1265	abk3532H	Human PON1(Paraoxonase 1) ELISA Kit
1266	abk4532H	Human HNP1-3(Neutrophil Peptide 1-3) ELISA Kit
1267	abk5532H	Human PFN1(Profilin 1) ELISA Kit
1268	abk6532H	Human ADAM22(A Disintegrin And Metalloprotease 22) ELISA Kit
1269	abk7532H	Human TP63(Tumor Protein p63) ELISA Kit
1270	abk8532H	Human TNF-β(Tumor Necrosis Factor Beta) ELISA Kit
1271	abk9532H	Human TL1A(Tumor Necrosis Factor Related Ligand 1A) ELISA Kit
1272	abk0632H	Human CCT8(T-complex protein 1 subunit theta) ELISA Kit
1273	abk1632H	Human CCT8L1(T-complex 1 protein subunit theta-like 1) ELISA Kit
1274	abk2632H	Human CCT5(T-complex 1 protein subunit epsilon) ELISA Kit
1275	abk3632H	Human CCT6A(T-complex 1 protein subunit zeta) ELISA Kit
1276	abk4632H	Human TSTA(Tumor Specific Transplantation Antigen) ELISA Kit
1277	abk5632H	Human RAG-1(Recombination Activating Gene 1) ELISA Kit
1278	abk6632H	Human RAG-2(Recombination Activating Gene 2) ELISA Kit
1279	abk7632H	Human CDK1(Cyclin Dependent Kinase 1) ELISA Kit
1280	abk8632H	Human CDK2(Cyclin Dependent Kinase 2) ELISA Kit
1281	abk9632H	Human CDK4(Cyclin Dependent Kinase 4) ELISA Kit
1282	abk0732H	Human CDK5(Cyclin Dependent Kinase 5) ELISA Kit
1283	abk1732H	Human CDK7(Cyclin Dependent Kinase 7) ELISA Kit
1284	abk2732H	Human CDK9(Cyclin Dependent Kinase 9) ELISA Kit
1285	abk3732H	Human CDKN2A(Cyclin Dependent Kinase Inhibitor 2A) ELISA Kit
1286	abk4732H	Human Ntn1(Netrin 1) ELISA Kit
1287	abk5732H	Human Ntn4(Netrin 4) ELISA Kit
1288	abk6732H	Human MHCC/HLA-C(Major Histocompatibility Complex Class I C) ELISA Kit
1289	abk7732H	Human MHCE/HLA-E(Major Histocompatibility Complex Class I E) ELISA Kit
1290	abk8732H	Human AMACR(Alpha-methylacyl-CoA racemase) ELISA Kit
1291	abk9732H	Human CCNA1(Cyclin-A1) ELISA Kit
1292	abk0832H	Human RELN(Reelin) ELISA Kit
1293	abk1832H	Human TGM2(Transglutaminase 2, Tissue) ELISA Kit
1294	abk2832H	Human TGM3(Transglutaminase 3, Epidermal) ELISA Kit
1295	abk3832H	Human IL-27(Interleukin 27) ELISA Kit
1296	abk4832H	Human EDN(Eosinophil-Derived Neurotoxin) ELISA Kit
1297	abk5832H	Human TTR(Transthyretin) ELISA Kit
1298	abk6832H	Human TCF20(Transcription Factor 20) ELISA Kit
1299	abk7832H	Human GAPDH(Glyceraldehyde-3-Phosphate Dehydrogenase) ELISA Kit
1300	abk8832H	Human TFR2(Transferrin Receptor 2) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
1301	abk9832H	Human TNPO1(Transportin 1) ELISA Kit
1302	abk0932H	Human Pax(Paxillin) ELISA Kit
1303	abk1932H	Human CD200R1(CD200 Receptor 1) ELISA Kit
1304	abk2932H	Human MAOB(Amine oxidase B) ELISA Kit
1305	abk3932H	Human FPN(Ferroportin) ELISA Kit
1306	abk4932H	Human HEPH(Hephaestin) ELISA Kit
1307	abk5932H	Human CALP(Calprotectin) ELISA Kit
1308	abk8932H	Human TPS(Tissue Polypeptide Specific Antigen) ELISA Kit
1309	abk9932H	Human SCARB1(Scavenger receptor class B member 1) ELISA Kit
1310	abk0042H	Human HPA(Heparanase) ELISA Kit
1311	abk1042H	Human TFPI2(Tissue Factor Pathway Inhibitor 2) ELISA Kit
1312	abk2042H	Human proPSA(precursor of Prostate Specific Antigen) ELISA Kit
1313	abk3042H	Human UBA52(Ubiquitin A 52 Residue Ribosomal Protein Fusion Product 1) ELISA Kit
1314	abk4042H	Human PON2(Paraoxonase 2) ELISA Kit
1315	abk5042H	Human NES1(Nesfatin 1) ELISA Kit
1316	abk6042H	Human PON3(Paraoxonase 3) ELISA Kit
1317	abk7042H	Human DDT(D-Dopachrome Tautomerase) ELISA Kit
1318	abk8042H	Human sC5b-9(Soluble Terminal Complement Complex C5b-9) ELISA Kit
1319	abk9042H	Human UCHL1(Ubiquitin Carboxyl Terminal Hydrolase L1) ELISA Kit
1320	abk0142H	Human UBR4(Ubiquitin Protein Ligase E3 Component N-Recognin 4) ELISA Kit
1321	abk1142H	Human UBE2A(Ubiquitin Conjugating Enzyme E2A) ELISA Kit
1322	abk3142H	Human BCAR1(Breast Cancer Anti Estrogen Resistance 1) ELISA Kit
1323	abk4142H	Human PLA2G2D(Phospholipase A2, Group II D) ELISA Kit
1324	abk5142H	Human UBAP2(Ubiquitin Associated Protein 2) ELISA Kit
1325	abk6142H	Human MASP1(Mannan Associated Serine Protease 1) ELISA Kit
1326	abk7142H	Human MASP2(Mannan Associated Serine Protease 2) ELISA Kit
1327	abk8142H	Human USP8(Ubiquitin Specific Peptidase 8) ELISA Kit
1328	abk9142H	Human DPP6(Dipeptidyl Peptidase VI) ELISA Kit
1329	abk0242H	Human MRC2(Mannose Receptor C Type 2) ELISA Kit
1330	abk1242H	Human MRC1(Mannose Receptor C Type 1) ELISA Kit
1331	abk2242H	Human GADL1(Glutamate Decarboxylase Like Protein 1) ELISA Kit
1332	abk3242H	Human GAD1(Glutamate Decarboxylase 1) ELISA Kit
1333	abk4242H	Human PPAR-δ(Peroxisome Proliferator Activated Receptor Delta) ELISA Kit
1334	abk5242H	Human NUP62(Nucleoporin 62kDa) ELISA Kit
1335	abk6242H	Human PCNA(Proliferating Cell Nuclear Antigen) ELISA Kit
1336	abk7242H	Human RANK(Receptor Activator of Nuclear Factor Kappa B) ELISA Kit
1337	abk9242H	Human MCAM(Melanoma Cell Adhesion Molecule) ELISA Kit
1338	abk0342H	Human MIA1(Melanoma Inhibitory Activity Protein 1) ELISA Kit
1339	abk1342H	Human MFI2(Melanotransferrin) ELISA Kit
1340	abk2342H	Human SPTβN4(Spectrin Beta, Non Erythrocytic 4) ELISA Kit
1341	abk3342H	Human IκBα(Inhibitory Subunit of NF Kappa B Alpha) ELISA Kit
1342	abk4342H	Human HNRPA2B1(Heterogeneous Nuclear Ribonucleoprotein A2/B1) ELISA Kit
1343	abk5342H	Human HNRPU(Heterogeneous Nuclear Ribonucleoprotein U) ELISA Kit
1344	abk6342H	Human RARα(Retinoic Acid Receptor Alpha) ELISA Kit
1345	abk7342H	Human RNASEN(Ribonuclease III, Nuclear) ELISA Kit
1346	abk8342H	Human MCSP(Melanoma Associated Chondroitin Sulfate Proteoglycan) ELISA Kit
1347	abk1442H	Human THRα(Thyroid Hormone Receptor Alpha) ELISA Kit
1348	abk2442H	Human NFRκB(Nuclear Factor Related To Kappa B Binding Protein) ELISA Kit
1349	abk3442H	Human TNKS1(Tankyrase 1) ELISA Kit
1350	abk4442H	Human TNKS2(Tankyrase 2) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
1351	abk5442H	Human PARP4(Poly ADP Ribose Polymerase 4) ELISA Kit
1352	abk6442H	Human ME20M(Melanoma Associated ME20) ELISA Kit
1353	abk7442H	Human PXR(Pregnane X Receptor) ELISA Kit
1354	abk8442H	Human ALOX15B(Arachidonate-15-Lipoxygenase, Type B) ELISA Kit
1355	abk9442H	Human RELB(V-Re1 Reticuloendotheliosis Viral Oncogene Homolog B) ELISA Kit
1356	abk0542H	Human IL-24(Interleukin 24) ELISA Kit
1357	abk1542H	Human AIM2(Absent In Melanoma 2) ELISA Kit
1358	abk2542H	Human TSPAN30/CD63(Tetraspanin 30) ELISA Kit
1359	abk3542H	Human BANP(BTG3 Associated Nuclear Protein) ELISA Kit
1360	abk4542H	Human PRKG1(Protein Kinase, cGMP Dependent Type I) ELISA Kit
1361	abk5542H	Human PDE4D(Phosphodiesterase 4D, cAMP Specific) ELISA Kit
1362	abk6542H	Human CKB(Creatine Kinase, Brain) ELISA Kit
1363	abk7542H	Human LOXL4(Lysyl Oxidase Like Protein 4) ELISA Kit
1364	abk8542H	Human TNFRSF1B(Tumor Necrosis Factor Receptor Superfamily, Member 1B) ELISA Kit
1365	abk9542H	Human LL-37(Antibacterial Protein LL-37) ELISA Kit
1366	abk0642H	Human Tα1(Thymosin-α1) ELISA Kit
1367	abk1642H	Human PACAP-38(Pituitary Adenylate Cyclase Activating Polypeptide 38) ELISA Kit
1368	abk2642H	Human PACAP-27(Pituitary Adenylate Cyclase Activating Polypeptide 27) ELISA Kit
1369	abk3642H	Human IL-35(Interleukin 35) ELISA Kit
1370	abk4642H	Human TK1(Thymidine Kinase 1, Soluble) ELISA Kit
1371	abk5642H	Human 12-LO(Arachidonate 12-Lipoxygenase) ELISA Kit
1372	abk6642H	Human VIL1(Villin 1) ELISA Kit
1373	abk7642H	Human PI3/Elafin(Peptidase Inhibitor 3, Skin Derived) ELISA Kit
1374	abk8642H	Human TTPα(Alpha Tocopherol Transfer Protein) ELISA Kit
1375	abk9642H	Human IL-21(Interleukin 21) ELISA Kit
1376	abk0742H	Human SSADH(Succinate-semialdehyde dehydrogenase, mitochondrial) ELISA Kit
1377	abk1742H	Human MMP-15(Matrix Metalloproteinase 15) ELISA Kit
1378	abk2742H	Human MMP-14(Matrix Metalloproteinase 14) ELISA Kit
1379	abk3742H	Human VLDL(Very Low Density Lipoprotein) ELISA Kit
1380	abk4742H	Human LCAT(Lecithin-Cholesterol Acyltransferase) ELISA Kit
1381	abk5742H	Human CYP2E1(Cytochrome P450, family 2, subfamily E, polypeptide 1) ELISA Kit
1382	abk6742H	Human MICB(Major Histocompatibility Complex Class I Chain-related gene B) ELISA Kit
1383	abk7742H	Human TNNC1(Troponin C Type 1) ELISA Kit
1384	abk8742H	Human GALT(Galactose-1-Phosphate Uridyltransferase) ELISA Kit
1385	abk9742H	Human FSHB(Follicotropin Subunit Beta) ELISA Kit
1386	abk0842H	Human LCE3B(Late Cornified Envelope Protein 3B) ELISA Kit
1387	abk1842H	Human CUBN (Cubilin) ELISA Kit
1388	abk2842H	Human CBG(Corticosteroid Binding Globulin) ELISA Kit
1389	abk3842H	Human PTK6(Protein Tyrosine Kinase 6) ELISA Kit
1390	abk4842H	Human TNFRSF5(Tumor Necrosis Factor Receptor Superfamily, Member 5) ELISA Kit
1391	abk5842H	Human HIF-2α(Hypoxia Inducible Factor 2 Alpha) ELISA Kit
1392	abk6842H	Human FCRL3(Fc receptor-like protein 3) ELISA Kit
1393	abk7842H	Human GPS1(G Protein Pathway Suppressor 1) ELISA Kit
1394	abk8842H	Human CYGB(Cytoglobin) ELISA Kit
1395	abk9842H	Human CDKN1A(Cyclin Dependent Kinase Inhibitor 1A) ELISA Kit
1396	abk0942H	Human CASP6(Caspase 6) ELISA Kit
1397	abk1942H	Human Itgam/CD11b(Integrin alpha-M) ELISA Kit
1398	abk2942H	Human CLSTN1(Calsyntenin 1) ELISA Kit
1399	abk3942H	Human SIRT2(Sirtuin 2) ELISA Kit
1400	abk4942H	Human TNFSF13(Tumor Necrosis Factor Ligand Superfamily, Member 13) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
1401	abk5942H	Human SKIL(SKI Like Oncogene ELISA Kit
1402	abk6942H	Human DSG3(Desmogleins 3) ELISA Kit
1403	abk8942H	Human DMBT1(Deleted In Malignant Brain Tumors 1) ELISA Kit
1404	abk9942H	Human PLXNA2(Plexin-A2) ELISA Kit
1405	abk0052H	Human TGM6(Transglutaminase-6) ELISA Kit
1406	abk1052H	Human MAP1B(Microtubule Associated Protein Light Chain 1B) ELISA Kit
1407	abk3052H	Human SI(Sucrase Isomaltase) ELISA Kit
1408	abk4052H	Human CART(Cocaine And Amphetamine Regulated Transcript) ELISA Kit
1409	abk5052H	Human DUOX2(Dual Oxidase 2) ELISA Kit
1410	abk6052H	Human GGT2(Gamma-glutamyltranspeptidase 2) ELISA Kit
1411	abk7052H	Human sβAPP(Soluble β Amyloid Precursor Protein) ELISA Kit
1412	abk8052H	Human sTNFRSF12A(Soluble Tumor Necrosis Factor Receptor Superfamily, Member 12A) ELISA Kit
1413	abk9052H	Human GLP1R(Glucagon-like Peptide 1 receptor) ELISA Kit
1414	abk0152H	Human AR(Aldose reductase) ELISA Kit
1415	abk1152H	Human MAGED2(Melanoma-associated antigen D2) ELISA Kit
1416	abk2152H	Human SYN1(Synapsin 1) ELISA Kit
1417	abk3152H	Human BTD(Biotinidase) ELISA Kit
1418	abk4152H	Human HABP2(Hyaluronan-binding protein 2) ELISA Kit
1419	abk5152H	Human SRGN(Serglycin) ELISA Kit
1420	abk6152H	Human SAT1(Spermidine/Spermine N1-Acetyltransferase 1) ELISA Kit
1421	abk7152H	Human RNLS(Renalase) ELISA Kit
1422	abk8152H	Human Cadm4(Cell adhesion molecule 4) ELISA Kit
1423	abk9152H	Human KEAP1(Kelch-like ECH-associated protein 1) ELISA Kit
1424	abk0252H	Human NEO1(Neogenin 1) ELISA Kit
1425	abk1252H	Human CD27BP(CD27 Binding Protein) ELISA Kit
1426	abk2252H	Human UCN3(Urocortin 3) ELISA Kit
1427	abk3252H	Human IRF1(Interferon Regulatory Factor 1) ELISA Kit
1428	abk4252H	Human ALPI(Intestinal alkaline phosphatase 1) ELISA Kit
1429	abk5252H	Human ADAM33 (A Disintegrin And Metalloprotease 33) ELISA Kit
1430	abk6252H	Human SRD5a1(Steroid 5 Alpha Reductase 1) ELISA Kit
1431	abk7252H	Human VASP(Vasodilator Stimulated Phosphoprotein) ELISA Kit
1432	abk8252H	Human DAG1(Dystrophin Associated Glycoprotein 1) ELISA Kit
1433	abk9252H	Human HABP1/C1QBP(Hyaluronan Binding Protein 1) ELISA Kit
1434	abk0352H	Human SOCS1(Suppressor of Cytokine Signaling 1) ELISA Kit
1435	abk1352H	Human SRD5a2(Steroid 5 Alpha Reductase 2) ELISA Kit
1436	abk3352H	Human IL-37(Interleukin 37) ELISA Kit
1437	abk4352H	Human CXCR2(CXC-Chemokine Receptor 2) ELISA Kit
1438	abk5352H	Human MEP1α(Meprin A Alpha) ELISA Kit
1439	abk6352H	Human SGPP1(Sphingosine-1-Phosphate Phosphatase 1) ELISA Kit
1440	abk7352H	Human SPHK1(Sphingosine Kinase 1) ELISA Kit
1441	abk8352H	Human suPAR(Soluble Urokinase-type Plasminogen Activator Receptor) ELISA Kit
1442	abk9352H	Human KLF4(Krueppel-Like Factor 4) ELISA Kit
1443	abk0452H	Human IL-11(Interleukin 11) ELISA Kit
1444	abk1452H	Human sdLDL(Small Dense Low Density Lipoprotein) ELISA Kit
1445	abk2452H	Human TSHR(Thyroid Stimulating Hormone Receptor) ELISA Kit
1446	abk3452H	Human TERF2IP(Telomeric Repeat Binding Factor 2 Interacting Protein) ELISA Kit
1447	abk4452H	Human DBT(Dihydrolipoamide branched chain transacylase E2) ELISA Kit
1448	abk5452H	Human PMEPA1(Prostate Transmembrane Protein androgen Induced 1) ELISA Kit
1449	abk6452H	Human RGC32(Response Gene To Complement 32) ELISA Kit
1450	abk7452H	Human KAL1(Anosmin-1) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
1451	abk8452H	Human AD(Adropin) ELISA Kit
1452	abk0552H	Human sENG(Soluble Endoglin) ELISA Kit
1453	abk1552H	Human c-MET/HGFR(hepatocyte growth factor receptor) ELISA Kit
1454	abk2552H	Human RECK(Reversion-inducing-cysteine-rich protein with kazal motifs) ELISA Kit
1455	abk3552H	Human TPH1(Tryptophan Hydroxylase 1) ELISA Kit
1456	abk4552H	Human pMAPT/pTAU(phosphorylated microtubule-associated protein tau) ELISA Kit
1457	abk5552H	Human HDAC1(Histone Deacetylase 1) ELISA Kit
1458	abk6552H	Human SCUBE1(Signal Peptide, CUB and EGF-like Domain-containing Protein 1) ELISA Kit
1459	abk7552H	Human GPX1(Glutathione Peroxidase 1) ELISA Kit
1460	abk8552H	Human PSMA(Prostate specific membrane antigen) ELISA Kit
1461	abk9552H	Human NID2(Nidogen 2) ELISA Kit
1462	abk0652H	Human MAC(Membrane Attack Complex) ELISA Kit
1463	abk2652H	Human GRIA1(Glutamate receptor 1) ELISA Kit
1464	abk3652H	Human MeCP2(Methyl CpG Binding Protein 2) ELISA Kit
1465	abk4652H	Human MBD1(Methyl CpG Binding Domain Protein 1) ELISA Kit
1466	abk5652H	Human MBD2(Methyl CpG Binding Domain Protein 2) ELISA Kit
1467	abk6652H	Human DNMT1(DNA Methyltransferase 1) ELISA Kit
1468	abk7652H	Human DNMT2(DNA Methyltransferase 2) ELISA Kit
1469	abk8652H	Human DNMT3A(DNA Methyltransferase 3A) ELISA Kit
1470	abk9652H	Human DNMT3B(DNA Methyltransferase 3B) ELISA Kit
1471	abk0752H	Human ROMO1(Reactive oxygen species modulator 1) ELISA Kit
1472	abk1752H	Human CTSS(Cathepsin S) ELISA Kit
1473	abk2752H	Human Ki67P(Ki-67 Protein) ELISA Kit
1474	abk3752H	Human HE4(epididymal protein 4) ELISA Kit
1475	abk4752H	Human FGF18(Fibroblast Growth Factor 18) ELISA Kit
1476	abk5752H	Human TACSTD2(Tumor-associated Calcium Signal Transducer 2) ELISA Kit
1477	abk6752H	Human PEBP4(Phosphatidylethanolamine Binding Protein 4) ELISA Kit
1478	abk7752H	Human MR-ProADM(Mid-regional pro-adrenomedullin) ELISA Kit
1479	abk8752H	Human SFRP4(secreted frizzled-related protein 4) ELISA Kit
1480	abk9752H	Human FCN3(Ficolin 3) ELISA Kit
1481	abk0852H	Human KL(Klotho) ELISA Kit
1482	abk1852H	Human ARL1(Aldose Reductase Like Protein 1) ELISA Kit
1483	abk2852H	Human PTP1B(protein tyrosine phosphatase 1B) ELISA Kit
1484	abk3852H	Human BIRC7(Baculoviral IAP repeat-containing protein 7) ELISA Kit
1485	abk4852H	Human C3d(Complement Fragment 3d) ELISA Kit
1486	abk5852H	Human FGF8(Fibroblast Growth Factor 8) ELISA Kit
1487	abk6852H	Human FGF10(Fibroblast Growth Factor 10) ELISA Kit
1488	abk7852H	Human ABCC5(ATP Binding Cassette Transporter C5) ELISA Kit
1489	abk8852H	Human LAG3(Lymphocyte Activation Gene 3) ELISA Kit
1490	abk9852H	Human HGS(Hepatocyte growth factor-regulated tyrosine kinase substrate) ELISA Kit
1491	abk0952H	Human LEFTY1(Left/Right Determination Factor 1) ELISA Kit
1492	abk1952H	Human SCGB2A2(Mammaglobin-A) ELISA Kit
1493	abk2952H	Human BLVRA(Biliverdin Reductase A) ELISA Kit
1494	abk4952H	Human SLC16A8(Solute Carrier Family 16, Member 8) ELISA Kit
1495	abk5952H	Human IER3(Immediate Early Response 3) ELISA Kit
1496	abk6952H	Human ADAMTS9(A Disintegrin And Metalloproteinase With Thrombospondin 9) ELISA Kit
1497	abk7952H	Human HES1(Transcription factor HES-1) ELISA Kit
1498	abk8952H	Human CSNK2A1(Casein Kinase 2, alpha 1 Polypeptide) ELISA Kit
1499	abk9952H	Human CSNK2B(Casein Kinase 2, beta) ELISA Kit
1500	abk0062H	Human SMRP(Soluble Mesothelin Related Peptide) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
1501	abk1062H	Human CCR4(Chemokine C-C-Motif Receptor 4) ELISA Kit
1502	abk2062H	Human Pro-ADM(Proadrenomedullin) ELISA Kit
1503	abk3062H	Human CECR1(Cat Eye Syndrome Chromosome Region, Candidate 1) ELISA Kit
1504	abk4062H	Human PROM1(Prominin-1) ELISA Kit
1505	abk5062H	Human Hpc25(Hepcidin 25) ELISA Kit
1506	abk6062H	Human HOXA10(Homeobox protein Hox-A10) ELISA Kit
1507	abk7062H	Human HOXB4(Homeobox protein Hox-B4) ELISA Kit
1508	abk8062H	Human MC5R(Melanocortin 5 Receptor) ELISA Kit
1509	abk9062H	Human QSOX1(Sulfhydryl oxidase 1) ELISA Kit
1510	abk0162H	Human ATRN(Attractin) ELISA Kit
1511	abk1162H	Human SCG3(Secretogranin III) ELISA Kit
1512	abk2162H	Human IL-28A(Interleukin 28A) ELISA Kit
1513	abk3162H	Human IL-28B(Interleukin 28B) ELISA Kit
1514	abk4162H	Human IL-29(Interleukin 29) ELISA Kit
1515	abk5162H	Human C4d(Complement Fragment 4d) ELISA Kit
1516	abk6162H	Human VNN1(Vanin 1) ELISA Kit
1517	abk7162H	Human ANXA1(Annexin A1) ELISA Kit
1518	abk8162H	Human BCL2A1(Bcl-2-related protein A1) ELISA Kit
1519	abk9162H	Human TIMD4(T-cell immunoglobulin and mucin domain-containing protein 4) ELISA Kit
1520	abk0262H	Human VSNL1(Visinin Like Protein 1) ELISA Kit
1521	abk2102H	Human IAP(Integrin Associated Protein) ELISA Kit
1522	abk3262H	Human ANXA3(Annexin A3) ELISA Kit
1523	abk4262H	Human AKAP12(A Kinase Anchor Protein 12) ELISA Kit
1524	abk5262H	Human C1QTNF2(C1q and Tumor Necrosis Factor Related Protein 2) ELISA Kit
1525	abk6262H	Human C1QTNF3(C1q and Tumor Necrosis Factor Related Protein 3) ELISA Kit
1526	abk7262H	Human C1QTNF9(C1q and Tumor Necrosis Factor Related Protein 9) ELISA Kit
1527	abk8262H	Human FAPα(Fibroblast Activation Protein Alpha) ELISA Kit
1528	abk9262H	Human GDN(Glia Derived Nexin) ELISA Kit
1529	abk0362H	Human PLOD3(Procollagen Lysine-2-Oxoglutarate-5-Dioxygenase 3) ELISA Kit
1530	abk1362H	Human CLIP(Corticotropin Like Intermediate Lobe Peptide) ELISA Kit
1531	abk2362H	Human SESN1(Sestrin 1) ELISA Kit
1532	abk3362H	Human LAMA4(Laminin Alpha 4) ELISA Kit
1533	abk4362H	Human CSNK2A2(Casein Kinase 2, alpha 2 Polypeptide) ELISA Kit
1534	abk5362H	Human STUB1(E3 ubiquitin-protein ligase CHIP) ELISA Kit
1535	abk6362H	Human pNF-H (Phosphorylated Neurofilament H) ELISA Kit
1536	abk7362H	Human DHFR(Dihydrofolate Reductase) ELISA Kit
1537	abk8362H	Human sFAS/sAPO-1(Soluble Factor-Related Apoptosis) ELISA Kit
1538	abk9362H	Human ADAMTS1(A Disintegrin and Metalloproteinase with Thrombospondin Motifs 1) ELISA Kit
1539	abk0462H	Human SPINK1(Serine Peptidase Inhibitor Kazal Type 1) ELISA Kit
1540	abk1462H	Human WISP1(WNT1 Inducible Signaling Pathway Protein 1) ELISA Kit
1541	abk2462H	Human SEMG1(Semenogelin I) ELISA Kit
1542	abk3462H	Human SFRP5(secreted frizzled-related protein 5) ELISA Kit
1543	abk4462H	Human TMPRSS6(Transmembrane Protease, Serine 6) ELISA Kit
1544	abk5462H	Human C1QTNF12(C1q and Tumor Necrosis Factor Related Protein 12) ELISA Kit
1545	abk6462H	Human SERPINA4(Kallistatin) ELISA Kit
1546	abk7462H	Human ADAMTS7(A Disintegrin And Metalloproteinase With Thrombospondin 7) ELISA Kit
1547	abk8462H	Human JAK1(Janus Kinase 1) ELISA Kit
1548	abk9462H	Human CAPN1(Calpain 1) ELISA Kit
1549	abk0562H	Human IRS1(Insulin Receptor Substrate 1) ELISA Kit
1550	abk1562H	Human CAMK2N2(Calcium/Calmodulin Dependent Protein Kinase II Inhibitor 2) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
1551	abk3562H	Human ADAMTS12(A Disintegrin And Metalloproteinase With Thrombospondin 12) ELISA Kit
1552	abk4562H	Human sRANKL(Soluble Receptor Activator of Nuclear factor-κB Ligand) ELISA Kit
1553	abk5562H	Human JAK3(Janus Kinase 3) ELISA Kit
1554	abk6562H	Human GLDL(Glycated Low Density Lipoprotein) ELISA Kit
1555	abk7562H	Human ITLN2(Intelectin 2) ELISA Kit
1556	abk8562H	Human TPMT(Thiopurine S-Methyltransferase) ELISA Kit
1557	abk9562H	Human CYP2A6(Cytochrome P450 2A6) ELISA Kit
1558	abk0662H	Human BCHE(Butyrylcholinesterase) ELISA Kit
1559	abk1662H	Human SCUBE3(Signal Peptide, CUB and EGF-like Domain-containing Protein 3) ELISA Kit
1560	abk2662H	Human CRY2(Cryptochrome 2) ELISA Kit
1561	abk3662H	Human RCN1(Reticulocalbin 1) ELISA Kit
1562	abk4662H	Human GLDC(Glycine Dehydrogenase) ELISA Kit
1563	abk5662H	Human CD5L(CD5 Antigen-like) ELISA Kit
1564	abk6662H	Human COX2(Cytochrome C Oxidase Subunit II) ELISA Kit
1565	abk7662H	Human PEPD(Peptidase D) ELISA Kit
1566	abk8662H	Human XBP1(X-Box Binding Protein 1) ELISA Kit
1567	abk9662H	Human XOD(Xanthine Oxidase) ELISA Kit
1568	abk0762H	Human pro-IL1β(pro interleukin 1 beta) ELISA Kit
1569	abk1762H	Human KRT1(Keratin 1) ELISA Kit
1570	abk2762H	Human ACC(Acetyl-CoA Carboxylase) ELISA Kit
1571	abk3762H	Human CRMP5(Collapsin Response Mediator Protein 5) ELISA Kit
1572	abk4762H	Human GRP78(Glucose Regulated Protein 78) ELISA Kit
1573	abk5762H	Human VCAM-1/CD106(Vascular Cell Adhesion Molecule 1) ELISA Kit
1574	abk6762H	Human HSL(Hormone Sensitive Lipase) ELISA Kit
1575	abk7762H	Human AITRL/TNFSF18(Activation-inducible TNF-related Ligand) ELISA Kit
1576	abk8762H	Human ADAMTS5(A Disintegrin And Metalloproteinase With Thrombospondin 5) ELISA Kit
1577	abk9762H	Human EAAT3(Excitatory Amino Acid Transporter 3) ELISA Kit
1578	abk0862H	Human STIP1(Stress Induced Phosphoprotein 1) ELISA Kit
1579	abk1862H	Human MIOX(Myoinositol Oxygenase) ELISA Kit
1580	abk2862H	Human FAAH(Fatty-acid amide hydrolase) ELISA Kit
1581	abk3862H	Human ATGL(Adipose Triglyceride Lipase) ELISA Kit
1582	abk4862H	Human ADAMTS8(A Disintegrin And Metalloproteinase With Thrombospondin 8) ELISA Kit
1583	abk5862H	Human NLRC4(NLR family CARD domain-containing protein 4) ELISA Kit
1584	abk6862H	Human NFATC2(Nuclear Factor Of Activated T-Cells, Cytoplasmic 2) ELISA Kit
1585	abk7862H	Human ELA1(Elastase 1, Pancreatic) ELISA Kit
1586	abk8862H	Human ADAMTS18(A Disintegrin And Metalloproteinase With Thrombospondin 18) ELISA Kit
1587	abk9862H	Human CST1(Cystatin-SN) ELISA Kit
1588	abk0962H	Human TBX2(T-Box Protein 2) ELISA Kit
1589	abk1962H	Human IREB2(Iron Responsive Element Binding Protein 2) ELISA Kit
1590	abk2962H	Human SPX(Spexin) ELISA Kit
1591	abk3962H	Human CTSA(Cathepsin A) ELISA Kit
1592	abk4962H	Human FGF3(Fibroblast Growth Factor 3) ELISA Kit
1593	abk5962H	Human STATH(Statherin) ELISA Kit
1594	abk6962H	Human UCP-3(Uncoupling Protein 3, Mitochondrial) ELISA Kit
1595	abk7962H	Human FABP12(Fatty Acid Binding Protein 12) ELISA Kit
1596	abk8962H	Human TGFBR1(TGF-beta receptor type 1) ELISA Kit
1597	abk9962H	Human ADORA2a(Adenosine A2a Receptor) ELISA Kit
1598	abk0072H	Human FGF14(Fibroblast Growth Factor 14) ELISA Kit
1599	abk2072H	Human HPX(Hemopexin) ELISA Kit
1600	abk3072H	Human CHIT1(Chitotriosidase-1) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
1601	abk4072H	Human ADP-HMW(High Molecular Weight Adiponectin) ELISA Kit
1602	abk5072H	Human EIF2αK3(Eukaryotic Translation Initiation Factor 2 Alpha K inase 3) ELISA Kit
1603	abk6072H	Human Smad2/MADH2(Mothers Against Decapentaplegic Homolog 2) ELISA Kit
1604	abk7072H	Human ADORA1(Adenosine A1 Receptor) ELISA Kit
1605	abk9072H	Human ART3(Ecto-ADP-ribosyltransferase 3) ELISA Kit
1606	abk0172H	Human TWIST1(Twist-related protein 1) ELISA Kit
1607	abk1172H	Human Smad6/MADH6(Mothers Against Decapentaplegic Homolog 6) ELISA Kit
1608	abk2172H	Human SH3YL1(SH3 domain-containing YSC84-like protein 1) ELISA Kit
1609	abk3172H	Human RXRα(Retinoid X Receptor Alpha) ELISA Kit
1610	abk4172H	Human CDH19(Cadherin 19) ELISA Kit
1611	abk5172H	Human ADORA3(Adenosine A3 Receptor) ELISA Kit
1612	abk6172H	Human TGFBR2(TGF-beta receptor type 2) ELISA Kit
1613	abk7172H	Human FGF22(Fibroblast Growth Factor 22) ELISA Kit
1614	abk8172H	Human CD3e(Cluster of Differentiation 3e) ELISA Kit
1615	abk9172H	Human EIF2a(Eukaryotic Translation Initiation Factor 2 Alpha) ELISA Kit
1616	abk0272H	Human SAA4(Serum Amyloid A4) ELISA Kit
1617	abk1272H	Human SSB(Sjogren Syndrome Antigen B) ELISA Kit
1618	abk2272H	Human PRKACA(cAMP-dependent protein kinase catalytic subunit alpha) ELISA Kit
1619	abk3272H	Human PLEC(Plectin) ELISA Kit
1620	abk4272H	Human SIK2(Serine/threonine protein kinase SIK2) ELISA Kit
1621	abk5272H	Human ANKRD12(Ankyrin Repeat Domain-containing Protein 12) ELISA Kit
1622	abk6272H	Human CD3g(Cluster of Differentiation 3g) ELISA Kit
1623	abk7272H	Human CMA1(Chymase 1, Mast Cell) ELISA Kit
1624	abk8272H	Human HDGF(Hepatoma Derived Growth Factor) ELISA Kit
1625	abk9272H	Human AgRP(Agouti Related Protein) ELISA Kit
1626	abk7611H	Human ANG(Angiogenin) ELISA Kit
1627	abk1372H	Human MCP-1(Monocyte Chemotactic Protein 1) ELISA Kit
1628	abk2372H	Human CFD(Complement Factor D) ELISA Kit
1629	abk3372H	Human IL-8(Interleukin 8) ELISA Kit
1630	abk4372H	Human CXCL16(Chemokine C-X-C-Motif Ligand 16) ELISA Kit
1631	abk5372H	Human ES(Endostatin) ELISA Kit
1632	abk6372H	Human HcpC(Hepcidin) ELISA Kit
1633	abk7372H	Human Gas6(Growth Arrest Specific Protein 6) ELISA Kit
1634	abk8372H	Human sgp130(Soluble Glycoprotein 130) ELISA Kit
1635	abk9372H	Human LEP(Leptin) ELISA Kit
1636	abk0472H	Human MAU(Microalbuminuria) ELISA Kit
1637	abk3472H	Human OxLDL(Oxidized Low Density Lipoprotein) ELISA Kit
1638	abk4472H	Human MMP-10(Matrix Metalloproteinase 10) ELISA Kit
1639	abk5472H	Human MMP-13(Matrix Metalloproteinase 13) ELISA Kit
1640	abk6472H	Human MYO(Myoglobin) ELISA Kit
1641	abk7472H	Human GLP-1(Glucagon Like Peptide 1) ELISA Kit
1642	abk8472H	Human C4(Complement Component 4) ELISA Kit
1643	abk9472H	Human TF(Transferrin) ELISA Kit
1644	abk0572H	Human KIM-1(Kidney Injury Molecule 1) ELISA Kit
1645	abk1572H	Human P I CP(Procollagen I C-Terminal Propeptide) ELISA Kit
1646	abk2572H	Human AChE(Acetylcholinesterase) ELISA Kit
1647	abk3572H	Human ALCAM(Activated Leukocyte Cell Adhesion Molecule) ELISA Kit
1648	abk4572H	Human ALDOA(Aldolase A, Fructose Bisphosphate) ELISA Kit
1649	abk5572H	Human ANGPTL2(Angiopoietin Like Protein 2) ELISA Kit
1650	abk6572H	Human FGFR4(Fibroblast Growth Factor Receptor 4) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
1651	abk7572H	Human CYP1A1(Cytochrome P450, family 1, subfamily A, polypeptide 1) ELISA Kit
1652	abk8572H	Human IGF-2(Insulin Like Growth Factor 2) ELISA Kit
1653	abk9572H	Human IGFBP-2(Insulin-like Growth Factor Binding Protein 2) ELISA Kit
1654	abk0672H	Human ApoA2(Apolipoprotein A2) ELISA Kit
1655	abk1672H	Human ApoF(Apolipoprotein F) ELISA Kit
1656	abk2672H	Human ApoA5(Apolipoprotein A5) ELISA Kit
1657	abk3672H	Human bFGF/FGF2(Basic Fibroblast Growth Factor) ELISA Kit
1658	abk4672H	Human FGF4(Fibroblast Growth Factor 4) ELISA Kit
1659	abk5672H	Human BMP-1(Bone Morphogenetic Protein 1) ELISA Kit
1660	abk6672H	Human BMP-15(Bone Morphogenetic Protein 15) ELISA Kit
1661	abk7672H	Human BMP-6(Bone Morphogenetic Protein 6) ELISA Kit
1662	abk9672H	Human CKLF(Chemokine Like Factor) ELISA Kit
1663	abk0772H	Human COL3(Collagen Type III) ELISA Kit
1664	abk1772H	Human COL10(Collagen Type X) ELISA Kit
1665	abk2772H	Human COLEC11(Collectin Sub-Family Member 11) ELISA Kit
1666	abk3772H	Human C1INH(Complement 1 Inhibitor) ELISA Kit
1667	abk4772H	Human C1q(Complement 1q) ELISA Kit
1668	abk5772H	Human C3(Complement Component 3) ELISA Kit
1669	abk6772H	Human C2(Complement Component 2) ELISA Kit
1670	abk7772H	Human CFB(Complement Factor B) ELISA Kit
1671	abk8772H	Human CFH(Complement Factor H) ELISA Kit
1672	abk9772H	Human CFI(Complement Factor I) ELISA Kit
1673	abk0872H	Human C4b(Complement Component 4b) ELISA Kit
1674	abk1872H	Human α2-M(Alpha-2 Macroglobulin) ELISA Kit
1675	abk2872H	Human GLβ(Galactosidase Beta) ELISA Kit
1676	abk3872H	Human Mly/CXCL9(Monocyte Interferon Gamma Inducing Factor) ELISA Kit
1677	abk4872H	Human PKM2(Pyruvate Kinase, Muscle) ELISA Kit
1678	abk5872H	Human VTN(Vitronectin) ELISA Kit
1679	abk6872H	Human DβH(Dopamine-β-Hydroxylase) ELISA Kit
1680	abk7872H	Human HIF-1α(Hypoxia Inducible Factor 1 Alpha) ELISA Kit
1681	abk8872H	Human LRG1(Leucine Rich Alpha-2-Glycoprotein 1) ELISA Kit
1682	abk9872H	Human GDH/GLDH(Glutamate dehydrogenase) ELISA Kit
1683	abk0972H	Human GAD2(Glutamate Decarboxylase 2) ELISA Kit
1684	abk1972H	Human GSTα1(Glutathione S Transferase Alpha 1) ELISA Kit
1685	abk3972H	Human TLN(Talin) ELISA Kit
1686	abk4972H	Human MMP-1(Matrix Metalloproteinase 1) ELISA Kit
1687	abk5972H	Human IgA2(Immunoglobulin A2) ELISA Kit
1688	abk7972H	Human I-PTH(intact Parathormone) ELISA Kit
1689	abk7453H	Human TG(Thyroglobulin) ELISA Kit
1690	abk9972H	Human HP(Haptoglobin) ELISA Kit
1691	abk0082H	Human MIP-3β/ELC/CCL19(Macrophage Inflammatory Protein 3β) ELISA Kit
1692	abk1082H	Human HBsAg(Hepatitis B surface antigen) ELISA Kit
1693	abk2082H	Human sST2(Soluble ST2) ELISA Kit
1694	abk3082H	Human EGFR2(Epidermal Growth Factor Receptor 2) ELISA Kit
1695	abk4082H	Human ApoB48(Apolipoprotein B48) ELISA Kit
1696	abk5082H	Human sTfR1(Soluble Transferrin Receptor 1) ELISA Kit
1697	abk6082H	Human ORM2(Orosomucoid 2) ELISA Kit
1698	abk7082H	Human IgG1(Immunoglobulin G1) ELISA Kit
1699	abk8082H	Human IgG3(Immunoglobulin G3) ELISA Kit
1700	abk9082H	Human IgG4(Immunoglobulin G4) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
1701	abk0182H	Human NT-proANP(N-Terminal Pro Atrial Natriuretic Peptide) ELISA Kit
1702	abk1182H	Human BGN(Biglycan) ELISA Kit
1703	abk2182H	Human NRG-1(Neuregulin 1) ELISA Kit
1704	abk4182H	Human NT-4(Neurotrophin 4) ELISA Kit
1705	abk5182H	Human ARG2(Arginase II) ELISA Kit
1706	abk6182H	Human CYPA(Cyclophilin A) ELISA Kit
1707	abk7182H	Human FETUB(Fetuin B) ELISA Kit
1708	abk8182H	Human KRT10(Keratin, Type I Cytoskeletal 10) ELISA Kit
1709	abk9182H	Human KISS1(Kisspeptin 1) ELISA Kit
1710	abk0282H	Human SHBG(Sex Hormone Binding Globulin) ELISA Kit
1711	abk1282H	Human TECK/CCL25(Thymus Expressed Chemokine) ELISA Kit
1712	abk2282H	Human OLFM4(Olfactomedin 4) ELISA Kit
1713	abk3282H	Human VE-Cadherin(Vascular Endothelial Cadherin) ELISA Kit
1714	abk6282H	Human FG(Fibrinogen) ELISA Kit
1715	abk7282H	Human LBP(Lipopolysaccharide Binding Protein) ELISA Kit
1716	abk8282H	Human DKK3(Dickkopf Related Protein 3) ELISA Kit
1717	abk9282H	Human iC3b(Complement iC3b) ELISA Kit
1718	abk0382H	Human SAA2(Serum Amyloid A2) ELISA Kit
1719	abk1382H	Human ICAM-1/CD54(intercellular adhesion molecule 1) ELISA Kit
1720	abk2382H	Human F8(Coagulation Factor VIII) ELISA Kit
1721	abk3382H	Human sVEGFR-1(soluble Vascular endothelial growth factor receptor 1) ELISA Kit
1722	abk4382H	Human ucOC(Undercarboxylated Osteocalcin) ELISA Kit
1723	abk5382H	Human ANG1(Angiopoietin 1) ELISA Kit
1724	abk6382H	Human Irisin ELISA Kit
1725	abk8382H	Human TLR4(Toll-Like Receptor 4) ELISA Kit
1726	abk9382H	Human IFN-α(Interferon Alpha) ELISA Kit
1727	abk1482H	Human D2D(D-Dimer) ELISA Kit
1728	abk2482H	Human DPP4(Dipeptidyl Peptidase IV) ELISA Kit
1729	abk3482H	Human sCD14(Soluble Cluster of Differentiation 14) ELISA Kit
1730	abk4482H	Human EMMPRIN/CD147(Extracellular Matrix Metalloproteinase Inducer) ELISA Kit
1731	abk5482H	Human CTSB(Cathepsin B) ELISA Kit
1732	abk6482H	Human FABP1(Fatty Acid Binding Protein 1, Liver) ELISA Kit
1733	abk7482H	Human AGP1(Alpha-1-Acid Glycoprotein 1) ELISA Kit
1734	abk8482H	Human BDNF(Brain Derived Neurotrophic Factor) ELISA Kit
1735	abk9482H	Human CASP1(Caspase 1) ELISA Kit
1736	abk0582H	Human CASP3(Caspase 3) ELISA Kit
1737	abk1582H	Human MIP-3α(Macrophage Inflammatory Protein 3 Alpha) ELISA Kit
1738	abk2582H	Human MDC/CCL22(Macrophage-Derived Chemokine) ELISA Kit
1739	abk3582H	Human CLU(Clusterin) ELISA Kit
1740	abk5582H	Human ENA-78/CXCL5(Epithelial Neutrophil Activating Peptide 78) ELISA Kit
1741	abk7582H	Human SDF-1/CXCL12(Stromal Cell Derived Factor 1) ELISA Kit
1742	abk9582H	Human APN(Aminopeptidase N) ELISA Kit
1743	abk0682H	Human Flt3L(FMS-like Tyrosine Kinase 3 Ligand) ELISA Kit
1744	abk1682H	Human G-CSF(Granulocyte Colony-stimulating Factor) ELISA Kit
1745	abk2682H	Human GM-CSF(Granulocyte-Macrophage Colony Stimulating Factor) ELISA Kit
1746	abk3682H	Human IGF-1(Insulin-like Growth Factor 1) ELISA Kit
1747	abk4682H	Human IGFBP-3(Insulin-like Growth Factor Binding Protein 3) ELISA Kit
1748	abk5682H	Human FGF7/KGF(Fibroblast Growth Factor 7) ELISA Kit
1749	abk6682H	Human IL-2(Interleukin 2) ELISA Kit
1750	abk8682H	Human ACAT2(Acetyl Coenzyme A Acetyltransferase 2) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
1751	abk0782H	Human IFABP/FABP2(Intestinal Fatty Acid Binding Protein) ELISA Kit
1752	abk3782H	Human IgG(Immunoglobulin G) ELISA Kit
1753	abk4782H	Human HCG(Chorionic Gonadotropin) ELISA Kit
1754	abk5782H	Human P I NP(Procollagen I N-Terminal Propeptide) ELISA Kit
1755	abk9782H	Human INHA(Inhibin Alpha Chain) ELISA Kit
1756	abk0882H	Human ACD(Adrenocortical Dysplasia Homolog) ELISA Kit
1757	abk2882H	Human ACO-2(Aconitase 2) ELISA Kit
1758	abk3882H	Human ACP1(Acid Phosphatase 1) ELISA Kit
1759	abk4882H	Human ADCY3(Adenylate Cyclase 3) ELISA Kit
1760	abk5882H	Human ADCY4(Adenylate Cyclase 4) ELISA Kit
1761	abk6882H	Human ADCYAP1/PACAP(Adenylate Cyclase Activating Polypeptide 1, Pituitary) ELISA Kit
1762	abk7882H	Human AFAP1(Actin Filament Associated Protein 1) ELISA Kit
1763	abk8882H	Human AIG1(Androgen Induced Protein 1) ELISA Kit
1764	abk9882H	Human AKAP1(A Kinase Anchor Protein 1) ELISA Kit
1765	abk0982H	Human ALK(Anaplastic Lymphoma Kinase) ELISA Kit
1766	abk3982H	Human ANGPTL3(Angiopoietin Like Protein 3) ELISA Kit
1767	abk5982H	Human RARS(Arginyl tRNA Synthetase) ELISA Kit
1768	abk0092H	Human iPLA2(Phospholipase A2, Calcium Independent) ELISA Kit
1769	abk1092H	Human GUK1(Guanylate Kinase 1) ELISA Kit
1770	abk2092H	Human HB(Hemoglobin) ELISA Kit
1771	abk3092H	Human HRG(Histidine Rich Glycoprotein) ELISA Kit
1772	abk4092H	Human HSPG1(Heparan Sulfate Proteoglycan 1) ELISA Kit
1773	abk5092H	Human IL-10Rβ(Interleukin 10 Receptor Beta) ELISA Kit
1774	abk6092H	Human AT-III(Antithrombin III) ELISA Kit
1775	abk9092H	Human IGFBP-5(Insulin-like Growth Factor Binding Protein 5) ELISA Kit
1776	abk0192H	Human IGFBP-4(Insulin-like Growth Factor Binding Protein 4) ELISA Kit
1777	abk1192H	Human IGFBP-6(Insulin-like Growth Factor Binding Protein 6) ELISA Kit
1778	abk4192H	Human ApoC3(Apolipoprotein C3) ELISA Kit
1779	abk5192H	Human ApoB100(Apolipoprotein B100) ELISA Kit
1780	abk6192H	Human ITIH4(Inter-alpha-trypsin Inhibitor Heavy Chain H4) ELISA Kit
1781	abk7192H	Human ARG1(Arginase 1) ELISA Kit
1782	abk8192H	Human ATXN2(Ataxin 2) ELISA Kit
1783	abk9192H	Human ATF5(Activating Transcription Factor 5) ELISA Kit
1784	abk0292H	Human ATF7(Activating Transcription Factor 7) ELISA Kit
1785	abk1292H	Human ABCG1(ATP Binding Cassette Transporter G1) ELISA Kit
1786	abk2292H	Human ABCA1(ATP Binding Cassette Transporter A1) ELISA Kit
1787	abk3292H	Human ABCG2(ATP Binding Cassette Transporter G2) ELISA Kit
1788	abk4292H	Human Aβ1-40(Amyloid Beta 1-40) ELISA Kit
1789	abk5292H	Human Aβ1-42(Amyloid Beta 1-42) ELISA Kit
1790	abk6292H	Human Bcl-3(B-cell Lymphoma/Leukemia 3) ELIS Kit ELISA Kit
1791	abk7292H	Human BLNK(B-Cell Linker Protein) ELISA Kit
1792	abk8292H	Human BTC(Betacellulin) ELISA Kit
1793	abk9292H	Human BMPR1A(Bone Morphogenetic Protein Receptor 1A) ELISA Kit
1794	abk0392H	Human IBSP(Bone Sialoprotein 2) ELISA Kit
1795	abk1392H	Human BFP(Brain Finger Protein) ELISA Kit
1796	abk2392H	Human C4BP(C4 Binding Protein) ELISA Kit
1797	abk3392H	Human CGRP1(Calcitonin Gene Related Peptide 1) ELISA Kit
1798	abk4392H	Human CaSR(Calcium Sensing Receptor) ELISA Kit
1799	abk5392H	Human CAP1(Adenylate Cyclase Associated Protein 1) ELISA Kit
1800	abk6392H	Human CA19-9(Carbohydrate antigen19-9) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
1801	abk7392H	Human CMLC-1(Cardiac Myosin Light Chain 1) ELISA Kit
1802	abk8392H	Human TNNI3/cTn-I(Troponin I Type 3, Cardiac) ELISA Kit
1803	abk9392H	Human CT1(Cardiotrophin 1) ELISA Kit
1804	abk0492H	Human CASP10(Caspase 10) ELISA Kit
1805	abk1492H	Human CASP12(Caspase 12) ELISA Kit
1806	abk2492H	Human Cav-1(Caveolin-1) ELISA Kit
1807	abk3492H	Human CCR1(CC-Chemokine Receptor 1) ELISA Kit
1808	abk4492H	Human CDC23(Cell Division Cycle Protein 23) ELISA Kit
1809	abk5492H	Human CENPE(Centromere Protein E) ELISA Kit
1810	abk6492H	Human CENPH(Centromere Protein H) ELISA Kit
1811	abk7492H	Human CENPI(Centromere Protein I) ELISA Kit
1812	abk8492H	Human CENPB(Centromere Protein B) ELISA Kit
1813	abk9492H	Human CEP110(Centrosomal Protein 110kDa) ELISA Kit
1814	abk0592H	Human CPN10(Chaperonin 10) ELISA Kit
1815	abk1592H	Human CCL1(Chemokine C-C-Motif Ligand 1) ELISA Kit
1816	abk2592H	Human CCL4L1(Chemokine C-C-Motif Ligand 4 Like Protein 1) ELISA Kit
1817	abk3592H	Human CCRL1(Chemokine C-C-Motif Receptor Like Protein 1) ELISA Kit
1818	abk4592H	Human CX3CR1(Chemokine C-X3-C-Motif Receptor 1) ELISA Kit
1819	abk5592H	Human CKLFSF7(Chemokine Like Factor Superfamily 7) ELISA Kit
1820	abk6592H	Human CKLFSF8(Chemokine Like Factor Superfamily 8) ELISA Kit
1821	abk7592H	Human CCKBR(Cholecystokinin B Receptor) ELISA Kit
1822	abk8592H	Human CHAc(Choline acetylase) ELISA Kit
1823	abk9592H	Human PC/CPG(Choline Phosphoglyceride) ELISA Kit
1824	abk0692H	Human CHD3(Chromodomain Helicase DNA Binding Protein 3) ELISA Kit
1825	abk1692H	Human CHGA(Chromogranin A) ELISA Kit
1826	abk2692H	Human CLDN14(Claudin 14) ELISA Kit
1827	abk3692H	Human CD109(Cluster of Differentiation 109) ELISA Kit
1828	abk5741H	Human MYC(Myc proto-oncogene protein) ELISA Kit
1829	abk5692H	Human CCSA-2(Colon Cancer Specific Antigen-2) ELISA Kit
1830	abk6692H	Human CCSA-3(Colon Cancer Specific Antigen-3) ELISA Kit
1831	abk7692H	Human CSF2Rα(Colony Stimulating Factor 2 Receptor Alpha) ELISA Kit
1832	abk8692H	Human P II NP(Procollagen II N-Terminal ProPeptide) ELISA Kit
1833	abk9692H	Human GCSFR(Colony Stimulating Factor Receptor, Granulocyte) ELISA Kit
1834	abk0792H	Human C1qR1(Complement Component 1q Receptor 1) ELISA Kit
1835	abk1792H	Human C1r(Complement Component 1R) ELISA Kit
1836	abk2792H	Human C5(Complement Component 5) ELISA Kit
1837	abk3792H	Human C8b(Complement Component 8b) ELISA Kit
1838	abk4792H	Human C8g(Complement Component 8g) ELISA Kit
1839	abk5792H	Human CTAPIII(Connective Tissue Activating Peptide III) ELISA Kit
1840	abk6792H	Human CX31(Connexin 31) ELISA Kit
1841	abk7792H	Human HELIX- II (Type II Collagen Helical Peptide) ELISA Kit
1842	abk8792H	Human MTHFR(Methylenetetrahydrofolate Reductase) ELISA Kit
1843	abk9792H	Human PRKAα1(Protein Kinase, AMP Activated Alpha 1) ELISA Kit
1844	abk0892H	Human IKKBK(Inhibitor Of Kappa-Light Polypeptide Gene Enhancer In B-Cells Kinase Beta) ELISA Kit
1845	abk1892H	Human PTPRB(Protein Tyrosine Phosphatase Receptor Type B) ELISA Kit
1846	abk2892H	Human CREB1(Cyclic AMP Response Element Binding Protein 1) ELISA Kit
1847	abk3892H	Human CXCR3(CXC-Chemokine Receptor 3) ELISA Kit
1848	abk4892H	Human CLEC2C(C-Type Lectin Domain Family 2, Member C) ELISA Kit
1849	abk5892H	Human SELE(E-Selectin) ELISA Kit
1850	abk6892H	Human FBXO2(F-Box Protein 2) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
1851	abk7892H	Human Flt3(FMS-like Tyrosine Kinase 3) ELISA Kit
1852	abk8892H	Human GATA1(GATA Binding Protein 1) ELISA Kit
1853	abk9892H	Human GATA2(GATA Binding Protein 2) ELISA Kit
1854	abk0992H	Human GATA5(GATA Binding Protein 5) ELISA Kit
1855	abk1992H	Human GDI1(GDP Dissociation Inhibitor 1) ELISA Kit
1856	abk2992H	Human GAP(Ras Gtpase Activating Protein) ELISA Kit
1857	abk3992H	Human RGS4(Regulator of G Protein Signaling 4) ELISA Kit
1858	abk4992H	Human NOD9(Nucleotide-binding oligomerization domain protein 9) ELISA Kit
1859	abk5992H	Human S100A1(S100 Calcium Binding Protein A1) ELISA Kit
1860	abk6992H	Human SHH(Sonic Hedgehog Protein) ELISA Kit
1861	abk7992H	Human TDP43(TAR DNA Binding Protein 43) ELISA Kit
1862	abk8992H	Human TAF1(Transcription initiation factor TFIID subunit 1) ELISA Kit
1863	abk9992H	Human TAF2(Transcription initiation factor TFIID subunit 2) ELISA Kit
1864	abk0003H	Human TAF13(Transcription initiation factor TFIID subunit 13) ELISA Kit
1865	abk1003H	Human TLR-2(Toll-like Receptor 2) ELISA Kit
1866	abk2003H	Human TU3A(TU3A Protein) ELISA Kit
1867	abk3003H	Human LAT(Linker for activation of T-cell) ELISA Kit
1868	abk4003H	Human TALLA-1/CD231(T-cell Acute Lymphoblastic Leukemia Antigen) ELISA Kit
1869	abk5003H	Human OCT1(Octamer-binding protein 1) ELISA Kit
1870	abk6003H	Human IL-1R1(Interleukin 1 Receptor Type I) ELISA Kit
1871	abk7003H	Human LCA/CD45(Leukocyte Common Antigen) ELISA Kit
1872	abk9003H	Human cPLA2(Cytosolic Phospholipase A2) ELISA Kit
1873	abk0103H	Human PDHB(Pyruvate dehydrogenase E1 component subunit Beta) ELISA Kit
1874	abk1103H	Human MlgA(Membrane Immunoglobulin A) ELISA Kit
1875	abk2103H	Human MlgD(Membrane Immunoglobulin D) ELISA Kit
1876	abk3103H	Human MlgM(Membrane Immunoglobulin M) ELISA Kit
1877	abk4103H	Human PDHA1(Pyruvate dehydrogenase E1 component subunit alpha) ELISA Kit
1878	abk5103H	Human VIM(Vimentin) ELISA Kit
1879	abk6103H	Human FOXP1(Forkhead Box Protein P1) ELISA Kit
1880	abk7103H	Human FGFBP2(Fibroblast Growth Factor Binding Protein 2) ELISA Kit
1881	abk8103H	Human FRS2(Fibroblast Growth Factor Receptor Substrate ELISA Kit
1882	abk9103H	Human VLA4(Very Late Appearing Antigen 4) ELISA Kit
1883	abk0203H	Human PRF1(Perforin 1) ELISA Kit
1884	abk3203H	Human PYGB(Glycogen Phosphorylase, Brain) ELISA Kit
1885	abk4203H	Human MAOA(Type A Monoamine Oxidase) ELISA Kit
1886	abk5203H	Human MCP-2/CCL8(Monocyte Chemotactic Protein 2) ELISA Kit
1887	abk6203H	Human SREBP(Sterol Regulatory Element Binding Protein) ELISA Kit
1888	abk7203H	Human EMILIN1(Elastin Microfibril Interface Located Protein 1) ELISA Kit
1889	abk8203H	Human PKB(Protein Kinase B) ELISA Kit
1890	abk9203H	Human PKD2(Protein Kinase D2) ELISA Kit
1891	abk0303H	Human LMP7/PSMB8(Low molecular mass protein 7/Proteasome subunit beta type-8) ELISA Kit
1892	abk2303H	Human PYY(Peptide YY) ELISA Kit
1893	abk3303H	Human DPYSL3(Dihydropyrimidinase Like Protein 3) ELISA Kit
1894	abk4303H	Human DFDT1(Farnesyl Diphosphate Farnesyltransferase 1) ELISA Kit
1895	abk5303H	Human SP-A1(Pulmonary surfactant-associated protein A1) ELISA Kit
1896	abk6303H	Human SP-B(Pulmonary surfactant-associated protein B) ELISA
1897	abk7303H	Human SP-D(Pulmonary surfactant-associated protein D) ELISA
1898	abk8303H	Human PRCP(Prolylcarboxypeptidase) ELISA Kit
1899	abk9303H	Human RPA1(Replication protein A1) ELISA Kit
1900	abk0403H	Human MBP-C/MBL2(Mannose Binding Protein C/Mannose Binding Lectin) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
1901	abk1403H	Human RNASEH1(Ribonuclease H1) ELISA Kit
1902	abk2403H	Human RNASET2(Ribonuclease T2) ELISA Kit
1903	abk3403H	Human NFkB-p105(Nuclear factor NF-kappa-B p105 subunit) ELISA Kit
1904	abk4403H	Human NFkB-p65(Nuclear factor NF-kappa-B p65 subunit) ELISA Kit
1905	abk5403H	Human βMSH(beta-Melanocyte Stimulating Hormone) ELISA Kit
1906	abk6403H	Human GRIK2(Glutamate Receptor, Ionotropic, Kainate 2) ELISA Kit
1907	abk8403H	Human CKMB(Creatine Kinase MB Isoenzyme) ELISA Kit
1908	abk0503H	Human CNTN4(Contactin 4) ELISA Kit
1909	abk1503H	Human UPP2(Uridine Phosphorylase 2) ELISA Kit
1910	abk2503H	Human MT1(Metallothionein 1) ELISA Kit
1911	abk3503H	Human SSFA2(Sperm Specific Antigen 2) ELISA Kit
1912	abk4503H	Human PD-1(Programmed Cell Death Protein 1) ELISA Kit
1913	abk5503H	Human PD-L1(Programmed Cell Death Protein 1 Ligand 1) ELISA Kit
1914	abk7503H	Human τPK1(Tau-Protein Kinase 1) ELISA Kit
1915	abk8503H	Human TRACP-5b(Tartrate Resistant Acid Phosphatase 5b) ELISA Kit
1916	abk9503H	Human HMGB-1(High Mobility Group Protein B1) ELISA Kit
1917	abk2603H	Human NFE2L2(Nuclear Factor, Erythroid Derived 2, Like 2) ELISA Kit
1918	abk4603H	Human PDGF-BB(Platelet Derived Growth Factor BB) ELISA Kit
1919	abk5603H	Human RBP4(Retinol Binding Protein 4) ELISA Kit
1920	abk6603H	Human TGF-β2(Transforming Growth Factor Beta 2) ELISA Kit
1921	abk7603H	Human TAP1(Transporter associated with Antigen Processing 1) ELISA Kit
1922	abk8603H	Human CD28(Cluster of Differentiation 28) ELISA Kit
1923	abk9603H	Human CD30L(Cluster of Differentiation 30 Ligand) ELISA Kit
1924	abk0703H	Human CD86(T-Lymphocyte Activation Antigen CD86) ELISA Kit
1925	abk1703H	Human sTNF(Soluble Tumor Necrosis Factor) ELISA Kit
1926	abk2703H	Human MHCG/HLA-G(Major Histocompatibility Complex Class I G) ELISA Kit
1927	abk3703H	Human PDGF-A(Platelet Derived Growth Factor Subunit A) ELISA Kit
1928	abk4703H	Human PDGF-B(Platelet Derived Growth Factor Subunit B) ELISA Kit
1929	abk5703H	Human CX37(Connexin 37) ELISA Kit
1930	abk6703H	Human LFA-3/CD58(Lymphocyte Function Associated Antigen 3) ELISA Kit
1931	abk7703H	Human pERK1/2(Phospho Extracellular Signal Regulated Kinase 1/2) ELISA Kit
1932	abk8703H	Human PLA2P(Phospholipase A2 Activating Protein) ELISA Kit
1933	abk9703H	Human PLA1(Phospholipase A1) ELISA Kit
1934	abk0803H	Human PLCβ1/PI-PLC(Phospholipase C Beta 1, Phosphoinositide Specific) ELISA Kit
1935	abk3803H	Human PGDS2(Prostaglandin D2 Synthase) ELISA Kit
1936	abk4803H	Human PIBF(Progesterone Induced Blocking Factor) ELISA Kit
1937	abk5803H	Human DEFα5(Defensin Alpha 5, Paneth Cell Specific) ELISA Kit
1938	abk6803H	Human IgM(Immunoglobulin M) ELISA Kit
1939	abk7803H	Human hs-Cys-C(high-sensitivity Cystatin C) ELISA Kit
1940	abk8803H	Human EP2(Prostaglandin E Receptor 2) ELISA Kit
1941	abk9803H	Human ACP3/PAP(Prostatic Acid Phosphatase) ELISA Kit
1942	abk0903H	Human DSG1(Desmogleins 1) ELISA Kit
1943	abk1903H	Human HJV(Hemojuvelin) ELISA Kit
1944	abk2903H	Human HSP-70/HSPA9(Heat Shock 70 kDa Protein 9) ELISA Kit
1945	abk3903H	Human HSF1(Heat Shock Factor 1) ELISA Kit
1946	abk4903H	Human WARS(Tryptophanyl tRNA Synthetase) ELISA Kit
1947	abk6903H	Human NT-3(Neurotrophin 3) ELISA Kit
1948	abk7903H	Human GROβ/CXCL2(Growth Regulated Oncogene Beta) ELISA Kit
1949	abk8903H	Human NE/ELA2(Neutrophil Elastase/Elastase-2) ELISA Kit
1950	abk9903H	Human STK39(Serine/Threonine Protein Kinase 39) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
1951	abk0013H	Human DAP6(Death Associated Protein 6) ELISA Kit
1952	abk1013H	Human MPO(Myeloperoxidase) ELISA Kit
1953	abk3013H	Human TPPP(Tubulin Polymerization Promoting Protein) ELISA Kit
1954	abk5013H	Human CTLA4(Cytotoxic T-Lymphocyte Associated Antigen 4) ELISA Kit
1955	abk6013H	Human CK-17/KRT17(Cytokeratin 17) ELISA Kit
1956	abk7013H	Human CK-18/KRT18(Cytokeratin 18) ELISA Kit
1957	abk8013H	Human CYP7A1(Cholesterol 7-Alpha-Hydroxylase) ELISA Kit
1958	abk9013H	Human CCND1(Cyclin D1) ELISA Kit
1959	abk0113H	Human CCND2(Cyclin D2) ELISA Kit
1960	abk1113H	Human FLNA(Filamin A, Alpha) ELISA Kit
1961	abk2113H	Human Plasmin ELISA Kit
1962	abk3113H	Human CKMT1A(Creatine kinase U-type, mitochondrial) ELISA Kit
1963	abk4113H	Human SNRPD1(Small Nuclear Ribonucleoprotein D1) ELISA Kit
1964	abk5113H	Human STAT5A(Signal Transducer and Activator of Transcription 5A) ELISA Kit
1965	abk6113H	Human SOX2(Sex Determining Region Y Box Protein 2) ELISA Kit
1966	abk8113H	Human HbH(Hemoglobin H) ELISA Kit
1967	abk9113H	Human SEPP1(Selenoprotein P) ELISA Kit
1968	abk1213H	Human CR1(Complement Receptor type 1) ELISA Kit
1969	abk3213H	Human PAF(Platelet Activating Factor) ELISA Kit
1970	abk4213H	Human PDGF-C(Platelet Derived Growth Factor C) ELISA Kit
1971	abk5213H	Human PDGF-D(Platelet Derived Growth Factor D) ELISA Kit
1972	abk7213H	Human NOX4(Nicotinamide Adenine Dinucleotide Phosphate Oxidase 4) ELISA Kit
1973	abk8213H	Human NOX5(Nicotinamide Adenine Dinucleotide Phosphate Oxidase 5) ELISA Kit
1974	abk1313H	Human FNDC5(Fibronectin type III domain-containing protein 5) ELISA Kit
1975	abk2313H	Human MEC/CCL28(Mucosae Associated Epithelia Chemokine) ELISA Kit
1976	abk4313H	Human TACE/ADAM17(TNF α Converting Enzyme) ELISA Kit
1977	abk5313H	Human TSGF(Tumor Specific Growth Factor/Tumor Supplied Group of Factor) ELISA Kit
1978	abk6313H	Human TFR1(Transferrin Receptor 1) ELISA Kit
1979	abk7313H	Human CD8(Cluster of Differentiation 8) ELISA Kit
1980	abk9313H	Human SEPCR/PROCR(Soluble Endothelial Protein C Receptor) ELISA Kit
1981	abk0413H	Human PLA2G2A(Phospholipase A2, Group II A) ELISA Kit
1982	abk1413H	Human PS1(Presenilin 1) ELISA Kit
1983	abk2413H	Human HMG-CoAR(3-Hydroxy-3-Methylglutaryl Coenzyme A Reductase) ELISA Kit
1984	abk3413H	Human IRF3(Interferon Regulatory Factor 3) ELISA Kit
1985	abk4413H	Human ZPI(Protein Z Dependent Protease Inhibitor) ELISA Kit
1986	abk5413H	Human CD74(HLA class II histocompatibility antigen gamma chain) ELISA Kit
1987	abk6413H	Human RTN4(Reticulon 4) ELISA Kit
1988	abk7413H	Human MAP2(Microtubule Associated Protein 2) ELISA Kit
1989	abk8413H	Human OPRM1(Mu-type opioid receptor) ELISA Kit
1990	abk9413H	Human sAGER(Soluble Advanced Glycosylation End Product Specific Receptor) ELISA Kit
1991	abk0513H	Human IL-26(Interleukin 26) ELISA Kit
1992	abk1513H	Human ITGa2B(Integrin Alpha 2B) ELISA Kit
1993	abk2513H	Human ALPL(Alkaline Phosphatase, Liver/Bone/Kidney) ELISA Kit
1994	abk4513H	Human TSH(Thyroid Stimulating Hormone) ELISA Kit
1995	abk5513H	Human EPO(Erythropoietin) ELISA Kit
1996	abk6513H	Human F-hCGβ(Free Chorionic Gonadotropin Beta subunit) ELISA Kit
1997	abk7513H	Human CGβ(Chorionic Gonadotrophin Beta) ELISA Kit
1998	abk8513H	Human TWEAK(Tumour Necrosis Factor Related Weak Inducer of Apoptosis) ELISA Kit
1999	abk9513H	Human hs-CRP(high-sensitivity C-Reactive Protein) ELISA Kit
2000	abk0613H	Human LTF/LF(Lactoferrin) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
2001	abk1613H	Human HSD17β10(17-Beta-Hydroxysteroid Dehydrogenase Type 10) ELISA Kit
2002	abk2613H	Human SORT1(Sortilin) ELISA Kit
2003	abk3613H	Human PD-L2(Programmed Cell Death Protein 1 Ligand 2) ELISA Kit
2004	abk4613H	Human PIVKA- II (Protein Induced by Vitamin K Absence or Antagonist- II) ELISA Kit
2005	abk5613H	Human N-ERC/Mesothelin ELISA Kit
2006	abk6613H	Human IL-31(Interleukin 31) ELISA Kit
2007	abk7613H	Human B3GALT1(Beta-1,3-glucosyltransferase) ELISA Kit
2008	abk8613H	Human CXCR4(Chemokine C-X-C-Motif Receptor 4) ELISA Kit
2009	abk9613H	Human C-ERC/Mesothelin ELISA Kit
2010	abk0713H	Human fFN(Fetal Fibronectin) ELISA Kit
2011	abk1713H	Human CF6(Coupling Factor 6) ELISA Kit
2012	abk3713H	Human γMSH(gamma-Melanocyte Stimulating Hormone) ELISA Kit
2013	abk4713H	Human Zonulin ELISA Kit
2014	abk5713H	Human OGDH(2-oxoglutarate dehydrogenase, mitochondrial) ELISA Kit
2015	abk6713H	Human AMPH(Amphiphysin) ELISA Kit
2016	abk7713H	Human ICAM-3/CD50(Intercellular Adhesion Molecule 3) ELISA Kit
2017	abk8713H	Human CD19(Cluster of Differentiation 19) ELISA Kit
2018	abk9713H	Human HNF1α(Hepatocyte Nuclear Factor 1 Alpha) ELISA Kit
2019	abk0813H	Human CoQ10B(Coenzyme Q10 Homolog B) ELISA Kit
2020	abk1813H	Human RANKL(Receptor Activator of Nuclear Factor Kappa B Ligand) ELISA Kit
2021	abk2813H	Human M I CA(major histocompatibility complex class I chain-related gene A) ELISA Kit
2022	abk3813H	Human CP(Ceruloplasmin) ELISA Kit
2023	abk4813H	Human NT-proBNP(N-Terminal Pro-Brain Natriuretic Peptide) ELISA Kit
2024	abk6813H	Human IL-6(Interleukin 6) ELISA Kit
2025	abk7813H	Human IL-10(Interleukin 10) ELISA Kit
2026	abk9813H	Human POSTN/OSF-2(Periostin) ELISA Kit
2027	abk1913H	Human GDF2(Growth Differentiation Factor 2) ELISA Kit
2028	abk2913H	Human NRP1(Neuropilin 1) ELISA Kit
2029	abk3913H	Human BNP(Brain Natriuretic Peptide) ELISA Kit
2030	abk4913H	Human DBP(Vitamin D-binding protein) ELISA Kit
2031	abk5913H	Human α-Glu(Alpha Glucosidase) ELISA Kit
2032	abk6913H	Human MIF(Macrophage Migration Inhibitory Factor) ELISA Kit
2033	abk7913H	Human ENG(Endoglin) ELISA Kit
2034	abk8913H	Human LEPR(Leptin Receptor) ELISA Kit
2035	abk9913H	Human VEGFR1/FLT1(Vascular Endothelial Growth Factor Receptor 1) ELISA Kit
2036	abk0023H	Human sTREM-1(soluble Triggering Receptor Expressed on Myeloid Cells-1) ELISA Kit
2037	abk1023H	Human ACE(Angiotensin I Converting Enzyme) ELISA Kit
2038	abk2023H	Human RANTES(Regulated On Activation, Normal T-Cell Expressed and Secreted) ELISA Kit
2039	abk3023H	Human SELP(P-Selectin) ELISA Kit
2040	abk4023H	Human IL17RA(Interleukin-17 receptor A) ELISA Kit
2041	abk6023H	Human PDGF-AA(Platelet Derived Growth Factor AA) ELISA Kit
2042	abk7023H	Human PF4(Platelet Factor 4) ELISA Kit
2043	abk0701H	Human ACTH(Adrenocorticotrophic Hormone) ELISA Kit
2044	abk7701H	Human α-CTx(Alpha Crosslaps) ELISA Kit
2045	abk0111H	Human 11-DH-TXB2(11-Dehydro-thromboxane B2) ELISA Kit
2046	abk5311H	Human ADH(Antidiuretic Hormone) ELISA Kit
2047	abk5611H	Human IAPP(Islet Amyloid Polypeptide) ELISA Kit
2048	abk3521H	Human APLN(Apelin) ELISA Kit
2049	abk1031H	Human ANP(Atrial Natriuretic Peptide) ELISA Kit
2050	abk9131H	Human β-EP(Beta-Endorphin) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
2051	abk8231H	Human BK(Bradykinin) ELISA Kit
2052	abk4331H	Human CNP(C-type Natriuretic Peptide) ELISA Kit
2053	abk0741H	Human CPP(Copeptin) ELISA Kit
2054	abk6841H	Human CTX-III(Cross Linked C-telopeptide of Type III Collagen) ELISA Kit
2055	abk2251H	Human Salusin α(Salusin Alpha) ELISA Kit
2056	abk0451H	Human β-CTx(Beta Crosslaps) ELISA Kit
2057	abk0751H	Human OXA(Orexin A) ELISA Kit
2058	abk0061H	Human NP(Nandrolone Phenylpropionate) ELISA Kit
2059	abk6261H	Human SOD1(Superoxide Dismutase 1, Soluble) ELISA Kit
2060	abk7261H	Human OGP(Osteogenic Growth Peptide) ELISA Kit
2061	abk1461H	Human TRH(Thyrotropin Releasing Hormone) ELISA Kit
2062	abk3461H	Human δSIP-α(Delta Sleep Inducing Peptide α) ELISA Kit
2063	abk5461H	Human SCT(Secretin) ELISA Kit
2064	abk4102H	Human MPA(Mycophenolic Acid) ELISA Kit
2065	abk9602H	Human PGF2α(Prostaglandin F2α) ELISA Kit
2066	abk5702H	Human Dyn(Big Dynorphin) ELISA Kit
2067	abk6902H	Human NKA(Neurokinin A) ELISA Kit
2068	abk7902H	Human NKB(Neurokinin B) ELISA Kit
2069	abk9902H	Human CL(Cardiolipin) ELISA Kit
2070	abk2012H	Human NPY(Neuropeptide Y) ELISA Kit
2071	abk1212H	Human GHRL(Ghrelin) ELISA Kit
2072	abk4212H	Human SST(Somatostatin) ELISA Kit
2073	abk0312H	Human OXB(Orexin B) ELISA Kit
2074	abk4812H	Human MT(Melatonin) ELISA Kit
2075	abk5812H	Human DHEA-S(Dehydroepiandrosterone Sulfate) ELISA Kit
2076	abk6812H	Human DPD(Deoxyypyridinoline) ELISA Kit
2077	abk9812H	Human SA(Sialic Acid) ELISA Kit
2078	abk9912H	Human UST2(Urotensin 2) ELISA Kit
2079	abk6022H	Human MTL(Motilin) ELISA Kit
2080	abk9022H	Human GT(Gastrin) ELISA Kit
2081	abk1622H	Human TMSβ10(Thymosin Beta 10) ELISA Kit
2082	abk2622H	Human TP5(Thymopentin) ELISA Kit
2083	abk4622H	Human VIP(Vasoactive Intestinal Peptide) ELISA Kit
2084	abk3822H	Human TXB2(Thromboxane B2) ELISA Kit
2085	abk4121H	Human GC(Glucagon) ELISA Kit
2086	abk0132H	Human GLP-2(Glucagon Like Peptide 2) ELISA Kit
2087	abk6932H	Human 11,12-DHETs(11,12-dihydroxyeicosatrienoic acids) ELISA Kit
2088	abk7932H	Human PCX(Podocalyxin) ELISA Kit
2089	abk2142H	Human SOD3(Superoxide Dismutase 3, Extracellular) ELISA Kit
2090	abk9342H	Human αMSH(Alpha-Melanocyte Stimulating Hormone) ELISA Kit
2091	abk0442H	Human EP(Erythrocyte Protoporphyrin) ELISA Kit
2092	abk2052H	Human DYNC2H1(Cytoplasmic Dynein 2 heavy Chain 1) ELISA Kit
2093	abk9452H	Human Salusin β(Salusin Beta) ELISA Kit
2094	abk3952H	Human BPn(biopyrrin) ELISA Kit
2095	abk2562H	Human BPDE(Benzo Pyrene Dihydrodiol Epoxide) ELISA Kit
2096	abk1072H	Human KISS-54(Kisspeptin-54) ELISA Kit
2097	abk8782H	Human AP17(Apelin 17) ELISA Kit
2098	abk1982H	Human Ang- II (Angiotensin II) ELISA Kit
2099	abk2982H	Human Ang- I (Angiotensin I) ELISA Kit
2100	abk4982H	Human Ang-III(Angiotensin III) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
2101	abk2192H	Human AP13(Apelin 13) ELISA Kit
2102	abk3192H	Human AP36(Apelin 36) ELISA Kit
2103	abk2203H	Human CRH(Corticotropin Releasing Hormone) ELISA Kit
2104	abk1803H	Human MS(Melatonin Sulfate) ELISA Kit
2105	abk2803H	Human EM2(Endomorphin 2) ELISA Kit
2106	abk5903H	Human NT(Neurotensin) ELISA Kit
2107	abk2013H	Human OB(Obestatin) ELISA Kit
2108	abk4013H	Human VIPR1(Vasoactive Intestinal Peptide Receptor 1) ELISA Kit
2109	abk7113H	Human TMSβ4(Thymosin Beta 4) ELISA Kit
2110	abk2213H	Human BMG/β2-MG(Beta-2-Microglobulin) ELISA Kit
2111	abk6213H	Human NADPH(Nicotinamide Adenine Dinucleotide Phosphate) ELISA Kit
2112	abk9213H	Human TAP(Trypsinogen Activation Peptide) ELISA Kit
2113	abk2713H	Human Ang1-7(Angiotensin 1-7) ELISA Kit
2114	abk0913H	Human IgG2(Immunoglobulin G2) ELISA Kit
2115	abk8023H	Human ACV-A(Activin A) ELISA Kit
2116	abk4422H	Human SOD2(Superoxide Dismutase 2, Mitochondrial) ELISA Kit
2117	abk0043H	Human TGF-β3(Transforming Growth Factor Beta 3) ELISA Kit
2118	abk4043H	Human βTG(Thromboglobulin, Beta) ELISA Kit
2119	abk4353H	Human RET(Ret Proto-Oncogene) ELISA Kit
2120	abk0213H	Human SAA(Serum Amyloid A) ELISA Kit
2121	abk2782H	Human FE(Ferritin) ELISA Kit
2122	abk7382H	Human ADP/Acrp30(Adiponectin) ELISA Kit
2123	abk7682H	Human IL-17A(Interleukin 17A) ELISA Kit
2124	abk6681H	Human PCT(Procalcitonin) ELISA Kit
2125	abk1701H	Human PRL(Prolactin) ELISA Kit
2126	abk7733H	Human Cys-C(Cystatin C) ELISA Kit
2127	abk4582H	Human CRP(C-Reactive Protein) ELISA Kit
2128	abk0002H	Human GPC4(Glypican 4) ELISA Kit
2129	abk1002H	Human TRUB1(Probable tRNA Pseudouridine Synthase 1) ELISA Kit
2130	abk1303H	Human VCAN(Versican/PG-M) ELISA Kit
2131	abk0482H	Human NGAL(Neutrophil Gelatinase Associated Lipocalin) ELISA Kit
2132	abk0271H	Human SEMA3A(Semaphorin 3A) ELISA Kit
2133	abk8531H	Human cTnT/TNNT2(Troponin T Type 2, Cardiac) ELISA Kit
2134	abk8313H	Human HS(Heparan Sulfate) ELISA Kit
2135	abk8802H	Human LZM(Lysozyme) ELISA Kit
2136	abk1241H	Human NOS3/eNOS(Nitric Oxide Synthase 3, Endothelial) ELISA Kit
2137	abk1203H	Human FSH(Follicle Stimulating Hormone) ELISA Kit
2138	abk1472H	Human LH(Luteinizing Hormone) ELISA Kit
2139	abk1782H	Human LP-a(Apolipoprotein a) ELISA Kit
2140	abk9682H	Human TnTf/TNNT3(Troponin T, fast skeletal muscle) ELISA Kit
2141	abk3182H	Human GFAP(Glial Fibrillary Acidic Protein) ELISA Kit
2142	abk1561H	Human PROC(Protein C) ELISA Kit
2143	abk8003H	Human NSE(Neuron Specific Enolase) ELISA Kit
2144	abk3313H	Human LpPLA2(Lipoprotein-associated Phospholipase A2) ELISA Kit
2145	abk7823H	Human PSA(Prostate Specific Antigen) ELISA Kit
2146	abk6972H	Human MMP-9(Matrix Metalloproteinase 9) ELISA Kit
2147	abk6733H	Human AMH(Anti-Mullerian Hormone) ELISA Kit
2148	abk5553H	Human CD81 antigen(CD81) ELISA kit
2149	abk3991H	Human Integrin beta-1(ITGB1) ELISA kit
2150	abk0653H	Human ZEB1(Zinc finger E-box-binding homeobox 1) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
2151	abk9553H	Human METRN(Meteorin) ELISA Kit
2152	abk8242H	Human IL-33(Interleukin 33) ELISA Kit
2153	abk8582H	Human αFP(Alpha-Fetoprotein) ELISA Kit
2154	abk0372H	Human CA125(Carbohydrate Antigen 125) ELISA Kit
2155	abk3481H	Human MMP-8(Matrix Metalloproteinase 8) ELISA Kit
2156	abk9531H	Human IL-7(Interleukin 7) ELISA Kit
2157	abk8672H	Human CEA(Carcinoembryonic Antigen) ELISA Kit
2158	abk9753H	Human GSDMD(Gasdermin D) ELISA Kit
2159	abk5031H	Human HBP/AZU(Heparin-binding protein/Azurocidin) ELISA Kit
2160	abk6201H	Human GROα/CXCL1(Growth Regulated Oncogene Alpha) ELISA Kit
2161	abk6601H	Human ApoA1(Apolipoprotein A1) ELISA Kit
2162	abk5813H	Human PTX 3/TSG-14(Pentraxin 3) ELISA Kit
2163	abk5023H	Human PCSK9(Proprotein Convertase Subtilisin/Kexin Type 9) ELISA Kit
2164	abk3513H	Human INS(Insulin) ELISA Kit
2165	abk7371H	Human S100B(S100 Calcium Binding Protein B) ELISA Kit
2166	abk4282H	Human IgE(Immunoglobulin E) ELISA Kit
2167	abk9501H	Human IFN-γ(Interferon Gamma) ELISA Kit
2168	abk7942H	Human CD200(Cluster of Differentiation 200) ELISA Kit
2169	abk4901H	Human TIMP-1(Tissue Inhibitors of Metalloproteinase 1) ELISA Kit
2170	abk8092H	Human IGFBP-1(Insulin-like Growth Factor Binding Protein 1) ELISA Kit
2171	abk8201H	Human IP-10/CXCL10(Interferon Gamma Induced Protein 10kDa) ELISA Kit
2172	abk2022H	Human PGA(Pepsinogen A) ELISA Kit
2173	abk3022H	Human PGC(Pepsinogen C) ELISA Kit
2174	abk0313H	Human CA242(Pancreatic Carcinoma Markers-CA242) ELISA Kit
2175	abk2472H	Human C-P(C-Peptide) ELISA Kit
2176	abk9201H	Human I-TAC/CXCL11(Interferon Inducible T-cell Alpha Chemoattractant) ELISA Kit
2177	abk2972H	Human IgA(Immunoglobulin A) ELISA Kit
2178	abk5002H	Human SCCA1(Squamous Cell Carcinoma Antigen 1) ELISA Kit
2179	abk5282H	Human ALB(Albumin) ELISA Kit
2180	abk0601H	Human TNF-α(Tumor Necrosis Factor Alpha) ELISA Kit
2181	abk4701H	Human IL-1β(Interleukin 1 Beta) ELISA Kit
2182	abk0054H	Human Amya (Amylase Alpha) ELISA Kit
2183	abk0563H	Human RIPK3(Receptor Interacting Serine Threonine Kinase 3)ELISA Kit
2184	abk2573H	Human IMA (Ischemia Modified Albumin) ELISA Kit
2185	abk2593H	Human PKBa (Protein Kinase B Alpha) ELISA Kit
2186	abk3073H	Human Asp(Asprosin) ELISA Kit
2187	abk3744H	Human LPO (Lactoperoxidase) ELISA Kit
2188	abk5820H	Human HSP gp96(Heat Shock Protein glycoprotein 96) ELISA Kit
2189	abk5883H	Human PIK3Ca (Phosphoinositide-3-Kinase Catalytic Alpha Polypeptide) ELISA Kit
2190	abk8453H	Human NLRP3(NLR Family, Pyrin Domain Containing Protein 3) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
1	abk3401M	Mouse FST(Follistatin) ELISA Kit
2	abk2603M	Mouse NFE2L2(Nuclear Factor, Erythroid Derived 2, Like 2) ELISA Kit
3	abk4701M	Mouse IL-1 β (Interleukin 1 Beta) ELISA Kit
4	abk1482M	Mouse D2D(D-Dimer) ELISA Kit
5	abk1561M	Mouse PROC(Protein C) ELISA Kit
6	abk0553M	Mouse ALT(Alanine Transaminase) ELISA Kit
7	abk8023M	Mouse ACV-A(Activin A) ELISA Kit
8	abk9023M	Mouse ANG3(Angiopoietin 3) ELISA Kit
9	abk2001M	Mouse BAFF/CD257(B-cell Activating Factor) ELISA Kit
10	abk6801M	Mouse LOX(Lysyl Oxidase) ELISA Kit
11	abk9001M	Mouse MIP-1 α (Macrophage Inflammatory Protein 1 Alpha) ELISA Kit
12	abk0101M	Mouse MIP-1 β (Macrophage Inflammatory Protein 1 Beta) ELISA Kit
13	abk1582M	Mouse MIP-3 α (Macrophage Inflammatory Protein 3 Alpha) ELISA Kit
14	abk2582M	Mouse MDC/CCL22(Macrophage-Derived Chemokine) ELISA Kit
15	abk7831M	Mouse CHEM(Chemerin) ELISA Kit
16	abk0123M	Mouse CHI3L1(Chitinase 3-like 1) ELISA Kit
17	abk6201M	Mouse GRO α /CXCL1(Growth Regulated Oncogene Alpha) ELISA Kit
18	abk7091M	Mouse ESM1(Endothelial Cell Specific Molecule 1) ELISA Kit
19	abk2401M	Mouse FGF21(Fibroblast Growth Factor 21) ELISA Kit
20	abk6401M	Mouse IFN- β (Interferon Beta) ELISA Kit
21	abk0501M	Mouse LIF(Leukemia Inhibitory Factor) ELISA Kit
22	abk7813M	Mouse IL-10(Interleukin 10) ELISA Kit
23	abk9201M	Mouse I-TAC/CXCL11(Interferon Inducible T-Cell Alpha Chemoattractant) ELISA Kit
24	abk7801M	Mouse GH(Growth Hormone) ELISA Kit
25	abk5631M	Mouse CASP8(Caspase 8) ELISA Kit
26	abk4472M	Mouse MMP-10(Matrix Metalloproteinase 10) ELISA Kit
27	abk8601M	Mouse MMP-12(Matrix Metalloproteinase 12) ELISA Kit
28	abk5472M	Mouse MMP-13(Matrix Metalloproteinase 13) ELISA Kit
29	abk7611M	Mouse ANG(Angiogenin) ELISA Kit
30	abk5671M	Mouse OPG(Osteoprotegerin) ELISA Kit
31	abk6472M	Mouse MYO(Myoglobin) ELISA Kit
32	abk2701M	Mouse PAPPA(Pregnancy Associated Plasma Protein A) ELISA Kit
33	abk5572M	Mouse ANGPTL2(Angiopoietin Like Protein 2) ELISA Kit
34	abk3711M	Mouse ANGPTL4(Angiopoietin Like Protein 4) ELISA Kit
35	abk6511M	Mouse AGT(Angiotensinogen) ELISA Kit
36	abk1001M	Mouse ANG2(Angiopoietin 2) ELISA Kit
37	abk6321M	Mouse ANXA5(Annexin A5) ELISA Kit
38	abk1521M	Mouse ANXA2(Annexin A2) ELISA Kit
39	abk2182M	Mouse NRG-1(Neuregulin 1) ELISA Kit
40	abk3732M	Mouse CDKN2A(Cyclin Dependent Kinase Inhibitor 2A) ELISA Kit
41	abk1671M	Mouse HDAC2(Histone Deacetylase 2) ELISA Kit
42	abk4941M	Mouse NRG-4(Neuregulin 4) ELISA Kit
43	abk5123M	Mouse ACF(Apobec 1 Complementation Factor) ELISA Kit
44	abk8521M	Mouse ApoC1(Apolipoprotein C1) ELISA Kit
45	abk0621M	Mouse ApoC4(Apolipoprotein C4) ELISA Kit
46	abk2621M	Mouse ApoE(Apolipoprotein E) ELISA Kit
47	abk1672M	Mouse ApoF(Apolipoprotein F) ELISA Kit
48	abk3621M	Mouse ApoH(Apolipoprotein H) ELISA Kit
49	abk6621M	Mouse AATF(Apoptosis Antagonizing Transcription Factor) ELISA Kit
50	abk6123M	Mouse AATK(Apoptosis Associated Tyrosine Kinase) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
51	abk8621M	Mouse AIF(Apoptosis Inducing Factor) ELISA Kit
52	abk1721M	Mouse ASK1(Apoptosis Signal Regulating Kinase 1) ELISA Kit
53	abk0721M	Mouse APAF1(Apoptosis Protease Activating Factor 1) ELISA Kit
54	abk4101M	Mouse SLC(Secondary Lymphoid Tissue Chemokine) ELISA Kit
55	abk2721M	Mouse AQP-1(Aquaporin 1) ELISA Kit
56	abk0112M	Mouse GRO γ /CXCL3(Growth Regulated Oncogene Gamma) ELISA Kit
57	abk3721M	Mouse AQP-2(Aquaporin 2) ELISA Kit
58	abk4721M	Mouse AQP-3(Aquaporin 3) ELISA Kit
59	abk5721M	Mouse AQP-4(Aquaporin 4) ELISA Kit
60	abk7721M	Mouse AVPR1A(Arginine Vasopressin Receptor 1A) ELISA Kit
61	abk5821M	Mouse AhR(Aryl Hydrocarbon Receptor) ELISA Kit
62	abk7123M	Mouse AIP(Aryl Hydrocarbon Receptor Interacting Protein) ELISA Kit
63	abk8123M	Mouse ARNT2(Aryl Hydrocarbon Receptor Nuclear Translocator 2) ELISA Kit
64	abk8821M	Mouse ASGR1(Asialoglycoprotein Receptor 1) ELISA Kit
65	abk1921M	Mouse AST(Aspartate Aminotransferase) ELISA Kit
66	abk4921M	Mouse ATM(Ataxia Telangiectasia Mutated) ELISA Kit
67	abk5921M	Mouse ATXN1(Ataxin 1) ELISA Kit
68	abk8192M	Mouse ATXN2(Ataxin 2) ELISA Kit
69	abk2292M	Mouse ABCA1(ATP Binding Cassette Transporter A1) ELISA Kit
70	abk1292M	Mouse ABCG1(ATP Binding Cassette Transporter G1) ELISA Kit
71	abk4031M	Mouse AXIN2(Axis Inhibition Protein 2) ELISA Kit
72	abk3672M	Mouse bFGF/FGF2(Basic Fibroblast Growth Factor) ELISA Kit
73	abk4672M	Mouse FGF4(Fibroblast Growth Factor 4) ELISA Kit
74	abk0131M	Mouse FGF6(Fibroblast Growth Factor 6) ELISA Kit
75	abk7031M	Mouse BAFFR(B-Cell Activation Factor Receptor) ELISA Kit
76	abk3131M	Mouse BAK1(BCL2 Antagonist/Killer 1) ELISA Kit
77	abk6131M	Mouse BMF(Bcl-2 Modifying Factor) ELISA Kit
78	abk8292M	Mouse BTC(Betacellulin) ELISA Kit
79	abk2111M	Mouse BACE1(Beta-Site APP Cleaving Enzyme 1) ELISA Kit
80	abk8413M	Mouse OPRM1(Mu-type opioid receptor) ELISA Kit
81	abk1231M	Mouse BID(BH3 Interacting Domain Death Agonist) ELISA Kit
82	abk8591M	Mouse IL-25(Interleukin 25) ELISA Kit
83	abk0301M	Mouse BLC(B-Lymphocyte Chemoattractant) ELISA Kit
84	abk5672M	Mouse BMP-1(Bone Morphogenetic Protein 1) ELISA Kit
85	abk3001M	Mouse BMP-2(Bone Morphogenetic Protein 2) ELISA Kit
86	abk4001M	Mouse BMP-4(Bone Morphogenetic Protein 4) ELISA Kit
87	abk5001M	Mouse BMP-7(Bone Morphogenetic Protein 7) ELISA Kit
88	abk9123M	Mouse BMP8B(Bone Morphogenetic Protein 8B) ELISA Kit
89	abk9292M	Mouse BMPR1A(Bone Morphogenetic Protein Receptor 1A) ELISA Kit
90	abk5231M	Mouse BMPR2(Bone Morphogenetic Protein Receptor II) ELISA Kit
91	abk9482M	Mouse CASP1(Caspase 1) ELISA Kit
92	abk8482M	Mouse BDNF(Brain Derived Neurotrophic Factor) ELISA Kit
93	abk7331M	Mouse CA724(Tumor Marker CA724) ELISA Kit
94	abk0223M	Mouse Ca-ATPase(Calcium ATPase) ELISA Kit
95	abk3282M	Mouse VE-Cadherin(Vascular Endothelial Cadherin) ELISA Kit
96	abk6001M	Mouse E-Cad(E-Cadherin) ELISA Kit
97	abk2011M	Mouse NCAD(Neural Cadherin) ELISA Kit
98	abk0431M	Mouse CABIN1(Calcineurin Binding Protein 1) ELISA Kit
99	abk1431M	Mouse CTR(Calcitonin Receptor) ELISA Kit
100	abk9601M	Mouse CT(Calcitonin) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
101	abk4392M	Mouse CaSR(Calcium Sensing Receptor) ELISA Kit
102	abk6431M	Mouse CNX(Calnexin) ELISA Kit
103	abk1223M	Mouse CAST(Calpastatin) ELISA Kit
104	abk7431M	Mouse CRT(Calreticulin) ELISA Kit
105	abk9431M	Mouse CALU(Calumenin) ELISA Kit
106	abk7001M	Mouse CA9(Carbonic Anhydrase IX) ELISA Kit
107	abk5531M	Mouse CA1(Carbonic Anhydrase I) ELISA Kit
108	abk0631M	Mouse CLS(Cardiolipin Synthase) ELISA Kit
109	abk1631M	Mouse CLCF1(Cardiotrophin Like Cytokine Factor 1) ELISA Kit
110	abk3631M	Mouse COMP(Cartilage Oligomeric Matrix Protein) ELISA Kit
111	abk0582M	Mouse CASP3(Caspase 3) ELISA Kit
112	abk6631M	Mouse CASP4(Caspase 4) ELISA Kit
113	abk7631M	Mouse CASP7(Caspase 7) ELISA Kit
114	abk9631M	Mouse CASP9(Caspase 9) ELISA Kit
115	abk0492M	Mouse CASP10(Caspase 10) ELISA Kit
116	abk4631M	Mouse CASP14(Caspase 14) ELISA Kit
117	abk2731M	Mouse CTSD(Cathepsin D) ELISA Kit
118	abk3731M	Mouse CTSK(Cathepsin K) ELISA Kit
119	abk4731M	Mouse CTSL(Cathepsin L) ELISA Kit
120	abk5731M	Mouse CDX2(Caudal Type Homeobox Transcription Factor 2) ELISA Kit
121	abk2223M	Mouse CEBPD(CCAAT/Enhancer Binding Protein Delta) ELISA Kit
122	abk2831M	Mouse CENPA(Centromere Protein A) ELISA Kit
123	abk8492M	Mouse CENPB(Centromere Protein B) ELISA Kit
124	abk5492M	Mouse CENPE(Centromere Protein E) ELISA Kit
125	abk6492M	Mouse CENPH(Centromere Protein H) ELISA Kit
126	abk9492M	Mouse CEP110(Centrosomal Protein 110kDa) ELISA Kit
127	abk3813M	Mouse CP(Ceruloplasmin) ELISA Kit
128	abk0931M	Mouse CCL3L1(Chemokine C-C-Motif Ligand 3 Like Protein 1) ELISA Kit
129	abk1931M	Mouse CCR8(Chemokine C-C-Motif Receptor 8) ELISA Kit
130	abk3592M	Mouse CCRL1(Chemokine C-C-Motif Receptor Like Protein 1) ELISA Kit
131	abk5201M	Mouse CX3CL1(Chemokine C-X3-C-Motif Ligand 1) ELISA Kit
132	abk4592M	Mouse CX3CR1(Chemokine C-X3-C-Motif Receptor 1) ELISA Kit
133	abk9672M	Mouse CKLF(Chemokine Like Factor) ELISA Kit
134	abk4931M	Mouse CKLFSF5(Chemokine Like Factor Superfamily 5) ELISA Kit
135	abk5931M	Mouse CKLFSF6(Chemokine Like Factor Superfamily 6) ELISA Kit
136	abk5592M	Mouse CKLFSF7(Chemokine Like Factor Superfamily 7) ELISA Kit
137	abk6592M	Mouse CKLFSF8(Chemokine Like Factor Superfamily 8) ELISA Kit
138	abk2041M	Mouse CCKAR(Cholecystokinin A Receptor) ELISA Kit
139	abk7592M	Mouse CCKBR(Cholecystokinin B Receptor) ELISA Kit
140	abk9592M	Mouse PC/CPG(choline phosphoglyceride) ELISA Kit
141	abk7372M	Mouse Gas6(Growth Arrest Specific Protein 6) ELISA Kit
142	abk7041M	Mouse CBX3(Chromobox Homolog 3) ELISA Kit
143	abk0692M	Mouse CHD3(Chromodomain Helicase DNA Binding Protein 3) ELISA Kit
144	abk8041M	Mouse CHD5(Chromodomain Helicase DNA Binding Protein 5) ELISA Kit
145	abk2201M	Mouse CNTF(Ciliary Neurotrophic Factor) ELISA Kit
146	abk3223M	Mouse CROP(Cisplatin Resistance Associated Overexpressed Protein) ELISA Kit
147	abk4223M	Mouse CC17(Clara Cell Protein) ELISA Kit
148	abk5223M	Mouse CD3(Cluster of differentiation 3) ELISA Kit
149	abk8141M	Mouse CD4(Cluster Of Differentiation 4) ELISA Kit
150	abk5141M	Mouse CD30(Cluster of differentiation 30) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
151	abk9603M	Mouse CD30L(Cluster of Differentiation 30 Ligand) ELISA Kit
152	abk6223M	Mouse CD40(Cluster of Differentiation 40) ELISA Kit
153	abk7223M	Mouse CD40L(Cluster of Differentiation 40 Ligand) ELISA Kit
154	abk6141M	Mouse CD97(Cluster of Differentiation 97) ELISA Kit
155	abk3692M	Mouse CD109(Cluster of Differentiation 109) ELISA Kit
156	abk5851M	Mouse GP4/CD36(Platelet Membrane Glycoprotein IV) ELISA Kit
157	abk7241M	Mouse F5(Coagulation Factor V) ELISA Kit
158	abk0341M	Mouse F13B(Coagulation Factor X III B Polypeptide) ELISA Kit
159	abk4831M	Mouse CFL1(Cofilin 1, Non-Muscle) ELISA Kit
160	abk1341M	Mouse CFL2(Cofilin 2, Muscle) ELISA Kit
161	abk0772M	Mouse COL3(Collagen Type III) ELISA Kit
162	abk8801M	Mouse COL4(Collagen Type IV) ELISA Kit
163	abk4341M	Mouse COL1α2(Collagen Type I Alpha 2) ELISA Kit
164	abk7341M	Mouse COL4α1(Collagen Type IV Alpha 1) ELISA Kit
165	abk1772M	Mouse COL10(Collagen Type X) ELISA Kit
166	abk8223M	Mouse CSF2Ra(Colony Stimulating Factor 2 Receptor Alpha) ELISA Kit
167	abk9223M	Mouse COL1(Collagen Type I) ELISA Kit
168	abk7911M	Mouse DCLK1(Doublecortin Like Kinase 1) ELISA Kit
169	abk9692M	Mouse GCSFR(Colony Stimulating Factor Receptor, Granulocyte) ELISA Kit
170	abk9991M	Mouse GPC3(Glypican 3) ELISA Kit
171	abk3772M	Mouse C1INH(Complement 1 Inhibitor) ELISA Kit
172	abk2792M	Mouse C5(Complement Component 5) ELISA Kit
173	abk7772M	Mouse CFB(Complement Factor B) ELISA Kit
174	abk2372M	Mouse CFD(Complement Factor D) ELISA Kit
175	abk3541M	Mouse C3a(Complement Component 3a) ELISA Kit
176	abk7901M	Mouse C4a(Complement Component 4a) ELISA Kit
177	abk9541M	Mouse CTGF(Connective Tissue Growth Factor) ELISA Kit
178	abk0641M	Mouse CX26(Connexin 26) ELISA Kit
179	abk6792M	Mouse CX31(Connexin 31) ELISA Kit
180	abk5703M	Mouse CX37(Connexin 37) ELISA Kit
181	abk0891M	Mouse CX40(Connexin 40) ELISA Kit
182	abk1891M	Mouse CX43(Connexin 43) ELISA Kit
183	abk0781M	Mouse CNTN1(Contactin 1) ELISA Kit
184	abk1781M	Mouse CNTN2(Contactin 2) ELISA Kit
185	abk2781M	Mouse CNTN3(Contactin 3) ELISA Kit
186	abk2091M	Mouse SIRT1(Sirtuin 1) ELISA Kit
187	abk2472M	Mouse C-P(C-Peptide) ELISA Kit
188	abk0381M	Mouse CKM(Creatine Kinase, Muscle) ELISA Kit
189	abk1232M	Mouse CRY1(Cryptochrome 1) ELISA Kit
190	abk4892M	Mouse CLEC2C(C-Type Lectin Domain Family 2, Member C) ELISA Kit
191	abk7101M	Mouse CTACK(Cutaneous T-Cell Attracting Chemokine) ELISA Kit
192	abk0323M	Mouse COL2(Collagen Type II) ELISA Kit
193	abk2732M	Mouse CDK9(Cyclin Dependent Kinase 9) ELISA Kit
194	abk1323M	Mouse CDKN1B(Cyclin Dependent Kinase Inhibitor 1B) ELISA Kit
195	abk6182M	Mouse CYP A(Cyclophilin A) ELISA Kit
196	abk2312M	Mouse CYPB(Cyclophilin B) ELISA Kit
197	abk3312M	Mouse CYPD(Cyclophilin D) ELISA Kit
198	abk7851M	Mouse CYSLTR2(Cysteiny Leukotriene Receptor 2) ELISA Kit
199	abk9122M	Mouse CYB561(Cytochrome b-561) ELISA Kit
200	abk1313M	Mouse FNDC5(Fibronectin type III domain-containing protein 5) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
201	abk7222M	Mouse CYTH1(Cytohesin 1) ELISA Kit
202	abk8222M	Mouse CYTIP(Cytohesin Interacting Protein) ELISA Kit
203	abk5013M	Mouse CTLA4(Cytotoxic T-Lymphocyte Associated Antigen 4) ELISA Kit
204	abk0013M	Mouse DAP6(Death Associated Protein 6) ELISA Kit
205	abk9132M	Mouse DCN(Decorin) ELISA Kit
206	abk2323M	Mouse DEFA5(Defensin Alpha 5, Paneth Cell Specific) ELISA Kit
207	abk6232M	Mouse VAP-1(Vascular Adhesion Protein 1) ELISA Kit
208	abk6781M	Mouse Des(Desmin) ELISA Kit
209	abk9961M	Mouse DAO(Diamine Oxidase) ELISA Kit
210	abk3323M	Mouse H3(Histone H3) ELISA Kit
211	abk9741M	Mouse DKK2(Dickkopf Related Protein 2) ELISA Kit
212	abk1071M	Mouse DLD(Dihydrolipoyl Dehydrogenase) ELISA Kit
213	abk2071M	Mouse DLAT(Dihydrolipoyl Transacetylase) ELISA Kit
214	abk0841M	Mouse D-LDH(D-Lactate Dehydrogenase) ELISA Kit
215	abk4323M	Mouse D1R(Dopamine D1 receptor) ELISA Kit
216	abk5323M	Mouse DRD3(Dopamine Receptor D3) ELISA Kit
217	abk6323M	Mouse DRD4(Dopamine Receptor D4) ELISA Kit
218	abk5961M	Mouse DAT(Dopamine Transporter) ELISA Kit
219	abk7832M	Mouse GAPDH(Glyceraldehyde-3-Phosphate Dehydrogenase) ELISA Kit
220	abk2381M	Mouse DMD(Dystrophin) ELISA Kit
221	abk9932M	Mouse SCARB1(Scavenger receptor class B member 1) ELISA Kit
222	abk8332M	Mouse EGR1(Early Growth Response Protein 1) ELISA Kit
223	abk9332M	Mouse EGR2(Early Growth Response Protein 2) ELISA Kit
224	abk0432M	Mouse EGR3(Early Growth Response Protein 3) ELISA Kit
225	abk1432M	Mouse EGR4(Early Growth Response Protein 4) ELISA Kit
226	abk8282M	Mouse DKK3(Dickkopf Related Protein 3) ELISA Kit
227	abk0561M	Mouse ELA4(Elastase 4) ELISA Kit
228	abk9461M	Mouse ELN(Elastin) ELISA Kit
229	abk7203M	Mouse EMILIN1(Elastin Microfibril Interface Located Protein 1) ELISA Kit
230	abk3032M	Mouse ELOA(Elongin A) ELISA Kit
231	abk9402M	Mouse EED(Embryonic Ectoderm Development) ELISA Kit
232	abk0302M	Mouse ERN1(Endoplasmic Reticulum To Nucleus Signalling 1) ELISA Kit
233	abk5372M	Mouse ES(Endostatin) ELISA Kit
234	abk7323M	Mouse EL(Endothelial lipase) ELISA Kit
235	abk1241M	Mouse NOS3/eNOS(Nitric Oxide Synthase 3, Endothelial) ELISA Kit
236	abk5202M	Mouse ECE1(Endothelin Converting Enzyme 1) ELISA Kit
237	abk4202M	Mouse ETRA(Endothelin Receptor A) ELISA Kit
238	abk5781M	Mouse ELMO2(Engulfment And Cell Motility 2) ELISA Kit
239	abk8323M	Mouse Ear 14(Eosinophil-associated, ribonuclease A family, member 14) ELISA Kit
240	abk7312M	Mouse ECFR(Eosinophil chemotactic factor receptor) ELISA Kit
241	abk5101M	Mouse MPIF2(Myeloid Progenitor Inhibitory Factor 2) ELISA Kit
242	abk6081M	Mouse EPOR(Erythropoietin receptor) ELISA Kit
243	abk1881M	Mouse UCP-2(Uncoupling Protein 2, Mitochondrial) ELISA Kit
244	abk6361M	Mouse ERα(Estrogen Receptor Alpha) ELISA Kit
245	abk6432M	Mouse EIF3A(Eukaryotic Translation Initiation Factor 3A) ELISA Kit
246	abk1222M	Mouse ECM1(Extracellular Matrix Protein 1) ELISA Kit
247	abk3222M	Mouse ERK1(Extracellular Signal Regulated Kinase 1) ELISA Kit
248	abk4222M	Mouse ERK2(Extracellular Signal Regulated Kinase 2) ELISA Kit
249	abk8301M	Mouse FAS/CD95(Factor Related Apoptosis) ELISA Kit
250	abk9181M	Mouse FSCN(Fascin) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
251	abk0161M	Mouse FABP5(Fatty Acid Binding Protein 5, Epidermal) ELISA Kit
252	abk6181M	Mouse FABP6(Fatty Acid Binding Protein 6, Ileal) ELISA Kit
253	abk6892M	Mouse FBXO2(F-Box Protein 2) ELISA Kit
254	abk5361M	Mouse ERβ(Estrogen Receptor Beta) ELISA Kit
255	abk2782M	Mouse FE(Ferritin) ELISA Kit
256	abk9322M	Mouse FBL(Fibrillarin) ELISA Kit
257	abk1332M	Mouse FBN1(Fibrillin 1) ELISA Kit
258	abk5822M	Mouse FDP(Fibrinogen Degradation Product) ELISA Kit
259	abk3422M	Mouse FGL2(Fibrinogen Like Protein 2) ELISA Kit
260	abk9261M	Mouse FGF9(Fibroblast Growth Factor 9) ELISA Kit
261	abk0361M	Mouse FGFBP1(Fibroblast Growth Factor Binding Protein 1) ELISA Kit
262	abk8103M	Mouse FRS2(Fibroblast Growth Factor Receptor Substrate 2) ELISA Kit
263	abk9323M	Mouse FRS3(Fibroblast Growth Factor Receptor Substrate 3) ELISA Kit
264	abk1322M	Mouse FMOD(Fibromodulin) ELISA Kit
265	abk2322M	Mouse FLRT2(Fibronectin Leucine Rich Transmembrane Protein 2) ELISA Kit
266	abk9801M	Mouse FN(Fibronectin) ELISA Kit
267	abk1113M	Mouse FLNa(Filamin A) ELISA Kit
268	abk7892M	Mouse Flt3(FMS-like Tyrosine Kinase 3) ELISA Kit
269	abk4302M	Mouse FAK(Focal Adhesion Kinase) ELISA Kit
270	abk9161M	Mouse FOXP3(Forkhead Box Protein P3) ELISA Kit
271	abk8581M	Mouse FPR2(Formyl Peptide Receptor 2) ELISA Kit
272	abk0423M	Mouse FLI1(Friend Leukemia Virus Integration 1) ELISA Kit
273	abk9411M	Mouse FUCA(Alpha-L-Fucosidase, Tissue) ELISA Kit
274	abk1423M	Mouse GNB1(G Protein Beta 1) ELISA Kit
275	abk2423M	Mouse GLB(Galactosidase, Beta) ELISA Kit
276	abk2471M	Mouse GALR1(Galanin Receptor 1) ELISA Kit
277	abk8851M	Mouse GAL1(Galectin 1) ELISA Kit
278	abk1951M	Mouse GAL2(Galectin 2) ELISA Kit
279	abk7581M	Mouse GAL3(Galectin 3) ELISA Kit
280	abk3951M	Mouse GAL4(Galectin 4) ELISA Kit
281	abk3423M	Mouse GAL6(Galectin 6) ELISA Kit
282	abk5951M	Mouse GAL7(Galectin 7) ELISA Kit
283	abk4951M	Mouse GAL8(Galectin 8) ELISA Kit
284	abk6951M	Mouse GAL9(Galectin 9) ELISA Kit
285	abk9851M	Mouse GAL12(Galectin 12) ELISA Kit
286	abk4423M	Mouse gGT1(Gamma Glutamyltransferase 1) ELISA Kit
287	abk5423M	Mouse GJB1(Gap Junction Protein Beta 1) ELISA Kit
288	abk0122M	Mouse GIP(Gastric Inhibitory Polypeptide) ELISA Kit
289	abk7022M	Mouse GT-34(Gastrin 34) ELISA Kit
290	abk5022M	Mouse GKN2(Gastrokine 2) ELISA Kit
291	abk8892M	Mouse GATA1(GATA Binding Protein 1) ELISA Kit
292	abk9892M	Mouse GATA2(GATA Binding Protein 2) ELISA Kit
293	abk9841M	Mouse GATA3(GATA Binding Protein 3) ELISA Kit
294	abk0941M	Mouse GATA4(GATA Binding Protein 4) ELISA Kit
295	abk0992M	Mouse GATA5(GATA Binding Protein 5) ELISA Kit
296	abk3453M	Mouse GS(Gelsolin) ELISA Kit
297	abk6423M	Mouse GTF II B(General Transcription Factor II B) ELISA Kit
298	abk7423M	Mouse GTF2H1(General Transcription Factor II H, Polypeptide 1) ELISA Kit
299	abk3182M	Mouse GFAP(glial fibrillary acidic protein) ELISA Kit
300	abk3671M	Mouse GR(glucocorticoid receptor) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
301	abk7502M	Mouse GCK(Glucokinase) ELISA Kit
302	abk5502M	Mouse G6PD(Glucose 6 Phosphate Dehydrogenase) ELISA Kit
303	abk1602M	Mouse GLUT3(Glucose Transporter 3) ELISA Kit
304	abk2602M	Mouse GLUT4(Glucose Transporter 4) ELISA Kit
305	abk6502M	Mouse GPI(Glucose 6 Phosphate Isomerase) ELISA Kit
306	abk7571M	Mouse GCLM(Glutamate Cysteine Ligase, Modifier Subunit) ELISA Kit
307	abk0972M	Mouse GAD2(Glutamate Decarboxylase 2, Acid) ELISA Kit
308	abk9872M	Mouse GDH/GLDH(Glutamate dehydrogenase) ELISA Kit
309	abk4453M	Mouse GR(Glutathione Reductase) ELISA Kit
310	abk8423M	Mouse GSTa1(Glutathione S Transferase Alpha 1) ELISA Kit
311	abk9423M	Mouse GSTo1(Glutathione S Transferase Omega 1) ELISA Kit
312	abk0523M	Mouse GSTt1(Glutathione S Transferase Theta 1) ELISA Kit
313	abk1523M	Mouse GSTt2(Glutathione S Transferase Theta 2) ELISA Kit
314	abk3203M	Mouse PYGB(Glycogen Phosphorylase, Brain) ELISA Kit
315	abk5471M	Mouse PYGL(Glycogen Phosphorylase, Liver) ELISA Kit
316	abk9281M	Mouse PYGM(Glycogen Phosphorylase, Muscle) ELISA Kit
317	abk7612M	Mouse GP5(Glycoprotein V, Platelet) ELISA Kit
318	abk0002M	Mouse GPC4(Glypican 4) ELISA Kit
319	abk7201M	Mouse GCP-2(Granulocyte Chemotactic Protein 2) ELISA Kit
320	abk3771M	Mouse GREM1(GremLin 1) ELISA Kit
321	abk1112M	Mouse GDF1(Growth Differentiation Factor 1) ELISA Kit
322	abk4112M	Mouse GDF3(Growth Differentiation Factor 3) ELISA Kit
323	abk2523M	Mouse GDF5(Growth Differentiation Factor 5) ELISA Kit
324	abk6112M	Mouse GDF7(Growth Differentiation Factor 7) ELISA Kit
325	abk7112M	Mouse GDF9(Growth Differentiation Factor 9) ELISA Kit
326	abk3112M	Mouse GDF11(Growth Differentiation Factor 11) ELISA Kit
327	abk4401M	Mouse GDF15(Growth Differentiation Factor 15) ELISA Kit
328	abk8212M	Mouse Grb2(Growth Factor Receptor Bound Protein 2) ELISA Kit
329	abk9212M	Mouse Grb7(Growth Factor Receptor Bound Protein 7) ELISA Kit
330	abk6212M	Mouse Grb10(Growth Factor Receptor Bound Protein 10) ELISA Kit
331	abk7212M	Mouse Grb14(Growth Factor Receptor Bound Protein 14) ELISA Kit
332	abk2461M	Mouse GHRH(Growth Hormone Releasing Hormone) ELISA Kit
333	abk5431M	Mouse CAM(Calmodulin) ELISA Kit
334	abk9972M	Mouse HP(Haptoglobin) ELISA Kit
335	abk3903M	Mouse HSF1(Heat Shock Factor 1) ELISA Kit
336	abk3523M	Mouse HSP-20(Heat Shock Protein 20) ELISA Kit
337	abk2802M	Mouse HSP-40(Heat Shock Protein 40) ELISA Kit
338	abk4523M	Mouse HSP-60(Heat Shock Protein 60) ELISA Kit
339	abk4802M	Mouse HSP-90(Heat Shock Protein 90) ELISA Kit
340	abk5802M	Mouse HSP gp96(Heat Shock Protein glycoprotein 96) ELISA Kit
341	abk6802M	Mouse HSF2(Heat Shock Transcription Factor 2) ELISA Kit
342	abk6992M	Mouse SHH(Sonic Hedgehog Protein) ELISA Kit
343	abk2941M	Mouse HHIP(Hedgehog Interacting Protein) ELISA Kit
344	abk2881M	Mouse HLTF(Helicase Like Transcription Factor) ELISA Kit
345	abk9381M	Mouse MMP-3(Matrix Metalloproteinase 3) ELISA Kit
346	abk9813M	Mouse POSTN/OSF-2(Periostin) ELISA Kit
347	abk5023M	Mouse PCSK9(Proprotein Convertase Subtilisin/Kexin Type 9) ELISA Kit
348	abk9371M	Mouse SCF(Stem Cell Factor) ELISA Kit
349	abk6941M	Mouse SELL(L-Selectin) ELISA Kit
350	abk3023M	Mouse SELP(P-Selectin) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
351	abk6453M	Mouse TPO(Thrombopoietin) ELISA Kit
352	abk6311M	Mouse ADIPOR1(Adiponectin Receptor 1) ELISA Kit
353	abk8542M	Mouse TNFRSF1B(Tumor Necrosis Factor Receptor Superfamily, Member 1B) ELISA Kit
354	abk1813M	Mouse RANKL(Receptor Activator of Nuclear Factor Kappa B Ligand) ELISA Kit
355	abk9291M	Mouse TSLP(Thymic Stromal Lymphopoietin) ELISA Kit
356	abk9913M	Mouse VEGFR1/FLT1(Vascular Endothelial Growth Factor Receptor 1) ELISA Kit
357	abk7242M	Mouse RANK(Receptor Activator of Nuclear Factor Kappa B) ELISA Kit
358	abk0042M	Mouse HPA(Heparanase) ELISA Kit
359	abk5523M	Mouse HC II(Heparin Cofactor II) ELISA Kit
360	abk6523M	Mouse HB-EGF(Heparin-binding Epidermal Growth Factor-like Growth Factor) ELISA Kit
361	abk7621M	Mouse LIPC(Lipase, Hepatic) ELISA Kit
362	abk7523M	Mouse HNF1b(Hepatocyte Nuclear Factor 1 Beta) ELISA Kit
363	abk8523M	Mouse HNF4a(Hepatocyte Nuclear Factor 4 Alpha) ELISA Kit
364	abk9523M	Mouse HNF6a(Hepatocyte Nuclear Factor 6 Alpha) ELISA Kit
365	abk0623M	Mouse HNF6b(Hepatocyte Nuclear Factor 6 Beta) ELISA Kit
366	abk8272M	Mouse HDGF(Hepatoma Derived Growth Factor) ELISA Kit
367	abk6372M	Mouse Hepc(Hepcidin) ELISA Kit
368	abk1623M	Mouse HNRPA1(Heterogeneous Nuclear Ribonucleoprotein A1) ELISA Kit
369	abk0012M	Mouse HK1(Hexokinase 1) ELISA Kit
370	abk5343M	Mouse histon-H2b ELISA Kit
371	abk2131M	Mouse Bcl-10(B-cell Lymphoma/Leukemia 10) ELIS Kit ELISA Kit
372	abk1922M	Mouse HAase(Hyaluronidase) ELISA Kit
373	abk7872M	Mouse HIF-1α(Hypoxia Inducible Factor 1 Alpha) ELISA Kit
374	abk3623M	Mouse ID(deiodinase) ELISA Kit
375	abk2972M	Mouse IgA(Immunoglobulin A) ELISA Kit
376	abk4623M	Mouse IgG2a(Immunoglobulin G2a) ELISA Kit
377	abk9141M	Mouse NOS2/iNOS(Nitric Oxide Synthase 2, Inducible) ELISA Kit
378	abk0611M	Mouse INHB(Inhibin B) ELISA Kit
379	abk5623M	Mouse INHbA(Inhibin Beta A) ELISA Kit
380	abk6623M	Mouse INHBP(Inhibin Binding Protein) ELISA Kit
381	abk7221M	Mouse INSL5(Insulin Like Protein 5) ELISA Kit
382	abk7972M	Mouse I-PTH(intact Parathormone) ELISA Kit
383	abk6161M	Mouse FOXO1(Forkhead Box Protein O1) ELISA Kit
384	abk0791M	Mouse UCP-1(Uncoupling Protein 1, Mitochondrial) ELISA Kit
385	abk8623M	Mouse FIL1h(Interleukin 1 Eta) ELISA Kit
386	abk9591M	Mouse IL-34(Interleukin 34) ELISA Kit
387	abk4501M	Mouse IL-3(Interleukin 3) ELISA Kit
388	abk9901M	Mouse IL-5(Interleukin 5) ELISA Kit
389	abk9531M	Mouse IL-7(Interleukin 7) ELISA Kit
390	abk0901M	Mouse IL-9(Interleukin 9) ELISA Kit
391	abk0452M	Mouse IL-11(Interleukin 11) ELISA Kit
392	abk6501M	Mouse IL-13(Interleukin 13) ELISA Kit
393	abk5111M	Mouse IL-15(Interleukin 15) ELISA Kit
394	abk3211M	Mouse IL-18(Interleukin 18) ELISA Kit
395	abk8501M	Mouse IL-23(Interleukin 23) ELISA Kit
396	abk3832M	Mouse IL-27(Interleukin 27) ELISA Kit
397	abk3642M	Mouse IL-35(Interleukin 35) ELISA Kit
398	abk0782M	Mouse IFABP/FABP2(Intestinal Fatty Acid Binding Protein) ELISA Kit
399	abk1791M	Mouse IL18BP(Interleukin 18 Binding Protein) ELISA Kit
400	abk6721M	Mouse AR(Amphiregulin) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
401	abk0613M	Mouse LTF/LF(Lactoferrin) ELISA Kit
402	abk1851M	Mouse LIFR(Leukemia Inhibitory Factor Receptor) ELISA Kit
403	abk9623M	Mouse PAL(L-Phenylalanine ammonia-lyase) ELISA Kit
404	abk7891M	Mouse Lptn/XCL1(Lymphotoctin) ELISA Kit
405	abk8802M	Mouse LZM(Lysozyme) ELISA Kit
406	abk3932M	Mouse FPN(Ferroportin) ELISA Kit
407	abk3922M	Mouse ITGB3(Integrin Beta 3) ELISA Kit
408	abk6913M	Mouse MIF(Macrophage Migration Inhibitory Factor) ELISA Kit
409	abk1353M	Mouse MHC I/H-2 I (Major Histocompatibility Complex I) ELISA Kit
410	abk2723M	Mouse MHC- II /H-2 II (Major Histocompatibility Complex- II) ELISA Kit
411	abk3723M	Mouse MHC-III/H-2III(Major Histocompatibility Complex-III) ELISA Kit
412	abk4972M	Mouse MMP-1(Matrix Metalloproteinase 1) ELISA Kit
413	abk8381M	Mouse MMP-2(Matrix Metalloproteinase 2) ELISA Kit
414	abk0481M	Mouse MMP-17(Matrix Metalloproteinase 17) ELISA Kit
415	abk1481M	Mouse MMP-24(Matrix Metalloproteinase 24) ELISA Kit
416	abk2481M	Mouse MMP-7(Matrix Metalloproteinase 7) ELISA Kit
417	abk3481M	Mouse MMP-8(Matrix Metalloproteinase 8) ELISA Kit
418	abk5842M	Mouse HIF-2α(Hypoxia Inducible Factor 2 Alpha) ELISA Kit
419	abk4723M	Mouse TRPV1(Transient Receptor Potential Cation Channel Subfamily V, Member 1) ELISA Kit
420	abk1101M	Mouse MCP-3(Monocyte Chemotactic Protein 3) ELISA Kit
421	abk3432M	Mouse MUC5AC(Mucin-5 subtype AC) ELISA Kit
422	abk4432M	Mouse MUC5B(Mucin-5 Subtype B) ELISA Kit
423	abk5723M	Mouse CHRM1(Cholinergic Receptor, Muscarinic 1) ELISA Kit
424	abk5041M	Mouse CHRM2(Cholinergic Receptor, Muscarinic 2) ELISA Kit
425	abk1801M	Mouse MBP(Myelin Basic Protein) ELISA Kit
426	abk4621M	Mouse ApoM(Apolipoprotein M) ELISA Kit
427	abk4051M	Mouse GAS(N-Acetylgalactosamine 6-Sulfatase) ELISA Kit
428	abk9941M	Mouse NAGase(N-Acetyl Beta-D-Glucosaminidase) ELISA Kit
429	abk8761M	Mouse NGF(Nerve growth factor) ELISA Kit
430	abk3012M	Mouse NCAM/CD56(Neural Cell Adhesion Molecule) ELISA Kit
431	abk2302M	Mouse NGB(Neuroglobin) ELISA Kit
432	abk4251M	Mouse Slit2(Slit Homolog 2) ELISA Kit
433	abk5012M	Mouse NTRK2(Neurotrophic Tyrosine Kinase Receptor Type 2) ELISA Kit
434	abk6012M	Mouse NTRK3(Neurotrophic Tyrosine Kinase Receptor Type 3) ELISA Kit
435	abk4182M	Mouse NT-4(Neurotrophin 4) ELISA Kit
436	abk0482M	Mouse NGAL(Neutrophil Gelatinase Associated Lipocalin) ELISA Kit
437	abk7213M	Mouse NOX4(Nicotinamide Adenine Dinucleotide Phosphate Oxidase 4) ELISA Kit
438	abk9751M	Mouse NFIL3(Nuclear Factor, Interleukin 3 Regulated) ELISA Kit
439	abk6971M	Mouse NRIP1(Nuclear Receptor Interacting Protein 1) ELISA Kit
440	abk3871M	Mouse NXF1(Nuclear RNA Export Factor 1) ELISA Kit
441	abk1971M	Mouse NUP35(Nucleoporin 35kDa) ELISA Kit
442	abk2971M	Mouse NUP37(Nucleoporin 37kDa) ELISA Kit
443	abk3971M	Mouse NUP50(Nucleoporin 50kDa) ELISA Kit
444	abk4971M	Mouse NUP85(Nucleoporin 85kDa) ELISA Kit
445	abk5971M	Mouse NUP98(Nucleoporin 98kDa) ELISA Kit
446	abk5871M	Mouse NUP107(Nucleoporin 107kDa) ELISA Kit
447	abk6871M	Mouse NUP133(Nucleoporin 133kDa) ELISA Kit
448	abk6723M	Mouse NUP153(Nucleoporin 153kDa) ELISA Kit
449	abk7871M	Mouse NUP155(Nucleoporin 155kDa) ELISA Kit
450	abk7723M	Mouse NUP160(Nucleoporin 160kDa) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
451	abk8871M	Mouse NUP188(Nucleoporin 188kDa) ELISA Kit
452	abk9871M	Mouse NUP205(Nucleoporin 205kDa) ELISA Kit
453	abk0971M	Mouse NUP214(Nucleoporin 214kDa) ELISA Kit
454	abk2282M	Mouse OLFM4(Olfactomedin 4) ELISA Kit
455	abk2912M	Mouse ITLN1(Intelectin 1/Omentin) ELISA Kit
456	abk2751M	Mouse OIT3(Oncoprotein Induced Transcript 3) ELISA Kit
457	abk6302M	Mouse ODC(Ornithine Decarboxylase) ELISA Kit
458	abk3342M	Mouse IκBα(Inhibitory Subunit of NF Kappa B Alpha) ELISA Kit
459	abk6671M	Mouse OC/BGP(Osteocalcin) ELISA Kit
460	abk4502M	Mouse OSCAR(Osteoclast Associated Receptor) ELISA Kit
461	abk4771M	Mouse ON(Osteonectin) ELISA Kit
462	abk0991M	Mouse pIκBα(Phospho-Inhibitory Subunit of NF Kappa B Alpha) ELISA Kit
463	abk0921M	Mouse MAPK14(Mitogen Activated Protein Kinase 14) ELISA Kit
464	abk2453M	Mouse PP(Pancreatic Polypeptide) ELISA Kit
465	abk9581M	Mouse PTHR2(Parathyroid Hormone Receptor 2) ELISA Kit
466	abk0932M	Mouse Pax(Paxillin) ELISA Kit
467	abk8051M	Mouse PDAP1(PDGFA Associated Protein 1) ELISA Kit
468	abk1022M	Mouse PP(Pepsin) ELISA Kit
469	abk2022M	Mouse PGA(Pepsinogen A) ELISA Kit
470	abk3022M	Mouse PGC(Pepsinogen C) ELISA Kit
471	abk8723M	Mouse PIN4(Peptidyl Prolyl Cis/Trans Isomerase NIMA Interacting 4 Protein) ELISA Kit
472	abk9771M	Mouse PMP3(Peroxisomal Membrane Protein 3) ELISA Kit
473	abk9703M	Mouse PLA1(Phospholipase A1) ELISA Kit
474	abk1861M	Mouse LDLR(Low Density Lipoprotein Receptor) ELISA Kit
475	abk9723M	Mouse PLCγ1(Phospholipase C Gamma 1) ELISA Kit
476	abk2902M	Mouse PEDF(Pigment Epithelium Derived Factor) ELISA Kit
477	abk9712M	Mouse PIR(Pirin) ELISA Kit
478	abk3361M	Mouse PACAPR1(Pituitary Adenylate Cyclase Activating Polypeptide Receptor 1) ELISA Kit
479	abk1823M	Mouse PACAP(Pituitary Adenylate Cyclase Activating Polypeptide) ELISA Kit
480	abk5091M	Mouse PGF(Placental Growth Factor) ELISA Kit
481	abk5922M	Mouse PKP2(Plakophilin 2) ELISA Kit
482	abk3322M	Mouse PAP(Plasmin-Antiplasmin Complex) ELISA Kit
483	abk7322M	Mouse tPA(Tissue-type Plasminogen Activator) ELISA Kit
484	abk3011M	Mouse PDGFRL(Platelet Derived Growth Factor Receptor Like Protein) ELISA Kit
485	abk2922M	Mouse GP II b/IIIa(Platelet Membrane Glycoprotein II b/IIIa complex) ELISA Kit
486	abk2823M	Mouse PDGFR-α(Platelet-Derived Growth Factor Soluble Receptor α) ELISA Kit
487	abk8961M	Mouse PTN(Pleiotrophin) ELISA Kit
488	abk3823M	Mouse PVRL1(Poliiovirus Receptor Related Protein 1) ELISA Kit
489	abk4823M	Mouse PVRL3(Poliiovirus Receptor Related Protein 3) ELISA Kit
490	abk7432M	Mouse PC4(activated RNA polymerase II transcription cofactor 4) ELISA Kit
491	abk5823M	Mouse Pro-ANP(Pro Atrial Natriuretic Peptide) ELISA Kit
492	abk5332M	Mouse PGR(Progesterone receptor) ELISA Kit
493	abk0232M	Mouse PHB(Prohibitin) ELISA Kit
494	abk6032M	Mouse PI(Proinsulin) ELISA Kit
495	abk5602M	Mouse PK2(Prokineticin 2) ELISA Kit
496	abk6602M	Mouse PKR1(Prokineticin Receptor 1) ELISA Kit
497	abk7552M	Mouse GPX1(Glutathione Peroxidase 1) ELISA Kit
498	abk6242M	Mouse PCNA(Proliferating Cell Nuclear Antigen) ELISA Kit
499	abk6823M	Mouse PGD2S(Prostaglandin D2 Synthase) ELISA Kit
500	abk8803M	Mouse EP2(Prostaglandin E Receptor 2) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
501	abk2702M	Mouse PTGES2(Prostaglandin E Synthase 2) ELISA Kit
502	abk7011M	Mouse PTGES(Prostaglandin E Synthase, Microsomal) ELISA Kit
503	abk4702M	Mouse PTGS2/COX-2(Prostaglandin Endoperoxide Synthase 2) ELISA Kit
504	abk7823M	Mouse PSA(Prostate Specific Antigen) ELISA Kit
505	abk4761M	Mouse PAR2(Protease Activated Receptor 2) ELISA Kit
506	abk1761M	Mouse PRSS1(Protease, Serine, 1) ELISA Kit
507	abk3561M	Mouse PDIA3(Protein Disulfide Isomerase A3) ELISA Kit
508	abk6761M	Mouse PDI(Protein Disulfide Isomerase) ELISA Kit
509	abk9251M	Mouse PIAS1(Protein Inhibitor Of Activated STAT 1) ELISA Kit
510	abk6561M	Mouse PICK1(Protein Interacting With Protein Kinase C, Alpha 1) ELISA Kit
511	abk5561M	Mouse PKBg(Protein Kinase B Gamma) ELISA Kit
512	abk8203M	Mouse PKB(Protein Kinase B) ELISA Kit
513	abk8823M	Mouse PKCb1(Protein Kinase C Beta 1) ELISA Kit
514	abk9823M	Mouse PKCd(Protein Kinase C Delta) ELISA Kit
515	abk0923M	Mouse PKCe(Protein Kinase C Epsilon) ELISA Kit
516	abk1923M	Mouse PKCg(Protein Kinase C Gamma) ELISA Kit
517	abk2923M	Mouse PKCi(Protein Kinase C Iota) ELISA Kit
518	abk3923M	Mouse PKCq(Protein Kinase C Theta) ELISA Kit
519	abk4923M	Mouse PKCz(Protein Kinase C Zeta) ELISA Kit
520	abk4661M	Mouse PKD1(Protein Kinase D1) ELISA Kit
521	abk9203M	Mouse PKD2(Protein Kinase D2) ELISA Kit
522	abk5923M	Mouse PKIa(Protein Kinase Inhibitor Alpha) ELISA Kit
523	abk6923M	Mouse PKIb(Protein Kinase Inhibitor Beta) ELISA Kit
524	abk7923M	Mouse PKIg(Protein Kinase Inhibitor Gamma) ELISA Kit
525	abk5661M	Mouse PKN1(Protein Kinase N1) ELISA Kit
526	abk8923M	Mouse PRKAa1(Protein Kinase, AMP Activated Alpha 1) ELISA Kit
527	abk5051M	Mouse PTPRO(Protein Tyrosine Phosphatase Receptor Type O) ELISA Kit
528	abk9232M	Mouse PCDH1(Protocadherin 1) ELISA Kit
529	abk2161M	Mouse PDP(Pyruvate Dehydrogenase Phosphatase) ELISA Kit
530	abk4471M	Mouse PK(Pyruvate Kinase, Liver And RBC) ELISA Kit
531	abk9923M	Mouse RAPGEF2(Rap Guanine Nucleotide Exchange Factor 2) ELISA Kit
532	abk0033M	Mouse RasGAP(Ras GTPase Activating Protein) ELISA Kit
533	abk4132M	Mouse REG4(Regenerating Islet Derived Protein 4) ELISA Kit
534	abk3992M	Mouse RGS4(Regulator Of G Protein Signaling 4) ELISA Kit
535	abk8412M	Mouse RXFP1(Relaxin/Insulin Like Family Peptide Receptor 1) ELISA Kit
536	abk1033M	Mouse RNASEA(Ribonuclease A) ELISA Kit
537	abk2033M	Mouse RNASEH(Ribonuclease H) ELISA Kit
538	abk2403M	Mouse RNASET2(Ribonuclease T2) ELISA Kit
539	abk0081M	Mouse RPL10(Ribosomal Protein L10) ELISA Kit
540	abk1081M	Mouse RPL26(Ribosomal Protein L26) ELISA Kit
541	abk5181M	Mouse RNF4(Ring Finger Protein 4) ELISA Kit
542	abk0251M	Mouse RUNX2(Runt Related Transcription Factor 2) ELISA Kit
543	abk2371M	Mouse S100A11(S100 Calcium Binding Protein A11) ELISA Kit
544	abk7371M	Mouse S100B(S100 Calcium Binding Protein B) ELISA Kit
545	abk5371M	Mouse S100A6(S100 Calcium Binding Protein A6) ELISA Kit
546	abk3902M	Mouse SDH(Sorbitol Dehydrogenase) ELISA Kit
547	abk1271M	Mouse sIgA(Secretory Immunoglobulin A) ELISA Kit
548	abk5061M	Mouse SEMA7A(Semaphorin 7A) ELISA Kit
549	abk2362M	Mouse SESN1(Sestrin 1) ELISA Kit
550	abk3033M	Mouse SESN2(Sestrin 2) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
551	abk4033M	Mouse SESN3(Sestrin 3) ELISA Kit
552	abk7522M	Mouse SOX3(Sex Determining Region Y Box Protein 3) ELISA Kit
553	abk1912M	Mouse SIGLEC3(Sialic Acid Binding Ig Like Lectin 3) ELISA Kit
554	abk0912M	Mouse SIGLEC12(Sialic Acid Binding Ig Like Lectin 12) ELISA Kit
555	abk5302M	Mouse OTC(Ornithine Transcarbamylase) ELISA Kit
556	abk6211M	Mouse ACTN2(Actinin Alpha 2) ELISA Kit
557	abk5033M	Mouse MYH7(Myosin Heavy Chain 7, Cardiac Muscle, Beta) ELISA Kit
558	abk4113M	Mouse SNRPD1(Small Nuclear Ribonucleoprotein D1) ELISA Kit
559	abk6033M	Mouse IL-1sR I (Soluble Interleukin-1 Receptor I) ELISA Kit
560	abk7033M	Mouse sMHC-1(soluble myosin heavy chain 1) ELISA Kit
561	abk8033M	Mouse sMHC-2(soluble myosin heavy chain 2) ELISA Kit
562	abk6291M	Mouse TRAIL/TNFSF10(Tumor Necrosis Factor Related Apoptosis Inducing Ligand) ELISA Kit
563	abk5212M	Mouse SSTR1(Somatostatin Receptor 1) ELISA Kit
564	abk4851M	Mouse SPOP(Speckle Type POZ Protein) ELISA Kit
565	abk3503M	Mouse SSFA2(Sperm Specific Antigen 2) ELISA Kit
566	abk6002M	Mouse SCCA2(Squamous Cell Carcinoma Antigen 2) ELISA Kit
567	abk0471M	Mouse SCFR(Stem Cell Growth Factor Receptor) ELISA Kit
568	abk5771M	Mouse SF1(Steroidogenic Factor 1) ELISA Kit
569	abk9033M	Mouse SDF-1β(Stromal Cell Derived Factor 1 Beta) ELISA Kit
570	abk2581M	Mouse SDF2L1(Stromal Cell Derived Factor 2 Like Protein 1) ELISA Kit
571	abk7781M	Mouse SSRP1(Structure Specific Recognition Protein 1) ELISA Kit
572	abk6222M	Mouse SOCS3(Suppressors Of Cytokine Signaling 3) ELISA Kit
573	abk8842M	Mouse CYGB(Cytoglobin) ELISA Kit
574	abk2812M	Mouse SYP(Synaptophysin) ELISA Kit
575	abk0133M	Mouse SYT1(Synaptotagmin 1) ELISA Kit
576	abk2891M	Mouse SDC4(Syndecan 4) ELISA Kit
577	abk9412M	Mouse TACR1(Tachykinin Receptor 1) ELISA Kit
578	abk0512M	Mouse TACR2(Tachykinin Receptor 2) ELISA Kit
579	abk3972M	Mouse TLN(Talin) ELISA Kit
580	abk4351M	Mouse THP(Tamm–Horsfall Glycoprotein) ELISA Kit
581	abk7992M	Mouse TDP43(TAR DNA Binding Protein 43) ELISA Kit
582	abk5351M	Mouse TBP(TATA Binding Protein) ELISA Kit
583	abk6351M	Mouse MAPτ(Microtubule Associated Protein Tau/Tau Protein) ELISA Kit
584	abk2181M	Mouse TXLNa(Taxilin Alpha) ELISA Kit
585	abk3181M	Mouse TXLNb(Taxilin Beta) ELISA Kit
586	abk1133M	Mouse TEK(Tyrosine Kinase, Endothelial) ELISA Kit
587	abk2801M	Mouse TE(telomerase) ELISA Kit
588	abk3281M	Mouse TNX(Tenascin X) ELISA Kit
589	abk2133M	Mouse TDG1(Teratocarcinoma Derived Growth Factor 1) ELISA Kit
590	abk8042M	Mouse sC5b-9(Soluble Terminal Complement Complex C5b-9) ELISA Kit
591	abk1102M	Mouse TBP2(Thioredoxin Binding Protein 2) ELISA Kit
592	abk0102M	Mouse TrxR1(Thioredoxin Reductase 1) ELISA Kit
593	abk9002M	Mouse Trx(Thioredoxin) ELISA Kit
594	abk8322M	Mouse TAFI(Thrombin activatable fibrinolysis inhibitor) ELISA Kit
595	abk3133M	Mouse TR(Thrombin Receptor) ELISA Kit
596	abk4133M	Mouse THBS2(Thrombospondin 2) ELISA Kit
597	abk6353M	Mouse TM(Thrombomodulin) ELISA Kit
598	abk2822M	Mouse TpP(Thrombus Precursor Protein) ELISA Kit
599	abk0832M	Mouse RELN(Reelin) ELISA Kit
600	abk7453M	Mouse TG(Thyroglobulin) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
601	abk4513M	Mouse TSH(Thyroid Stimulating Hormone) ELISA Kit
602	abk1942M	Mouse Itgam/CD11b(Integrin alpha-M) ELISA Kit
603	abk0461M	Mouse TEF(Thyrotroph Embryonic Factor) ELISA Kit
604	abk3681M	Mouse TBG(Thyroxine Binding Globulin) ELISA Kit
605	abk3801M	Mouse TFPI(Tissue Factor Pathway Inhibitor) ELISA Kit
606	abk5453M	Mouse TF(Tissue factor) ELISA Kit
607	abk5481M	Mouse TIMP-2(Tissue Inhibitors of Metalloproteinase 2) ELISA Kit
608	abk6481M	Mouse TIMP-3(Tissue Inhibitors of Metalloproteinase 3) ELISA Kit
609	abk7481M	Mouse TIMP-4(Tissue Inhibitors of Metalloproteinase 4) ELISA Kit
610	abk8351M	Mouse TLR9(Toll-like Receptor 9) ELISA Kit
611	abk6133M	Mouse FBXO32(F-Box Protein 32) ELISA Kit
612	abk7313M	Mouse CD8(Cluster of Differentiation 8) ELISA Kit
613	abk7133M	Mouse MuRF1(muscle-specific RING-finger protein 1) ELISA Kit
614	abk6832M	Mouse TCF20(Transcription Factor 20) ELISA Kit
615	abk9442M	Mouse RELB(V-Rel Reticuloendotheliosis Viral Oncogene Homolog B) ELISA Kit
616	abk8133M	Mouse TMED1(Transmembrane Emp24 Domain Containing Protein 1) ELISA Kit
617	abk9832M	Mouse TNPO1(Transportin 1) ELISA Kit
618	abk5832M	Mouse TTR(Transthyretin) ELISA Kit
619	abk1902M	Mouse TFF2(Trefoil Factor 2) ELISA Kit
620	abk2261M	Mouse TFF3(Trefoil Factor 3, Intestinal) ELISA Kit
621	abk8392M	Mouse TNNI3/cTn-I(Troponin I Type 3, Cardiac) ELISA Kit
622	abk9133M	Mouse TUBb(Tubulin Beta) ELISA Kit
623	abk6071M	Mouse Ub(Ubiquitin) ELISA Kit
624	abk9302M	Mouse PLAU/uPA(Urokinase-Type Plasminogen Activator) ELISA Kit
625	abk0722M	Mouse VEGF-B(Vascular Endothelial cell Growth Factor B) ELISA Kit
626	abk1722M	Mouse VEGFR3/Flt4(Vascular Endothelial Cell Growth Factor Receptor 3) ELISA Kit
627	abk5762M	Mouse VCAM-1/CD106(Vascular cell adhesion molecule 1) ELISA Kit
628	abk5831M	Mouse FOS(V-Fos FBJ Murine Osteosarcoma Viral Oncogene Homolog) ELISA Kit
629	abk8202M	Mouse VASPIN(Visceral Adipose Specific Serine Protease Inhibitor) ELISA Kit
630	abk9202M	Mouse VF(Visfatin) ELISA Kit
631	abk3722M	Mouse VWF(Von Willebrand Factor) ELISA Kit
632	abk0233M	Mouse α1-MG(α1-Microglobulin) ELISA Kit
633	abk6451M	Mouse α2-AP(α2-Antiplasmin) ELISA Kit
634	abk1233M	Mouse α2-PI(α2-plasmin inhibitor) ELISA Kit
635	abk2353M	Mouse Gal(α-Galactosyl Epitopes) ELISA Kit
636	abk2233M	Mouse IDUA(Alpha-L-iduronidase) ELISA Kit
637	abk3233M	Mouse αNAG(αN-acetylglucosaminidase) ELISA Kit
638	abk0651M	Mouse β Manase(β mannosidase) ELISA Kit
639	abk4233M	Mouse β-Glu(Beta Glucosidase) ELISA Kit
640	abk5712M	Mouse GLTP(Glycolipid Transfer Protein) ELISA Kit
641	abk9642M	Mouse IL-21(Interleukin 21) ELISA Kit
642	abk0391M	Mouse PLAUR/uPAR(Plasminogen Activator, Urokinase Receptor) ELISA Kit
643	abk9502M	Mouse GLUT1(Glucose Transporter 1) ELISA Kit
644	abk6111M	Mouse 5-LO(Arachidonate 5-Lipoxygenase) ELISA Kit
645	abk5642M	Mouse 12-LO(Arachidonate 12-Lipoxygenase) ELISA Kit
646	abk1111M	Mouse 15-LO(Arachidonate 15-Lipoxygenase) ELISA Kit
647	abk1141M	Mouse NOS1/nNOS(Nitric Oxide Synthase 1, Neuronal) ELISA Kit
648	abk2402M	Mouse GS(Glutamine synthetase) ELISA Kit
649	abk0031M	Mouse ATF6(Activating Transcription Factor 6) ELISA Kit
650	abk5233M	Mouse SEMA3C(Semaphorin-3C) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
651	abk2712M	Mouse GSK3β(Glycogen Synthase Kinase 3 Beta) ELISA Kit
652	abk4552M	Mouse pMAPT/pTAU(phosphorylated microtubule-associated protein tau) ELISA Kit
653	abk2152M	Mouse SYN1(Synapsin 1) ELISA Kit
654	abk5811M	Mouse CPE(Carboxypeptidase E) ELISA Kit
655	abk9732M	Mouse CCNA1(Cyclin-A1) ELISA Kit
656	abk6233M	Mouse ANTXR2(Anthrax Toxin Receptor 2) ELISA Kit
657	abk7233M	Mouse CCNE1(Cyclin-E1) ELISA Kit
658	abk8233M	Mouse RARβ(Retinoic Acid Receptor Beta) ELISA Kit
659	abk4011M	Mouse LUM(Lumican) ELISA Kit
660	abk9233M	Mouse CCNE2(Cyclin-E2) ELISA Kit
661	abk0333M	Mouse PI16(Peptidase Inhibitor 16) ELISA Kit
662	abk1333M	Mouse TOP1(Topoisomerase I) ELISA Kit
663	abk1392M	Mouse BFP(Brain Finger Protein) ELISA Kit
664	abk2333M	Mouse S100A10(S100 Calcium-binding Protein A10) ELISA Kit
665	abk4142M	Mouse PLA2G2D(Phospholipase A2, Group II D) ELISA Kit
666	abk0092M	Mouse iPLA2(Phospholipase A2, Calcium Independent) ELISA Kit
667	abk0382M	Mouse SAA2(Serum Amyloid A2) ELISA Kit
668	abk6413M	Mouse RTN4(Reticulon 4) ELISA Kit
669	abk1413M	Mouse PS1(Presenilin 1) ELISA Kit
670	abk9272M	Mouse AgRP(Agouti Related Protein) ELISA Kit
671	abk3333M	Mouse G6PC(Glucose-6-Phosphatase, Catalytic) ELISA Kit
672	abk0861M	Mouse LDL(Low Density Lipoprotein) ELISA Kit
673	abk3021M	Mouse FASN(Fatty Acid Synthase) ELISA Kit
674	abk2191M	Mouse HDL(High Density Lipoprotein) ELISA Kit
675	abk2303M	Mouse PYY(Peptide YY) ELISA Kit
676	abk4333M	Mouse IgG2b(Immunoglobulin G2b) ELISA Kit
677	abk5333M	Mouse IgG2c(Immunoglobulin G2c) ELISA Kit
678	abk8082M	Mouse IgG3(Immunoglobulin G3) ELISA Kit
679	abk8281M	Mouse FABP3(Fatty Acid Binding Protein 3, Muscle and Heart) ELISA Kit
680	abk7411M	Mouse FABP4(Fatty Acid Binding Protein 4, Adipocyte) ELISA Kit
681	abk8582M	Mouse αFP(Alpha-Fetoprotein) ELISA Kit
682	abk3701M	Mouse CSTA(Cystatin A) ELISA Kit
683	abk1151M	Mouse P-Cadherin(Cadherin, Placental) ELISA Kit
684	abk2451M	Mouse PCDH15(Protocadherin 15) ELISA Kit
685	abk2502M	Mouse PDPN(Podoplanin) ELISA Kit
686	abk3791M	Mouse NFκB2(Nuclear Factor Kappa B2) ELISA Kit
687	abk8012M	Mouse NPHN(Nephrin) ELISA Kit
688	abk0401M	Mouse AFGF/FGF1(Acidic Fibroblast Growth Factor 1) ELISA Kit
689	abk8261M	Mouse FGF23(Fibroblast Growth Factor 23) ELISA Kit
690	abk0811M	Mouse FGFR1(Fibroblast Growth Factor Receptor 1) ELISA Kit
691	abk8382M	Mouse TLR4(Toll-Like Receptor 4) ELISA Kit
692	abk6681M	Mouse PCT(Procalcitonin) ELISA Kit
693	abk3691M	Mouse TH(Tyrosine Hydroxylase) ELISA Kit
694	abk8311M	Mouse ADM(Adrenomedullin) ELISA Kit
695	abk5482M	Mouse CTSB(Cathepsin B) ELISA Kit
696	abk6381M	Mouse MMP-11(Matrix Metalloproteinase 11) ELISA Kit
697	abk2601M	Mouse EREG(Epiregulin) ELISA Kit
698	abk4451M	Mouse α1-AT(Alpha 1-Antitrypsin) ELISA Kit
699	abk2742M	Mouse MMP-14(Matrix Metalloproteinase 14) ELISA Kit
700	abk3482M	Mouse sCD14(Soluble Cluster of Differentiation 14) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
701	abk8101M	Mouse FcεR II /CD23(Receptor II for the Fc Region of Immunoglobulin E) ELISA Kit
702	abk2991M	Mouse PCK1(Phosphoenolpyruvate Carboxykinase 1) ELISA Kit
703	abk6333M	Mouse ACTa2(Actin Alpha 2, Smooth Muscle) ELISA Kit
704	abk0091M	Mouse SOST(Sclerostin) ELISA Kit
705	abk4742M	Mouse LCAT(Lecithin-Cholesterol Acyltransferase) ELISA Kit
706	abk3901M	Mouse PIIINP(N-Terminal Procollagen III Propeptide) ELISA Kit
707	abk7333M	Mouse CD163(Cluster of Differentiation 163) ELISA Kit
708	abk2482M	Mouse DPP4(Dipeptidyl Peptidase IV) ELISA Kit
709	abk5301M	Mouse EG-VEGF(Endocrine Gland Derived Vascular Endothelial Growth Factor) ELISA Kit
710	abk9401M	Mouse IL-1R2(Interleukin 1 Receptor Type II) ELISA Kit
711	abk1501M	Mouse TNFSF14(Tumor Necrosis Factor Ligand Superfamily, Member 14) ELISA Kit
712	abk2501M	Mouse MCSF(Macrophage Colony Stimulating Factor 1) ELISA Kit
713	abk7501M	Mouse IL-22(Interleukin 22) ELISA Kit
714	abk1453M	Mouse PL(Pancreatic Lipase) ELISA Kit
715	abk1091M	Mouse IRF5(Interferon Regulatory Factor 5) ELISA Kit
716	abk7591M	Mouse IL-12 p35(Interleukin 12 p35) ELISA Kit
717	abk8712M	Mouse FTH(Ferritin, Heavy Polypeptide) ELISA Kit
718	abk0011M	Mouse IL-6R(Interleukin 6 Receptor) ELISA Kit
719	abk4211M	Mouse IL-19(Interleukin 19) ELISA Kit
720	abk1411M	Mouse IL-20(Interleukin 20) ELISA Kit
721	abk3572M	Mouse ALCAM(Activated Leukocyte Cell Adhesion Molecule) ELISA Kit
722	abk8371M	Mouse SDC1(Syndecan 1) ELISA Kit
723	abk4192M	Mouse ApoC3(Apolipoprotein C3) ELISA Kit
724	abk5742M	Mouse CYP2E1(Cytochrome P450, family 2, subfamily E, polypeptide 1) ELISA Kit
725	abk6231M	Mouse MOG(Myelin Oligodendrocyte Glycoprotein) ELISA Kit
726	abk8981M	Mouse AZGP1(Zinc-alpha-2-glycoprotein) ELISA Kit
727	abk9331M	Mouse MUC1(Mucin 1) ELISA Kit
728	abk8333M	Mouse IL8Rα(Interleukin 8 Receptor Alpha) ELISA Kit
729	abk9333M	Mouse IL8Rβ(Interleukin 8 Receptor Beta) ELISA Kit
730	abk9042M	Mouse UCHL1(Ubiquitin Carboxyl Terminal Hydrolase L1) ELISA Kit
731	abk2932M	Mouse MAOB(Amine oxidase B) ELISA Kit
732	abk0433M	Mouse Slc18a2(Solute Carrier Family 18, Member 2) ELISA Kit
733	abk1433M	Mouse Slc18a1(Solute Carrier Family 18, Member 1) ELISA Kit
734	abk2433M	Mouse S1PR1(Sphingosine 1 Phosphate Receptor 1) ELISA Kit
735	abk3433M	Mouse S1PR2(Sphingosine 1 Phosphate Receptor 2) ELISA Kit
736	abk4433M	Mouse S1PR3(Sphingosine 1 Phosphate Receptor 3) ELISA Kit
737	abk5433M	Mouse S1PR4(Sphingosine 1 Phosphate Receptor 4) ELISA Kit
738	abk6433M	Mouse S1PR5(Sphingosine 1 Phosphate Receptor 5) ELISA Kit
739	abk7433M	Mouse DEFb14(beta-defensin 14) ELISA Kit
740	abk8433M	Mouse AK1(Adenylate Kinase 1) ELISA Kit
741	abk6281M	Mouse MLCK/MYLK(Myosin light chain kinase) ELISA Kit
742	abk3713M	Mouse γMSH(gamma-Melanocyte Stimulating Hormone) ELISA Kit
743	abk2572M	Mouse AChE(Acetylcholinesterase) ELISA Kit
744	abk2503M	Mouse MT1(Metallothionein 1) ELISA Kit
745	abk2443M	Mouse AMT(Aminomethyltransferase) ELISA Kit
746	abk7151M	Mouse RHOA(Ras Homolog Gene Family, Member A) ELISA Kit
747	abk1042M	Mouse TFPI2(Tissue Factor Pathway Inhibitor 2) ELISA Kit
748	abk8242M	Mouse IL-33(Interleukin 33) ELISA Kit
749	abk9262M	Mouse GDN(Glia Derived Nexin) ELISA Kit
750	abk8372M	Mouse sgp130(Soluble Glucoprotein 130) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
751	abk8122M	Mouse CVL/EZR(Cytovillin/Ezrin) ELISA Kit
752	abk0533M	Mouse sCD146(Soluble Cluster of Differentiation 146) ELISA Kit
753	abk5362M	Mouse STUB1(E3 ubiquitin-protein ligase CHIP) ELISA Kit
754	abk1533M	Mouse TITF1(Thyroid transcription factor 1) ELISA Kit
755	abk4531M	Mouse CNR1(Cannabinoid Receptor 1, Brain) ELISA Kit
756	abk2533M	Mouse CNR2(Cannabinoid Receptor 2) ELISA Kit
757	abk3353M	Mouse KARSL(Lysyl tRNA Synthetase) ELISA Kit
758	abk5982M	Mouse RARS(Arginyl tRNA Synthetase) ELISA Kit
759	abk4533M	Mouse IARS(Isoleucyl tRNA Synthetase) ELISA Kit
760	abk9182M	Mouse KISS1(Kisspeptin 1) ELISA Kit
761	abk5533M	Mouse ACV-AB(Activin AB) ELISA Kit
762	abk6533M	Mouse DNMT3L(DNA[cytosine-5]-methyltransferase 3-like) ELISA Kit
763	abk7272M	Mouse CMA1(Chymase 1, Mast Cell) ELISA Kit
764	abk7533M	Mouse LAMP1(Lysosomal Associated Membrane Protein 1) ELISA Kit
765	abk8533M	Mouse MMP-19(Matrix Metalloproteinase 19) ELISA Kit
766	abk3082M	Mouse EGFR2(Epidermal Growth Factor Receptor 2) ELISA Kit
767	abk9533M	Mouse AP-1(Transcription Factor AP-1) ELISA Kit
768	abk0633M	Mouse APC(Activated Protein C) ELISA Kit
769	abk0503M	Mouse CNTN4(Contactin 4) ELISA Kit
770	abk5271M	Mouse GPR120(G Protein Coupled Receptor 120) ELISA Kit
771	abk1633M	Mouse TBK1(TANK Binding Kinase 1) ELISA Kit
772	abk2633M	Mouse sAPPα(soluble amyloid precursor protein α) ELISA Kit
773	abk3221M	Mouse TARS(Threonyl tRNA Synthetase, cytoplasmic) ELISA Kit
774	abk3633M	Mouse CCL6(Chemokine C-C-Motif Ligand 6) ELISA Kit
775	abk9041M	Mouse NEF3(Neurofilament 3) ELISA Kit
776	abk5042M	Mouse NES1(Nesfatin 1) ELISA Kit
777	abk4633M	Mouse ABHD6(Monoacylglycerol Lipase ABHD6) ELISA Kit
778	abk5633M	Mouse RAB35(Ras-related protein Rab-35) ELISA Kit
779	abk5511M	Mouse AGRN(Agrin) ELISA Kit
780	abk6633M	Mouse LIPA(Lipase A, Lysosomal Acid) ELISA Kit
781	abk3512M	Mouse MyD88(Myeloid Differentiation Factor 88) ELISA Kit
782	abk0172M	Mouse TWIST1(Twist-related protein 1) ELISA Kit
783	abk2172M	Mouse SH3YL1(SH3 domain-containing YSC84-like protein 1) ELISA Kit
784	abk7633M	Mouse Nrf1(Nuclear Respiratory Factor 1) ELISA Kit
785	abk8633M	Mouse CTSG(Cathepsin G) ELISA Kit
786	abk9633M	Mouse TLR7(Toll Like Receptor 7) ELISA Kit
787	abk4762M	Mouse GRP78(Glucose Regulated Protein 78) ELISA Kit
788	abk0733M	Mouse RARγ(Retinoic Acid Receptor Gamma) ELISA Kit
789	abk1733M	Mouse ACTN1(Actinin Alpha 1) ELISA Kit
790	abk2733M	Mouse TFAM(Transcription Factor A, Mitochondrial) ELISA Kit
791	abk0871M	Mouse PPARγC1α(Peroxisome Proliferator Activated Receptor Gamma Coactivator 1 Alpha) ELISA Kit
792	abk3733M	Mouse NEO(Neogenin) ELISA Kit
793	abk6301M	Mouse ET-1(Endothelin 1) ELISA Kit
794	abk4733M	Mouse PPARγC1β(Peroxisome Proliferator Activated Receptor Gamma Coactivator 1 beta) ELISA Kit
795	abk6682M	Mouse IL-2(Interleukin 2) ELISA Kit
796	abk5501M	Mouse IL-4(Interleukin 4) ELISA Kit
797	abk1372M	Mouse MCP-1(Monocyte Chemotactic Protein 1) ELISA Kit
798	abk5733M	Mouse CCL12/MCP-5(Monocyte Chemotactic Protein 5) ELISA Kit
799	abk3101M	Mouse TARC(Thymus Activation Regulated Chemokine) ELISA Kit
800	abk3682M	Mouse IGF-1(Insulin-like Growth Factor 1) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
801	abk4682M	Mouse IGFBP-3(Insulin-like Growth Factor Binding Protein 3) ELISA Kit
802	abk9372M	Mouse LEP(Leptin) ELISA Kit
803	abk4292M	Mouse Aβ1-40(Amyloid Beta 1-40) ELISA Kit
804	abk5292M	Mouse Aβ1-42(Amyloid Beta 1-42) ELISA Kit
805	abk1701M	Mouse PRL(Prolactin) ELISA Kit
806	abk7472M	Mouse GLP-1(Glucagon Like Peptide 1) ELISA Kit
807	abk1023M	Mouse ACE(Angiotensin I Converting Enzyme) ELISA Kit
808	abk5382M	Mouse ANG1(Angiopietin 1) ELISA Kit
809	abk6733M	Mouse AMH(Anti-Mullerian Hormone) ELISA Kit
810	abk6601M	Mouse ApoA1(Apolipoprotein A1) ELISA Kit
811	abk7521M	Mouse ApoB(Apolipoprotein B) ELISA Kit
812	abk2511M	Mouse AGER(Total Advanced Glycosylation End Product Specific Receptor) ELISA Kit
813	abk7672M	Mouse BMP-6(Bone Morphogenetic Protein 6) ELISA Kit
814	abk8472M	Mouse C4(Complement Component 4) ELISA Kit
815	abk8772M	Mouse CFH(Complement Factor H) ELISA Kit
816	abk4641M	Mouse NTX I (Cross Linked N-Telopeptide of Type I Collagen) ELISA Kit
817	abk3641M	Mouse CTX I (Cross Linked C-telopeptide of Type I Collagen) ELISA Kit
818	abk7733M	Mouse Cys-C(Cystatin C) ELISA Kit
819	abk8903M	Mouse NE/ELA2(Neutrophil Elastase/Elastase-2) ELISA Kit
820	abk2101M	Mouse ECF/CCL11(Eosinophil Chemotactic Factor) ELISA Kit
821	abk5892M	Mouse SELE(E-selectin) ELISA Kit
822	abk8681M	Mouse GDNF(Glial Cell Line Derived Neurotrophic Factor) ELISA Kit
823	abk2092M	Mouse HB(Hemoglobin) ELISA Kit
824	abk0771M	Mouse OPN(Osteopontin) ELISA Kit
825	abk6722M	Mouse HO1(Heme Oxygenase 1) ELISA Kit
826	abk5282M	Mouse ALB(Albumin) ELISA Kit
827	abk5401M	Mouse HGF(Hepatocyte Growth Factor) ELISA Kit
828	abk4282M	Mouse IgE(Immunoglobulin E) ELISA Kit
829	abk7082M	Mouse IgG1(Immunoglobulin G1) ELISA Kit
830	abk6803M	Mouse IgM(Immunoglobulin M) ELISA Kit
831	abk0572M	Mouse KIM-1(Kidney Injury Molecule 1) ELISA Kit
832	abk8913M	Mouse LEPR(Leptin Receptor) ELISA Kit
833	abk5322M	Mouse PAI1(Plasminogen Activator Inhibitor 1) ELISA Kit
834	abk5603M	Mouse RBP4(Retinol Binding Protein 4) ELISA Kit
835	abk9113M	Mouse SEPP1(Selenoprotein P) ELISA Kit
836	abk0213M	Mouse SAA(Serum Amyloid A) ELISA Kit
837	abk7582M	Mouse SDF-1/CXCL12(Stromal Cell Derived Factor 1) ELISA Kit
838	abk1281M	Mouse TNC(Tenascin C) ELISA Kit
839	abk9271M	Mouse S100A8(S100 Calcium Binding Protein A8) ELISA Kit
840	abk0371M	Mouse S100A9(S100 Calcium Binding Protein A9) ELISA Kit
841	abk6482M	Mouse FABP1(Fatty Acid Binding Protein 1, Liver) ELISA Kit
842	abk8733M	Mouse Klotho(Klotho) ELISA Kit
843	abk1472M	Mouse LH(Luteinizing Hormone) ELISA Kit
844	abk9382M	Mouse IFNα(Interferon Alpha) ELISA Kit
845	abk3861M	Mouse RETN(Resistin) ELISA Kit
846	abk7382M	Mouse ADP/Acrp30(Adiponectin) ELISA Kit
847	abk2023M	Mouse RANTES(Regulated On Activation, Normal T-Cell Expressed and Secreted) ELISA Kit
848	abk8201M	Mouse IP-10/CXCL10(Interferon Gamma Induced Protein 10kDa) ELISA Kit
849	abk1301M	Mouse Cyt-C(Cytochrome C) ELISA Kit
850	abk2301M	Mouse DKK1(Dickkopf Related Protein 1) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
851	abk3301M	Mouse EGF(Epidermal Growth Factor) ELISA Kit
852	abk9301M	Mouse FASL/TNFSF6(Factor Related Apoptosis Ligand) ELISA Kit
853	abk0682M	Mouse Flt3L(FMS-like Tyrosine Kinase 3 Ligand) ELISA Kit
854	abk1682M	Mouse G-CSF(Granulocyte Colony-stimulating Factor) ELISA Kit
855	abk8401M	Mouse IL-1ra/IL-1F3(Interleukin 1 Receptor Antagonist) ELISA Kit
856	abk6813M	Mouse IL-6(Interleukin 6) ELISA Kit
857	abk0542M	Mouse IL-24(Interleukin 24) ELISA Kit
858	abk4601M	Mouse REN(Renin) ELISA Kit
859	abk3472M	Mouse OxLDL(Oxidized Low Density Lipoprotein) ELISA Kit
860	abk9682M	Mouse TnTf/TNNT3(Troponin T, fast skeletal muscle) ELISA Kit
861	abk8003M	Mouse NSE(Neuron Specific Enolase) ELISA Kit
862	abk5711M	Mouse ANG-R-Tie1(Angiopietin Receptor Tie1) ELISA Kit
863	abk9842M	Mouse CDKN1A(Cyclin Dependent Kinase Inhibitor 1A) ELISA Kit
864	abk6092M	Mouse AT-III(Antithrombin III) ELISA Kit
865	abk5192M	Mouse ApoB100(Apolipoprotein B100) ELISA Kit
866	abk9733M	Mouse AQP-0(Aquaporin 0) ELISA Kit
867	abk7192M	Mouse ARG1(Arginase I) ELISA Kit
868	abk3601M	Mouse Bcl-2(B-cell Lymphoma/Leukemia 2) ELISA Kit
869	abk6292M	Mouse Bcl-3(B-cell Lymphoma/Leukemia 3) ELIS Kit ELISA Kit
870	abk5131M	Mouse BAX(Bcl-2 Associated X Protein) ELISA Kit
871	abk4131M	Mouse BAD(Bcl2 Associated Death Promoter) ELISA Kit
872	abk0331M	Mouse BRAK/CXCL14(Breast and Kidney Expressed Chemokine) ELISA Kit
873	abk3392M	Mouse CGRP1(Calcitonin Gene Related Peptide 1) ELISA Kit
874	abk2431M	Mouse CACYBP(Calcylin Binding Protein) ELISA Kit
875	abk8431M	Mouse CR(Calretinin) ELISA Kit
876	abk6531M	Mouse CA2(Carbonic Anhydrase II) ELISA Kit
877	abk1572M	Mouse P I CP(Procollagen I C-Terminal Propeptide) ELISA Kit
878	abk8672M	Mouse CEA(Carcinoembryonic Antigen) ELISA Kit
879	abk5782M	Mouse P I NP(Procollagen I N-Terminal Propeptide) ELISA Kit
880	abk9392M	Mouse CT1(Cardiotrophln 1) ELISA Kit
881	abk8631M	Mouse CAD(Caspase Activated Deoxyribonuclease) ELISA Kit
882	abk0731M	Mouse CTNNα1(Catenin, Alpha 1) ELISA Kit
883	abk1731M	Mouse CTNNβ1(Catenin, Beta 1) ELISA Kit
884	abk2492M	Mouse Cav-1(Caveolin-1) ELISA Kit
885	abk7492M	Mouse CENPI(Centromere Protein I) ELISA Kit
886	abk0833M	Mouse CXCL15(Chemokine C-X-C-Motif Ligand 15) ELISA Kit
887	abk4372M	Mouse CXCL16(Chemokine C-X-C-Motif Ligand 16) ELISA Kit
888	abk3931M	Mouse CKLFSF4(Chemokine Like Factor Superfamily 4) ELISA Kit
889	abk3041M	Mouse CETP(Cholesteryl Ester Transfer Protein) ELISA Kit
890	abk8713M	Mouse CD19(Cluster of Differentiation 19) ELISA Kit
891	abk2232M	Mouse SEMA4D(Semaphorin 4D) ELISA Kit
892	abk1833M	Mouse F2(Coagulation Factor II) ELISA Kit
893	abk2382M	Mouse F8(Coagulation Factor VIII) ELISA Kit
894	abk3241M	Mouse F9(Coagulation Factor IX) ELISA Kit
895	abk4241M	Mouse F10(coagulation factor X) ELISA Kit
896	abk6241M	Mouse F12(Coagulation Factor XII) ELISA Kit
897	abk5772M	Mouse C3(Complement Component 3) ELISA Kit
898	abk0792M	Mouse C1qR1(Complement Component 1q Receptor 1) ELISA Kit
899	abk8901M	Mouse C5a(Complement Component 5a) ELISA Kit
900	abk8403M	Mouse CKMB(Creatine Kinase MB Isoenzyme) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
901	abk7881M	Mouse CRYαB(Crystallin Alpha B) ELISA Kit
902	abk8641M	Mouse PIII(Procollagen III C-Terminal ProPeptide) ELISA Kit
903	abk3892M	Mouse CXCR3(CXC-Chemokine Receptor 3) ELISA Kit
904	abk2892M	Mouse CREB1(Cyclic AMP Response Element Binding Protein 1) ELISA Kit
905	abk5741M	Mouse MYC(Myc proto-oncogene protein) ELISA Kit
906	abk0732M	Mouse CDK5(Cyclin Dependent Kinase 5) ELISA Kit
907	abk9013M	Mouse CCND1(Cyclin D1) ELISA Kit
908	abk0113M	Mouse CCND2(Cyclin D2) ELISA Kit
909	abk2833M	Mouse CCND3(Cyclin D3) ELISA Kit
910	abk8013M	Mouse CYP7A1(Cholesterol 7-Alpha-Hydroxylase) ELISA Kit
911	abk3833M	Mouse DNASE1(Deoxyribonuclease 1) ELISA Kit
912	abk4071M	Mouse DNM1(Dynamin 1) ELISA Kit
913	abk4833M	Mouse EPCAM(Epithelial Cell Adhesion Molecule) ELISA Kit
914	abk5833M	Mouse CXCL5(Chemokine C-X-C-Motif Ligand 5) ELISA Kit
915	abk1422M	Mouse FGα(Fibrinogen Alpha) ELISA Kit
916	abk2422M	Mouse FGβ(Fibrinogen Beta) ELISA Kit
917	abk6282M	Mouse FG(Fibrinogen) ELISA Kit
918	abk5682M	Mouse FGF7/KGF(Fibroblast Growth Factor 7) ELISA Kit
919	abk1203M	Mouse FSH(Follicle Stimulating Hormone) ELISA Kit
920	abk8902M	Mouse GMFB(Glia Maturation Factor, Beta) ELISA Kit
921	abk0712M	Mouse GRα(Glucocorticoid Receptor Alpha) ELISA Kit
922	abk1712M	Mouse GRβ(Glucocorticoid Receptor Beta) ELISA Kit
923	abk0602M	Mouse GLUT2(Glucose Transporter 2) ELISA Kit
924	abk3651M	Mouse GUSβ(β-glucuronidase) ELISA Kit
925	abk3712M	Mouse GSK3α(Glycogen Synthase Kinase 3 Alpha) ELISA Kit
926	abk8391M	Mouse GzmA(Granzyme A) ELISA Kit
927	abk9391M	Mouse GzmB(Granzyme B) ELISA Kit
928	abk1913M	Mouse GDF2(Growth Differentiation Factor 2) ELISA Kit
929	abk2112M	Mouse GDF10(Growth Differentiation Factor 10) ELISA Kit
930	abk1802M	Mouse HSP-27/HSPB1(Heat Shock Protein 27) ELISA Kit
931	abk2903M	Mouse HSP-70/HSPA9(Heat Shock Protein 70) ELISA Kit
932	abk6191M	Mouse PDGF-AB(Platelet Derived Growth Factor AB) ELISA Kit
933	abk4603M	Mouse PDGF-BB(Platelet Derived Growth Factor BB) ELISA Kit
934	abk8291M	Mouse TREM-1(Triggering Receptor Expressed on Myeloid Cells-1) ELISA Kit
935	abk3391M	Mouse VEGFR-2/KDR(Vascular Endothelial Growth Factor Receptor 2) ELISA Kit
936	abk7722M	Mouse HO2(Heme Oxygenase 2) ELISA Kit
937	abk1903M	Mouse HJV(Hemojuvelin) ELISA Kit
938	abk9713M	Mouse HNF1α(Hepatocyte Nuclear Factor 1 Alpha) ELISA Kit
939	abk9503M	Mouse HMGB-1(High Mobility Group Protein B1) ELISA Kit
940	abk5552M	Mouse HDAC1(Histone Deacetylase 1) ELISA Kit
941	abk6833M	Mouse HABP(Hyaluronate binding protein) ELISA Kit
942	abk9252M	Mouse HABP1/C1QBP(Hyaluronan Binding Protein 1) ELISA Kit
943	abk9102M	Mouse Fcy(Fc Fragment of IgG) ELISA Kit
944	abk3782M	Mouse IgG(Immunoglobulin G) ELISA Kit
945	abk9782M	Mouse INHA(Inhibin Alpha Chain) ELISA Kit
946	abk2521M	Mouse INSR(Insulin Receptor) ELISA Kit
947	abk8092M	Mouse IGFBP-1(Insulin-like Growth Factor Binding Protein 1) ELISA Kit
948	abk9572M	Mouse IGFBP-2(Insulin-like Growth Factor Binding Protein 2) ELISA Kit
949	abk0192M	Mouse IGFBP-4(Insulin-like Growth Factor Binding Protein 4) ELISA Kit
950	abk7833M	Mouse ICAM-2/CD102(Intercellular Adhesion Molecule 2) ELIS Kit

SR. NO.	AbChek Cat No.	Product Name
951	abk8833M	Mouse FIL1d(Interleukin 1 Delta) ELISA Kit
952	abk3003M	Mouse LAT(Linker for activation of T-cell) ELISA Kit
953	abk3313M	Mouse LpPLA2(Lipoprotein-associated Phospholipase A2) ELISA
954	abk0392M	Mouse IBSP(Bone Sialoprotein 2) ELISA Kit
955	abk0082M	Mouse MIP-3β/ELC/CCL19(Macrophage Inflammatory Protein 3β) ELISA Kit
956	abk9833M	Mouse HAT1(Histone Acetyltransferase 1) ELISA Kit
957	abk1103M	Mouse MlgA(Membrane Immunoglobulin A) ELISA Kit
958	abk0472M	Mouse MAU(Microalbuminuria) ELISA Kit
959	abk6251M	Mouse Smad4/MADH4(Mothers Against Decapentaplegic Homolog 4) ELISA Kit
960	abk0933M	Mouse ATPase(ATPase, Na+/K+ Transporting) ELISA Kit
961	abk0141M	Mouse NEFL(Neurofilament, Light Polypeptide) ELISA Kit
962	abk2032M	Mouse NOX1(Nicotinamide Adenine Dinucleotide Phosphate Oxidase 1) ELISA Kit
963	abk4813M	Mouse NT-proBNP(N-terminal pro-Brain Natriuretic Peptide) ELISA Kit
964	abk3403M	Mouse NFKB-p105(Nuclear factor NF-kappa-B p105 subunit) ELISA Kit
965	abk4403M	Mouse NFKB-p65(Nuclear Factor Kappa B p65) ELISA Kit
966	abk5003M	Mouse OCT1(Octamer-binding protein 1) ELISA Kit
967	abk4751M	Mouse OCT2(Octamer-binding protein 2) ELISA Kit
968	abk8132M	Mouse OSM(Oncostatin M) ELISA Kit
969	abk0203M	Mouse PRF1(Perforin 1) ELISA Kit
970	abk2871M	Mouse PPAR-γ(Peroxisome Proliferator Activated Receptor Gamma) ELISA Kit
971	abk1871M	Mouse PPAR-α(Peroxisome Proliferator Activated Receptor Alpha) ELISA Kit
972	abk7703M	Mouse pERK1/2(Phospho Extracellular Signal Regulated Kinase 1/2) ELISA Kit
973	abk6991M	Mouse PI3K2α(Phosphoinositide-3-Kinase Class-2-Alpha Polypeptide) ELISA Kit
974	abk8703M	Mouse PLA2P(Phospholipase A2 Activating Protein) ELISA Kit
975	abk6322M	Mouse PAI2(Plasminogen Activator Inhibitor 2) ELISA Kit
976	abk4213M	Mouse PDGF-C(Platelet Derived Growth Factor C) ELISA Kit
977	abk5213M	Mouse PDGF-D(Platelet Derived Growth Factor D) ELISA Kit
978	abk3703M	Mouse PDGF-A(Platelet Derived Growth Factor Subunit A) ELISA
979	abk4703M	Mouse PDGF-B(Platelet Derived Growth Factor Subunit B) ELISA
980	abk3591M	Mouse PECAM1/CD31(Platelet/Endothelial Cell Adhesion Molecule 1) ELISA Kit
981	abk9822M	Mouse PHLDA2(Pleckstrin Homology Like Domain Family A, Member 2) ELISA Kit
982	abk7932M	Mouse PCX(Podocalyxin) ELISA Kit
983	abk1933M	Mouse PVRL2(Poliiovirus Receptor Related Protein 2) ELISA Kit
984	abk9651M	Mouse POMC(Pro-Opiomelanocortin) ELISA Kit
985	abk3702M	Mouse PTGS1(Prostaglandin Endoperoxide Synthase 1) ELISA
986	abk9803M	Mouse ACP3/PAP(Prostatic Acid Phosphatase) ELISA Kit
987	abk3051M	Mouse PTPRN(Protein Tyrosine Phosphatase Receptor Type N) ELISA Kit
988	abk5303M	Mouse SP-A1(Pulmonary surfactant-associated protein A1) ELISA
989	abk7303M	Mouse SP-D(Pulmonary surfactant-associated protein D) ELISA
990	abk1161M	Mouse PDK4(Pyruvate Dehydrogenase Kinase Isozyme 4) ELISA
991	abk4103M	Mouse PDHA1(Pyruvate dehydrogenase E1 component subunit alpha) ELISA Kit
992	abk4872M	Mouse PKM2(Pyruvate Kinase, Muscle) ELISA Kit
993	abk2933M	Mouse SPLA2(Secreted Phospholipase A2) ELISA Kit
994	abk9903M	Mouse STK39(Serine/Threonine Protein Kinase 39) ELISA Kit
995	abk0282M	Mouse SHBG(Sex Hormone Binding Globulin) ELISA Kit
996	abk9422M	Mouse SRP9(Signal Recognition Particle 9kDa) ELISA Kit
997	abk1522M	Mouse STAT1(Signal Transducer and Activator of Transcription 1)
998	abk2522M	Mouse STAT2(Signal Transducer and Activator of Transcription 2) ELISA Kit
999	abk3522M	Mouse STAT4(Signal Transducer and Activator of Transcription 4) ELISA Kit
1000	abk5113M	Mouse STAT5A(Signal Transducer and Activator of Transcription 5A) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
1001	abk4522M	Mouse STAT6(Signal Transducer and Activator of Transcription 6) ELISA Kit
1002	abk5251M	Mouse Smad1/MADH1(Mothers Against Decapentaplegic Homolog 1) ELISA Kit
1003	abk7251M	Mouse Smad7/MADH7(Mothers Against Decapentaplegic Homolog 7) ELISA Kit
1004	abk3933M	Mouse sCD86(Soluble Cluster of Differentiation 86) ELISA Kit
1005	abk9313M	Mouse SEPCR/PROCR(Soluble Endothelial Protein C Receptor) ELISA Kit
1006	abk4933M	Mouse IL-2R α /CD25(Interleukin-2 Receptor alpha chain) ELISA
1007	abk5933M	Mouse IL-2R β /CD122(Interleukin-2 Receptor beta chain) ELISA Kit
1008	abk5992M	Mouse S100A1(S100 Calcium Binding Protein A1) ELISA Kit
1009	abk6203M	Mouse SREBP(Sterol Regulatory Element Binding Protein) ELISA Kit
1010	abk1581M	Mouse SDF-2(Stromal Cell Derived Factor 2) ELISA Kit
1011	abk3581M	Mouse SDF-4(Stromal Cell Derived Factor 4) ELISA Kit
1012	abk3202M	Mouse SMN1(Survival Of Motor Neuron 1, Telomeric) ELISA Kit
1013	abk6933M	Mouse SNAP-25(Synaptosome Associated Protein 25) ELISA Kit
1014	abk2551M	Mouse SNC α (Synuclein Alpha) ELISA Kit
1015	abk8503M	Mouse TRACP-5b(Tartrate Resistant Acid Phosphatase 5b) ELISA Kit
1016	abk8992M	Mouse TAF1(Transcription initiation factor TFIID subunit 1) ELISA Kit
1017	abk9992M	Mouse TAF2(Transcription initiation factor TFIID subunit 2) ELISA Kit
1018	abk7933M	Mouse TAF4(Transcription initiation factor TFIID subunit 4) ELISA Kit
1019	abk7351M	Mouse TIEG1(TGF Beta Inducible Early Response Gene 1) ELISA Kit
1020	abk3291M	Mouse TSP-1(Thrombospondin-1) ELISA Kit
1021	abk4402M	Mouse TAT(Thrombin/Antithrombin Complex) ELISA Kit
1022	abk5932M	Mouse CALP(Calprotectin) ELISA Kit
1023	abk1282M	Mouse TECK/CCL25(Thymus Expressed Chemokine) ELISA Kit
1024	abk4881M	Mouse TJP1(Tight junction protein ZO-1) ELISA Kit
1025	abk9472M	Mouse TF(Transferrin) ELISA Kit
1026	abk8832M	Mouse TFR2(Transferrin Receptor 2) ELISA Kit
1027	abk6313M	Mouse TFR1(Transferrin Receptor 1) ELISA Kit
1028	abk6691M	Mouse TGF β I(Transforming Growth Factor Beta Induced) ELISA Kit
1029	abk9933M	Mouse TGF- α (Transforming Growth Factor Alpha) ELISA Kit
1030	abk6603M	Mouse TGF- β 2(Transforming Growth Factor Beta 2) ELISA Kit
1031	abk0043M	Mouse TGF- β 3(Transforming Growth Factor Beta 3) ELISA Kit
1032	abk1832M	Mouse TGM2(Transglutaminase 2, Tissue) ELISA Kit
1033	abk1043M	Mouse TUBa3C(Tubulin Alpha 3c) ELISA Kit
1034	abk3013M	Mouse TPPP(Tubulin Polymerization Promoting Protein) ELISA Kit
1035	abk8532M	Mouse TNF- β (Tumor Necrosis Factor Beta) ELISA Kit
1036	abk7051M	Mouse TP53(Tumor Protein p53) ELISA Kit
1037	abk5313M	Mouse TSGF(Tumor Specific Growth Factor/Tumor Supplied Group of Factor) ELISA Kit
1038	abk2791M	Mouse TYR(Tyrosinase) ELISA Kit
1039	abk1391M	Mouse VEGF-C(Vascular Endothelial Growth Factor C) ELISA Kit
1040	abk2391M	Mouse VEGF-D(Vascular Endothelial Growth Factor D) ELISA Kit
1041	abk7482M	Mouse AGP1(Alpha-1-Acid Glycoprotein 1) ELISA Kit
1042	abk5913M	Mouse α -Glu(Alpha Glucosidase) ELISA Kit
1043	abk3043M	Mouse β -Hex A(β -hexosaminidase A) ELISA Kit
1044	abk4043M	Mouse β TG/PBP/CXCL7/NAP2(Thromboglobulin, Beta) ELISA Kit
1045	abk5112M	Mouse GDF6(Growth Differentiation Factor 6) ELISA Kit
1046	abk9632M	Mouse CDK4(Cyclin Dependent Kinase 4) ELISA Kit
1047	abk5043M	Mouse IFN- γ R1(Interferon Gamma Receptor 1) ELISA Kit
1048	abk9171M	Mouse SPARCL1(SPARC Like Protein 1) ELISA Kit
1049	abk1601M	Mouse VEGF-A(Vascular Endothelial Cell Growth Factor A) ELISA Kit
1050	abk9481M	Mouse TNNI2(Troponin I Type 2, Fast Skeletal) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
1051	abk6043M	MouseMAP1LC3A(Microtubule-associated proteins 1A/1B light chain 3A) ELISA Kit
1052	abk7043M	Mouse MAP1LC3B(Microtubule-associated proteins 1A/1B light chain 3B) ELISA Kit
1053	abk9003M	Mouse cPLA2(Cytosolic Phospholipase A2) ELISA Kit
1054	abk5403M	Mouse β MSH(beta-Melanocyte Stimulating Hormone) ELISA Kit
1055	abk0642M	Mouse T α 1(Thymosin- α 1) ELISA Kit
1056	abk1642M	Mouse PACAP-38(Pituitary Adenylate Cyclase Activating Polypeptide 38) ELISA Kit
1057	abk2642M	Mouse PACAP-27(Pituitary Adenylate Cyclase Activating Polypeptide 27) ELISA Kit
1058	abk3513M	Mouse INS(Insulin) ELISA Kit
1059	abk8531M	Mouse cTnT/TNNT2(Troponin T Type 2, Cardiac) ELISA Kit
1060	abk1003M	Mouse TLR-2(Toll-like Receptor 2) ELISA Kit
1061	abk6903M	Mouse NT-3(Neurotrophin 3) ELISA Kit
1062	abk3982M	Mouse ANGPTL3(Angiopoietin Like Protein 3) ELISA Kit
1063	abk9092M	Mouse IGFBP-5(Insulin-like Growth Factor Binding Protein 5) ELISA Kit
1064	abk8043M	Mouse LRP-1(Low-Density Lipoprotein-Receptor-Related Protein 1) ELISA Kit
1065	abk4203M	Mouse MAOA(Monoamine Oxidase Type A) ELISA Kit
1066	abk1692M	Mouse CHGA(Chromogranin A) ELISA Kit
1067	abk4803M	Mouse PIBF(Progesterone Induced Blocking Factor) ELISA Kit
1068	abk5911M	Mouse CYP1A2(Cytochrome P450, family 1, subfamily A, polypeptide 2) ELISA Kit
1069	abk2241M	Mouse F7(Coagulation Factor VII) ELISA Kit
1070	abk7282M	Mouse LBP(Lipopolysaccharide Binding Protein) ELISA Kit
1071	abk0143M	Mouse IFN α 4(Interferon Alpha 4) ELISA Kit
1072	abk2513M	Mouse ALPL(Alkaline Phosphatase, Liver/Bone/Kidney) ELISA Kit
1073	abk0672M	Mouse ApoA2(Apolipoprotein A2) ELISA Kit
1074	abk6382M	Mouse Irisin ELISA Kit
1075	abk5701M	Mouse IL-12(Interleukin 12) ELISA Kit
1076	abk5513M	Mouse EPO(Erythropoietin) ELISA Kit
1077	abk7401M	Mouse IL-1 α (Interleukin 1 Alpha) ELISA Kit
1078	abk0422M	Mouse GHR(Growth Hormone Receptor) ELISA Kit
1079	abk1592M	Mouse CCL1(Chemokine C-C-Motif Ligand 1) ELISA Kit
1080	abk2143M	Mouse Pro-MMP-9(Pro-Matrix Metalloproteinase-9) ELISA Kit
1081	abk1192M	Mouse IGFBP-6(Insulin-like Growth Factor Binding Protein 6) ELISA Kit
1082	abk4301M	Mouse EGFR(Epidermal Growth Factor Receptor) ELISA Kit
1083	abk1121M	Mouse FETUA(Fetuin A) ELISA Kit
1084	abk3582M	Mouse CLU(Clusterin) ELISA Kit
1085	abk3143M	Mouse gp130(Glucoprotein 130) ELISA Kit
1086	abk5203M	Mouse MCP-2/CCL8(Monocyte Chemotactic Protein 2) ELISA Kit
1087	abk3872M	Mouse Mly/CXCL9(Monocyte Interferon Gamma Inducing Factor) ELISA Kit
1088	abk8572M	Mouse IGF-2(Insulin Like Growth Factor 2) ELISA Kit
1089	abk6023M	Mouse PDGF-AA(Platelet Derived Growth Factor AA) ELISA Kit
1090	abk7023M	Mouse PF4(Platelet Factor 4) ELISA Kit
1091	abk0701M	Mouse ACTH(Adrenocorticotrophic Hormone) ELISA Kit
1092	abk5311M	Mouse ADH(Antidiuretic Hormone) ELISA Kit
1093	abk1031M	Mouse ANP(Atrial Natriuretic Peptide) ELISA Kit
1094	abk9131M	Mouse β -EP(Beta-Endorphin) ELISA Kit
1095	abk8231M	Mouse BK(Bradykinin) ELISA Kit
1096	abk3913M	Mouse BNP(Brain Natriuretic Peptide) ELISA Kit
1097	abk2203M	Mouse CRH(Corticotropin Releasing Hormone) ELISA Kit
1098	abk7701M	Mouse α -CTx(Alpha Crosslaps) ELISA Kit
1099	abk4331M	Mouse CNP(C-type natriuretic peptide) ELISA Kit
1100	abk0451M	Mouse β -CTx(Beta Crosslaps) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
1101	abk5812M	Mouse DHEA-S(Dehydroepiandrosterone sulfate) ELISA Kit
1102	abk6812M	Mouse DPD(Deoxypyridinoline) ELISA Kit
1103	abk1212M	Mouse GHRL(Ghrelin) ELISA Kit
1104	abk4121M	Mouse GC(Glucagon) ELISA Kit
1105	abk0312M	Mouse OXB(Orexin B) ELISA Kit
1106	abk0723M	Mouse LHRH(Luteinizing Hormone-Releasing Hormone) ELISA Kit
1107	abk4812M	Mouse MT(Melatonin) ELISA Kit
1108	abk2012M	Mouse NPY(Neuropeptide Y) ELISA Kit
1109	abk6213M	Mouse NADPH(Nicotinamide Adenine Dinucleotide Phosphate) ELISA Kit
1110	abk0751M	Mouse OXA(Orexin A) ELISA Kit
1111	abk6571M	Mouse OFQ/N(Orphanin FQ/Nociceptin) ELISA Kit
1112	abk5461M	Mouse SCT(Secretin) ELISA Kit
1113	abk4212M	Mouse SST(Somatostatin) ELISA Kit
1114	abk3822M	Mouse TXB2(Thromboxane B2) ELISA Kit
1115	abk9602M	Mouse PGF2α(Prostaglandin F2 Alpha) ELISA Kit
1116	abk6261M	Mouse SOD1(Superoxide Dismutase 1, Soluble) ELISA Kit
1117	abk3521M	Mouse APLN(Apelin) ELISA Kit
1118	abk4422M	Mouse SOD2(Superoxide Dismutase 2, Mitochondrial) ELISA Kit
1119	abk2142M	Mouse SOD3(Superoxide Dismutase 3, Extracellular) ELISA Kit
1120	abk5611M	Mouse IAPP(Islet Amyloid Polypeptide) ELISA Kit
1121	abk0132M	Mouse GLP-2(Glucagon Like Peptide 2) ELISA Kit
1122	abk9022M	Mouse GT(Gastrin) ELISA Kit
1123	abk1461M	Mouse TRH(Thyrotropin Releasing Hormone) ELISA Kit
1124	abk9213M	Mouse TAP(Trypsinogen Activation Peptide) ELISA Kit
1125	abk2043M	Mouse U-T4(Ultrasensitivity Thyroxine) ELISA Kit
1126	abk4622M	Mouse VIP(Vasoactive Intestinal Peptide) ELISA Kit
1127	abk2982M	Mouse Ang- I (Angiotensin I) ELISA Kit
1128	abk1982M	Mouse Ang- II (Angiotensin II) ELISA Kit
1129	abk2713M	Mouse Ang1-7(Angiotensin 1-7) ELISA Kit
1130	abk9043M	Mouse Ang1-9(Angiotensin 1-9) ELISA Kit
1131	abk9342M	Mouse αMSH(Alpha-Melanocyte Stimulating Hormone) ELISA Kit
1132	abk2213M	Mouse BMG/β2-MG(Beta-2-Microglobulin) ELISA Kit
1133	abk4582M	Mouse CRP(C-Reactive-Protein) ELISA Kit
1134	abk4692M	Mouse GRN(Granulin) ELISA Kit
1135	abk9513M	Mouse hs-CRP(high-sensitivity C-Reactive Protein) ELISA Kit
1136	abk8622M	Mouse IDO1 (Indoleamine 2,3-dioxygenase 1) ELISA kit
1137	abk1013M	Mouse MPO(Myeloperoxidase) ELISA Kit
1138	abk2682M	Mouse GM-CSF(Granulocyte-Macrophage Colony Stimulating Factor) ELISA Kit
1139	abk0601M	Mouse TNF-α(Tumor Necrosis Factor Alpha) ELISA Kit
1140	abk7682M	Mouse IL-17A(Interleukin 17A) ELISA Kit
1141	abk9501M	Mouse IFN-γ(Interferon Gamma) ELISA Kit
1142	abk4901M	Mouse TIMP-1(Tissue Inhibitors of Metalloproteinase 1) ELISA Kit
1143	abk9111M	Mouse IL-16(Interleukin 16) ELISA Kit
1144	abk6701M	Mouse IL-12 p40(Interleukin 12 p40) ELISA Kit
1145	abk6003M	Mouse IL-1R1(Interleukin 1 Receptor Type I) ELISA Kit
1146	abk7191M	Mouse PGRN(progranulin) ELISA Kit
1147	abk6972M	Mouse MMP-9(MatrixMetalloproteinase 9) ELISA Kit
1148	abk4111M	Mouse TNFRSF1A(Tumor Necrosis Factor Receptor Superfamily, Member 1A) ELISA Kit
1149	abk1382M	Mouse ICAM-1/CD54(Intercellular Adhesion Molecule 1) ELISA Kit
1150	abk7913M	Mouse ENG(Endoglin) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
1	abk5932R	Rat CALP(Calprotectin) ELISA Kit
2	abk6472R	Rat MYO(Myoglobin) ELISA Kit
3	abk9943R	Rat ICOSLG(Inducible T-Cell Co Stimulator Ligand) ELISA Kit
4	abk8453R	Rat NLRP3 (NACHT, LRR and PYD domains-containing protein 3) ELISA kit
5	abk1482R	Rat D2D(D-Dimer) ELISA Kit
6	abk1561R	Rat PROC(Protein C) ELISA Kit
7	abk8023R	Rat ACV-A(Activin A) ELISA Kit
8	abk3001R	Rat BMP-2(Bone Morphogenetic Protein 2) ELISA Kit
9	abk6201R	Rat GROα/CXCL1(Growth Regulated Oncogene Alpha) ELISA Kit
10	abk9461R	Rat ELN(Elastin) ELISA Kit
11	abk1301R	Rat Cyt-C(Cytochrome C) ELISA Kit
12	abk5513R	Rat EPO(Erythropoietin) ELISA Kit
13	abk2682R	Rat GM-CSF(Granulocyte-Macrophage Colony Stimulating Factor) ELISA Kit
14	abk7401R	Rat IL-1α(Interleukin 1 Alpha) ELISA Kit
15	abk5501R	Rat IL-4(Interleukin 4) ELISA Kit
16	abk0472R	Rat MAU(Microalbuminuria) ELISA Kit
17	abk1472R	Rat LH(Luteinizing Hormone) ELISA Kit
18	abk4601R	Rat REN(Renin) ELISA Kit
19	abk6801R	Rat LOX(Lysyl Oxidase) ELISA Kit
20	abk2841R	Rat COL1α1(Collagen Type I Alpha 1) ELISA Kit
21	abk4472R	Rat MMP-10(Matrix Metalloproteinase 10) ELISA Kit
22	abk8601R	Rat MMP-12(Matrix Metalloproteinase 12) ELISA Kit
23	abk5472R	Rat MMP-13(Matrix Metalloproteinase 13) ELISA Kit
24	abk9601R	Rat CT(Calcitonin) ELISA Kit
25	abk2701R	Rat PAPPA(Pregnancy Associated Plasma Protein A) ELISA Kit
26	abk8003R	Rat NSE(Neuron Specific Enolase) ELISA Kit
27	abk6701R	Rat IL-12 p40(Interleukin 12 p40) ELISA Kit
28	abk5701R	Rat IL-12(Interleukin 12) ELISA Kit
29	abk4721R	Rat AQP-3(Aquaporin 3) ELISA Kit
30	abk5721R	Rat AQP-4(Aquaporin 4) ELISA Kit
31	abk7721R	Rat AVPR1A(Arginine Vasopressin Receptor 1A) ELISA Kit
32	abk5821R	Rat AhR(Aryl Hydrocarbon Receptor) ELISA Kit
33	abk7123R	Rat AIP(Aryl Hydrocarbon Receptor Interacting Protein) ELISA Kit
34	abk8123R	Rat ARNT2(Aryl Hydrocarbon Receptor Nuclear Translocator 2) ELISA Kit
35	abk8821R	Rat ASGR1(Asialoglycoprotein Receptor 1) ELISA Kit
36	abk1921R	Rat AST(Aspartate Aminotransferase) ELISA Kit
37	abk0962R	Rat TBX2(T-Box Protein 2) ELISA Kit
38	abk4921R	Rat ATM(Ataxia Telangiectasia Mutated) ELISA Kit
39	abk5921R	Rat ATXN1(Ataxin 1) ELISA Kit
40	abk8192R	Rat ATXN2(Ataxin 2) ELISA Kit
41	abk4031R	Rat AXIN2(Axis Inhibition Protein 2) ELISA Kit
42	abk0131R	Rat FGF6(Fibroblast Growth Factor 6) ELISA Kit
43	abk9261R	Rat FGF9(Fibroblast Growth Factor 9) ELISA Kit
44	abk3672R	Rat bFGF/FGF2(Basic Fibroblast Growth Factor) ELISA Kit
45	abk7031R	Rat BAFFR(B-Cell Activation Factor Receptor) ELISA Kit
46	abk2131R	Rat Bcl-10(B-cell Lymphoma/Leukemia 10) ELIS Kit ELISA Kit
47	abk6292R	Rat Bcl-3(B-cell Lymphoma/Leukemia 3) ELISA Kit
48	abk6131R	Rat BMF(BCL-2 Modifying Factor) ELISA Kit
49	abk4131R	Rat BAD(Bcl2 Associated Death Promoter) ELISA Kit
50	abk8292R	Rat BTC(Betacellulin) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
51	abk2111R	Rat BACE1(Beta-Site APP Cleaving Enzyme 1) ELISA Kit
52	abk0231R	Rat BACE2(Beta-Site APP Cleaving Enzyme 2) ELISA Kit
53	abk8413R	Rat OPRM1(Mu-type opioid receptor) ELISA Kit
54	abk4143R	Rat hs-INS(high-sensitivity Insulin) ELISA Kit
55	abk1231R	Rat BID(BH3 Interacting Domain Death Agonist) ELISA Kit
56	abk0301R	Rat BLC(B-Lymphocyte Chemoattractant) ELISA Kit
57	abk5672R	Rat BMP-1(Bone Morphogenetic Protein 1) ELISA Kit
58	abk4411R	Rat BMP-3(Bone Morphogenetic Protein 3) ELISA Kit
59	abk4001R	Rat BMP-4(Bone Morphogenetic Protein 4) ELISA Kit
60	abk6591R	Rat BMP-5(Bone Morphogenetic Protein 5) ELISA Kit
61	abk7672R	Rat BMP-6(Bone Morphogenetic Protein 6) ELISA Kit
62	abk5001R	Rat BMP-7(Bone Morphogenetic Protein 7) ELISA Kit
63	abk5231R	Rat BMPR2(Bone Morphogenetic Protein Receptor II) ELISA Kit
64	abk1001R	Rat ANG2(Angiopoietin 2) ELISA Kit
65	abk3282R	Rat VE-Cadherin(Vascular Endothelial Cadherin) ELISA Kit
66	abk2011R	Rat NCAD(Neural Cadherin) ELISA Kit
67	abk0431R	Rat CABIN1(Calcineurin Binding Protein 1) ELISA Kit
68	abk4392R	Rat CASR(Calcium Sensing Receptor) ELISA Kit
69	abk5143R	Rat CAMK2(Calcium/Calmodulin-Dependent Protein Kinase- II) ELISA Kit
70	abk4252R	Rat ALPI(Intestinal alkaline phosphatase 1) ELISA Kit
71	abk5431R	Rat CAM(Calmodulin) ELISA Kit
72	abk6431R	Rat CNX(Calnexin) ELISA Kit
73	abk7431R	Rat CRT(Calreticulin) ELISA Kit
74	abk8431R	Rat CR(Calretinin) ELISA Kit
75	abk9431R	Rat CALU(Calumenin) ELISA Kit
76	abk5531R	Rat CA1(Carbonic Anhydrase I) ELISA Kit
77	abk6531R	Rat CA2(Carbonic Anhydrase II) ELISA Kit
78	abk8672R	Rat CEA(Carcinoembryonic Antigen) ELISA Kit
79	abk8582R	Rat αFP(Alpha-Fetoprotein) ELISA Kit
80	abk0631R	Rat CLS(Cardiolipin Synthase) ELISA Kit
81	abk9392R	Rat CT1(Cardiotrophin 1) ELISA Kit
82	abk1631R	Rat CLCF1(Cardiotrophin Like Cytokine Factor 1) ELISA Kit
83	abk6143R	Rat CACT(Carnitine-Acylcarnitine Translocase) ELISA Kit
84	abk3631R	Rat COMP(Cartilage Oligomeric Matrix Protein) ELISA Kit
85	abk0582R	Rat CASP3(Caspase 3) ELISA Kit
86	abk6631R	Rat CASP4(Caspase 4) ELISA Kit
87	abk7631R	Rat CASP7(Caspase 7) ELISA Kit
88	abk9631R	Rat CASP9(Caspase 9) ELISA Kit
89	abk1492R	Rat CASP12(Caspase 12) ELISA Kit
90	abk4631R	Rat CASP14(Caspase 14) ELISA Kit
91	abk2731R	Rat CTSD(Cathepsin D) ELISA Kit
92	abk3731R	Rat CTSK(Cathepsin K) ELISA Kit
93	abk4731R	Rat CTSL(Cathepsin L) ELISA Kit
94	abk5731R	Rat CDX2(Caudal Type Homeobox Transcription Factor 2) ELISA Kit
95	abk2492R	Rat CAV-1(Caveolin-1) ELISA Kit
96	abk2831R	Rat CENPA(Centromere Protein A) ELISA Kit
97	abk8492R	Rat CENPB(Centromere Protein B) ELISA Kit
98	abk5492R	Rat CENPE(Centromere Protein E) ELISA Kit
99	abk9492R	Rat CEP110(Centrosomal Protein 110kDa) ELISA Kit
100	abk3813R	Rat CP(Ceruloplasmin) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
101	abk5831R	Rat FOS(V-Fos FBJ Murine Osteosarcoma Viral Oncogene Homolog) ELISA Kit
102	abk1931R	Rat CCR8(Chemokine C-C-Motif Receptor 8) ELISA Kit
103	abk3592R	Rat CCRL1(Chemokine C-C-Motif Receptor Like Protein 1) ELISA Kit
104	abk7143R	Rat SCD(Stearoyl Coenzyme A Desaturase) ELISA Kit
105	abk4592R	Rat CX3CR1(Chemokine C-X3-C-Motif Receptor 1) ELISA Kit
106	abk4372R	Rat CXCL16(Chemokine C-X-C-Motif Ligand 16) ELISA Kit
107	abk9672R	Rat CKLF(Chemokine Like Factor) ELISA Kit
108	abk3931R	Rat CKLFSF4(Chemokine Like Factor Superfamily 4) ELISA Kit
109	abk4931R	Rat CKLFSF5(Chemokine Like Factor Superfamily 5) ELISA Kit
110	abk5931R	Rat CKLFSF6(Chemokine Like Factor Superfamily 6) ELISA Kit
111	abk5592R	Rat CKLFSF7(Chemokine Like Factor Superfamily 7) ELISA Kit
112	abk2041R	Rat CCKAR(Cholecystokinin A Receptor) ELISA Kit
113	abk6511R	Rat AGT(Angiotensinogen) ELISA Kit
114	abk5041R	Rat CHRM2(Cholinergic Receptor, Muscarinic 2) ELISA Kit
115	abk0692R	Rat CHD3(Chromodomain Helicase DNA Binding Protein 3) ELISA Kit
116	abk8041R	Rat CHD5(Chromodomain Helicase DNA Binding Protein 5) ELISA Kit
117	abk2201R	Rat CNTF(Ciliary Neurotrophic Factor) ELISA Kit
118	abk3223R	Rat CROP(Cisplatin Resistance Associated Overexpressed Protein) ELISA Kit
119	abk5312R	Rat CS(Citrate Synthase) ELISA Kit
120	abk2141R	Rat C-jun(Proto-oncogene c-Jun) ELISA Kit
121	abk4223R	Rat CC17(Clara Cell Protein) ELISA Kit
122	abk3692R	Rat CD109(Cluster Of Differentiation 109) ELISA Kit
123	abk5223R	Rat CD3(Cluster of Differentiation 3) ELISA Kit
124	abk5141R	Rat CD30(Cluster of differentiation 30) ELISA Kit
125	abk7223R	Rat CD40L(Cluster Of Differentiation 40 Ligand) ELISA Kit
126	abk8143R	Rat CD6(Cluster of Differentiation 6) ELISA Kit
127	abk7313R	Rat CD8(Cluster of Differentiation 8) ELISA Kit
128	abk6141R	Rat CD97(Cluster of Differentiation 97) ELISA Kit
129	abk3582R	Rat CLU(Clusterin) ELISA Kit
130	abk4831R	Rat CFL1(Cofilin 1, Non-Muscle) ELISA Kit
131	abk1341R	Rat CFL2(Cofilin 2, Muscle) ELISA Kit
132	abk0323R	Rat COL2(Collagen Type II) ELISA Kit
133	abk0772R	Rat COL3(Collagen Type III) ELISA Kit
134	abk4341R	Rat COL1α2(Collagen Type I Alpha 2) ELISA Kit
135	abk7341R	Rat COL4α1(Collagen Type IV Alpha 1) ELISA Kit
136	abk9143R	Rat COL9α1(Collagen Type IX Alpha 1) ELISA Kit
137	abk1772R	Rat COL10(Collagen Type X) ELISA Kit
138	abk0243R	Rat Collagenase I ELISA Kit
139	abk6671R	Rat OC/BGP(Osteocalcin) ELISA Kit
140	abk9692R	Rat GCSFR(Colony Stimulating Factor Receptor, Granulocyte) ELISA Kit
141	abk3772R	Rat C1INH(Complement 1 Inhibitor) ELISA Kit
142	abk8472R	Rat C4(Complement Component 4) ELISA Kit
143	abk1792R	Rat C1R(Complement Component 1Q Receptor) ELISA Kit
144	abk1541R	Rat C9(Complement Component 9) ELISA Kit
145	abk7772R	Rat CFB(Complement Factor B) ELISA Kit
146	abk3541R	Rat C3a(Complement Component 3a) ELISA Kit
147	abk7901R	Rat C4a(Complement Component 4a) ELISA Kit
148	abk9541R	Rat CTGF(Connective Tissue Growth Factor) ELISA Kit
149	abk0641R	Rat CX26(Connexin 26) ELISA Kit
150	abk6792R	Rat CX31(Connexin 31) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
151	abk5703R	Rat CX37(Connexin 37) ELISA Kit
152	abk0891R	Rat CX40(Connexin 40) ELISA Kit
153	abk1891R	Rat CX43(Connexin 43) ELISA Kit
154	abk0781R	Rat CNTN1(Contactin 1) ELISA Kit
155	abk1781R	Rat CNTN2(Contactin 2) ELISA Kit
156	abk2781R	Rat CNTN3(Contactin 3) ELISA Kit
157	abk0381R	Rat CKM(Creatine Kinase, Muscle) ELISA Kit
158	abk4641R	Rat NTX I (Cross Linked N-Telopeptide of Type I Collagen) ELISA Kit
159	abk5631R	Rat CASP8(Caspase 8) ELISA Kit
160	abk1232R	Rat CRY1(Cryptochrome 1) ELISA Kit
161	abk7101R	Rat CTACK(Cutaneous T-Cell Attracting Chemokine) ELISA Kit
162	abk2732R	Rat CDK9(Cyclin Dependent Kinase 9) ELISA Kit
163	abk3732R	Rat CDKN2A(Cyclin Dependent Kinase Inhibitor 2A) ELISA Kit
164	abk9013R	Rat CCND1(Cyclin D1) ELISA Kit
165	abk0113R	Rat CCND2(Cyclin D2) ELISA Kit
166	abk6182R	Rat CYP(A) (Cyclophilin A) ELISA Kit
167	abk2312R	Rat CYP(B) (Cyclophilin B) ELISA Kit
168	abk3312R	Rat CYP(D) (Cyclophilin D) ELISA Kit
169	abk7733R	Rat Cys-C(Cystatin C) ELISA Kit
170	abk7851R	Rat CYSLTR2(Cysteiny Leukotriene Receptor 2) ELISA Kit
171	abk9122R	Rat CYB561(Cytochrome b-561) ELISA Kit
172	abk5742R	Rat CYP2E1(Cytochrome P450, family 2, subfamily E, polypeptide1) ELISA Kit
173	abk7222R	Rat CYTH1(Cytohesin 1) ELISA Kit
174	abk8222R	Rat CYTIP(Cytohesin Interacting Protein) ELISA Kit
175	abk1243R	Rat CVL(Cytovillin) ELISA Kit
176	abk5412R	Rat DAF(Decay Accelerating Factor) ELISA Kit
177	abk9132R	Rat DCN(Decorin) ELISA Kit
178	abk2243R	Rat D-5(Defensin-5) ELISA Kit
179	abk7001R	Rat CA9(Carbonic Anhydrase IX) ELISA Kit
180	abk3833R	Rat DNASE1(Deoxyribonuclease I) ELISA Kit
181	abk6781R	Rat Des(Desmin) ELISA Kit
182	abk2301R	Rat DKK1(Dickkopf Related Protein 1) ELISA Kit
183	abk9741R	Rat DKK2(Dickkopf Related Protein 2) ELISA Kit
184	abk1071R	Rat DLD(Dihydrolipoyl Dehydrogenase) ELISA Kit
185	abk2071R	Rat DLAT(Dihydrolipoyl Transacetylase) ELISA Kit
186	abk0841R	Rat D-LDH(D-Lactate Dehydrogenase) ELISA Kit
187	abk3243R	Rat DRD1(Dopamine Receptor D1) ELISA Kit
188	abk4243R	Rat DRD2(Dopamine Receptor D2) ELISA Kit
189	abk5323R	Rat DRD3(Dopamine Receptor D3) ELISA Kit
190	abk6323R	Rat DRD4(Dopamine Receptor D4) ELISA Kit
191	abk5961R	Rat DAT(Dopamine Transporter) ELISA Kit
192	abk4801R	Rat TM(Thrombin) ELISA Kit
193	abk2381R	Rat DMD(Dystrophin) ELISA Kit
194	abk8332R	Rat EGR1(Early Growth Response Protein 1) ELISA Kit
195	abk9332R	Rat EGR2(Early Growth Response Protein 2) ELISA Kit
196	abk0432R	Rat EGR3(Early Growth Response Protein 3) ELISA Kit
197	abk9502R	Rat GLUT1(Glucose Transporter 1) ELISA Kit
198	abk2572R	Rat AChE(Acetylcholinesterase) ELISA Kit
199	abk5101R	Rat MPlF2(Myeloid Progenitor Inhibitory Factor 2) ELISA Kit
200	abk0561R	Rat ELA4(Elastase 4) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
201	abk3032R	Rat ELOA(Elongin A) ELISA Kit
202	abk9402R	Rat EED(Embryonic Ectoderm Development) ELISA Kit
203	abk5282R	Rat ALB(Albumin) ELISA Kit
204	abk0302R	Rat ERN1(Endoplasmic Reticulum To Nucleus Signalling 1) ELISA Kit
205	abk5372R	Rat ES(Endostatin) ELISA Kit
206	abk7323R	Rat EL(Endothelial lipase) ELISA Kit
207	abk1241R	Rat NOS3/eNOS(Nitric Oxide Synthase 3, Endothelial) ELISA Kit
208	abk2101R	Rat ECF/CCL11(Eosinophil Chemotactic Factor) ELISA Kit
209	abk3301R	Rat EGF(Epidermal Growth Factor) ELISA Kit
210	abk9482R	Rat CASP1(Caspase 1) ELISA Kit
211	abk8301R	Rat FAS/CD95(Factor Related Apoptosis) ELISA Kit
212	abk5243R	Rat ACS(Fatty Acyl-CoA synthetase) ELISA Kit
213	abk5822R	Rat FDP(Fibrinogen Degradation Product) ELISA Kit
214	abk9323R	Rat FRS3(Fibroblast Growth Factor Receptor Substrate 3) ELISA Kit
215	abk2322R	Rat FLRT2(Fibronectin Leucine Rich Transmembrane Protein 2) ELISA Kit
216	abk4302R	Rat FAK(Focal Adhesion Kinase) ELISA Kit
217	abk3401R	Rat FST(Follistatin) ELISA Kit
218	abk5201R	Rat CX3CL1(Chemokine C-X3-C-Motif Ligand 1) ELISA Kit
219	abk3862R	Rat ATGL(Adipose Triglyceride Lipase) ELISA Kit
220	abk9411R	Rat FUCA(Alpha-L-Fucosidase, Tissue) ELISA Kit
221	abk1423R	Rat GNB1(G Protein Beta 1) ELISA Kit
222	abk3021R	Rat FASN(Fatty Acid Synthase) ELISA Kit
223	abk1471R	Rat GAL(Galanin) ELISA Kit
224	abk1951R	Rat GAL2(Galectin 2) ELISA Kit
225	abk7581R	Rat GAL3(Galectin 3) ELISA Kit
226	abk3951R	Rat GAL4(Galectin 4) ELISA Kit
227	abk5951R	Rat GAL7(Galectin 7) ELISA Kit
228	abk4951R	Rat GAL8(Galectin 8) ELISA Kit
229	abk6951R	Rat GAL9(Galectin 9) ELISA Kit
230	abk8651R	Rat γGT1(Gamma Glutamyltransferase 1) ELISA Kit
231	abk5423R	Rat GJb1(Gap Junction Protein Beta 1) ELISA Kit
232	abk0122R	Rat GIP(Gastric Inhibitory Polypeptide) ELISA Kit
233	abk7022R	Rat GT-34(Gastrin 34) ELISA Kit
234	abk8892R	Rat GATA1(GATA Binding Protein 1) ELISA Kit
235	abk9892R	Rat GATA2(GATA Binding Protein 2) ELISA Kit
236	abk9841R	Rat GATA3(GATA Binding Protein 3) ELISA Kit
237	abk0941R	Rat GATA4(GATA Binding Protein 4) ELISA Kit
238	abk0992R	Rat GATA5(GATA Binding Protein 5) ELISA Kit
239	abk3453R	Rat GS(Gelsolin) ELISA Kit
240	abk6243R	Rat GTF II B(General Transcription Factor II B) ELISA Kit
241	abk7243R	Rat G-6-Pase/G6PT/GSD(Glucose 6 Phosphate) ELISA Kit
242	abk8902R	Rat GMFB(Glia Maturation Factor, Beta) ELISA Kit
243	abk8681R	Rat GDNF(Glial Cell Line Derived Neurotrophic Factor) ELISA Kit
244	abk3671R	Rat GR(Glucocorticoid Receptor) ELISA Kit
245	abk7502R	Rat GCK(Glucokinase) ELISA Kit
246	abk3143R	Rat gp130(Glucoprotein 130) ELISA Kit
247	abk5502R	Rat G6PD(Glucose 6 Phosphate Dehydrogenase) ELISA Kit
248	abk1602R	Rat GLUT3(Glucose Transporter 3) ELISA Kit
249	abk2602R	Rat GLUT4(Glucose Transporter 4) ELISA Kit
250	abk7571R	Rat GCLM(Glutamate Cysteine Ligase, Modifier Subunit) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
251	abk0972R	Rat GAD2(Glutamate Decarboxylase 2, Acid) ELISA Kit
252	abk9872R	Rat GDH/GLDH(Glutamate Dehydrogenase) ELISA Kit
253	abk6403R	Rat GRIK2(Glutamate Receptor, Ionotropic, Kainate 2) ELISA Kit
254	abk8423R	Rat GSTA1(Glutathione S Transferase Alpha 1) ELISA Kit
255	abk9423R	Rat GSTO1(Glutathione S Transferase Omega 1) ELISA Kit
256	abk8243R	Rat GSTT1(Glutathione S Transferase Theta 1) ELISA Kit
257	abk9243R	Rat GSTT2(Glutathione S Transferase Theta 2) ELISA Kit
258	abk3203R	Rat PYGB(Glycogen Phosphorylase, Brain) ELISA Kit
259	abk9281R	Rat PYGM(Glycogen Phosphorylase, Muscle) ELISA Kit
260	abk7612R	Rat GP5(Glycoprotein V, Platelet) ELISA Kit
261	abk0343R	Rat GnRHR(Gonadotropin-Releasing Hormone Receptor) ELISA Kit
262	abk1682R	Rat G-CSF(Granulocyte Colony Stimulating Factor) ELISA Kit
263	abk3771R	Rat GREM1(GremLin 1) ELISA Kit
264	abk1112R	Rat GDF1(Growth Differentiation Factor 1) ELISA Kit
265	abk1913R	Rat GDF2(Growth Differentiation Factor 2) ELISA Kit
266	abk4112R	Rat GDF3(Growth Differentiation Factor 3) ELISA Kit
267	abk2523R	Rat GDF5(Growth Differentiation Factor 5) ELISA Kit
268	abk7112R	Rat GDF9(Growth Differentiation Factor 9) ELISA Kit
269	abk3112R	Rat GDF11(Growth Differentiation Factor 11) ELISA Kit
270	abk4401R	Rat GDF15(Growth Differentiation Factor 15) ELISA Kit
271	abk6212R	Rat GRB10(Growth Factor Receptor Bound Protein 10) ELISA Kit
272	abk7212R	Rat GRB14(Growth Factor Receptor Bound Protein 14) ELISA Kit
273	abk8212R	Rat GRB2(Growth Factor Receptor Bound Protein 2) ELISA Kit
274	abk2461R	Rat GHRH(Growth Hormone Releasing Hormone) ELISA Kit
275	abk9972R	Rat HP(Haptoglobin) ELISA Kit
276	abk1343R	Rat HSPA4(Heat Shock 70kDa Protein 4) ELISA Kit
277	abk3903R	Rat HSF1(Heat Shock Factor 1) ELISA Kit
278	abk2802R	Rat HSP-40(Heat Shock Protein 40) ELISA Kit
279	abk4523R	Rat HSP-60(Heat Shock Protein 60) ELISA Kit
280	abk4802R	Rat HSP-90(Heat Shock Protein 90) ELISA Kit
281	abk5802R	Rat HSP GP96(Heat Shock Protein Glycoprotein 96) ELISA Kit
282	abk2343R	Rat HSP47(Heat Shock Protein 47) ELISA Kit
283	abk6802R	Rat HSF2(Heat Shock Transcription Factor 2) ELISA Kit
284	abk6992R	Rat SHH(Sonic Hedgehog Protein) ELISA Kit
285	abk2941R	Rat HHIP(Hedgehog Interacting Protein) ELISA Kit
286	abk2881R	Rat HLTF(Helicase Like Transcription Factor) ELISA Kit
287	abk0042R	Rat HPA(Heparanase) ELISA Kit
288	abk5523R	Rat HC II (Heparin Cofactor II) ELISA Kit
289	abk7621R	Rat LIPC(Lipase, Hepatic) ELISA Kit
290	abk5401R	Rat HGF(Hepatocyte Growth Factor) ELISA Kit
291	abk3343R	Rat HNF1a(Hepatocyte Nuclear Factor 1 Alpha) ELISA Kit
292	abk8272R	Rat HDGF(Hepatoma Derived Growth Factor) ELISA Kit
293	abk6372R	Rat Hepc(Hepcidin) ELISA Kit
294	abk1623R	Rat HNRPA1(Heterogeneous Nuclear Ribonucleoprotein A1) ELISA Kit
295	abk0012R	Rat HK1(Hexokinase 1) ELISA Kit
296	abk2191R	Rat HDL(High Density Lipoprotein) ELISA Kit
297	abk9503R	Rat HMGB-1(High Mobility Group Protein B1) ELISA Kit
298	abk4343R	Rat HAT1(Histone Acetyltransferase 1) ELISA Kit
299	abk5343R	Rat Histon-H2b ELISA Kit
300	abk7872R	Rat HIF-1α(Hypoxia Inducible Factor 1 Alpha) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
301	abk6343R	Rat IDS(Iduronate Sulfatase) ELISA Kit
302	abk2413R	Rat HMG-CoAR(3-Hydroxy-3-Methylglutaryl Coenzyme A Reductase) ELISA Kit
303	abk4282R	Rat IgE(Immunoglobulin E) ELISA Kit
304	abk9141R	Rat NOS2/iNOS(Nitric Oxide Synthase 2, Inducible) ELISA Kit
305	abk6623R	Rat INHBP(Inhibin Binding Protein) ELISA Kit
306	abk7343R	Rat IKBα(Inhibitory Subunit of NF-κBα) ELISA Kit
307	abk3642R	Rat IL-35(Interleukin 35) ELISA Kit
308	abk8572R	Rat IGF-2(Insulin Like Growth Factor 2) ELISA Kit
309	abk7972R	Rat I-PTH(Intact Parathormone) ELISA Kit
310	abk7833R	Rat ICAM-2/CD102(Intercellular Adhesion Molecule 2) ELISA Kit
311	abk6401R	Rat IFN-β(Interferon Beta) ELISA Kit
312	abk8343R	Rat IFI16/p16(Interferon-Inducible Protein 16) ELISA Kit
313	abk9343R	Rat PEPCK(Phosphoenolpyruvate Carboxykinase) ELISA Kit
314	abk8833R	Rat FIL1d(Interleukin 1 Delta) ELISA Kit
315	abk8623R	Rat FIL1h(Interleukin 1 Eta) ELISA Kit
316	abk4742R	Rat LCAT(Lecithin-Cholesterol Acyltransferase) ELISA Kit
317	abk4501R	Rat IL-3(Interleukin 3) ELISA Kit
318	abk9901R	Rat IL-5(Interleukin 5) ELISA Kit
319	abk9531R	Rat IL-7(Interleukin 7) ELISA Kit
320	abk0452R	Rat IL-11(Interleukin 11) ELISA Kit
321	abk6501R	Rat IL-13(Interleukin 13) ELISA Kit
322	abk5111R	Rat IL-15(Interleukin 15) ELISA Kit
323	abk9111R	Rat IL-16(Interleukin 16) ELISA Kit
324	abk3211R	Rat IL-18(Interleukin 18) ELISA Kit
325	abk9642R	Rat IL-21(Interleukin 21) ELISA Kit
326	abk8501R	Rat IL-23(Interleukin 23) ELISA Kit
327	abk3832R	Rat IL-27(Interleukin 27) ELISA Kit
328	abk0782R	Rat IFABP/FABP2(Intestinal Fatty Acid Binding Protein) ELISA Kit
329	abk2531R	Rat MUC2(Mucin 2) ELISA Kit
330	abk6721R	Rat AR(Amphiregulin) ELISA Kit
331	abk9801R	Rat FN(Fibronectin) ELISA Kit
332	abk0861R	Rat LDL(Low Density Lipoprotein) ELISA Kit
333	abk9372R	Rat LEP(Leptin) ELISA Kit
334	abk0501R	Rat LIF(Leukemia Inhibitory Factor) ELISA Kit
335	abk9491R	Rat LOX-1(Lectin Like Oxidized Low Density Lipoprotein Receptor 1) ELISA Kit
336	abk6941R	Rat SELL(L-Selectin) ELISA Kit
337	abk7891R	Rat Lptn/XCL1(Lymphotactin) ELISA Kit
338	abk9001R	Rat MIP-1α(Macrophage Inflammatory Protein 1 Alpha) ELISA Kit
339	abk0101R	Rat MIP-1β(Macrophage Inflammatory Protein 1 Beta) ELISA Kit
340	abk5471R	Rat PYGL(Glycogen Phosphorylase, Liver) ELISA Kit
341	abk1582R	Rat MIP-3α(Macrophage Inflammatory Protein 3 Alpha) ELISA Kit
342	abk6913R	Rat MIF(Macrophage Migration Inhibitory Factor) ELISA Kit
343	abk2582R	Rat MDC/CCL22(Macrophage-Derived Chemokine) ELISA Kit
344	abk0443R	Rat MHC II /AgB II /H-1 II (Major Histocompatibility Complex II) ELISA Kit
345	abk1443R	Rat MHCIII/AgBIII/H-1III(Major Histocompatibility ComplexIII) ELISA Kit
346	abk3861R	Rat RETN(Resistin) ELISA Kit
347	abk2443R	Rat MAT(Matrilysin) ELISA Kit
348	abk4972R	Rat MMP-1(Matrix Metalloproteinase 1) ELISA Kit
349	abk8381R	Rat MMP-2(Matrix Metalloproteinase 2) ELISA Kit
350	abk9381R	Rat MMP-3(Matrix Metalloproteinase 3) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
351	abk0481R	Rat MMP-17(Matrix Metalloproteinase 17) ELISA Kit
352	abk1481R	Rat MMP-24(Matrix Metalloproteinase 24) ELISA Kit
353	abk2481R	Rat MMP-7(Matrix Metalloproteinase 7) ELISA Kit
354	abk7611R	Rat ANG(Angiogenin) ELISA Kit
355	abk5382R	Rat ANG1(Angiopoietin 1) ELISA Kit
356	abk4101R	Rat SLC(Secondary Lymphoid Tissue Chemokine) ELISA Kit
357	abk2503R	Rat MT1(Metallothionein 1) ELISA Kit
358	abk1372R	Rat MCP-1(Monocyte Chemotactic Protein 1) ELISA Kit
359	abk1101R	Rat MCP-3(Monocyte Chemotactic Protein 3) ELISA Kit
360	abk1801R	Rat MBP(Myelin Basic Protein) ELISA Kit
361	abk2511R	Rat AGER(Total Advanced Glycosylation End Product Specific Receptor) ELISA Kit
362	abk3443R	Rat PTEN(Phosphatase and Tensin Homolog) ELISA Kit
363	abk9941R	Rat NAGase(N-Acetyl Beta-D-Glucosaminidase) ELISA Kit
364	abk4051R	Rat GAS(N-Acetylgalactosamine 6-Sulfatase) ELISA Kit
365	abk8761R	Rat NGF(Nerve Growth Factor) ELISA Kit
366	abk2302R	Rat NGB(Neuroglobin) ELISA Kit
367	abk5012R	Rat NTRK2(Neurotrophic Tyrosine Kinase Receptor Type 2) ELISA Kit
368	abk6012R	Rat NTRK3(Neurotrophic Tyrosine Kinase Receptor Type 3) ELISA Kit
369	abk6903R	Rat NT-3(Neurotrophin 3) ELISA Kit
370	abk4182R	Rat NT-4(Neurotrophin 4) ELISA Kit
371	abk0482R	Rat NGAL(Neutrophil Gelatinase Associated Lipocalin) ELISA Kit
372	abk4443R	Rat NOX1(Nicotinamide Adenine Dinucleotide Phosphate Oxidase 1) ELISA Kit
373	abk5443R	Rat NOX4(Nicotinamide Adenine Dinucleotide Phosphate Oxidase 4) ELISA Kit
374	abk9751R	Rat NFIL3(Nuclear Factor, Interleukin 3 Regulated) ELISA Kit
375	abk6971R	Rat NRIP1(Nuclear Receptor Interacting Protein 1) ELISA Kit
376	abk3871R	Rat NXF1(Nuclear RNA Export Factor 1) ELISA Kit
377	abk5871R	Rat NUP107(Nucleoporin 107kDa) ELISA Kit
378	abk6723R	Rat NUP153(Nucleoporin 153kDa) ELISA Kit
379	abk7871R	Rat NUP155(Nucleoporin 155kDa) ELISA Kit
380	abk7723R	Rat NUP160(Nucleoporin 160kDa) ELISA Kit
381	abk8871R	Rat NUP188(Nucleoporin 188kDa) ELISA Kit
382	abk9871R	Rat NUP205(Nucleoporin 205kDa) ELISA Kit
383	abk0971R	Rat NUP214(Nucleoporin 214kDa) ELISA Kit
384	abk1971R	Rat NUP35(Nucleoporin 35kDa) ELISA Kit
385	abk2971R	Rat NUP37(Nucleoporin 37kDa) ELISA Kit
386	abk4971R	Rat NUP85(Nucleoporin 85kDa) ELISA Kit
387	abk6443R	Rat NUP198(Nucleoporin 198kDa) ELISA Kit
388	abk2912R	Rat ITLN1(Intelectin 1/Omentin) ELISA Kit
389	abk2751R	Rat OIT3(Oncoprotein Induced Transcript 3) ELISA Kit
390	abk1532R	Rat OSMR(Oncostatin M Receptor) ELISA Kit
391	abk6302R	Rat ODC(Ornithine Decarboxylase) ELISA Kit
392	abk4502R	Rat OSCAR(Osteoclast Associated Receptor) ELISA Kit
393	abk3502R	Rat OSF(Osteoclast Stimulating Factor 1) ELISA Kit
394	abk4771R	Rat ON(Osteonectin) ELISA Kit
395	abk0771R	Rat OPN(Osteopontin) ELISA Kit
396	abk0921R	Rat MAPK14(Mitogen Activated Protein Kinase 14) ELISA Kit
397	abk4381R	Rat PUMA(P53 Upregulated Modulator of Apoptosis) ELISA Kit
398	abk9581R	Rat PTHR2(Parathyroid Hormone Receptor 2) ELISA Kit
399	abk8051R	Rat PDAP1(Pdgfa Associated Protein 1) ELISA Kit
400	abk1022R	Rat PP(Pepsin) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
401	abk2303R	Rat PYY(Peptide YY) ELISA Kit
402	abk0203R	Rat PRF1(Perforin 1) ELISA Kit
403	abk7443R	Rat PMP3(Peroxisomal Membrane Protein 3) ELISA Kit
404	abk6991R	Rat PI3K2α(Phosphoinositide 3 Kinase 2 Alpha) ELISA Kit
405	abk9703R	Rat PLA1(Phospholipase A1) ELISA Kit
406	abk9723R	Rat PLCγ1(Phospholipase C Gamma 1) ELISA Kit
407	abk8443R	Rat pACCCase(Phosphorylated Acetyl-CoA Carboxylase) ELISA Kit
408	abk1572R	Rat P I CP(Procollagen I C-terminal Propeptide) ELISA Kit
409	abk5091R	Rat PGF(Placental Growth Factor) ELISA Kit
410	abk3322R	Rat PAP(Plasmin-Antiplasmin Complex) ELISA Kit
411	abk7322R	Rat tPA(Tissue-type Plasminogen Activator) ELISA Kit
412	abk4322R	Rat Plg(Plasminogen) ELISA Kit
413	abk4213R	Rat PDGF-C(Platelet Derived Growth Factor C) ELISA Kit
414	abk5213R	Rat PDGF-D(Platelet Derived Growth Factor D) ELISA Kit
415	abk3011R	Rat PDGFR(Platelet Derived Growth Factor Receptor Like Protein) ELISA Kit
416	abk2823R	Rat PDGFR-α(Platelet-Derived Growth Factor Soluble Receptor α) ELISA Kit
417	abk8961R	Rat PTN(Pleiotrophin) ELISA Kit
418	abk7932R	Rat PCX(Podocalyxin) ELISA Kit
419	abk3823R	Rat PVRL1(Poliavirus Receptor Related Protein 1) ELISA Kit
420	abk1933R	Rat PVRL2(Poliavirus Receptor Related Protein 2) ELISA Kit
421	abk7432R	Rat PC4(Positive Cofactor 4) ELISA Kit
422	abk5332R	Rat PGR(Progesterone Receptor) ELISA Kit
423	abk5801R	Rat FSTL1(Follistatin Like Protein 1) ELISA Kit
424	abk0232R	Rat PHB(Prohibitin) ELISA Kit
425	abk6032R	Rat PI(Proinsulin) ELISA Kit
426	abk5602R	Rat PK2(Prokineticin 2) ELISA Kit
427	abk6602R	Rat PKR1(Prokineticin Receptor 1) ELISA Kit
428	abk6242R	Rat PCNA(Proliferating Cell Nuclear Antigen) ELISA Kit
429	abk9221R	Rat GZMK(Granzyme K) ELISA Kit
430	abk8303R	Rat PRCP(Prolylcarboxypeptidase) ELISA Kit
431	abk8803R	Rat EP2(Prostaglandin E Receptor 2) ELISA Kit
432	abk2702R	Rat PTGES2(Prostaglandin E Synthase 2) ELISA Kit
433	abk7011R	Rat PTGES(Prostaglandin E Synthase, Microsomal) ELISA Kit
434	abk2182R	Rat NRG-1(Neuregulin 1) ELISA Kit
435	abk3702R	Rat PTGS1(Prostaglandin Endoperoxide Synthase 1) ELISA Kit
436	abk4702R	Rat PTGS2/COX-2(Prostaglandin Endoperoxide Synthase 2) ELISA Kit
437	abk0702R	Rat FP(Prostaglandin F Receptor) ELISA Kit
438	abk7823R	Rat PSA(Prostate Specific Antigen) ELISA Kit
439	abk4761R	Rat PAR2(Protease Activated Receptor 2) ELISA Kit
440	abk1761R	Rat PRSS1(Protease, Serine, 1) ELISA Kit
441	abk3561R	Rat PDIA3(Protein Disulfide Isomerase A3) ELISA Kit
442	abk6761R	Rat PDI(Protein Disulfide Isomerase) ELISA Kit
443	abk9251R	Rat PIAS1(Protein Inhibitor Of Activated STAT 1) ELISA Kit
444	abk9443R	Rat PICK1(Protein Interacting With Protein Kinase C, Alpha 1) ELISA Kit
445	abk5561R	Rat PKBG(Protein Kinase B Gamma) ELISA Kit
446	abk8203R	Rat PKB(Protein Kinase B) ELISA Kit
447	abk8823R	Rat PKCB1(Protein Kinase C Beta 1) ELISA Kit
448	abk9823R	Rat PKCD(Protein Kinase C Delta) ELISA Kit
449	abk0923R	Rat PKCE(Protein Kinase C Epsilon) ELISA Kit
450	abk1923R	Rat PKCG(Protein Kinase C Gamma) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
451	abk2923R	Rat PKCι(Protein Kinase C Iota) ELISA Kit
452	abk3923R	Rat PKCq(Protein Kinase C Theta) ELISA Kit
453	abk4923R	Rat PKCZ(Protein Kinase C Zeta) ELISA Kit
454	abk0543R	Rat PKD1(Protein Kinase D1) ELISA Kit
455	abk1543R	Rat PKD2(Protein Kinase D2) ELISA Kit
456	abk5923R	Rat PKIA(Protein Kinase Inhibitor Alpha) ELISA Kit
457	abk6923R	Rat PKIB(Protein Kinase Inhibitor Beta) ELISA Kit
458	abk7923R	Rat PKIG(Protein Kinase Inhibitor Gamma) ELISA Kit
459	abk2543R	Rat PKN1(Protein Kinase N1) ELISA Kit
460	abk8923R	Rat PRKAA1(Protein Kinase, AMP Activated Alpha 1) ELISA Kit
461	abk3051R	Rat PTPRN(Protein Tyrosine Phosphatase Receptor Type N) ELISA Kit
462	abk5051R	Rat PTPRO(Protein Tyrosine Phosphatase Receptor Type O) ELISA Kit
463	abk9232R	Rat PCDH1(Protocadherin 1) ELISA Kit
464	abk3543R	Rat GRAP2(GRB2-Related Adapter protein 2) ELISA Kit
465	abk5171R	Rat SPC(Pulmonary Surfactant Associated Protein C) ELISA Kit
466	abk4543R	Rat PDK4(Pyruvate Dehydrogenase Kinase Isozyme 4) ELISA Kit
467	abk2161R	Rat PDP(Pyruvate Dehydrogenase Phosphatase) ELISA Kit
468	abk0051R	Rat NRG-2(Neuregulin 2) ELISA Kit
469	abk1813R	Rat RANKL(Receptor Activator of Nuclear Factor Kappa B Ligand) ELISA Kit
470	abk4132R	Rat REG4(Regenerating Islet Derived Protein 4) ELISA Kit
471	abk2023R	Rat RANTES(Regulated On Activation, Normal T-Cell Expressed and Secreted) ELISA Kit
472	abk3992R	Rat RGS4(Regulator Of G Protein Signaling 4) ELISA Kit
473	abk3872R	Rat Mly/CXCL9(Monocyte Interferon Gamma Inducing Factor) ELISA Kit
474	abk1033R	Rat RNASEA(Ribonuclease A) ELISA Kit
475	abk2403R	Rat RNASET2(Ribonuclease T2) ELISA Kit
476	abk0081R	Rat RPL10(Ribosomal Protein L10) ELISA Kit
477	abk5181R	Rat RNF4(Ring Finger Protein 4) ELISA Kit
478	abk0251R	Rat RUNX2(Runt Related Transcription Factor 2) ELISA Kit
479	abk2371R	Rat S100A11(S100 Calcium Binding Protein A11) ELISA Kit
480	abk7371R	Rat S100B(S100 Calcium Binding Protein B) ELISA Kit
481	abk8281R	Rat FABP3(Fatty Acid Binding Protein 3, Muscle and Heart) ELISA Kit
482	abk1271R	Rat sIgA(Secretory Immunoglobulin A) ELISA Kit
483	abk9113R	Rat SEPP1(Selenoprotein P) ELISA Kit
484	abk5061R	Rat SEMA7A(Semaphorin 7A) ELISA Kit
485	abk6113R	Rat SOX2(Sex Determining Region Y Box Protein 2) ELISA Kit
486	abk5543R	Rat STAT5B(Signal Transducer And Activator Of Transcription 5B) ELISA Kit
487	abk6211R	Rat ACTN2(Actinin Alpha 2) ELISA Kit
488	abk4251R	Rat SLIT2(Slit Homolog 2) ELISA Kit
489	abk4113R	Rat SNRPD1(Small Nuclear Ribonucleoprotein D1) ELISA Kit
490	abk3482R	Rat sCD14(Soluble Cluster of Differentiation 14) ELISA Kit
491	abk7301R	Rat EPCR(Endothelial Protein C Receptor) ELISA Kit
492	abk8362R	Rat sFAS/sAPO-1(Soluble Factor-Related Apoptosis) ELISA Kit
493	abk0011R	Rat IL-6R(Interleukin 6 Receptor) ELISA Kit
494	abk9913R	Rat VEGFR1/FLT1(Vascular Endothelial Growth Factor Receptor 1) ELISA Kit
495	abk6543R	Rat SSTR1(Somatostatin Receptor 1) ELISA Kit
496	abk3503R	Rat SSFA2(Sperm Specific Antigen 2) ELISA Kit
497	abk0471R	Rat SCFR(Stem Cell Growth Factor Receptor) ELISA Kit
498	abk9371R	Rat SCF(Stem Cell Factor) ELISA Kit
499	abk5771R	Rat SF1(Steroidogenic Factor 1) ELISA Kit
500	abk2581R	Rat SDF2L1(Stromal Cell Derived Factor 2 Like Protein 1) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
501	abk7781R	Rat SSRP1(Structure Specific Recognition Protein 1) ELISA Kit
502	abk6222R	Rat SOCS3(Suppressors Of Cytokine Signaling 3) ELISA Kit
503	abk3202R	Rat SMN1(Survival of Motor Neuron 1, Telomeric) ELISA Kit
504	abk2812R	Rat SYP(Synaptophysin) ELISA Kit
505	abk6933R	Rat SNAP-25(Synaptosome Associated Protein 25) ELISA Kit
506	abk2891R	Rat SDC4(Syndecan 4) ELISA Kit
507	abk7543R	Rat ChAT(Choline Acetyltransferase) ELISA Kit
508	abk9412R	Rat TACR1(Tachykinin Receptor 1) ELISA Kit
509	abk0512R	Rat TACR2(Tachykinin Receptor 2) ELISA Kit
510	abk4351R	Rat THP(Tamm-Horsfall Glycoprotein) ELISA Kit
511	abk5351R	Rat TBP(Tata Binding Protein) ELISA Kit
512	abk2181R	Rat TXLNA(Taxilin Alpha) ELISA Kit
513	abk3181R	Rat TXLNB(Taxilin Beta) ELISA Kit
514	abk1133R	Rat TEK(Tyrosine Kinase, Endothelial) ELISA Kit
515	abk2801R	Rat TE(Telomerase) ELISA Kit
516	abk1281R	Rat TNC(Tenascin C) ELISA Kit
517	abk3281R	Rat TNX(Tenascin X) ELISA Kit
518	abk4581R	Rat TDGF1(Teratocarcinoma Derived Growth Factor 1) ELISA Kit
519	abk8042R	Rat sC5b-9(Soluble Terminal Complement Complex C5b-9) ELISA Kit
520	abk0102R	Rat TrxR1(Thioredoxin Reductase 1) ELISA Kit
521	abk9002R	Rat TRX(Thioredoxin) ELISA Kit
522	abk8322R	Rat TAFI(Thrombin Activatable Fibrinolysis Inhibitor) ELISA Kit
523	abk3133R	Rat TR(Thrombin Receptor) ELISA Kit
524	abk8543R	Rat MSRA(Methionine Sulfoxide Reductase A) ELISA Kit
525	abk9543R	Rat MSRB1(Methionine Sulfoxide Reductase B1) ELISA Kit
526	abk6353R	Rat TM(Thrombomodulin) ELISA Kit
527	abk6453R	Rat TPO(Thrombopoietin) ELISA Kit
528	abk2822R	Rat TPP(Thrombus Precursor Protein) ELISA Kit
529	abk9291R	Rat TSLP(Thymic Stromal Lymphopoietin) ELISA Kit
530	abk3101R	Rat TARC(Thymus Activation Regulated Chemokine) ELISA Kit
531	abk7453R	Rat TG(Thyroglobulin) ELISA Kit
532	abk1671R	Rat HDAC2(Histone Deacetylase 2) ELISA Kit
533	abk9991R	Rat GPC3(Glypican 3) ELISA Kit
534	abk5552R	Rat HDAC1(Histone Deacetylase 1) ELISA Kit
535	abk4513R	Rat TSH(Thyroid Stimulating Hormone) ELISA Kit
536	abk6341R	Rat COL3a1(Collagen Type III Alpha 1) ELISA Kit
537	abk0461R	Rat TEF(Thyrotroph Embryonic Factor) ELISA Kit
538	abk3801R	Rat TFPI(Tissue Factor Pathway Inhibitor) ELISA Kit
539	abk5453R	Rat TF(Tissue Factor) ELISA Kit
540	abk5481R	Rat TIMP-2(Tissue Inhibitors of Metalloproteinase 2) ELISA Kit
541	abk6481R	Rat TIMP-3(Tissue Inhibitors of Metalloproteinase 3) ELISA Kit
542	abk7481R	Rat TIMP-4(Tissue Inhibitors of Metalloproteinase 4) ELISA Kit
543	abk0643R	Rat SERP1(Stress-associated Endoplasmic Reticulum Protein 1) ELISA Kit
544	abk8382R	Rat TLR4(Toll-Like Receptor 4) ELISA Kit
545	abk8351R	Rat TLR9(Toll-Like Receptor 9) ELISA Kit
546	abk8371R	Rat SDC1(Syndecan 1) ELISA Kit
547	abk6111R	Rat 5-LO(Arachidonate 5-Lipoxygenase) ELISA Kit
548	abk5642R	Rat 12-LO(Arachidonate 12-Lipoxygenase) ELISA Kit
549	abk7013R	Rat CK-18/KRT18(Cytokeratin 18) ELISA Kit
550	abk6832R	Rat TCF20(Transcription factor 20) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
551	abk1643R	Rat ERRα(Estrogen Related Receptor Alpha) ELISA Kit
552	abk6451R	Rat α2-AP(α2-Antiplasmin) ELISA Kit
553	abk0542R	Rat IL-24(Interleukin 24) ELISA Kit
554	abk9933R	Rat TGF-α(Transforming Growth Factor Alpha) ELISA Kit
555	abk5851R	Rat GP4/CD36(Platelet Membrane Glycoprotein IV) ELISA Kit
556	abk7691R	Rat TSC22(Transforming Growth Factor Beta Stimulated Protein Clone 22) ELISA Kit
557	abk8972R	Rat TG(Transglutaminase) ELISA Kit
558	abk8133R	Rat TMED1(Transmembrane Emp24 Domain Containing Protein 1) ELISA Kit
559	abk9832R	Rat TNPO1(Transportin 1) ELISA Kit
560	abk5832R	Rat TTR(Transthyretin) ELISA Kit
561	abk1902R	Rat TFF2(Trefoil Factor 2) ELISA Kit
562	abk2261R	Rat TFF3(Trefoil Factor 3, Intestinal) ELISA Kit
563	abk0611R	Rat INHB(Inhibin B) ELISA Kit
564	abk9133R	Rat TUBB(Tubulin Beta) ELISA Kit
565	abk3013R	Rat PPP(Tubulin Polymerization Promoting Protein) ELISA Kit
566	abk8542R	Rat TNFRSF1B(Tumor Necrosis Factor Receptor Superfamily, Member 1B) ELISA Kit
567	abk6291R	Rat TRAIL/TNFSF10(Tumor Necrosis Factor Related Apoptosis Inducing Ligand) ELISA Kit
568	abk5313R	Rat TSGF(Tumor Specific Growth Factor/Tumor Supplied Group of Factor) ELISA Kit
569	abk2791R	Rat TYR(Tyrosinase) ELISA Kit
570	abk9302R	Rat PLAU/uPA(Urokinase-Type Plasminogen Activator) ELISA Kit
571	abk0391R	Rat PLAUR/uPAR(Plasminogen Activator, Urokinase Receptor) ELISA Kit
572	abk0722R	Rat VEGF-B(Vascular Endothelial Cell Growth Factor B) ELISA Kit
573	abk2643R	Rat HRAS(V-Ha-Ras Harvey Rat Sarcoma Viral Oncogene Homolog) ELISA Kit
574	abk8202R	Rat VASPIN(Visceral Adipose Specific Serine Protease Inhibitor) ELISA Kit
575	abk9202R	Rat VF(Visfatin) ELISA Kit
576	abk4913R	Rat DBP(Vitamin D Binding Protein) ELISA Kit
577	abk3722R	Rat vWF(Von Willebrand Factor) ELISA Kit
578	abk3643R	Rat XIAP(X-Linked Inhibitor Of Apoptosis Protein) ELISA Kit
579	abk0233R	Rat α1-MG(α1-Microglobulin) ELISA Kit
580	abk5712R	Rat GLTP(Glycolipid Transfer Protein) ELISA Kit
581	abk4111R	Rat TNFRSF1A(Tumor Necrosis Factor Receptor Superfamily, Member 1A) ELISA Kit
582	abk4643R	Rat MDH1(Malate Dehydrogenase 1) ELISA Kit
583	abk8913R	Rat LEPR(Leptin Receptor) ELISA Kit
584	abk2901R	Rat PCIII(Procollagen III) ELISA Kit
585	abk5643R	Rat FLAP(5-Lipoxygenase Activating Protein) ELISA Kit
586	abk3532R	Rat PON1(Paraoxonase 1) ELISA Kit
587	abk6652R	Rat DNMT1(DNA Methyltransferase 1) ELISA Kit
588	abk2091R	Rat SIRT1(Sirtuin 1) ELISA Kit
589	abk6552R	Rat SCUBE1(Signal Peptide, CUB and EGF-like Domain-containing Protein 1) ELISA Kit
590	abk1313R	Rat FNDC5(Fibronectin type III domain-containing protein 5) ELISA Kit
591	abk6643R	Rat GNMT(Glycine N-methyltransferase) ELISA Kit
592	abk9742R	Rat FSHB(Follicle Stimulating Subunit Beta) ELISA Kit
593	abk7643R	Rat NTSR2(Neurotensin Receptor Type 2) ELISA Kit
594	abk2513R	Rat ALPL(Alkaline Phosphatase, Liver/Bone/Kidney) ELISA Kit
595	abk7511R	Rat AIF1(Allograft Inflammatory Factor 1) ELISA Kit
596	abk0562R	Rat IRS1(Insulin Receptor Substrate 1) ELISA Kit
597	abk8842R	Rat CYGB(Cytoglobin) ELISA Kit
598	abk9353R	Rat GLO1(Glyoxalase I) ELISA Kit
599	abk9632R	Rat CDK4(Cyclin Dependent Kinase 4) ELISA Kit
600	abk5092R	Rat IL-10Rβ(Interleukin 10 Receptor Beta) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
601	abk2521R	Rat INSR(Insulin Receptor) ELISA Kit
602	abk7131R	Rat BECN1(Beclin 1) ELISA Kit
603	abk8521R	Rat ApoC1(Apolipoprotein C1) ELISA Kit
604	abk8643R	Rat HSD11β1(11-Beta-Hydroxysteroid Dehydrogenase Type 1) ELISA Kit
605	abk6282R	Rat FG(Fibrinogen) ELISA Kit
606	abk7233R	Rat CCNE1(Cyclin-E1) ELISA Kit
607	abk4453R	Rat GR(Glutathione Reductase) ELISA Kit
608	abk9643R	Rat RGN(Regucalcin) ELISA Kit
609	abk2932R	Rat MAOB(Amine oxidaseB) ELISA Kit
610	abk1743R	Rat IRF4(Interferon Regulatory Factor 4) ELISA Kit
611	abk1872R	Rat α2-M(Alpha-2 Macroglobulin) ELISA Kit
612	abk1822R	Rat SERT(Serotonin Transporter) ELISA Kit
613	abk6772R	Rat C2(Complement Component 2) ELISA Kit
614	abk9233R	Rat CCNE2(Cyclin-E2) ELISA Kit
615	abk2743R	Rat UGT1A1(UDP Glucuronosyltransferase 1 Family, Polypeptide A1) ELISA Kit
616	abk3881R	Rat MT2(Metalllothionein 2) ELISA Kit
617	abk0791R	Rat UCP-1(Uncoupling Protein 1, Mitochondrial) ELISA Kit
618	abk9521R	Rat ApoC2(Apolipoprotein C2) ELISA Kit
619	abk1242R	Rat MRC1(Mannose Receptor C Type 1) ELISA Kit
620	abk3743R	Rat PEMT(Phosphatidylethanolamine N-methyltransferase) ELISA Kit
621	abk4743R	Rat FECH(Ferrochelatase) ELISA Kit
622	abk5743R	Rat BCL2L1(Bcl-2-Like Protein 1) ELISA Kit
623	abk3742R	Rat VLDL(Very Low Density Lipoprotein) ELISA Kit
624	abk6743R	Rat MC4R(Melanocortin Receptor 4) ELISA Kit
625	abk6402R	Rat F1+2(Prothrombin Fragment 1+2) ELISA Kit
626	abk7252R	Rat VASP(Vasodilator Stimulated Phosphoprotein) ELISA Kit
627	abk8652R	Rat DNMT3A(DNA Methyltransferase 3A) ELISA Kit
628	abk6232R	Rat VAP-1(Vascular Adhesion Protein 1) ELISA Kit
629	abk7743R	Rat ICD(Isocitrate Dehydrogenase) ELISA Kit
630	abk0672R	Rat ApoA2(Apolipoprotein A2) ELISA Kit
631	abk7521R	Rat ApoB(Apolipoprotein B) ELISA Kit
632	abk8732R	Rat AMACR(Alpha-methylacyl-CoA racemase) ELISA Kit
633	abk9732R	Rat CCNA1(Cyclin-A1) ELISA Kit
634	abk0832R	Rat RELN(Reelin) ELISA Kit
635	abk8743R	Rat GPIb(Platelet Membrane Glycoprotein Ib) ELISA Kit
636	abk1132R	Rat JAK2(Janus Kinase 2) ELISA Kit
637	abk3932R	Rat FPN(Ferroportin) ELISA Kit
638	abk2621R	Rat ApoE(Apolipoprotein E) ELISA Kit
639	abk9652R	Rat DNMT3B(DNA Methyltransferase 3B) ELISA Kit
640	abk8482R	Rat BDNF(Brain Derived Neurotrophic Factor) ELISA Kit
641	abk5292R	Rat Aβ1-42(Amyloid Beta 1-42) ELISA Kit
642	abk8013R	Rat CYP7A1(Cholesterol 7-Alpha-Hydroxylase) ELISA Kit
643	abk2471R	Rat GALR1(Galanin Receptor 1) ELISA Kit
644	abk9743R	Rat GALR2(Galanin Receptor2) ELISA Kit
645	abk0843R	Rat GALR3(Galanin Receptor3) ELISA Kit
646	abk1832R	Rat TGM2(Transglutaminase 2, Tissue) ELISA Kit
647	abk3182R	Rat GFAP(Glial Fibrillary Acidic Protein) ELISA Kit
648	abk1413R	Rat PS1(Presenilin 1) ELISA Kit
649	abk9272R	Rat AgRP(Agouti Related Protein) ELISA Kit
650	abk3691R	Rat TH(Tyrosine Hydroxylase) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
651	abk1141R	Rat NOS1/nNOS(Nitric Oxide Synthase 1, Neuronal) ELISA Kit
652	abk3641R	Rat CTX I (Cross Linked C-telopeptide of Type I Collagen) ELISA Kit
653	abk6301R	Rat ET-1(Endothelin 1) ELISA Kit
654	abk6413R	Rat RTN4(Reticulon 4) ELISA Kit
655	abk3721R	Rat AQP-2(Aquaporin 2) ELISA Kit
656	abk4192R	Rat ApoC3(Apolipoprotein C3) ELISA Kit
657	abk5511R	Rat AGRN(Agrin) ELISA Kit
658	abk1823R	Rat PACAP(Pituitary Adenylate Cyclase Activating Polypeptide) ELISA Kit
659	abk6321R	Rat ANXA5(Annexin A5) ELISA Kit
660	abk2843R	Rat GPIb-IX-V(Platelet Membrane Glycoprotein GPIb-IX-V Complex) ELISA Kit
661	abk9932R	Rat SCARB1(Scavenger receptor class B member 1) ELISA Kit
662	abk8591R	Rat IL-25(Interleukin 25) ELISA Kit
663	abk9813R	Rat POSTN/OSF-2(Periostin) ELISA Kit
664	abk2001R	Rat BAFF/CD257(B-Cell Activating Factor) ELISA Kit
665	abk3791R	Rat NFκB2(Nuclear Factor Kappa B2) ELISA Kit
666	abk9442R	Rat RELB(V-Rel Reticuloendotheliosis Viral Oncogene Homolog B) ELISA Kit
667	abk2502R	Rat PDPN(Podoplanin) ELISA Kit
668	abk6681R	Rat PCT(Procalcitonin) ELISA Kit
669	abk1023R	Rat ACE(Angiotensin I Converting Enzyme) ELISA Kit
670	abk1322R	Rat FMOD(Fibromodulin) ELISA Kit
671	abk1171R	Rat NNE(Non-Neuronal Enolase) ELISA Kit
672	abk3491R	Rat cADPRH/CD38(Cyclic ADP Ribose Hydrolase) ELISA Kit
673	abk8311R	Rat ADM(Adrenomedullin) ELISA Kit
674	abk8012R	Rat NPHN(Nephrin) ELISA Kit
675	abk0401R	Rat AFGF/FGF1(Acidic Fibroblast Growth Factor 1) ELISA Kit
676	abk2401R	Rat FGF21(Fibroblast Growth Factor 21) ELISA Kit
677	abk1401R	Rat FGF19(Fibroblast Growth Factor 19) ELISA Kit
678	abk0811R	Rat FGFR1(Fibroblast Growth Factor Receptor 1) ELISA Kit
679	abk2902R	Rat PEDF(Pigment Epithelium Derived Factor) ELISA Kit
680	abk5482R	Rat CTSB(Cathepsin B) ELISA Kit
681	abk4623R	Rat IgG2a(Immunoglobulin G2a) ELISA Kit
682	abk4333R	Rat IgG2b(Immunoglobulin G2b) ELISA Kit
683	abk5333R	Rat IgG2c(Immunoglobulin G2c) ELISA Kit
684	abk6381R	Rat MMP-11(Matrix Metalloproteinase 11) ELISA Kit
685	abk2601R	Rat EREG(Epiregulin) ELISA Kit
686	abk4451R	Rat α1-AT(Alpha 1-Antitrypsin) ELISA Kit
687	abk2742R	Rat MMP-14(Matrix Metalloproteinase 14) ELISA Kit
688	abk8101R	Rat FcεR II/CD23(Receptor II for the Fc Region of Immunoglobulin E) ELISA Kit
689	abk9591R	Rat IL-34(Interleukin 34) ELISA Kit
690	abk2991R	Rat PCK1(Phosphoenolpyruvate Carboxykinase 1) ELISA Kit
691	abk6333R	Rat ACTα2(Actin Alpha 2, Smooth Muscle) ELISA Kit
692	abk1791R	Rat IL18BP(Interleukin 18 Binding Protein) ELISA Kit
693	abk9161R	Rat FOXP3(Forkhead Box Protein P3) ELISA Kit
694	abk5341R	Rat COL2α1(Collagen Type II Alpha 1) ELISA Kit
695	abk7333R	Rat CD163(Cluster of Differentiation 163) ELISA Kit
696	abk2372R	Rat CFD(Complement Factor D) ELISA Kit
697	abk9201R	Rat I-TAC/CXCL11(Interferon Inducible T-Cell Alpha Chemoattractant) ELISA Kit
698	abk7913R	Rat ENG(Endoglin) ELISA Kit
699	abk7411R	Rat FABP4(Fatty Acid Binding Protein 4, Adipocyte) ELISA Kit
700	abk6003R	Rat IL-1R1(Interleukin 1 Receptor Type I) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
701	abk7501R	Rat IL-22(Interleukin 22) ELISA Kit
702	abk1453R	Rat PL(Pancreatic Lipase) ELISA Kit
703	abk7591R	Rat IL-12 p35(Interleukin 12 p35) ELISA Kit
704	abk3431R	Rat CALD(Caldesmon) ELISA Kit
705	abk8712R	Rat FTH(Ferritin, Heavy Polypeptide) ELISA Kit
706	abk4211R	Rat IL-19(Interleukin 19) ELISA Kit
707	abk3572R	Rat ALCAM(Activated Leukocyte Cell Adhesion Molecule) ELISA Kit
708	abk3982R	Rat ANGPTL3(Angiopietin Like Protein 3) ELISA Kit
709	abk8711R	Rat MAG(Myelin Associated Glycoprotein) ELISA Kit
710	abk1121R	Rat FETUA(Fetuin A) ELISA Kit
711	abk8282R	Rat DKK3(Dickkopf Related Protein 3) ELISA Kit
712	abk3411R	Rat ACE2(Angiotensin I Converting Enzyme 2) ELISA Kit
713	abk1421R	Rat IDE(Insulin Degrading Enzyme) ELISA Kit
714	abk8141R	Rat CD4(Cluster of Differentiation 4) ELISA Kit
715	abk2652R	Rat GRIA1(Glutamate receptor 1) ELISA Kit
716	abk3902R	Rat SDH(Sorbitol Dehydrogenase) ELISA Kit
717	abk0532R	Rat GYS2(Glycogen Synthase 2, Liver) ELISA Kit
718	abk7572R	Rat CYP1A1(Cytochrome P450, family 1, subfamily A, polypeptide 1) ELISA Kit
719	abk5911R	Rat CYP1A2(Cytochrome P450, family 1, subfamily A, polypeptide 2) ELISA Kit
720	abk0121R	Rat CYP1B1(Cytochrome P450, family 1, subfamily B, polypeptide 1) ELISA Kit
721	abk3843R	Rat CYP2C9(Cytochrome P450, family 2, subfamily C, polypeptide 9) ELISA Kit
722	abk4843R	Rat CYP2C19(Cytochrome P450, family 2, subfamily C, polypeptide 19) ELISA Kit
723	abk5843R	Rat ALAD(Aminolevulinate Delta Dehydratase) ELISA Kit
724	abk6843R	Rat LIPE(Lipase, Hormone Sensitive) ELISA Kit
725	abk4382R	Rat ucOC(Undercarboxylated Osteocalcin) ELISA Kit
726	abk9042R	Rat UCHL1(Ubiquitin Carboxyl Terminal Hydrolase L1) ELISA Kit
727	abk7992R	Rat TDP43(TAR DNA Binding Protein 43) ELISA Kit
728	abk7843R	Rat sCD80(Soluble Cluster of Differentiation 80) ELISA Kit
729	abk7091R	Rat ESM1(Endothelial Cell Specific Molecule 1) ELISA Kit
730	abk2062R	Rat Pro-ADM(Proadrenomedullin) ELISA Kit
731	abk6861R	Rat APP(Amyloid Precursor Protein) ELISA Kit
732	abk7552R	Rat GPX1(Glutathione Peroxidase 1) ELISA Kit
733	abk7831R	Rat CHEM(Chemerin) ELISA Kit
734	abk8843R	Rat PKCα(Protein Kinase Cα) ELISA Kit
735	abk2402R	Rat GS(Glutamine synthetase) ELISA Kit
736	abk9843R	Rat SDF-1α(Stromal Cell Derived Factor 1α) ELISA Kit
737	abk0943R	Rat dGHRH(Desacyl Ghrelin) ELISA Kit
738	abk8611R	Rat ANG II R1(Angiotensin II Receptor 1) ELISA Kit
739	abk1943R	Rat ADRβ1(Adrenergic Receptor Beta 1) ELISA Kit
740	abk2943R	Rat ADRβK1(Adrenergic Receptor Beta Kinase 1) ELISA Kit
741	abk4061R	Rat OCLN(Occludin) ELISA Kit
742	abk3943R	Rat F-actin ELISA Kit
743	abk0442R	Rat EP(Erythrocyte Protoporphyrin) ELISA Kit
744	abk7072R	Rat ADORA1(Adenosine A1 Receptor) ELISA Kit
745	abk9962R	Rat ADORA2a(Adenosine A2a Receptor) ELISA Kit
746	abk4943R	Rat ADORA2b(Adenosine A2b Receptor) ELISA Kit
747	abk5172R	Rat ADORA3(Adenosine A3 Receptor) ELISA Kit
748	abk2762R	Rat ACC(Acetyl-CoA Carboxylase) ELISA Kit
749	abk5943R	Rat PC(Pyruvate Carboxylase) ELISA Kit
750	abk0413R	Rat PLA2G2A(Phospholipase A2, Group II A) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
751	abk6943R	Rat CAPN2(Calpain-2 catalytic subunit) ELISA Kit
752	abk3043R	Rat β-Hex A(β-hexosaminidase A) ELISA Kit
753	abk7943R	Rat β-Hex B(β-hexosaminidase B) ELISA Kit
754	abk8603R	Rat CD28(Cluster of Differentiation 28) ELISA Kit
755	abk8943R	Rat ICOS(Inducible T-cell costimulator) ELISA Kit
756	abk2291R	Rat TPO(Thyroid Peroxidase) ELISA Kit
757	abk1713R	Rat CF6(Coupling Factor 6) ELISA Kit
758	abk4222R	Rat ERK2(Extracellular Signal Regulated Kinase 2) ELISA Kit
759	abk0053R	Rat ERK1/2(Extracellular Signal Regulated Kinase 1/2) ELISA Kit
760	abk9182R	Rat KISS1(Kisspeptin 1) ELISA Kit
761	abk8692R	Rat P II NP(Procollagen II N-Terminal Propeptide) ELISA Kit
762	abk9321R	Rat HYOU1(Hypoxia Up Regulated 1) ELISA Kit
763	abk9262R	Rat GDN(Glia Derived Nexin) ELISA Kit
764	abk0141R	Rat NEFL(Neurofilament, Light Polypeptide) ELISA Kit
765	abk9462R	Rat CAPN1(calpain 1) ELISA Kit
766	abk1053R	Rat AZU/HBP(Azurocidin/Heparin Binding Protein) ELISA Kit
767	abk8372R	Rat sgp130(soluble Glucoprotein 130) ELISA Kit
768	abk2053R	Rat CD68(Cluster of Differentiation 68) ELISA Kit
769	abk3053R	Rat sE-selectin(soluble E-selectin) ELISA Kit
770	abk4611R	Rat AMY1(Amylase Alpha 1, Salivary) ELISA Kit
771	abk7032R	Rat AMY2(Amylase Alpha 2, Pancreatic) ELISA Kit
772	abk1131R	Rat LDHA(Lactate Dehydrogenase A) ELISA Kit
773	abk6641R	Rat LDHB(Lactate Dehydrogenase B) ELISA Kit
774	abk7741R	Rat LDHC(Lactate Dehydrogenase C) ELISA Kit
775	abk9781R	Rat TNNT1(Troponin T Type 1, Slow Skeletal) ELISA Kit
776	abk5641R	Rat CTX II (Cross Linked C-telopeptide of Type II Collagen) ELISA Kit
777	abk5053R	Rat GSTP(Glutathione S-transferase P) ELISA Kit
778	abk6053R	Rat HRH4(Histamine Receptor H4) ELISA Kit
779	abk4551R	Rat P II CP(Procollagen II C-Terminal ProPeptide) ELISA Kit
780	abk9041R	Rat NEF3(Neurofilament 3) ELISA Kit
781	abk9952R	Rat CSNK2B(Casein Kinase 2, beta) ELISA Kit
782	abk8452R	Rat AD(Adropin) ELISA Kit
783	abk9762R	Rat EAAT3(Excitatory Amino Acid Transporter 3) ELISA Kit
784	abk7053R	Rat RFTN1(Raft Linking Protein) ELISA Kit
785	abk0031R	Rat ATF6(Activating Transcription Factor 6) ELISA Kit
786	abk0371R	Rat S100A9(S100 Calcium Binding Protein A9) ELISA Kit
787	abk8053R	Rat LXRβ(liver X Receptor Beta) ELISA Kit
788	abk9053R	Rat APLNR(Apelin Receptor) ELISA Kit
789	abk0153R	Rat NEFH(Neurofilament, Heavy Polypeptide) ELISA Kit
790	abk1153R	Rat WNT10B(Wingless-Type MMTV Integration Site Family, Member 10B) ELISA Kit
791	abk3212R	Rat GAP43(Growth Associated Protein 43) ELISA Kit
792	abk0852R	Rat KL(Klotho) ELISA Kit
793	abk3492R	Rat CCR1(CC-Chemokine Receptor 1) ELISA Kit
794	abk2153R	Rat ACTα1(Actin Alpha 1, Skeletal Muscle) ELISA Kit
795	abk1861R	Rat LDLR(Low Density Lipoprotein Receptor) ELISA Kit
796	abk9662R	Rat XOD(Xanthine Oxidase) ELISA Kit
797	abk6521R	Rat ApoA4(Apolipoprotein A4) ELISA Kit
798	abk7151R	Rat RHOA(Ras Homolog Gene Family, Member A) ELISA Kit
799	abk1862R	Rat MIOX(Myoinositol Oxygenase) ELISA Kit
800	abk3153R	Rat MC3R(Melanocortin Receptor 3) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
801	abk2672R	Rat ApoA5(Apolipoprotein A5) ELISA Kit
802	abk3682R	Rat IGF-1(Insulin-like Growth Factor 1) ELISA Kit
803	abk7801R	Rat GH(Growth Hormone) ELISA Kit
804	abk2472R	Rat C-P(C-peptide) ELISA Kit
805	abk5671R	Rat OPG(Osteoprotegerin) ELISA Kit
806	abk1701R	Rat PRL(Prolactin) ELISA Kit
807	abk7472R	Rat GLP-1(Glucagon Like Peptide 1) ELISA Kit
808	abk9223R	Rat COL1(Collagen Type I) ELISA Kit
809	abk8801R	Rat COL4(Collagen Type IV) ELISA Kit
810	abk8772R	Rat CFH(Complement Factor H) ELISA Kit
811	abk8901R	Rat C5a(Complement Component 5a) ELISA Kit
812	abk7382R	Rat ADP/Acrp30(Adiponectin) ELISA Kit
813	abk9961R	Rat DAO(Diamine Oxidase) ELISA Kit
814	abk2972R	Rat IgA(Immunoglobulin A) ELISA Kit
815	abk6803R	Rat IgM(Immunoglobulin M) ELISA Kit
816	abk4682R	Rat IGFBP-3(Insulin-Like Growth Factor Binding Protein 3) ELISA Kit
817	abk2782R	Rat FE(Ferritin) ELISA Kit
818	abk0572R	Rat KIM-1(Kidney Injury Molecule 1) ELISA Kit
819	abk3481R	Rat MMP-8(Matrix Metalloproteinase 8) ELISA Kit
820	abk6733R	Rat AMH(Anti-Mullerian Hormone) ELISA Kit
821	abk4813R	Rat NT-proBNP(N-Terminal Pro-Brain Natriuretic Peptide) ELISA Kit
822	abk3901R	Rat PIIINP(N-Terminal Procollagen III Propeptide) ELISA Kit
823	abk5322R	Rat PAI1(Plasminogen Activator Inhibitor 1) ELISA Kit
824	abk0213R	Rat SAA(Serum Amyloid A) ELISA Kit
825	abk7582R	Rat SDF-1/CXCL12(Stromal Cell Derived Factor 1) ELISA Kit
826	abk9472R	Rat TF(Transferrin) ELISA Kit
827	abk6601R	Rat ApoA1(Apolipoprotein A1) ELISA Kit
828	abk4292R	Rat Aβ1-40(Amyloid Beta 1-40) ELISA Kit
829	abk8261R	Rat FGF23(Fibroblast Growth Factor 23) ELISA Kit
830	abk0091R	Rat SOST(Sclerostin) ELISA Kit
831	abk6482R	Rat FABP1(Fatty Acid Binding Protein 1, Liver) ELISA Kit
832	abk3513R	Rat INS(Insulin) ELISA Kit
833	abk9382R	Rat IFN-α(Interferon Alpha) ELISA Kit
834	abk5192R	Rat ApoB100(Apolipoprotein B100) ELISA Kit
835	abk5303R	Rat SP-A1(Pulmonary surfactant-associated protein A1) ELISA Kit
836	abk7251R	Rat Smad7/MADH7(Mothers Against Decapentaplegic Homolog 7) ELISA Kit
837	abk8903R	Rat NE/ELA2(Neutrophil Elastase/Elastase-2) ELISA Kit
838	abk0403R	Rat MBP-C/MBL2(Mannose Binding Protein C/Mannose Binding Lectin) ELISA Kit
839	abk7192R	Rat ARG1(Arginase I) ELISA Kit
840	abk6303R	Rat SP-B(Pulmonary surfactant-associated protein B) ELISA Kit
841	abk2292R	Rat ABCA1(ATP Binding Cassette Transporter A1) ELISA Kit
842	abk1292R	Rat ABCG1(ATP Binding Cassette Transporter G1) ELISA Kit
843	abk3601R	Rat Bcl-2(B-cell Lymphoma/Leukemia 2) ELISA Kit
844	abk3131R	Rat BAK1(BCL2 Antagonist/Killer 1) ELISA Kit
845	abk5131R	Rat BAX(Bcl-2 Associated X Protein) ELISA Kit
846	abk9292R	Rat BMPR1A(Bone Morphogenetic Protein Receptor 1A) ELISA Kit
847	abk0331R	Rat BRAK/CXCL14(Breast and Kidney Expressed Chemokine) ELISA Kit
848	abk3392R	Rat CGRP1(Calcitonin Gene Related Peptide 1) ELISA Kit
849	abk8531R	Rat cTnT/TNNT2(Troponin T Type 2, Cardiac) ELISA Kit
850	abk0731R	Rat CTNNα1(Catenin, Alpha 1) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
851	abk1731R	Rat CTNNβ1(Catenin, Beta 1) ELISA Kit
852	abk8731R	Rat C/EBPδ(CCAAT/Enhancer Binding Protein, Delta) ELISA Kit
853	abk7492R	Rat CENPI(Centromere Protein I) ELISA Kit
854	abk3041R	Rat CETP(Cholesteryl Ester Transfer Protein) ELISA Kit
855	abk1692R	Rat CHGA(Chromogranin A) ELISA Kit
856	abk4153R	Rat JNK1(C-Jun N-Terminal Kinases 1) ELISA Kit
857	abk5153R	Rat CD44(Cluster of Differentiation 44) ELISA Kit
858	abk1833R	Rat F2(Coagulation Factor II) ELISA Kit
859	abk6153R	Rat F3(Coagulation Factor III) ELISA Kit
860	abk4241R	Rat F10(Coagulation Factor X) ELISA Kit
861	abk5241R	Rat F13(Coagulation Factor XIII) ELISA Kit
862	abk6241R	Rat F12(Coagulation Factor XII) ELISA Kit
863	abk3241R	Rat F9(Coagulation Factor IX) ELISA Kit
864	abk7241R	Rat F5(Coagulation Factor V) ELISA Kit
865	abk0392R	Rat IBSP(Bone Sialoprotein 2) ELISA Kit
866	abk5772R	Rat C3(Complement Component 3) ELISA Kit
867	abk2792R	Rat C5(Complement Component 5) ELISA Kit
868	abk0503R	Rat CNTN4(Contactin 4) ELISA Kit
869	abk7153R	Rat CTX-III(Cross Linked C-telopeptide of Type III Collagen) ELISA Kit
870	abk4892R	Rat CLEC2C(C-Type Lectin Domain Family 2, Member C) ELISA Kit
871	abk3892R	Rat CXCR3(CXC-Chemokine Receptor 3) ELISA Kit
872	abk2892R	Rat CREB1(Cyclic AMP Response Element Binding Protein 1) ELISA Kit
873	abk0732R	Rat CDK5(Cyclin Dependent Kinase 5) ELISA Kit
874	abk2833R	Rat CCND3(Cyclin D3) ELISA Kit
875	abk5682R	Rat FGF7/KGF(Fibroblast Growth Factor 7) ELISA Kit
876	abk9842R	Rat CDKN1A(Cyclin Dependent Kinase Inhibitor 1A) ELISA Kit
877	abk9003R	Rat cPLA2(Cytosolic Phospholipase A2) ELISA Kit
878	abk5013R	Rat CTLA4(Cytotoxic T-Lymphocyte Associated Antigen 4) ELISA Kit
879	abk0013R	Rat DAP6(Death Associated Protein 6) ELISA Kit
880	abk2482R	Rat DPP4(Dipeptidyl Peptidase IV) ELISA Kit
881	abk4071R	Rat DNM1(Dynamin 1) ELISA Kit
882	abk0602R	Rat GLUT2(Glucose Transporter 2) ELISA Kit
883	abk7203R	Rat EMILIN1(Elastin Microfibril Interface Located Protein 1) ELISA Kit
884	abk5301R	Rat EG-VEGF(Endocrine Gland Derived Vascular Endothelial Growth Factor) ELISA Kit
885	abk3222R	Rat ERK1(Extracellular Signal Regulated Kinase 1) ELISA Kit
886	abk1422R	Rat FGα(Fibrinogen Alpha) ELISA Kit
887	abk6852R	Rat FGF10(Fibroblast Growth Factor 10) ELISA Kit
888	abk1203R	Rat FSH(Follicle Stimulating Hormone) ELISA Kit
889	abk2872R	Rat GLβ(Galactosidase Beta) ELISA Kit
890	abk5022R	Rat GKN2(Gastrokinase 2) ELISA Kit
891	abk0712R	Rat GRα(Glucocorticoid Receptor Alpha) ELISA Kit
892	abk1712R	Rat GRβ(Glucocorticoid Receptor Beta) ELISA Kit
893	abk3651R	Rat GUSβ(β-glucuronidase) ELISA Kit
894	abk3712R	Rat GSK3α(Glycogen Synthase Kinase 3 Alpha) ELISA Kit
895	abk8391R	Rat GzmA(Granzyme A) ELISA Kit
896	abk9391R	Rat GzmB(Granzyme B) ELISA Kit
897	abk1802R	Rat HSP-27/HSPB1(Heat Shock Protein 27) ELISA Kit
898	abk2903R	Rat HSP-70/HSPA9(Heat Shock Protein 70) ELISA Kit
899	abk6722R	Rat HO1(Heme Oxygenase 1) ELISA Kit
900	abk7722R	Rat HO2(Heme Oxygenase 2) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
901	abk8153R	Rat HSPG(Heparan Sulfate Proteoglycan) ELISA Kit
902	abk6523R	Rat HB-EGF(Heparin-binding Epidermal Growth Factor-like Growth Factor) ELISA Kit
903	abk8523R	Rat HNF4a(Hepatocyte Nuclear Factor 4 Alpha) ELISA Kit
904	abk4582R	Rat CRP(C-Reactive Protein) ELISA Kit
905	abk9252R	Rat HABP-1/C1QBP(Hyaluronan Binding Protein 1) ELISA Kit
906	abk3782R	Rat IgG(Immunoglobulin G) ELISA Kit
907	abk9782R	Rat INHA(Inhibin Alpha Chain) ELISA Kit
908	abk9153R	Rat INHbC(Inhibin Beta C) ELISA Kit
909	abk9092R	Rat IGFBP-5(Insulin-like Growth Factor Binding Protein 5) ELISA
910	abk1192R	Rat IGFBP-6(Insulin-like Growth Factor Binding Protein 6) ELISA
911	abk8092R	Rat IGFBP-1(Insulin-like Growth Factor Binding Protein 1) ELISA
912	abk9572R	Rat IGFBP-2(Insulin-like Growth Factor Binding Protein 2) ELISA
913	abk0192R	Rat IGFBP-4(Insulin-like Growth Factor Binding Protein 4)ELISA
914	abk4603R	Rat PDGF-BB(Platelet Derived Growth Factor BB) ELISA Kit
915	abk4901R	Rat TIMP-1(Tissue Inhibitors of Metalloproteinase 1) ELISA Kit
916	abk8201R	Rat IP-10/CXCL10(Interferon Gamma Induced Protein 10kDa) ELISA Kit
917	abk8401R	Rat IL-1ra/IL-1F3(Interleukin 1 Receptor Antagonist) ELISA Kit
918	abk7682R	Rat IL-17A(Interleukin 17A) ELISA Kit
919	abk0613R	Rat LTF/LF(Lactoferrin) ELISA Kit
920	abk8043R	Rat LRP-1(Low-Density Lipoprotein-Receptor-Related Protein 1) ELISA Kit
921	abk7282R	Rat LBP(Lipopolysaccharide Binding Protein) ELISA Kit
922	abk3313R	Rat LpPLA2(Lipoprotein-associated Phospholipase A2) ELISA Kit
923	abk2501R	Rat MCSF(Macrophage Colony Stimulating Factor 1) ELISA Kit
924	abk0253R	Rat MHC I /AgB I /H-1 I (Major Histocompatibility Complex I) ELISA Kit
925	abk9331R	Rat MUC1(Mucin 1) ELISA Kit
926	abk5251R	Rat Smad1/MADH1(Mothers Against Decapentaplegic Homolog 1) ELISA Kit
927	abk6251R	Rat Smad4/MADH4(Mothers Against Decapentaplegic Homolog 4) ELISA Kit
928	abk3012R	Rat NCAM/CD56(Neural Cell Adhesion Molecule) ELISA Kit
929	abk3403R	Rat NFkB-p105(Nuclear factor NF-kappa-B p105 subunit) ELISA Kit
930	abk4403R	Rat NFkB-p65(Nuclear Factor Kappa B p65) ELISA Kit
931	abk5751R	Rat OCT4(Octamer-binding protein 4) ELISA Kit
932	abk8132R	Rat OSM(Oncostatin M) ELISA Kit
933	abk7903R	Rat GROβ/CXCL2(Growth Regulated Oncogene Beta) ELISA Kit
934	abk0112R	Rat GROγ/CXCL3(Growth Regulated Oncogene Gamma) ELISA Kit
935	abk3552R	Rat TPH1(Tryptophan Hydroxylase 1) ELISA Kit
936	abk3472R	Rat OxLDL(Oxidized Low Density Lipoprotein) ELISA Kit
937	abk5571R	Rat Smad3/MADH3(Mothers Against Decapentaplegic Homolog 3) ELISA Kit
938	abk0681R	Rat PTHRP(Parathyroid Hormone Related Protein) ELISA Kit
939	abk2871R	Rat PPAR-γ(Peroxisome Proliferator Activated Receptor Gamma) ELISA Kit
940	abk1871R	Rat PPAR-α(Peroxisome Proliferator Activated Receptor Alpha) ELISA Kit
941	abk7703R	Rat pERK1/2(Phospho Extracellular Signal Regulated Kinase 1/2) ELISA Kit
942	abk0991R	Rat plkBα(Phospho-Inhibitory Subunit of NF-KBα) ELISA Kit
943	abk8703R	Rat PLA2P(Phospholipase A2 Activating Protein) ELISA Kit
944	abk6322R	Rat PAI2(Plasminogen Activator Inhibitor 2) ELISA Kit
945	abk3703R	Rat PDGF-A(Platelet Derived Growth Factor Subunit A) ELISA Kit
946	abk4703R	Rat PDGF-B(Platelet Derived Growth Factor Subunit B) ELISA Kit
947	abk7023R	Rat PF4(Platelet Factor 4) ELISA Kit
948	abk2922R	Rat GP II b/IIIa(Platelet Membrane Glycoprotein II b/IIIa complex) ELISA Kit
949	abk3591R	Rat PECAM1/CD31(Platelet/Endothelial Cell Adhesion Molecule 1) ELISA Kit
950	abk9822R	Rat PHLDA2(Pleckstrin Homology Like Domain Family A, Member 2) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
951	abk8241R	Rat F11(Coagulation Factor XI) ELISA Kit
952	abk5823R	Rat Pro-ANP(Pro Atrial Natriuretic Peptide) ELISA Kit
953	abk6823R	Rat PGD2S(Prostaglandin D2 Synthase) ELISA Kit
954	abk9803R	Rat ACP3/PAP(Prostatic Acid Phosphatase) ELISA Kit
955	abk7303R	Rat SP-D(Pulmonary surfactant-associated protein D) ELISA Kit
956	abk4103R	Rat PDHA1(Pyruvate dehydrogenase E1 component subunit alpha) ELISA Kit
957	abk4471R	Rat PK(Pyruvate Kinase, Liver And RBC) ELISA Kit
958	abk4872R	Rat PKM2(Pyruvate Kinase, Muscle) ELISA Kit
959	abk2992R	Rat GAP(Ras Gtpase Activating Protein) ELISA Kit
960	abk1253R	Rat REG(Regenerating Islet Derived Protein) ELISA Kit
961	abk5603R	Rat RBP4(Retinol Binding Protein 4) ELISA Kit
962	abk1403R	Rat RNASEH1(Ribonuclease H1) ELISA Kit
963	abk1081R	Rat RPL26(Ribosomal Protein L26) ELISA Kit
964	abk4313R	Rat TACE/ADAM17(TNF α Converting Enzyme) ELISA Kit
965	abk5371R	Rat S100A6(S100 Calcium Binding Protein A6) ELISA Kit
966	abk2933R	Rat SPLA2(Secreted Phospholipase A2) ELISA Kit
967	abk2271R	Rat SLPI(Secretory Leukocyte Protease Inhibitor) ELISA Kit
968	abk9903R	Rat STK39(Serine/Threonine Protein Kinase 39) ELISA Kit
969	abk0282R	Rat SHBG(Sex Hormone Binding Globulin) ELISA Kit
970	abk4933R	Rat IL-2Rα/CD25(Interleukin-2 Receptor alpha chain) ELISA Kit
971	abk5992R	Rat S100A1(S100 Calcium Binding Protein A1) ELISA Kit
972	abk1003R	Rat TLR-2(Toll-like Receptor 2) ELISA Kit
973	abk2253R	Rat TLR6(Toll-Like Receptor 6) ELISA Kit
974	abk6203R	Rat SREBP(Sterol Regulatory Element Binding Protein) ELISA Kit
975	abk1581R	Rat SDF-2(Stromal Cell Derived Factor 2) ELISA Kit
976	abk3581R	Rat SDF-4(Stromal Cell Derived Factor 4) ELISA Kit
977	abk8503R	Rat TRACP-5b(Tartrate Resistant Acid Phosphatase 5b) ELISA Kit
978	abk3253R	Rat TAF1A(Tata Box Binding Protein Associated Factor 1A) ELISA Kit
979	abk6351R	Rat MAPτ(Microtubule Associated Protein Tau/Tau Protein) ELISA Kit
980	abk7351R	Rat TIEG1(TGF Beta Inducible Early Response Gene 1) ELISA Kit
981	abk4402R	Rat TAT(Thrombin/Antithrombin Complex) ELISA Kit
982	abk3291R	Rat TSP-1(Thrombospondin-1) ELISA Kit
983	abk4253R	Rat TXB3(Thromboxanes B3) ELISA Kit
984	abk1282R	Rat TECK/CCL25(Thymus Expressed Chemokine) ELISA Kit
985	abk4881R	Rat TJP1(Tight junction protein ZO-1) ELISA Kit
986	abk5741R	Rat MYC(Myc proto-oncogene protein) ELISA Kit
987	abk1111R	Rat 15-LO(Arachidonate 15-Lipoxygenase) ELISA Kit
988	abk5253R	Rat PTX-3/TSG-14(Pentraxin 3) ELISA Kit
989	abk8832R	Rat TFR2(Transferrin Receptor 2) ELISA Kit
990	abk6313R	Rat TFR1(Transferrin Receptor 1) ELISA Kit
991	abk6253R	Rat TGF-βRII(Transforming Growth Factor Receptor II) ELISA Kit
992	abk8291R	Rat TREM-1(Triggering Receptor Expressed on Myeloid Cells-1) ELISA Kit
993	abk8532R	Rat TNF-β(Tumor Necrosis Factor Beta) ELISA Kit
994	abk7051R	Rat TP53(Tumor Protein p53) ELISA Kit
995	abk2603R	Rat NFE2L2(Nuclear Factor, Erythroid Derived 2, Like 2) ELISA Kit
996	abk1391R	Rat VEGF-C(Vascular Endothelial Growth Factor C) ELISA Kit
997	abk2391R	Rat VEGF-D(Vascular Endothelial Growth Factor D) ELISA Kit
998	abk5762R	Rat VCAM-1/CD106(Vascular Cell Adhesion Molecule 1) ELISA Kit
999	abk9103R	Rat VLA4(Very Late Appearing Antigen 4) ELISA Kit
1000	abk7253R	Rat ADAMTS13/VWFCP(Von Willebrand Factor Cleaving Protease) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
1001	abk5913R	Rat α -Glu(Alpha Glucosidase) ELISA Kit
1002	abk5112R	Rat GDF6(Growth Differentiation Factor 6) ELISA Kit
1003	abk4552R	Rat pMAPT/pTAU(phosphorylated microtubule-associated protein tau) ELISA Kit
1004	abk0642R	Rat Ta1(Thymosin- α 1) ELISA Kit
1005	abk2842R	Rat CBG(Corticosteroid Binding Globulin) ELISA Kit
1006	abk4772R	Rat C1q(Complement 1q) ELISA Kit
1007	abk7392R	Rat CMLC-1(Cardiac Myosin Light Chain 1) ELISA Kit
1008	abk9821R	Rat ASM(Acid Sphingomyelinase) ELISA Kit
1009	abk5713R	Rat OGDH(2-oxoglutarate dehydrogenase, mitochondrial) ELISA Kit
1010	abk8253R	Rat PFK-A(6-phosphofructokinase Type A) ELISA Kit
1011	abk2551R	Rat SNC α (Synuclein Alpha) ELISA Kit
1012	abk8392R	Rat TNNT3/cTn-I(Troponin I Type 3, Cardiac) ELISA Kit
1013	abk8403R	Rat CKMB(Creatine Kinase MB Isoenzyme) ELISA Kit
1014	abk8751R	Rat IL-2R β (Interleukin 2 Receptor Beta) ELISA Kit
1015	abk3022R	Rat PGC(Pepsinogen C) ELISA Kit
1016	abk0182R	Rat NT-proANP(N-Terminal Pro Atrial Natriuretic Peptide) ELISA Kit
1017	abk5782R	Rat P I NP(Procollagen I N-Terminal Propeptide) ELISA Kit
1018	abk5403R	Rat β MSH(beta-Melanocyte Stimulating Hormone) ELISA Kit
1019	abk9301R	Rat FASL/TNFSF6(Factor Related Apoptosis Ligand) ELISA Kit
1020	abk2241R	Rat F7(Coagulation Factor VII) ELISA Kit
1021	abk2382R	Rat F8(Coagulation Factor VIII) ELISA Kit
1022	abk7191R	Rat PGRN(progranulin) ELISA Kit
1023	abk4491R	Rat CLDN5(Claudin 5) ELISA Kit
1024	abk5042R	Rat NES1(Nesfatin 1) ELISA Kit
1025	abk0271R	Rat SEMA3A(Semaphorin 3A) ELISA Kit
1026	abk3613R	Rat PD-L2(Programmed Cell Death Protein 1 Ligand 2) ELISA Kit
1027	abk9682R	Rat TnTf/TNNT3(Troponin T Type 3, Fast Skeletal) ELISA Kit
1028	abk8411R	Rat AFP-L3(Alpha-Fetoprotein Lens Culinaris Agglutinin 3) ELISA Kit
1029	abk4203R	Rat MAOA(Monoamine Oxidase Type A) ELISA Kit
1030	abk9702R	Rat LXR α (Liver X Receptor Alpha) ELISA Kit
1031	abk1922R	Rat HAase(Hyaluronidase) ELISA Kit
1032	abk6072R	Rat Smad2/MADH2(Mothers Against Decapentaplegic Homolog 2) ELISA Kit
1033	abk1903R	Rat HJV(Hemojuvelin) ELISA Kit
1034	abk1601R	Rat VEGF-A(Vascular Endothelial Cell Growth Factor A) ELISA Kit
1035	abk6382R	Rat Irisin ELISA Kit
1036	abk9513R	Rat hs-CRP(high-sensitivity C-Reactive Protein) ELISA Kit
1037	abk9253R	Rat Notch 2(Notch Homolog 2) ELISA Kit
1038	abk2913R	Rat NRP1(Neuropilin 1) ELISA Kit
1039	abk8851R	Rat GAL1(Galectin 1) ELISA Kit
1040	abk5833R	Rat CXCL5(Chemokine C-X-C-Motif Ligand 5) ELISA Kit
1041	abk3023R	Rat SELP(P-Selectin) ELISA Kit
1042	abk4043R	Rat β TG/PBP/CXCL7/NAP2(Thromboglobulin, Beta) ELISA Kit
1043	abk2072R	Rat HPX(Hemopexin) ELISA Kit
1044	abk3391R	Rat VEGFR-2/KDR(Vascular Endothelial Growth Factor Receptor 2) ELISA Kit
1045	abk4312R	Rat NRP2(Neuropilin 2) ELISA Kit
1046	abk3271R	Rat SAP(Serum Amyloid P Component) ELISA Kit
1047	abk0353R	Rat GFRA1(Glial Cell Line Derived Neurotrophic Factor Receptor Alpha 1) ELISA Kit
1048	abk1031R	Rat ANP(Atrial Natriuretic Peptide) ELISA Kit
1049	abk4812R	Rat MT(Melatonin) ELISA Kit
1050	abk0701R	Rat ACTH(Adrenocorticotrophic Hormone) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
1051	abk9131R	Rat β -EP(Beta-Endorphin) ELISA Kit
1052	abk8231R	Rat BK(Bradykinin) ELISA Kit
1053	abk3913R	Rat BNP(Brain Natriuretic Peptide) ELISA Kit
1054	abk7113R	Rat TMS β 4(Thymosin Beta 4) ELISA Kit
1055	abk2203R	Rat CRH(Corticotropin Releasing Hormone) ELISA Kit
1056	abk4331R	Rat CNP(C-Type Natriuretic Peptide) ELISA Kit
1057	abk5812R	Rat DHEA-S(Dehydroepiandrosterone Sulfate) ELISA Kit
1058	abk6812R	Rat DPD(Deoxypyridinoline) ELISA Kit
1059	abk4121R	Rat GC(Glucagon) ELISA Kit
1060	abk9022R	Rat GT(Gastrin) ELISA Kit
1061	abk5311R	Rat ADH(Antidiuretic Hormone) ELISA Kit
1062	abk0723R	Rat LHRH(Luteinizing Hormone-Releasing Hormone) ELISA Kit
1063	abk6022R	Rat MTL(Motilin) ELISA Kit
1064	abk0312R	Rat OXB(Orexin B) ELISA Kit
1065	abk2012R	Rat NPY(Neuropeptide Y) ELISA Kit
1066	abk6213R	Rat NADPH(Nicotinamide Adenine Dinucleotide Phosphate) ELISA Kit
1067	abk0751R	Rat OXA(Orexin A) ELISA Kit
1068	abk9602R	Rat PGF2 α (Prostaglandin F2 Alpha) ELISA Kit
1069	abk1212R	Rat GHRL(Ghrelin) ELISA Kit
1070	abk5461R	Rat SCT(Secretin) ELISA Kit
1071	abk4212R	Rat SST(Somatostatin) ELISA Kit
1072	abk3822R	Rat TXB2(Thromboxane B2) ELISA Kit
1073	abk2043R	Rat U-T4(Ultrasensitivity Thyroxine) ELISA Kit
1074	abk0743R	Rat OXM(Oxyntomodulin) ELISA Kit
1075	abk7701R	Rat α -CTx(Alpha Crosslaps) ELISA Kit
1076	abk0451R	Rat β -CTx(Beta Crosslaps) ELISA Kit
1077	abk0111R	Rat 11-DH-TXB2(11-Dehydrothromboxane B2) ELISA Kit
1078	abk6261R	Rat SOD1(Superoxide Dismutase 1, Soluble) ELISA Kit
1079	abk4422R	Rat SOD2(Superoxide Dismutase 2, Mitochondrial) ELISA Kit
1080	abk2142R	Rat SOD3(Superoxide Dismutase 3, Extracellular) ELISA Kit
1081	abk1843R	Rat SOD4(Superoxide Dismutase 4, Copper Chaperone) ELISA Kit
1082	abk1642R	Rat PACAP-38(Pituitary Adenylate Cyclase Activating Polypeptide 38) ELISA Kit
1083	abk2642R	Rat PACAP-27(Pituitary Adenylate Cyclase Activating Polypeptide 27) ELISA Kit
1084	abk3521R	Rat APLN(Apelin) ELISA Kit
1085	abk5611R	Rat IAPP(Islet Amyloid Polypeptide) ELISA Kit
1086	abk4053R	Rat δ SIP(Delta Sleep Inducing Peptide) ELISA Kit
1087	abk9912R	Rat UST2(Urotensin 2) ELISA Kit
1088	abk2192R	Rat AP13(Apelin 13) ELISA Kit
1089	abk2013R	Rat OB(Obestatin) ELISA Kit
1090	abk3192R	Rat AP36(Apelin 36) ELISA Kit
1091	abk6571R	Rat OFQ/N(Orphanin FQ/Nociceptin) ELISA Kit
1092	abk1461R	Rat TRH(Thyrotropin Releasing Hormone) ELISA Kit
1093	abk9213R	Rat TAP(Trypsinogen Activation Peptide) ELISA Kit
1094	abk4622R	Rat VIP(Vasoactive Intestinal Peptide) ELISA Kit
1095	abk2982R	Rat Ang- I (Angiotensin I) ELISA Kit
1096	abk1982R	Rat Ang- II (Angiotensin II) ELISA Kit
1097	abk2713R	Rat Ang1-7(Angiotensin 1-7) ELISA Kit
1098	abk6603R	Rat TGF- β 2(Transforming Growth Factor Beta 2) ELISA Kit
1099	abk6813R	Rat IL-6(Interleukin 6) ELISA Kit
1100	abk9501R	Rat IFN- γ (Interferon Gamma) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
1	abk57101	Testosterone ELISA Kit
2	abk90101	Cortisol ELISA Kit
3	abk81101	Pg(Progesterone) ELISA Kit
4	abk31001	TGF-β1(Transforming Growth Factor Beta 1) ELISA Kit
5	abk66001	TMAO(Trimethylamine-N-oxide) ELISA Kit
6	abk75611	KYNA(Kynurenic Acid)ELISA Kit
7	abk09001	HCy(Homocysteine) ELISA Kit
8	abk20001	PGE1(Prostaglandin E1) ELISA Kit
9	abk11101	VD(Vitamin D) ELISA Kit
10	abk14101	VD2(Vitamin D2) ELISA Kit
11	abk79001	VD3(Vitamin D3) ELISA Kit
12	abk72101	25-HVD3(25-Hydroxy Vitamin D3) ELISA Kit
13	abk99001	DHVD3(1,25-Dihydroxyvitamin D3) ELISA Kit
14	abk68001	VK1(Vitamin K1) ELISA Kit
15	abk01001	MEK(Methionine-Enkephalin) ELISA Kit
16	abk59001	PGI2(Prostacyclin) ELISA Kit
17	abk55101	GSH(Glutathione) ELISA Kit
18	abk17101	OT(Oxytocin) ELISA Kit
19	abk21101	DHT(Dihydrotestosterone) ELISA Kit
20	abk24101	HA(Histamine) ELISA Kit
21	abk66101	5-HT(5-Hydroxytryptamine) ELISA Kit
22	abk13101	PGE2(Prostaglandin E2) ELISA Kit
23	abk40001	F-E3(Free Estriol) ELISA Kit
24	abk34101	3-NT(3-Nitrotyrosine) ELISA Kit
25	abk70101	8-epi-PGF2a(8-Epi Prostaglandin F2 Alpha) ELISA Kit
26	abk65101	ADMA(Asymmetrical Dimethylarginine) ELISA Kit
27	abk43101	EPI(Epinephrine) ELISA Kit
28	abk53101	DA(Dopamine) ELISA Kit
29	abk23101	NE(Noradrenaline) ELISA Kit
30	abk84001	VB2(Vitamin B2) ELISA Kit
31	abk85101	AA(Arachidonic Acid) ELISA Kit
32	abk15101	LXA4(Lipoxin A4) ELISA Kit
33	abk45101	TXA2(Thromboxane A2) ELISA Kit
34	abk97001	IP3(Inositol Triphosphate) ELISA Kit
35	abk56101	MDA(Malondialdehyde) ELISA Kit
36	abk03101	LTB4(Leukotriene B4) ELISA Kit
37	abk26101	E1(Estrone) ELISA Kit
38	abk00101	PGD2(Prostaglandin D2) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
39	abk90001	SP(Substance P) ELISA Kit
40	abk73101	ASD(Androstenedione) ELISA Kit
41	abk33101	ALD(Aldosterone) ELISA Kit
42	abk61101	ADT(Androsterone) ELISA Kit
43	abk69001	Bb(Bilirubin) ELISA Kit
44	abk91101	HVA(Homovanillic Acid) ELISA Kit
45	abk92101	ACH(Acetylcholine) ELISA Kit
46	abk28001	cGMP(Cyclic Guanosine Monophosphate) ELISA Kit
47	abk41101	E1S(Estrone Sulfate) ELISA Kit
48	abk30101	LTA4(Leukotriene A4) ELISA Kit
49	abk51101	PTD(Pentosidine) ELISA Kit
50	abk19001	LTC4(Leukotriene C4) ELISA Kit
51	abk72001	ACTb(Actin Beta) ELISA Kit
52	abk23001	LTD4(Leukotriene D4) ELISA Kit
53	abk11001	LTE4(Leukotriene E4) ELISA Kit
54	abk83101	E3(Estriol) ELISA Kit
55	abk50001	F-TESTO(Free Testosterone) ELISA Kit
56	abk29001	cAMP(Cyclic Adenosine Monophosphate) ELISA Kit
57	abk82101	VB6(Vitamin B6) ELISA Kit
58	abk46101	FA(Folic Acid) ELISA Kit
59	abk93101	17-OHP(17-Hydroxyprogesterone) ELISA Kit
60	abk42001	CTN(Creatinine) ELISA Kit
61	abk10101	Car(Carnosine) ELISA Kit
62	abk80001	GnRH(Gonadotropin Releasing Hormone) ELISA Kit
63	abk16101	5-HIAA(5-Hydroxyindoleacetic Acid) ELISA Kit
64	abk96101	fT3(Free Triiodothyronine) ELISA Kit
65	abk28101	P5(Pregnenolone) ELISA Kit
66	abk42101	rT3(Reverse Triiodothyronine) ELISA Kit
67	abk86101	fT4(Free Thyroxine) ELISA Kit
68	abk41001	AcCoA(Acetyl CoA) ELISA Kit
69	abk70001	4-HNE(4-Hydroxynonenal) ELISA Kit
70	abk65001	VMA(Vanillylmandelic Acid) ELISA Kit
71	abk60001	VA(Vitamin A) ELISA Kit
72	abk07001	SPH(Sphingomyelin) ELISA Kit
73	abk58101	8-OHdG(8-Hydroxydeoxyguanosine) ELISA Kit
74	abk77101	E2(Estradiol) ELISA Kit
75	abk31101	Cort(Corticosterone) ELISA Kit
76	abk27101	gABA(Gamma-Aminobutyric Acid) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
1	abk6813HF	MF-Human IL-6(Interleukin 6) ELISA Kit
2	abk9382HF	MF-Human IFN-α(Interferon Alpha) ELISA Kit
3	abk6813MF	MF-Mouse IL-6(Interleukin 6) ELISA Kit
4	abk0601MF	MF-Mouse TNF-α(Tumor Necrosis Factor Alpha) ELISA Kit
5	abk6201MF	MF-Mouse GROα/CXCL1(Growth Regulated Oncogene Alpha) ELISA Kit
6	abk9501MF	MF-Mouse IFN-γ(Interferon Gamma) ELISA Kit
7	abk2682HF	MF-Human GM-CSF(Granulocyte-Macrophage Colony Stimulating Factor) ELISA Kit
8	abk10013F	MF-TGF-β1(Transforming Growth Factor Beta 1) ELISA Kit
9	abk0901MF	MF-Mouse IL-9(Interleukin 9) ELISA Kit
10	abk7682MF	MF-Mouse IL-17A(Interleukin 17A) ELISA Kit
11	abk5501HF	HF-Human IL-4(Interleukin 4) ELISA Kit
12	abk6682RF	MF-Rat IL-2(Interleukin 2) ELISA Kit
13	abk0601HF	MF-Human TNF-α(Tumor Necrosis Factor Alpha) ELISA Kit
14	abk9501HF	MF-Human IFN-γ(Interferon Gamma) ELISA Kit
15	abk8532HF	MF-Human TNF-β(Tumor Necrosis Factor Beta) ELISA Kit
16	abk5671HF	MF-Human OPG(Osteoprotegerin) ELISA Kit
17	abk9501RF	MF-Rat IFN-γ(Interferon Gamma) ELISA Kit
18	abk4401HF	MF-Human GDF15(Growth Differentiation Factor 15) ELISA Kit
19	abk9101HF	MF-Human sCD40L(Soluble Cluster of Differentiation 40 Ligand) ELISA Kit
20	abk0901HF	MF-Human IL-9(Interleukin 9) ELISA Kit
21	abk5892HF	MF-Human SELE(E-Selectin) ELISA Kit
22	abk5892RF	MF-Rat SELE(E-Selectin) ELISA Kit
23	abk5892MF	MF-Mouse SELE(E-Selectin) ELISA Kit
24	abk9372HF	MF-Human LEP(Leptin) ELISA Kit
25	abk9372RF	MF-Rat LEP(Leptin) ELISA Kit
26	abk9372MF	MF-Mouse LEP(Leptin) ELISA Kit
27	abk4601HF	MF-Human REN(Renin) ELISA Kit
28	abk4601RF	MF-Rat REN(Renin) ELISA Kit
29	abk4601MF	MF-Mouse REN(Renin) ELISA Kit
30	abk3861HF	MF-Human RETN(Resistin) ELISA Kit
31	abk3861RF	MF-Rat RETN(Resistin) ELISA Kit
32	abk3861MF	MF-Mouse RETN(Resistin) ELISA Kit
33	abk6831RF	MF-Rat IL-6(Interleukin 6) ELISA Kit
34	abk9382RF	MF-Rat IFN-α(Interferon Alpha) ELISA Kit
35	abk9382MF	MF-Mouse IFN-α(Interferon Alpha) ELISA Kit
36	abk0601RF	MF-Rat TNF-α(Tumor Necrosis Factor Alpha) ELISA Kit
37	abk6201HF	MF-Human GROα/CXCL1(Growth Regulated Oncogene Alpha) ELISA Kit
38	abk6201RF	MF-Rat GROα/CXCL1(Growth Regulated Oncogene Alpha) ELISA Kit
39	abk5501RF	MF-Rat IL-4(Interleukin 4) ELISA Kit
40	abk5501MF	MF-Mouse IL-4(Interleukin 4) ELISA Kit
41	abk2682RF	MF-Rat GM-CSF(Granulocyte-Macrophage Colony Stimulating Factor) ELISA Kit
42	abk7682RF	MF-Rat IL-17A(Interleukin 17A) ELISA Kit
43	abk8532RF	MF-Rat TNF-β(Tumor Necrosis Factor Beta) ELISA Kit
44	abk8532MF	MF-Mouse TNF-β(Tumor Necrosis Factor Beta) ELISA Kit
45	abk6682HF	MF-Human IL-2(Interleukin 2) ELISA Kit
46	abk6682MF	MF-Mouse IL-2(Interleukin 2) ELISA Kit
47	abk4401RF	MF-Rat GDF15(Growth Differentiation Factor 15) ELISA Kit
48	abk4401MF	MF-Mouse GDF15(Growth Differentiation Factor 15) ELISA Kit
49	abk5671RF	MF-Rat OPG(Osteoprotegerin) ELISA Kit
50	abk5671MF	MF-Mouse OPG(Osteoprotegerin) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
51	abk4701HF	MF-Human IL-1β(Interleukin 1 Beta) ELISA Kit
52	abk4701RF	MF-Rat IL-1β(Interleukin 1 Beta) ELISA Kit
53	abk7401HF	MF-Human IL-1α(Interleukin 1 Alpha) ELISA Kit
54	abk7401RF	MF-Rat IL-1α(Interleukin 1 Alpha) ELISA Kit
55	abk7401MF	MF-Mouse IL-1α(Interleukin 1 Alpha) ELISA Kit
56	abk7813HF	MF-Human IL-10(Interleukin 10) ELISA Kit
57	abk7813RF	MF-Rat IL-10(Interleukin 10) ELISA Kit
58	abk7813MF	MF-Mouse IL-10(Interleukin 10) ELISA Kit
59	abk2482HF	MF-Human DPP4(Dipeptidyl Peptidase IV) ELISA Kit
60	abk2482RF	MF-Rat DPP4(Dipeptidyl Peptidase IV) ELISA Kit
61	abk2482MF	MF-Mouse DPP4(Dipeptidyl Peptidase IV) ELISA Kit
62	abk5701HF	MF-Human IL-12(Interleukin 12) ELISA Kit
63	abk5701RF	MF-Rat IL-12(Interleukin 12) ELISA Kit
64	abk5701MF	MF-Mouse IL-12(Interleukin 12) ELISA Kit
65	abk6701HF	MF-Human IL-12 p40(Interleukin 12 p40) ELISA Kit
66	abk6701MF	MF-Mouse IL-12 p40(Interleukin 12 p40) ELISA Kit
67	abk7501MF	MF-Mouse IL-22(Interleukin 22) ELISA Kit
68	abk1191HF	MF-Human IGFBP-6(Insulin-like Growth Factor Binding Protein 6) ELISA Kit
69	abk1191RF	MF-Rat IGFBP-6(Insulin-like Growth Factor Binding Protein 6) ELISA Kit
70	abk1191MF	MF-Mouse IGFBP-6(Insulin-like Growth Factor Binding Protein 6) ELISA Kit
71	abk7913HF	MF-Human ENG(Endoglin) ELISA Kit
72	abk7913MF	MF-Mouse ENG(Endoglin) ELISA Kit
73	abk1682HF	MF-Human G-CSF(Granulocyte Colony-stimulating Factor) ELISA Kit
74	abk1682MF	MF-Mouse G-CSF(Granulocyte Colony-stimulating Factor) ELISA Kit
75	abk1382HF	MF-Human ICAM-1/CD54(intercellular adhesion molecule 1) ELISA Kit
76	abk1382RF	MF-Rat ICAM-1/CD54(intercellular adhesion molecule 1) ELISA Kit
77	abk1382MF	MF-Mouse ICAM-1/CD54(intercellular adhesion molecule 1) ELISA Kit
78	abk3481HF	MF-Human MMP-8(Matrix Metalloproteinase 8) ELISA Kit
79	abk3481RF	MF-Rat MMP-8(Matrix Metalloproteinase 8) ELISA Kit
80	abk3481MF	MF-Mouse MMP-8(Matrix Metalloproteinase 8) ELISA Kit
81	abk5101HF	MF-Human MPIF2(Myeloid Progenitor Inhibitory Factor 2) ELISA Kit
82	abk5101MF	MF-Mouse MPIF2(Myeloid Progenitor Inhibitory Factor 2) ELISA Kit
83	abk5401HF	MF-Human HGF(Hepatocyte Growth Factor) ELISA Kit
84	abk5401RF	MF-Rat HGF(Hepatocyte Growth Factor) ELISA Kit
85	abk5401MF	MF-Mouse HGF(Hepatocyte Growth Factor) ELISA Kit
86	abk3682HF	MF-Human IGF-1(Insulin-like Growth Factor 1) ELISA Kit
87	abk3682RF	MF-Rat IGF-1(Insulin-like Growth Factor 1) ELISA Kit
88	abk3682MF	MF-Mouse IGF-1(Insulin-like Growth Factor 1) ELISA Kit
89	abk1372HF	MF-Human MCP-1(Monocyte Chemotactic Protein 1) ELISA Kit
90	abk1371RF	MF-Rat MCP-1(Monocyte Chemotactic Protein 1) ELISA Kit
91	abk5733MF	MF-Mouse CCL12/MCP-5(Monocyte Chemotactic Protein 5) ELISA Kit
92	abk0572HF	MF-Human KIM-1(Kidney Injury Molecule 1) ELISA Kit
93	abk0572RF	MF-Rat KIM-1(Kidney Injury Molecule 1) ELISA Kit
94	abk0572MF	MF-Mouse KIM-1(Kidney Injury Molecule 1) ELISA Kit
95	abk1601RF	MF-Rat VEGF-A(Vascular Endothelial Cell Growth Factor A) ELISA Kit
96	abk1601MF	MF-Mouse VEGF-A(Vascular Endothelial Cell Growth Factor A) ELISA Kit
97	abk3211RF	MF-Rat IL-18(Interleukin 18) ELISA Kit
98	abk3211MF	MF-Mouse IL-18(Interleukin 18) ELISA Kit
99	abk6401HF	MF-Human IFN-β(Interferon Beta) ELISA Kit
100	abk6401RF	MF-Rat IFN-β(Interferon Beta) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
101	abk6401MF	MF-Mouse IFN- β (Interferon Beta) ELISA Kit
102	abk0422HF	MF-Human GHR(Growth Hormone Receptor) ELISA Kit
103	abk0422MF	MF-Mouse GHR(Growth Hormone Receptor) ELISA Kit
104	abk8913HF	MF-Human LEPR(Leptin Receptor) ELISA Kit
105	abk8913RF	MF-Rat LEPR(Leptin Receptor) ELISA Kit
106	abk8913MF	MF-Mouse LEPR(Leptin Receptor) ELISA Kit
107	abk1752HF	MF-Human CTSS(Cathepsin S) ELISA Kit
108	abk8542HF	MF-Human TNFRSF1B(Tumor Necrosis Factor Receptor Superfamily, Member 1B) ELISA Kit
109	abk8542RF	MF-Rat TNFRSF1B(Tumor Necrosis Factor Receptor Superfamily, Member 1B) ELISA Kit
110	abk8542MF	MF-Mouse TNFRSF1B(Tumor Necrosis Factor Receptor Superfamily, Member 1B) ELISA Kit
111	abk0452HF	MF-Human IL-11(Interleukin 11) ELISA Kit
112	abk0452RF	MF-Rat IL-11(Interleukin 11) ELISA Kit
113	abk0452MF	MF-Mouse IL-11(Interleukin 11) ELISA Kit
114	abk2501HF	MF-Human MCSF(Macrophage Colony Stimulating Factor 1) ELISA Kit
115	abk2501RF	MF-Rat MCSF(Macrophage Colony Stimulating Factor 1) ELISA Kit
116	abk2501MF	MF-Mouse MCSF(Macrophage Colony Stimulating Factor 1) ELISA Kit
117	abk4501HF	MF-Human IL-3(Interleukin 3) ELISA Kit
118	abk4501RF	MF-Rat IL-3(Interleukin 3) ELISA Kit
119	abk4501MF	MF-Mouse IL-3(Interleukin 3) ELISA Kit
120	abk9901HF	MF-Human IL-5(Interleukin 5) ELISA Kit
121	abk9901RF	MF-Rat IL-5(Interleukin 5) ELISA Kit
122	abk9901MF	MF-Mouse IL-5(Interleukin 5) ELISA Kit
123	abk3352HF	MF-Human IL-37(Interleukin 37) ELISA Kit
124	abk2511HF	MF-Human AGER(Total Advanced Glycosylation End Product Specific Receptor) ELISA Kit
125	abk2511RF	MF-Rat AGER(Total Advanced Glycosylation End Product Specific Receptor) ELISA Kit
126	abk2511MF	MF-Mouse AGER(Total Advanced Glycosylation End Product Specific Receptor) ELISA Kit
127	abk8482HF	MF-Human BDNF(Brain Derived Neurotrophic Factor) ELISA Kit
128	abk8482RF	MF-Rat BDNF(Brain Derived Neurotrophic Factor) ELISA Kit
129	abk8482MF	MF-Mouse BDNF(Brain Derived Neurotrophic Factor) ELISA Kit
130	abk4292HF	MF-Human A β 1-40(Amyloid Beta 1-40) ELISA Kit
131	abk4292RF	MF-Rat A β 1-40(Amyloid Beta 1-40) ELISA Kit
132	abk4292MF	MF-Mouse A β 1-40(Amyloid Beta 1-40) ELISA Kit
133	abk5292HF	MF-Human A β 1-42(Amyloid Beta 1-42) ELISA Kit
134	abk5292RF	MF-Rat A β 1-42(Amyloid Beta 1-42) ELISA Kit
135	abk5292MF	MF-Mouse A β 1-42(Amyloid Beta 1-42) ELISA Kit
136	abk6301HF	MF-Human ET-1(Endothelin 1) ELISA Kit
137	abk6301RF	MF-Rat ET-1(Endothelin 1) ELISA Kit
138	abk6301MF	MF-Mouse ET-1(Endothelin 1) ELISA Kit
139	abk8901HF	MF-Human C5a(Complement Component 5a) ELISA Kit
140	abk8901RF	MF-Rat C5a(Complement Component 5a) ELISA Kit
141	abk8901MF	MF-Mouse C5a(Complement Component 5a) ELISA Kit
142	abk8572HF	MF-Human IGF-2(Insulin Like Growth Factor 2) ELISA Kit
143	abk8572RF	MF-Rat IGF-2(Insulin Like Growth Factor 2) ELISA Kit
144	abk8572MF	MF-Mouse IGF-2(Insulin Like Growth Factor 2) ELISA Kit
145	abk2401HF	MF-Human FGF21(Fibroblast Growth Factor 21) ELISA Kit
146	abk2401RF	MF-Rat FGF21(Fibroblast Growth Factor 21) ELISA Kit
147	abk2401MF	MF-Mouse FGF21(Fibroblast Growth Factor 21) ELISA Kit
148	abk8501RF	MF-Rat IL-23(Interleukin 23) ELISA Kit
149	abk8501MF	MF-Mouse IL-23(Interleukin 23) ELISA Kit
150	abk8242MF	MF-Mouse IL-33(Interleukin 33) ELISA Kit

SR. NO.	AbChek Cat No.	Product Name
151	abk3301RF	MF-Rat EGF(Epidermal Growth Factor) ELISA Kit
152	abk3301MF	MF-Mouse EGF(Epidermal Growth Factor) ELISA Kit
153	abk1203HF	MF-Human FSH(Follicle Stimulating Hormone) ELISA Kit
154	abk1203RF	MF-Rat FSH(Follicle Stimulating Hormone) ELISA Kit
155	abk1203MF	MF-Mouse FSH(Follicle Stimulating Hormone) ELISA Kit
156	abk1472HF	MF-Human LH(Luteinizing Hormone) ELISA Kit
157	abk1472RF	MF-Rat LH(Luteinizing Hormone) ELISA Kit
158	abk1472MF	MF-Mouse LH(Luteinizing Hormone) ELISA Kit
159	abk1701HF	MF-Human PRL(Prolactin) ELISA Kit
160	abk1701RF	MF-Rat PRL(Prolactin) ELISA Kit
161	abk1701MF	MF-Mouse PRL(Prolactin) ELISA Kit
162	abk8392HF	MF-Human TNNI3/cTn-I(Troponin I Type 3, Cardiac) ELISA Kit
163	abk8392RF	MF-Rat TNNI3/cTn-I(Troponin I Type 3, Cardiac) ELISA Kit
164	abk8392MF	MF-Mouse TNNI3/cTn-I(Troponin I Type 3, Cardiac) ELISA Kit
165	abk8531HF	MF-Human cTnT/TNNT2(Troponin T Type 2, Cardiac) ELISA Kit
166	abk8531RF	MF-Rat cTnT/TNNT2(Troponin T Type 2, Cardiac) ELISA Kit
167	abk8531MF	MF-Mouse cTnT/TNNT2(Troponin T Type 2, Cardiac) ELISA Kit
168	abk2101HF	MF-Human ECF/CCL11(Eosinophil Chemotactic Factor) ELISA Kit
169	abk2101RF	MF-Rat ECF/CCL11(Eosinophil Chemotactic Factor) ELISA Kit
170	abk2101MF	MF-Mouse ECF/CCL11(Eosinophil Chemotactic Factor) ELISA Kit
171	abk3392HF	MF-Human CGRP1(Calcitonin Gene Related Peptide 1) ELISA Kit
172	abk3392RF	MF-Rat CGRP1(Calcitonin Gene Related Peptide 1) ELISA Kit
173	abk3392MF	MF-Mouse CGRP1(Calcitonin Gene Related Peptide 1) ELISA Kit
174	abk9141HF	MF-Human NOS2/iNOS(Nitric Oxide Synthase 2, Inducible) ELISA Kit
175	abk9141RF	MF-Rat NOS2/iNOS(Nitric Oxide Synthase 2, Inducible) ELISA Kit
176	abk9141MF	MF-Mouse NOS2/iNOS(Nitric Oxide Synthase 2, Inducible) ELISA Kit
177	abk7371HF	MF-Human S100B(S100 Calcium Binding Protein B) ELISA Kit
178	abk7371RF	MF-Rat S100B(S100 Calcium Binding Protein B) ELISA Kit
179	abk7371MF	MF-Mouse S100B(S100 Calcium Binding Protein B) ELISA Kit
180	abk4603HF	MF-Human PDGF-BB(Platelet Derived Growth Factor BB) ELISA Kit
181	abk4603RF	MF-Rat PDGF-BB(Platelet Derived Growth Factor BB) ELISA Kit
182	abk4603MF	MF-Mouse PDGF-BB(Platelet Derived Growth Factor BB) ELISA Kit
183	abk5322HF	MF-Human PAI-1(Plasminogen Activator Inhibitor 1) ELISA Kit
184	abk5322RF	MF-Rat PAI-1(Plasminogen Activator Inhibitor 1) ELISA Kit
185	abk5322MF	MF-Mouse PAI-1(Plasminogen Activator Inhibitor 1) ELISA Kit
186	abk8401HF	MF-Human IL-1ra/IL-1F3(Interleukin 1 Receptor Antagonist) ELISA Kit
187	abk8401RF	MF-Rat IL-1ra/IL-1F3(Interleukin 1 Receptor Antagonist) ELISA Kit
188	abk8401MF	MF-Mouse IL-1ra/IL-1F3(Interleukin 1 Receptor Antagonist) ELISA Kit
189	abk1752MF	MF-Mouse CTSS(Cathepsin S) ELISA Kit
190	abk1752RF	MF-Rat CTSS(Cathepsin S) ELISA Kit
191	abk1601HF	MF-Human VEGF-A(Vascular Endothelial Cell Growth Factor A) ELISA Kit
192	abk3082HF	MF-Human EGFR2(Epidermal Growth Factor Receptor 2) ELISA Kit
193	abk7682HF	MF-Human IL-17A(Interleukin 17A) ELISA Kit
194	abk4701MF	MF-Mouse IL-1 β (Interleukin 1 Beta) ELISA Kit
195	abk7501HF	MF-Human IL-22(Interleukin 22) ELISA Kit
196	abk3211HF	MF-Human IL-18(Interleukin 18) ELISA Kit
197	abk3372HF	MF-Human IL-8(Interleukin 8) ELISA Kit
198	abk8501HF	MF-Human IL-23(Interleukin 23) ELISA Kit
199	abk8242HF	MF-Human IL-33(Interleukin 33) ELISA Kit
200	abk6972MF	MF-Mouse MMP-9(MatrixMetalloproteinase 9) ELISA Kit

What is an ELISA?

ELISA (Enzyme Linked Immunosorbent Assays) are a type of immunoassay that are commonly used to quantify levels of a specific target within a sample. Samples routinely used in ELISA include serum, plasma, cell culture supernates, cell lysates, saliva, tissue lysates, and urine.


ELISA are usually run in 96-well microplates coated with a capture antibody specific for the analyte of interest. Upon incubation with experimental samples, standards, or controls, the target analyte is captured by this antibody. A conjugated detection antibody binds to a different epitope on the target analyte. A substrate solution is subsequently added to produce a signal that is proportional to the amount of analyte bound. ELISA can have different formats. Descriptions and diagrams of these can be found in the next section.

ELISA Types

The four main types of ELISA are indirect, direct, sandwich, and competitive. Each type of ELISA has its own advantages & disadvantages.

Direct ELISA

In a direct ELISA, an antigen or sample is immobilized directly on the plate and a conjugated detection antibody binds to the target protein. Substrate is then added, producing a signal that is proportional to the amount of analyte in the sample. Since only one antibody is used in a direct ELISA, they are less specific than a sandwich ELISA.

**When to Use**

Assessing antibody affinity and specificity. Investigating blocking/ inhibitory interactions

Advantages


- Fast and simple protocol

Disadvantages

- Less specific since you are only using 1 antibody.
- Potential for high background if all proteins from a sample are immobilized in well.

Indirect ELISA

An indirect ELISA is similar to a direct ELISA in that an antigen is immobilized on a plate, but it includes an additional amplification detection step. First, an unconjugated primary detection antibody is added and binds to the specific antigen. A conjugated secondary antibody directed against the host species of the primary antibody is then added. Substrate then produces a signal proportional to the amount of antigen bound in the well.

**When to Use**

Measuring endogenous antibodies

Advantages

- Amplification using a secondary antibody

Disadvantages

- Potential for cross-reactivity caused by secondary antibody

Sandwich ELISA

An indirect ELISA is similar to a direct ELISA in that an antigen is immobilized on a plate, but it includes an additional amplification detection step. First, an unconjugated primary detection antibody is added and binds to the specific antigen. A conjugated secondary antibody directed against the host species of the primary antibody is then added. Substrate then produces a signal proportional to the amount of antigen bound in the well.

**When to Use**

Determining analyte concentration in a biological sample

Advantages


- Highest specificity and sensitivity
- Compatible with complex sample matrices

Disadvantages

- Longer protocol
- Challenging to develop

Competitive ELISA

An indirect ELISA is similar to a direct ELISA in that an antigen is immobilized on a plate, but it includes an additional amplification detection step. First, an unconjugated primary detection antibody is added and binds to the specific antigen. A conjugated secondary antibody directed against the host species of the primary antibody is then added. Substrate then produces a signal proportional to the amount of antigen bound in the well.

**When to Use**

Determining concentrations of small molecules & hormones

Advantages

- Ability to quantitate small molecules

Disadvantages

- Less specific since you are only using 1 antibody
- Requires a conjugated antigen

Key Features

- ★ Low sample volume
- ★ Faster detection times
- ★ Optimized Performance
- ★ High reproducibility



Save Samples



High Sensitivity



Faster Results

Sample loading: 25 μ L

Sensitivity: Higher

Detection Time: Reduced by 30 minutes

Operation steps: 3 steps

Storage: 4 $^{\circ}$ C

Versus

Sample loading: 100 μ L

Sensitivity: High

Detection Time: more than 3 hours

Operation steps: 4 steps

Storage: 4 $^{\circ}$ C / -20 $^{\circ}$ C

AbChek ELISA Kits

Conventional ELISA Kits



Manufactured by-

NuLife, F-4, Sector-8, Noida,
U.P., 201301- INDIA

+91 9899464625

012-04139569
011-46112691



www.abchek.com
www.nulifecare.in



info@nulifecare.in
technical@nulife.in